



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Motability Retail Prescription Form

MAY 2024
Customer Ref:

Customer Ref\*, Contact No\*, Delivery Address\*, Account No\*, Order No\*, Order Qty\*, Order date\*

AVIVA RX20 Modulite M21 707041 - Weekly Price £45.77

Standard Features

- Max user weight 136kg (21 Stone)
Crash tested ISO 7176 -Part 19
Modulite seating system
Flex3 adjustable Backrest Plates with Laguna Backrest cushion
Standard T-Armrests, detachable with Long armpad
Manually settable Seat Tilt 0 - 15°
Manual Settable Recliner (in 3.75° steps) 90° - 120°
80° hangers with fixed footplate
Remote LiNX DLX-REM110 with swing away bracket
8" Castor, puncture proof solid tyre
Solid inlay black wheels
6 km/h 4mph motors
60Ah batteries
Anti tippers
4 Point tie down system
Lapbelt with autolock (adjustable both sides)

Standard Base price: £5,977



Please make your selection with a cross in the box [ ] Free of charge [O] standard feature [◆] POA special request

Customer weight = KG

Table with columns: Code, Description, Price. Rows include BASE MODEL\* (JBC0000), CHASSIS WIDTH\* (JBL2810), SEATING TYPE\* (JBC2100), SEAT FRAME\* (JBC2111).

<b>SEAT WIDTH*</b>			
JBC2121	<input type="checkbox"/>	380 mm Seat Plate	<input type="radio"/>
JBC2122	<input type="checkbox"/>	430 mm Seat Plate	<input type="radio"/>
JBC2123	<input type="checkbox"/>	480 mm Seat Plate	<input type="radio"/>
JBC2124	<input type="checkbox"/>	530 mm Seat Plate	<input type="radio"/>
<b>SEAT DEPTH*</b>			
JBC2132	<input type="checkbox"/>	Seat Depth 410 mm	<input type="radio"/>
JBC2134	<input type="checkbox"/>	Seat Depth 460 mm	<input type="radio"/>
JBC2136	<input type="checkbox"/>	Seat Depth 510 mm	<input type="radio"/>
<b>SEAT CUSHION*</b>			
JBC2142	<input type="checkbox"/>	Matrx Flo-tech Cushion Flo-tech Lite	<input checked="" type="radio"/>
<b>BACKREST SUPPORT*</b>			
JBC2162	<input type="checkbox"/>	Backposts tall (540 mm)	<input checked="" type="radio"/>
<b>PUSH HANDLES FOR BACKPOSTS</b>			
JBC2108	<input type="checkbox"/>	Backposts without push handles	<input checked="" type="radio"/>
<b>BACKREST UPHOLSTERY / CUSHION*</b>			
JBC2173	<input type="checkbox"/>	Flex3 adjustable Backrest Plates with Laguna Backrest cushion	<input checked="" type="radio"/>
<b>UPHOLSTERY TYPE*</b>			
JBC2152	<input type="checkbox"/>	Black Black multi stretch PU coated fabric Upholstery Upholstery	<input type="radio"/>
JBC2155	<input type="checkbox"/>	Black Spacer Fabric Upholstery	<input type="radio"/>
<b>ARMREST (Pair)*</b>			
JBC2191	<input type="checkbox"/>	T-Arm detachable, Standard	<input checked="" type="radio"/>
<b>ARMPAD (Pair)*</b>			
JBC2196	<input type="checkbox"/>	Long Armpad (375 mm)	<input checked="" type="radio"/>
<b>SEAT LIFTER / SEAT TILT*</b>			
JBC1091	<input type="checkbox"/>	Manually settable Seat Tilt 0 - 15°	<input checked="" type="radio"/>
<b>SEAT TO GROUND HEIGHT*</b>			
JBC2857	<input type="checkbox"/>	Seat-to-Floor Height 435 mm	<input type="radio"/>
JBC2858	<input type="checkbox"/>	Seat-to-Floor Height 460 mm	<input type="radio"/>
JBC2859	<input type="checkbox"/>	Seat-to-Floor Height 485 mm	<input type="radio"/>
<b>BACKREST RECLINER*</b>			
JBC0975	<input type="checkbox"/>	Manual Settable Recliner (in 3.75° steps) 90° - 120°	<input checked="" type="radio"/>
<b>FOOTREST HANGER TYPE*</b>			
JBC0526	<input type="checkbox"/>	Standard hanger 80° pair Lower leg length: 350 - 450 mm	<input checked="" type="radio"/>
<b>FOOTPLATE TYPE*</b>			
JBC0512	<input type="checkbox"/>	Fixed Footplate pair	<input checked="" type="radio"/>
<b>CONTROLS*</b>			
JBC2201	<input type="checkbox"/>	Remote LiNX DLX-REM110 = Drive only 1 Drive function with Drive only power module	<input checked="" type="radio"/>
<b>CONTROL BRACKETRY*</b>			
JBC1707	<input type="checkbox"/>	Swing Away Bracket RH	<input type="radio"/>
JBC1708	<input type="checkbox"/>	Swing Away Bracket LH	<input type="radio"/>

<b>SHROUD COLOUR*</b>			
JBC2871	<input type="checkbox"/>	Mid range - Black/Anthracite Shroud	<input type="radio"/>
JBC2872	<input type="checkbox"/>	Mid range - Silver Shroud	<input type="radio"/>
JBC2873	<input type="checkbox"/>	Mid range - Bengal Red Shroud	<input type="radio"/>
JBC2874	<input type="checkbox"/>	Mid range - Mercury Blue Shroud	<input type="radio"/>
JBC2875	<input type="checkbox"/>	Mid range - Midnight Purple Shroud	<input type="radio"/>
<b>RIM COLOUR*</b>			
JBC2882	<input type="checkbox"/>	Mid range - Black Rim Inserts	<input type="radio"/>
JBC2883	<input type="checkbox"/>	Mid range - Silver Rim Inserts	<input type="radio"/>
JBC2884	<input type="checkbox"/>	Mid range - Bengal Red Rim Inserts	<input type="radio"/>
JBC2885	<input type="checkbox"/>	Mid range - Mercury Blue Rim Inserts	<input type="radio"/>
JBC2886	<input type="checkbox"/>	Mid range - Midnight Purple Rim Inserts	<input type="radio"/>
<b>FORKS*</b>			
JBC2632	<input type="checkbox"/>	Double sided fork Without castor rim insert (enlarges wide base to 656 mm)	<input checked="" type="checkbox"/>
<b>FRONT WHEEL *</b>			
JBC0828	<input type="checkbox"/>	8" (200 mm) x 2.50-4 (65 mm) Castor, puncture proof solid tyre	<input checked="" type="checkbox"/>
<b>DRIVE WHEEL *</b>			
JBC0705	<input type="checkbox"/>	DW 14" (355 mm) x 3" (75 mm), Solid inlay Dual density foam insert	<input checked="" type="checkbox"/>
<b>C.T.C SUSPENSION*</b>			
JBC2870	<input type="checkbox"/>	CTC Suspension Control - Traction - Comfort	<input checked="" type="checkbox"/>
<b>SUSPENSION SETTING*</b>			
JBC1910	<input type="checkbox"/>	Soft (suggested for up to 70 kg User Weight)	<input type="radio"/>
JBC1912	<input type="checkbox"/>	Medium (suggested for up to 90 kg User Weight)	<input type="radio"/>
JBC1911	<input type="checkbox"/>	Firm (suggested for up to 110 kg User Weight)	<input type="radio"/>
JBC1915	<input type="checkbox"/>	Very firm (suggested for up to 136 kg User Weight)	<input type="radio"/>
<b>MOTOR TYPE*</b>			
JBC1282	<input type="checkbox"/>	6 km/h	<input checked="" type="checkbox"/>
<b>BATTERY TYPE*</b>			
JBC1259	<input type="checkbox"/>	60 Ah Gel	<input checked="" type="checkbox"/>
<b>CHARGER*</b>			
JBC1343	<input type="checkbox"/>	Opti Charge Charger 8A (IP44)	<input checked="" type="checkbox"/>
<b>ANTI-TIPPERS*</b>			
JBC1205	<input type="checkbox"/>	Anti-tipper	<input checked="" type="checkbox"/>
<b>TRANSPORTATION KIT*</b>			
JBC1494	<input type="checkbox"/>	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19) Adequate headrest is compulsory	<input checked="" type="checkbox"/>
<b>LAP BELT*</b>			
JBC2113	<input type="checkbox"/>	Lap belt with Auto Lock Flexible mount in side rails	<input checked="" type="checkbox"/>

<b>TOTAL PRICE</b>	
--------------------	--

**Non Conforming Options:**



For a more complex customer requirement please refer to our Invacare Unique special order form. email [specialsuk@invacare.com](mailto:specialsuk@invacare.com) or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Description of special request (please include relevant option code or part number):