



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

MAY 2024
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

E-Fix & Action 3NG M82 707362 - Weekly Price £46.32

Standard Features

- Action 3NG + E-Fix
- Max user weight 125kg (19.7 stone)
- Crash tested ISO- 7176-Part 19
- Seat widths 380mm to 505mm (15" -20")
- Seat depths 400mm or 450mm (16" or 18")
- Seat heights 460mm to 510mm (18" -20")
- Nylon with hook and loop strips seat upholstery
- Fixed 7°Backrest 2 heights 400mm/430mm
- Lumbar support padded backrest with shaped cover
- 80°swing in/swing out hangers
- Flip up footplates
- Removable full arms
- Basic 50mm foam cushion
- Adjustable calf strap
- Choice of 8 standard colours
- 8" Black solid front wheels, with 24" Black solid tyre spoked rim
- Aluminium Pre-anodized handrim
- Quick release axle
- Push to lock brakes
- Transportation kit (including buckle belt)
- Alber Efix E35 1566293 (24" wheels)
- Brackets
- Alber anti-tippers
- Adapter Kit
- Mobility Plus package
- Fitting at Invacare



Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

Code	Description	
BASE PRODUCT*		
CEE0003	<input type="checkbox"/> ACTION 3NG SELF PROPEL	(NC 1195) ◆
PROPULSION*		
CEE9259	<input type="checkbox"/> E-Fix E36	£6,086

SEAT WIDTH*			
CEE0180		15" / 380 MM SEAT WIDTH (lapbelt adjustment 420mm - 1520mm)	○
CEE0190		16" / 405 MM SEAT WIDTH (lapbelt adjustment 420mm - 1520mm)	○
CEE0200		17" / 430 MM SEAT WIDTH (lapbelt adjustment 420mm - 1520mm)	◆
CEE0210		18" / 455 MM SEAT WIDTH (lapbelt adjustment 420mm - 1880mm)	○
CEE0220		19" / 480 MM SEAT WIDTH (lapbelt adjustment 420mm - 1880mm)	○
CEE0230		20" / 505 MM SEAT WIDTH (lapbelt adjustment 420mm - 1880mm)	○
SEAT DEPTH*			
CEE0240		16" / 400 MM SEAT DEPTH (NC 0879)	○
CEE0250		18" / 450 MM SEAT DEPTH (NC 0878)	◆
SEAT HEIGHT*			
CEE0260		20" / 510 MM SEAT HEIGHT	○
CEE0270		19" / 485 MM SEAT HEIGHT	○
CEE0280		18" / 460 MM SEAT HEIGHT	◆
SEAT UPHOLSTERY*			
CEE0330		NYLON W/O HOOK AND LOOP STRIPS	◆
BACKREST*			
CEE0410		FIXED 7° BACKREST 2 HEIGHTS 16" / 400 or 17" / 430 MM FIX PUSH HANDLES	◆
BACKREST UPHOLSTERY*			
CEE0474		LUMBAR SUPPORT BI MATERIAL PADDED COVER	◆
LEGREST HANGER RH*			
CEE0510		80° RH SWING IN/SWING OUT	◆
LEGREST HANGER LH*			
CEE0560		80° LH SWING IN/SWING OUT	◆
FOOTREST RH*			
CEE0600		FLIP UP RH FOOTPLATE	◆
FOOTREST LH*			
CEE0620		FLIP UP LH FOOTPLATE	◆
ARMREST RH*			
CEE0640		REMOVABLE RH (N°1) DESK ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED (NC 0760, 0839)	○
CEE0650		REMOVABLE RH (N°1) FULL LENGTH ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED (NC 0760, 0839)	○
CEE0660		SWING AWAY DETACHABLE RH (N°2) DESK 240 MM HEIGHT (NC 0740, 0750, 0839)	◆
ARMREST LH*			
CEE0740		REMOVABLE LH (N°1) DESK ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED	○
CEE0750		REMOVABLE LH (N°1) FULL LENGTH ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED	○
CEE0760		SWING AWAY DETACHABLE LH (N°2) DESK 240 MM HEIGHT	◆
CEE0839		TUBULAR ARMREST LH, SWING AWAY PADDED TUBE WITH REMOVABLE MUDGUARD LH (N°6) (CUSHION IS RECOMMENDED FOR SEAT HEIGHT 410 AND 435 MM) (NC 1200, 1345)	○

SEAT CUSHION			
CEE0849	<input type="checkbox"/>	BASIC 50 MM FOAM CUSHION NYLON AND ANTISLIP COVER	◆
SEAT CUSHION DEPTH			
CEE0878	<input type="checkbox"/>	FOR 400 MM SEAT DEPTH	○
CEE0879	<input type="checkbox"/>	FOR 450 MM SEAT DEPTH	○
COLOUR*			
CEE0961	<input type="checkbox"/>	SLATE GREY (RAL7024 FINE TEXTURE)	○
CEE0964	<input type="checkbox"/>	PEARL WHITE	○
CEE0965	<input type="checkbox"/>	ELECTRIC GREEN	○
CEE0966	<input type="checkbox"/>	MATT BLACK	○
CEE0974	<input type="checkbox"/>	BLUE	○
CEE0976	<input type="checkbox"/>	DARK PURPLE	○
CEE0969	<input type="checkbox"/>	RED MATT	○
CEE0972	<input type="checkbox"/>	WATER BLUE	○
CEE0980	<input type="checkbox"/>	ANTHRACITE BLACK	○
CASTOR HOUSING*			
CEE0995	<input type="checkbox"/>	STANDARD	◆
FRONT CASTORS*			
CEE1010	<input type="checkbox"/>	8" X 1.1/4" BLACK SOLID TYRE	◆
CEE1030	<input type="checkbox"/>	8" X 2" BLACK SOLID TYRE	○
REAR WHEEL SIZE*			
CEE1068	<input type="checkbox"/>	24"	◆
REAR WHEELS*			
CEE1071	<input type="checkbox"/>	STANDARD (BASIC) (BLACK HUB, COLORLESS RIM AND ZINC PLATED SPOKES)	◆
TYRES*			
CEE1082	<input type="checkbox"/>	SOLID BLACK 1.3/8" / 35MM	◆
CEE1088	<input type="checkbox"/>	PNEUMATIC LOW PROFILE BLACK 1" / 25MM	○
HANDRIM*			
CEE1140	<input type="checkbox"/>	ALUMINIUM PREANODIZED HANDRIM	◆
REAR WHEEL AXLE*			
CEE1180	<input type="checkbox"/>	Q/R AXLE	◆
REAR WHEEL POSITION*			
CEE1190	<input type="checkbox"/>	STANDARD REAR WHEELS POSITION	◆
CEE1195	<input type="checkbox"/>	MID WHEEL PROPULSION REAR WHEEL POSITION	○
CEE1200	<input type="checkbox"/>	PASSIVE REAR WHEELS POSITION	○
BRAKES*			
CEE1270	<input type="checkbox"/>	PUSH TO LOCK BRAKES	◆
TRANSPORTATION KIT*			
CEE1345	<input type="checkbox"/>	TRANSPORTATION BRACKETS	◆
SAFETY OPTIONS*			
CEE1360	<input type="checkbox"/>	POSTURE BELT	◆

TOTAL PRICE	<input type="text"/>
--------------------	----------------------

Non Conforming Options:

CEE0740	REMOVABLE LH (N°1) DESK ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED
CEE0750	REMOVABLE LH (N°1) FULL LENGTH ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED
CEE0760	SWING AWAY DETACHABLE LH (N°2) DESK 240 MM HEIGHT
CEE0839	TUBULAR ARMREST LH, SWING AWAY PADDED TUBE WITH REMOVABLE MUDGUARD LH (N°6) (CUSHION IS RECOMMENDED FOR SEAT HEIGHT 410 AND 435 MM)
CEE0878	FOR 400 MM SEAT DEPTH
CEE0879	FOR 450 MM SEAT DEPTH
CEE1195	MID WHEEL PROPULSION REAR WHEEL POSITION
CEE1200	PASSIVE REAR WHEELS POSITION
CEE1345	TRANSPORTATION BRACKETS



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):