



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Motability Retail Prescription Form

MAY 2024
Customer Ref:

Customer Ref\*, Contact No\*, Delivery Address\*, Account No\*, Order No\*, Order Qty\*, Order date\*

Storm4 Modulite M30 708645 - Weekly Price £52.38

Standard Features

- Max user weight 150Kg (23.5 stone)
Crash tested ISO 7176 -Part 19
Modulite seating system
Matrix Flo-tech Lite Cushion
Tension Adjustable Backrest
Standard T-Armrests, detachable with long armpad
Manual settable seat Tilt 0° - 20°
Manual settable Recliner (0° - 30° in 3.75° steps)
VARI F hangers with Fixed Footplate
LiNX Control REM216 withswing away bracket
Pneumatic Castors
Pneumatic Drive Wheels
73.5 Ah Batteries
12km/h (8mph) Motors
LED light system with indicators and Mirror
4 Point tie down system
Lapbelt with autolock (adjustable both sides)
Anti-Tippers



Standard Base price: £7,361

Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

Table with columns: Code, Description, Price. Rows include CHAIR TYPE\* (Storm 4 Modulite), SEATING TYPE\* (Modulite Seating Unit), and SEATING OPTIONS\* (Telescopic Seat Frame).

<b>SEAT WIDTH*</b>			
ATBXXXX	<input type="checkbox"/>	Seat width 330mm (13") seat width achieved on a 380mm seat frame (2121) with inrigged armrests (2183)	<input type="radio"/>
ATB2121	<input type="checkbox"/>	Set at 380mm (adjustable from 380mm to 430mm (15" - 17") with stepless adjustment on armrests, legrests and backrest)	<input type="radio"/>
ATB2122	<input type="checkbox"/>	Set at 430mm (adjustable from 430mm to 480mm (17" - 19) with stepless adjustment on armrests, legrests and backrest)	<input type="radio"/>
ATB2123	<input type="checkbox"/>	Set at 480mm (adjustable from 480mm to 530mm (19" - 21") with stepless adjustment on armrests, legrests and backrest)	<input type="radio"/>
ATB2124	<input type="checkbox"/>	Set at 530mm (adjustable from 530mm to 580mm (21" - 23") with stepless adjustment on armrests, legrests and backrest)	<input type="radio"/>
<b>SEAT DEPTH*</b>			
ATB2132	<input type="checkbox"/>	410mm (16") (adjustable 410mm - 510mm stepless adjustment) (NC 2196)	<input type="radio"/>
ATB2134	<input type="checkbox"/>	460mm (18") (adjustable 410mm - 510mm stepless adjustment)	<input type="radio"/>
ATB2136	<input type="checkbox"/>	510mm (20") (adjustable 410mm - 510mm stepless adjustment)	<input type="radio"/>
<b>SEAT CUSHION*</b>			
ATB2142	<input type="checkbox"/>	Matrx Flo-tech Lite	<input checked="" type="radio"/>
<b>BACKREST SUPPORT*</b>			
ATB2162	<input type="checkbox"/>	Backposts 540mm (Tall) must be selected for Flex3 Laguna backrest	<input type="radio"/>
ATB2163	<input type="checkbox"/>	Backposts (Short)	<input type="radio"/>
<b>PUSH HANDLES FOR BACKPOSTS*</b>			
ATB2107	<input type="checkbox"/>	Backposts WITH push handles	<input checked="" type="radio"/>
<b>BACKREST UPHOLSTERY / CUSHION*</b>			
ATB2171	<input type="checkbox"/>	Basic Tension Adjustable Backrest (540mm or 480mm upholstery height) Straps easily accessible and adjustable from the rear, Black fabric only	<input checked="" type="radio"/>
<b>UPHOLSTREY TYPE*</b>			
ATB2151	<input type="checkbox"/>	Black fabric	<input checked="" type="radio"/>
<b>ARMREST (PAIR)*</b>			
ATB2191	<input type="checkbox"/>	T-Armrest Standard, Range of adjustment from seat plate 245mm to 315mm	<input checked="" type="radio"/>
<b>ARMPAD (PAIR)*</b>			
ATB2196	<input type="checkbox"/>	Long Armpad (375mm) (NC 2132)	<input checked="" type="radio"/>
<b>SEAT HEIGHT AND SEAT RISER*</b>			
ATB1023	<input type="checkbox"/>	Fixed Seat Height 450mm (17.75")	<input type="radio"/>
ATB1025	<input type="checkbox"/>	Fixed Seat Height 480mm (19")	<input type="radio"/>
<b>SEAT TILT*</b>			
ATB1002	<input type="checkbox"/>	Manual Seat Tilt by spindle	<input checked="" type="radio"/>
<b>BACKREST RECLINER*</b>			
ATB0975	<input type="checkbox"/>	Manual Settable Recliner (0° - 30°, in 3.75° steps) Foldable	<input checked="" type="radio"/>
<b>FOOTREST HANGER TYPE (PAIR)*</b>			
ATB0503	<input type="checkbox"/>	Vari F : Angle adjustable with tools 0° - 70°, length 290mm - 460mm	<input checked="" type="radio"/>
<b>FOOTPLATE TYPE*</b>			
ATB0512	<input type="checkbox"/>	Fixed Footplate	<input checked="" type="radio"/>
<b>CONTROLS*</b>			
ATB2204	<input type="checkbox"/>	REM216 Drive, Actuators & Lights	<input checked="" type="radio"/>

<b>CONTROL BRACKETRY*</b>			
ATB1707	<input type="checkbox"/>	Swing Away RH Bracket	<input type="radio"/>
ATB1708	<input type="checkbox"/>	Swing Away LH Bracket	<input type="radio"/>
<b>SHROUD COLOUR*</b>			
ATB0202	<input type="checkbox"/>	Blue	<input type="radio"/>
ATB0228	<input type="checkbox"/>	Black	<input type="radio"/>
ATB0220	<input type="checkbox"/>	Champagne	<input type="radio"/>
ATB0204	<input type="checkbox"/>	Transparent Cherry Red	<input type="radio"/>
<b>FRONT WHEEL*</b>			
ATB0851	<input type="checkbox"/>	3.00 - 6 Grey Pneumatic	<input checked="" type="radio"/>
<b>REAR WHEEL*</b>			
ATB0771	<input type="checkbox"/>	3.00 - 8 Grey Pneumatic	<input checked="" type="radio"/>
<b>MOTOR TYPE*</b>			
ATB1909	<input type="checkbox"/>	12 km/h (8mph) (includes LED light system with indicators)	<input checked="" type="radio"/>
<b>BATTERY*</b>			
ATB1364	<input type="checkbox"/>	70 Ah Gel	<input checked="" type="radio"/>
<b>SUSPENSION*</b>			
ATB1240	<input type="checkbox"/>	Suspension low (up to 80kg user weight)	<input type="radio"/>
ATB1241	<input type="checkbox"/>	Suspension middle (from 70kg up to 100kg user weight)	<input type="radio"/>
ATB1242	<input type="checkbox"/>	Suspension hard (from 90kg up to 150kg user weight)	<input type="radio"/>
<b>CHARGER*</b>			
ATB1343	<input type="checkbox"/>	OptiCharge Charger 8A (200 - 250V)	<input checked="" type="radio"/>
<b>ANTI TIPPER*</b>			
ATB1205	<input type="checkbox"/>	Anti-tipper (standard version)	<input checked="" type="radio"/>
<b>TRANSPORTATION KIT*</b>			
ATB1494	<input type="checkbox"/>	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory, see options below)	<input checked="" type="radio"/>
<b>FIXATION BELT*</b>			
ATB1329	<input type="checkbox"/>	Lap Belt with Auto Lock	<input checked="" type="radio"/>
<b>REAR VIEW MIRROR*</b>			
ATB1415	<input type="checkbox"/>	Rear view Mirror RH - Must be selected with 10 km/h Motor	<input type="radio"/>
ATB1339	<input type="checkbox"/>	Rear view Mirror LH - Must be selected with 10 km/h Motor	<input type="radio"/>
<b>TOTAL PRICE</b>			<input type="text"/>

**Non Conforming Options:**

ATB2132            410mm (16") (adjustable 410mm - 510mm stepless adjustment)

ATB2196            Long Armpad (375mm)



For a more complex customer requirement please refer to our Invacare Unique special order form. email [specialsuk@invacare.com](mailto:specialsuk@invacare.com) or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Description of special request (please include relevant option code or part number):