



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Motability Retail Prescription Form

MAY 2024
Customer Ref:

Customer Ref\*, Contact No\*, Delivery Address\*, Account No\*, Order No\*, Order Qty\*, Order date\*

Storm4 Modulite M31 1001116 - Weekly
Price £56.74

Standard Features

- Max user weight 150Kg (23.5 stone)
Crash tested ISO 7176 -Part 19
Modulite seating system
Flex 3 Laguna Backrest
Standard T-Armrests, detachable with Long armpad
Electrical Tilt 0° to 25°
Manual settable Recliner (0° - 30° in 3.75° steps)
VARI F with Fixed Footplate
LinX Control REM216 with swing away bracket
Pneumatic Castors
Pneumatic Drive Wheels
12 km/h (8mph) Motors, LED light system with indicators
4 Point tie down system
Lapbelt with autolock (adjustable both sides)
Anti-Tippers



Standard Base price: £8,334

Please make your selection with a cross in the box
Free of charge
standard feature
POA
special request

Customer weight = KG

Table with columns: Code, Description, Price. Rows include CHAIR TYPE\* (Storm 4 Modulite), SEATING\* (Modulite Seating Unit), and SEATING OPTION\* (Telescopic Seat Frame).

| <b>SEAT WIDTH*</b>  |                          |   |                                  |
|---|--------------------------|---|----------------------------------|
| ATBXXXX   | <input type="checkbox"/> | Seat width 330mm (13") seat width achieved on a 380mm seat frame (2121) with inrigged armrests (2183)                 | <input type="radio"/>            |
| ATB2121   | <input type="checkbox"/> | Set at 380mm (adjustable from 380mm to 430mm (15" - 17") with stepless adjustment on armrests, legrests and backrest) | <input type="radio"/>            |
| ATB2122   | <input type="checkbox"/> | Set at 430mm (adjustable from 430mm to 480mm (17" - 19) with stepless adjustment on armrests, legrests and backrest)  | <input type="radio"/>            |
| ATB2123   | <input type="checkbox"/> | Set at 480mm (adjustable from 480mm to 530mm (19" - 21") with stepless adjustment on armrests, legrests and backrest) | <input type="radio"/>            |
| ATB2124   | <input type="checkbox"/> | Set at 530mm (adjustable from 530mm to 580mm (21" - 23") with stepless adjustment on armrests, legrests and backrest) | <input type="radio"/>            |
| <b>SEAT DEPTH*</b>  |                          |   |                                  |
| ATB2132   | <input type="checkbox"/> | 410mm (16") (adjustable 410mm - 510mm stepless adjustment) (NC 2196)  | <input type="radio"/>            |
| ATB2134   | <input type="checkbox"/> | 460mm (18") (adjustable 410mm - 510mm stepless adjustment)  | <input type="radio"/>            |
| ATB2136   | <input type="checkbox"/> | 510mm (20") (adjustable 410mm - 510mm stepless adjustment)  | <input type="radio"/>            |
| <b>SEAT CUSHION*</b>  |                          |   |                                  |
| ATB2142   | <input type="checkbox"/> | Matrx Flo-tech Lite   | <input checked="" type="radio"/> |
| <b>BACKREST SUPPORT*</b>  |                          |   |                                  |
| ATB2162   | <input type="checkbox"/> | Backposts 540mm (Tall) must be selected for Flex3 Laguna backrest   | <input checked="" type="radio"/> |
| <b>PUSH HANDLES FOR BACKPOSTS*</b>                                      |                          |   |                                  |
| ATB2107   | <input type="checkbox"/> | Backposts WITH push handles   | <input checked="" type="radio"/> |
| <b>BACKREST UPHOLSTERY / CUSHION*</b>                                   |                          |   |                                  |
| ATB2173   | <input type="checkbox"/> | Flex3 adjustable Backrest Plates with Laguna Backrest cushion   | <input checked="" type="radio"/> |
| <b>UPHOLSTERY TYPE*</b>   |                          |   |                                  |
| ATB2152   | <input type="checkbox"/> | Black multi stretch PU coated fabric  | <input type="radio"/>            |
| ATB2155   | <input type="checkbox"/> | Black Spacer fabric   | <input type="radio"/>            |
| <b>ARMREST (PAIR)*</b>  |                          |   |                                  |
| ATB2191   | <input type="checkbox"/> | T-Armrest Standard, Range of adjustment from seat plate 245mm to 315mm  | <input checked="" type="radio"/> |
| <b>ARMPAD (PAIR)*</b>   |                          |   |                                  |
| ATB2196   | <input type="checkbox"/> | Long Armpad (375mm) (NC 2132)   | <input checked="" type="radio"/> |
| <b>HEADREST / NECKREST - Adequate headrest option must be selected*</b> |                          |   |                                  |
| ATB1350   | <input type="checkbox"/> | neckrest: height and depth adjustable   | <input checked="" type="radio"/> |
| <b>SEAT HEIGHT AND SEAT RISER*</b>                                      |                          |   |                                  |
| ATB1023   | <input type="checkbox"/> | Fixed Seat Height 450mm (17.75")  | <input type="radio"/>            |
| ATB1025   | <input type="checkbox"/> | Fixed Seat Height 480mm (19")   | <input type="radio"/>            |
| <b>SEAT TILT*</b>   |                          |   |                                  |
| ATB1014   | <input type="checkbox"/> | Electrical Tilt 0° to 25°   | <input checked="" type="radio"/> |
| <b>BACKREST RECLINER*</b>   |                          |   |                                  |
| ATB0975   | <input type="checkbox"/> | Manual Settable Recliner (0° - 30°, in 3.75° steps) Foldable  | <input checked="" type="radio"/> |
| <b>FOOTREST HANGER TYPE (PAIR)*</b>                                     |                          |   |                                  |
| ATB0503   | <input type="checkbox"/> | Vari F : Angle adjustable with tools 0° - 70°, length 290mm - 460mm   | <input type="radio"/>            |
| <b>FOOTPLATE TYPE*</b>  |                          |   |                                  |
| ATB0512   | <input type="checkbox"/> | Fixed Footplate   | <input type="radio"/>            |

| <b>CONTROLS*</b>           |                          |   |             |
|----------------------------|--------------------------|---|-------------|
| ATB2204                    | <input type="checkbox"/> | REM216 Drive, Actuators & Lights  | ◆           |
| <b>CONTROL BRACKETRY*</b>  |                          |   |             |
| ATB1707                    | <input type="checkbox"/> | Swing Away RH Bracket   | ○           |
| ATB1708                    | <input type="checkbox"/> | Swing Away LH Bracket   | ○           |
| <b>SHROUD COLOUR*</b>      |                          |   |             |
| ATB0202                    | <input type="checkbox"/> | Blue  | ○           |
| ATB0228                    | <input type="checkbox"/> | Black   | ○           |
| ATB0220                    | <input type="checkbox"/> | Champagne   | ○           |
| ATB0204                    | <input type="checkbox"/> | Transparent Cherry Red  | ○           |
| <b>FRONT WHEEL*</b>        |                          |   |             |
| ATB0851                    | <input type="checkbox"/> | 3.00 - 6 Grey Pneumatic   | ◆           |
| <b>REAR WHEEL*</b>         |                          |   |             |
| ATB0771                    | <input type="checkbox"/> | 3.00 - 8 Grey Pneumatic   | ◆           |
| <b>MOTOR TYPE*</b>         |                          |   |             |
| ATB1909                    | <input type="checkbox"/> | 12 km/h (8mph) (includes LED light system with indicators)  | (NC 1437) ◆ |
| <b>BATTERY*</b>            |                          |   |             |
| ATB1364                    | <input type="checkbox"/> | 70 Ah Gel   | ◆           |
| <b>SUSPENSION*</b>         |                          |   |             |
| ATB1240                    | <input type="checkbox"/> | Suspension low (up to 80kg user weight)   | ○           |
| ATB1241                    | <input type="checkbox"/> | Suspension middle (from 70kg up to 100kg user weight)   | ○           |
| ATB1242                    | <input type="checkbox"/> | Suspension hard (from 90kg up to 150kg user weight)   | ○           |
| <b>CHARGER*</b>            |                          |   |             |
| ATB1343                    | <input type="checkbox"/> | OptiCharge Charger 8A (200 - 250V)  | ◆           |
| <b>LIGHTS*</b>             |                          |   |             |
| ATB1437                    | <input type="checkbox"/> | LED Light system with indicators  | (NC 1909) ◆ |
| <b>ANTI TIPPER*</b>        |                          |   |             |
| ATB1205                    | <input type="checkbox"/> | Anti-tipper (standard version)  | ◆           |
| <b>TRANSPORTATION KIT*</b> |                          |   |             |
| ATB1494                    | <input type="checkbox"/> | 4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory, see options below) | ◆           |
| <b>FIXATION BELT*</b>      |                          |   |             |
| ATB2113                    | <input type="checkbox"/> | Lap belt with auto lock (flexible mount in side rails)  | ◆           |
| <b>REAR VIEW MIRROR</b>    |                          |   |             |
| ATB1415                    | <input type="checkbox"/> | Rear view Mirror RH - Must be selected with 10 km/h Motor   | ○           |
| ATB1339                    | <input type="checkbox"/> | Rear view Mirror LH - Must be selected with 10 km/h Motor   | ○           |

|                    |  |
|--------------------|--|
| <b>TOTAL PRICE</b> |  |
|--------------------|--|

**Non Conforming Options:**

ATB2132 410mm (16") (adjustable 410mm - 510mm stepless adjustment)

ATB2196 Long Armpad (375mm)

ATB1909 12 km/h (8mph) (includes LED light system with indicators)

ATB1437 LED Light system with indicators



For a more complex customer requirement please refer to our Invacare Unique special order form. email [specialsuk@invacare.com](mailto:specialsuk@invacare.com) or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Description of special request (please include relevant option code or part number):