



Motability Retail Prescription Form

Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

MAY 2024
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Storm4 X-Plore Modulite M36 1001119 - Weekly Price £59.30

Standard Features

- Max user weight 150Kg (23.5 stone)
- Crash tested ISO7176-Part 19
- Modulite seating system
- Flex3 Laguna Backrest
- Standard T-Armrests, detachable with Long Armpads
- Electrical Seat Riser with Manual settable seat Tilt 0° - 20°
- Manual settable Recliner (0° - 30° in 3.75° steps)
- VARI F hangers with Fixed Footplate
- LiNX control REM216 with swing away bracket
- Puncture proof castors
- Puncture proof Drive Wheels
- 12 km/h (8mph) Motors, LED light system with indicators
- 4 Point-tie down system
- Lapbelt with autolock (adjustable both sides)
- Anti-Tippers
- X-plore suspension 50mm all round



Standard Base price: £9,113

Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

Code	Description	
CHAIR TYPE*		
ATE0017	<input type="checkbox"/> Storm 4 X-plore	£9,113
SEATING TYPE*		
ATE2100	<input type="checkbox"/> Invacare Modulite 'Flex3'	◆
SEAT FRAME*		
ATE2111	<input type="checkbox"/> Telescopic Seat Frame with metal seat plate	◆

SEAT WIDTH*			
ATEXXXX	<input type="checkbox"/>	Seat width 330mm (13")	<input type="radio"/>
ATE2121	<input type="checkbox"/>	Set at 380mm	<input type="radio"/>
ATE2122	<input type="checkbox"/>	Set at 430mm	<input type="radio"/>
ATE2123	<input type="checkbox"/>	Set at 480mm	<input type="radio"/>
ATE2124	<input type="checkbox"/>	Set at 530mm	<input type="radio"/>
SEAT DEPTH*			
ATE2132	<input type="checkbox"/>	410mm (16") (NC 2196)	<input type="radio"/>
ATE2134	<input type="checkbox"/>	460mm (18")	<input type="radio"/>
ATE2136	<input type="checkbox"/>	510mm (20")	<input type="radio"/>
SEAT CUSHION*			
ATE2142	<input type="checkbox"/>	Matrx Flo-tech Lite	<input checked="" type="radio"/>
BACKREST SUPPORT*			
ATE2162	<input type="checkbox"/>	Backposts 540mm (Tall)	<input type="radio"/>
PUSH HANDLES FOR BACKPOSTS*			
ATE2107	<input type="checkbox"/>	Backposts WITH push handles	<input checked="" type="radio"/>
BACKREST UPHOLSTREY / CUSHION*			
ATE2173	<input type="checkbox"/>	Flex3 adjustable Backrest Plates with Laguna Backrest cushion	<input checked="" type="radio"/>
UPHOLSTERY COLOUR*			
ATE2152	<input type="checkbox"/>	Black multi stretch PU coated fabric	<input type="radio"/>
ATE2155	<input type="checkbox"/>	Black Spacer fabric	<input type="radio"/>
ARMREST*			
ATE2191	<input type="checkbox"/>	T-Armrest Standard	<input checked="" type="radio"/>
ARMPAD*			
ATE2196	<input type="checkbox"/>	Long Armpad (375mm)	<input checked="" type="radio"/>
HEADREST / NECKREST - Adequate headrest option must be selected*			
ATE1350	<input type="checkbox"/>	Neckrest: height and depth Adjustable	<input checked="" type="radio"/>
SEAT HEIGHT SETTING & SEAT RISER*			
ATE1040	<input type="checkbox"/>	Electrical Seat Riser (440mm - 690mm)	<input checked="" type="radio"/>
BACKREST RECLINER*			
ATE0975	<input type="checkbox"/>	Manual Settable Recliner (0° - 30°, in 3.75° steps) Foldable	<input checked="" type="radio"/>
FOOTREST HANGER TYPE (PAIR)*			
ATE0503	<input type="checkbox"/>	Vari F : Angle adjustable with tools 0° - 70°, length 290mm - 460mm	<input checked="" type="radio"/>
FOOTPLATE TYPE*			
ATE0512	<input type="checkbox"/>	Fixed Footplate	<input checked="" type="radio"/>
CONTROLS*			
ATE2204	<input type="checkbox"/>	REM216 Drive, Actuators & Lights	<input checked="" type="radio"/>
CONTROL BRACKETRY*			
ATE1707	<input type="checkbox"/>	Swing Away RH Bracket	<input type="radio"/>
ATE1708	<input type="checkbox"/>	Swing Away LH Bracket	<input type="radio"/>
SHROUD COLOUR*			
ATE0202	<input type="checkbox"/>	Blue	<input type="radio"/>
ATE0228	<input type="checkbox"/>	Black	<input type="radio"/>
ATE0220	<input type="checkbox"/>	Champagne	<input type="radio"/>
ATE0204	<input type="checkbox"/>	Transparent Cherry Red	<input type="radio"/>

CASTOR WHEEL*		
ATE0853	<input type="checkbox"/>	3.00 - 6 Grey Puncture Proof (sealed) ◆
DRIVE WHEEL*		
ATE0773	<input type="checkbox"/>	3.00 - 8 Grey Puncture Proof (sealed) ◆
MOTOR TYPE*		
ATE1909	<input type="checkbox"/>	12 km/h (8mph) includes LED light system with indicators ◆
BATTERY*		
ATE1364	<input type="checkbox"/>	73.5 Ah GEL ◆
REAR SUSPENSION*		
ATE1233	<input type="checkbox"/>	Hitec spring X-plore ○
FRONT SUSPENSION*		
ATE1238	<input type="checkbox"/>	Front suspension (suggested for up to 100kg) ○
ATE1239	<input type="checkbox"/>	Front suspension (suggested for above 100kg) ○
CHARGER*		
ATE1343	<input type="checkbox"/>	OptiCharge Charger 8A (200 - 250V) ◆
LIGHTS*		
ATZ1006	<input type="checkbox"/>	Lights Included with Motor Package 10 & 12kmh ◆
ANTI-TIPPERS*		
ATE1205	<input type="checkbox"/>	Anti-tipper (standard version) ◆
TRANSPORTATION*		
ATE1494	<input type="checkbox"/>	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory). ◆
FIXATION BELT*		
ATE2113	<input type="checkbox"/>	Lapbelt (adjustable both sides) adjustable mounting on seat rail extrusion ◆
REAR VIEW MIRROR		
ATE1415	<input type="checkbox"/>	Rear view Mirror RH - Must be selected with 10 km/h Motor ○
ATE1339	<input type="checkbox"/>	Rear view Mirror LH - Must be selected with 10 km/h Motor ○

TOTAL PRICE	
--------------------	--

Non Conforming Options:

ATE2196 Long Armpad (375mm)



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):