



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Motability Retail Prescription Form

MAY 2024
Customer Ref:

Customer Ref\*, Contact No\*, Delivery Address\*, Account No\*, Order No\*, Order Qty\*, Order date\*

Storm4 X-Plore Modulite M37 1001120 - Weekly Price £59.76

Standard Features

- Max user weight 150Kg (23.5 stone)
Crash tested ISO7176-Part 19
Modulite seating system
Flex3 Laguna Backrest
Standard T-Armrests, detachable with Long Armpads
Neckrest
Electrical Seat Riser, Electrical Tilt from 0° to 25°
Manual settable Recliner (0° - 30° in 3.75° steps)
VARI F hangers with Fixed Footplate
LiNX control REM216 with swing away bracket
Puncture proof castors
Puncture proof Drive Wheels
12 km/h (8mph) Motors, LED light system with indicators
4 Point-tie down system
Lapbelt with autolock (adjustable both sides)
Anti-Tippers
X-plore suspension 50mm all round

Standard Base price: £9,653



Please make your selection with a cross in the box [ ] ( ) Free of charge [ ] [ ] standard feature [ ] POA special request [ ]

Customer weight = KG

Table with columns: Code, Description, Price. Rows include CHAIR TYPE\* (ATE0017 Storm 4 X-plore £9,653), SEATING TYPE\* (ATE2100 Invacare Modulite 'Flex3'), SEAT FRAME\* (ATE2111 Telescopic Seat Frame with metal seat plate).

<b>SEAT WIDTH*</b>			
ATEXXXX	<input type="checkbox"/>	Seat width 330mm (13")	<input type="radio"/>
ATE2121	<input type="checkbox"/>	Set at 380mm	<input type="radio"/>
ATE2122	<input type="checkbox"/>	Set at 430mm	<input type="radio"/>
ATE2123	<input type="checkbox"/>	Set at 480mm	<input type="radio"/>
ATE2124	<input type="checkbox"/>	Set at 530mm	<input type="radio"/>
<b>SEAT DEPTH*</b>			
ATE2132	<input type="checkbox"/>	410mm (16") (NC 2196)	<input type="radio"/>
ATE2134	<input type="checkbox"/>	460mm (18")	<input type="radio"/>
ATE2136	<input type="checkbox"/>	510mm (20")	<input type="radio"/>
<b>SEAT CUSHION*</b>			
ATE2142	<input type="checkbox"/>	Matrx Flo-tech Lite	<input checked="" type="checkbox"/>
<b>BACKREST SUPPORT*</b>			
ATE2162	<input type="checkbox"/>	Backposts 540mm (Tall)	<input checked="" type="checkbox"/>
<b>PUSH HANDLES FOR BACKPOSTS*</b>			
ATE2107	<input type="checkbox"/>	Backposts WITH push handles	<input checked="" type="checkbox"/>
<b>BACKREST UPHOLSTREY / CUSHION*</b>			
ATE2173	<input type="checkbox"/>	Flex3 adjustable Backrest Plates with Laguna Backrest cushion	<input checked="" type="checkbox"/>
<b>UPHOLSTERY COLOUR*</b>			
ATE2152	<input type="checkbox"/>	Black multi stretch PU coated fabric	<input type="radio"/>
ATE2155	<input type="checkbox"/>	Black Spacer fabric	<input type="radio"/>
<b>ARMREST*</b>			
ATE2191	<input type="checkbox"/>	T-Armrest Standard	<input checked="" type="checkbox"/>
<b>ARMPAD*</b>			
ATE2196	<input type="checkbox"/>	Long Armpad (375mm)	<input checked="" type="checkbox"/>
<b>HEADREST / NECKREST - Adequate headrest option must be selected*</b>			
ATE1350	<input type="checkbox"/>	Neckrest: height and depth Adjustable	<input checked="" type="checkbox"/>
<b>SEAT HEIGHT SETTING &amp; SEAT RISER*</b>			
ATE1040	<input type="checkbox"/>	Electrical Seat Riser (440mm - 690mm)	<input checked="" type="checkbox"/>
<b>SEAT TILT*</b>			
ATE1014	<input type="checkbox"/>	Electrical Tilt from 0° to 25°	<input checked="" type="checkbox"/>
<b>BACKREST RECLINER*</b>			
ATE0975	<input type="checkbox"/>	Manual Settable Recliner (0° - 30°, in 3.75° steps) Foldable	<input checked="" type="checkbox"/>
<b>FOOTREST HANGER TYPE (PAIR)*</b>			
ATE0503	<input type="checkbox"/>	Vari F : Angle adjustable with tools 0° - 70°, length 290mm - 460mm	<input checked="" type="checkbox"/>
<b>FOOTPLATE TYPE*</b>			
ATE0512	<input type="checkbox"/>	Fixed Footplate	<input checked="" type="checkbox"/>
<b>CONTROLS*</b>			
ATE2204	<input type="checkbox"/>	REM216 Drive, Actuators & Lights	<input checked="" type="checkbox"/>
<b>CONTROL BRACKETRY*</b>			
ATE1707	<input type="checkbox"/>	Swing Away RH Bracket	<input type="radio"/>
ATE1708	<input type="checkbox"/>	Swing Away LH Bracket	<input type="radio"/>

<b>SHROUD COLOUR*</b>			
ATE0202	<input type="checkbox"/>	Blue	<input type="radio"/>
ATE0228	<input type="checkbox"/>	Black	<input type="radio"/>
ATE0220	<input type="checkbox"/>	Champagne	<input type="radio"/>
ATE0204	<input type="checkbox"/>	Transparent Cherry Red	<input type="radio"/>
<b>CASTOR WHEEL*</b>			
ATE0853	<input type="checkbox"/>	3.00 - 6 Grey Puncture Proof (sealed)	<input type="checkbox"/>
<b>DRIVE WHEEL*</b>			
ATE0773	<input type="checkbox"/>	3.00 - 8 Grey Puncture Proof (sealed)	<input type="checkbox"/>
<b>MOTOR TYPE*</b>			
ATE1909	<input type="checkbox"/>	12 km/h (8mph) includes LED light system with indicators	<input type="checkbox"/>
<b>BATTERY*</b>			
ATE1364	<input type="checkbox"/>	73.5 Ah GEL	<input type="checkbox"/>
<b>REAR SUSPENSION*</b>			
ATE1233	<input type="checkbox"/>	Hitec spring X-plore	<input type="radio"/>
<b>FRONT SUSPENSION*</b>			
ATE1238	<input type="checkbox"/>	Front suspension (suggested for up to 100kg)	<input type="radio"/>
ATE1239	<input type="checkbox"/>	Front suspension (suggested for above 100kg)	<input type="radio"/>
<b>CHARGER*</b>			
ATE1343	<input type="checkbox"/>	OptiCharge Charger 8A (200 - 250V)	<input type="checkbox"/>
<b>LIGHTS*</b>			
ATZ1006	<input type="checkbox"/>	Lights Included with Motor Package 10 & 12kmh	<input type="checkbox"/>
<b>ANTI-TIPPERS*</b>			
ATE1205	<input type="checkbox"/>	Anti-tipper (standard version)	<input type="checkbox"/>
<b>TRANSPORTATION*</b>			
ATE1494	<input type="checkbox"/>	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory).	<input type="checkbox"/>
<b>FIXATION BELT*</b>			
ATE2113	<input type="checkbox"/>	Lapbelt (adjustable both sides) adjustable mounting on seat rail extrusion	<input type="checkbox"/>
<b>REAR VIEW MIRROR</b>			
ATE1415	<input type="checkbox"/>	Rear view Mirror RH - Must be selected with 10 km/h Motor	<input type="radio"/>
ATE1339	<input type="checkbox"/>	Rear view Mirror LH - Must be selected with 10 km/h Motor	<input type="radio"/>

<b>TOTAL PRICE</b>	<input type="text"/>
--------------------	----------------------

**Non Conforming Options:**

ATE2196          Long Armpad (375mm)

---



For a more complex customer requirement please refer to our Invacare Unique special order form. email [specialsuk@invacare.com](mailto:specialsuk@invacare.com) or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Description of special request (please include relevant option code or part number):