



Motability Retail Prescription Form

Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

MAY 2024
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

**TDX SP2 HD M17 1001588 - Weekly Price
£59.81**

Standard Features

- Max user weight 160kg (28 stone)
- Crash tested ISO 7176 -Part 19
- Modulite Seating System
- Matrix Flo-tech Cushion
- Flex 3 adjustable Backrest plates with Laguna backrest cushion
- Standard T-armrests, detachable with long armpad
- Elan Headrest
- Seat lifter with electrical tilt (300 mm, 30°) with CoG shift
- Manual settable backrest Recline 0° - 30° in 3.75° steps
- Centre mounted Manual legrest
- LiNX control REM216 swing away bracket
- Solid Black tyres
- 12km/h (8mph) 4 pole motors
- 4-Point tie-down system (Q'Straint QLK-150 compatible)
- indicators, LED headlights and Mirror
- Lapbelt with autolock (adjustable both sides)



Standard Base price: £8,939

Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

Code	Description	
CHAIR TYPE*		
ARJ0000	<input type="checkbox"/> TDX SP2 Heavy Duty	◆
User Weight*		
ARJ1918	<input type="checkbox"/> User weight less 100 kg	○
ARJ1919	<input type="checkbox"/> User weight 101 to 125 kg	○
ARJ1921	<input type="checkbox"/> User weight 126 to 136 kg	○
ARJ1922	<input type="checkbox"/> User weight 137 to 150 kg	○
BASE TYPE*		
ARJ2100	<input type="checkbox"/> Modulite Seating unit	£8,939

SEAT WIDTH*			
ARJ2121		Set at 380mm (15") (adjustable from 380mm to 430mm (15" - 17") with stepless adjustment on armrests, legrests and backrest) (NC 2138 + 2142)	○
ARJ2122		Set at 430mm (17") (adjustable from 430mm to 480mm (17" - 19") with stepless adjustment on armrests, legrests and backrest) (NC 2138 + 2142)	○
ARJ2123		Set at 480mm (19") (adjustable from 480mm to 530mm (19" - 21") with stepless adjustment on armrests, legrests and backrest) (NC 2138 + 2142)	○
ARJ2124		Set at 530mm (21") (adjustable from 530mm to 580mm (21" - 23") with stepless adjustment on armrests, legrests and backrest)	○
SEAT DEPTH RAIL LENGTH*			
ARJ2129		Standard Depth Package (16" - 20", 410mm - 510mm simple stepless adjustment from the side) (NC 2138, 2196)	○
ARJ2130		Extended Depth Package (20" - 22", 510mm - 560mm simple stepless adjustment from the side) (NC 2132, 2134)	○
SEAT DEPTH SETTING*			
ARJ2132		410mm (16") (410mm - 510mm simple stepless adjustment from the side) (NC 2196)	○
ARJ2134		460mm (18") (410mm - 510mm simple stepless adjustment from the side)	○
ARJ2136		510mm (20") (410mm - 510mm simple stepless adjustment from the side)	○
ARJ2138		560mm (22") (510mm - 560mm simple stepless adjustment from the side) (Cushion (option 2142) is only available with Seat Width (option 2124) 530mm)	○
SEAT CUSHION*			
ARJ2142		Matrx Flo-tech Lite	◆
BACKREST SUPPORT*			
ARJ2162		Backposts 540mm (Tall) must be selected for Flex3 adjustable backrest plates	◆
PUSH HANDLES FOR BACKPOSTS*			
ARJ2107		Backposts WITH push handles	◆
BACKREST UPHOLSTERY / CUSHION*			
ARJ2173		Flex3 adjustable Backrest Plates with Laguna Backrest cushion	◆
UPHOLSTERY COLOUR*			
ARJ2152		Black Black multi stretch PU coated fabric Upholstery	○
ARJ2155		Black Spacer fabric	○
ARMREST*			
ARJ2191		T-Armrest Standard, Range of adjustment from seat plate 245mm to 315mm	◆
ARMPAD*			
ARJ2196		Long Armpad (375mm) (NC 2132)	◆
HEADREST / NECKREST - Adequate headrest option must be selected*			
ARJ2324		Elan Standard Pad 355mm (10") wide with standard cover complete with Elan Headrest Hardware	◆
SEAT LIFTER / SEAT TILT*			
ARJ1094		Seat lifter with Electrical Tilt (300mm, 30°) with CoG shift Seat Height 450mm - 480mm (Also choose option 1265)	◆
SEAT HEIGHT SETTING FOR 20 DEGREES SEAT TILT			
ARZ1003		Seat height set at 460mm (18") (from floor to front of seat plate at 0° tilt)	◆
BACKREST RECLINER*			
ARJ0975		Manual Settable Recliner (0° - 30°, in 3.75° steps) Foldable	◆

FOOTREST HANGER TYPE*			
ARJ0556	<input type="checkbox"/>	Centre Mounted: Manual : Hanger, centre mounted, angle settable (0° - 90°) : knee to heel length 390 - 510mm, (including calf pad and angle adjustable footplate)	◆
FOOTPLATE TYPE*			
ARJ0559	<input type="checkbox"/>	Angle Adjustable Footplate for Centre Mounted Hanger	◆
CONTROLS*			
ARJ2204	<input type="checkbox"/>	REM216 Drive, Actuators & Lights	◆
CONTROL BRACKETRY*			
ARJ1707	<input type="checkbox"/>	Swing Away RH Bracket	○
ARJ1708	<input type="checkbox"/>	Swing Away LH Bracket	○
SHROUD COLOUR*			
ARJ0228	<input type="checkbox"/>	Wet Black	○
ARJ2601	<input type="checkbox"/>	SpaceStation Silver	○
ARJ2602	<input type="checkbox"/>	Island Blue	○
ARJ2603	<input type="checkbox"/>	Grape JellyBean (Purple)	○
ARJ2607	<input type="checkbox"/>	Rockstar Red	○
ARJ2608	<input type="checkbox"/>	Lights Out Black - Matt	○
ARJ2609	<input type="checkbox"/>	Invacare Blue	○
RIM COLOUR*			
ARJ2611	<input type="checkbox"/>	SpaceStation Silver	○
ARJ2612	<input type="checkbox"/>	Island Blue	○
ARJ2613	<input type="checkbox"/>	Grape JellyBean (Purple)	○
ARJ2617	<input type="checkbox"/>	Rockstar Red	○
ARJ2618	<input type="checkbox"/>	Lights Out Black - Matt	○
ARJ2619	<input type="checkbox"/>	Invacare Blue	○
ARJ2620	<input type="checkbox"/>	Wet Black	○
ARJ2621	<input type="checkbox"/>	Without Rim insert	○
FORKS*			
ARJ2632	<input type="checkbox"/>	Double sided fork	◆
CENTRE WHEELS (DRIVING WHEELS)*			
ARJ0733	<input type="checkbox"/>	355mm x 76mm (14" x 3") Black Tyres with Gel Foam Inserts	◆
MOTOR TYPE*			
ARJ1224	<input type="checkbox"/>	12 km/h (8mph) (includes light system and indicators) G-Trac recommended	◆
BATTERY*			
ARJ1260	<input type="checkbox"/>	73.5 Ah (Gel)	◆
CHARGER*			
ARJ1343	<input type="checkbox"/>	OptiCharge Charger 8A (200 - 250V)	◆
TRANSPORTATION KIT*			
ARJ1494	<input type="checkbox"/>	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory, see options below)	◆
FIXATION BELT*			
ARJ2113	<input type="checkbox"/>	Lapbelt (adjustable both sides) adjustable mounting on seat rail extrusion	◆

REAR VIEW MIRROR*

ARJ1415	<input type="checkbox"/>	Rear view Mirror RH - Must be selected with 10 km/h Motor	<input type="radio"/>
ARJ1339	<input type="checkbox"/>	Rear view Mirror LH - Must be selected with 10 km/h Motor	<input type="radio"/>

TOTAL PRICE			
--------------------	--	--	--

Non Conforming Options:

ARJ2132	410mm (16") (410mm - 510mm simple stepless adjustment from the side)
ARJ2134	460mm (18") (410mm - 510mm simple stepless adjustment from the side)
ARJ2138	560mm (22") (510mm - 560mm simple stepless adjustment from the side) (Cushion (option 2142) is only available with Seat Width (option 2124) 530mm)
ARJ2196	Long Armpad (375mm)



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):