



# Motability Retail Prescription Form

**Yes, you can:**

Tel: 01656 776222 Fax: 01656 776220  
Email: ordersuk@invacare.com  
Online Spares available at: www.invacare.co.uk

MAY 2024  
Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

**TDX SP2 HD M12 735799 - Weekly Price  
£59.21**

### Standard Features

- Max user weight 150kg (23.5 stone)
- Crash tested ISO 7176 -Part 19
- Modulite Seating System
- Matrix Flo-tech Cushion
- Flex 3 adjustable Backrest plates with Laguna backrest cushion
- Standard T-armrests, detachable with short armpad
- Elan headrest
- Seat lifter with Electrical Tilt (300mm, 30°) with CoG shift
- Manual Centre mount legrest
- LiNX control REM216 swing away bracket
- Solid Black tyres
- 10km/h (6mph) 4 pole motors
- 4-Point tie-down system (Q'Straint QLK-150 compatible)
- indicators, LED headlights and Mirror
- Lapbelt with autolock (adjustable both sides)



**Standard Base price: £8,939**

Please make your selection with a cross in the box     Free of charge     standard feature    POA    special request

Customer weight =                      KG

Code	Description	
<b>CHAIR TYPE*</b>		
ARJ0000	<input type="checkbox"/> TDX SP2 Heavy Duty	◆
<b>User Weight*</b>		
ARJ1918	<input type="checkbox"/> User weight less 100 kg	○
ARJ1919	<input type="checkbox"/> User weight 101 to 125 kg	○
ARJ1921	<input type="checkbox"/> User weight 126 to 136 kg	○
ARJ1922	<input type="checkbox"/> User weight 137 to 150 kg	○
<b>BASE TYPE*</b>		
ARJ2100	<input type="checkbox"/> Modulite Seating unit	£8,939

SEAT WIDTH*			
ARJ2121	<input type="checkbox"/>	Set at 380mm (15") (adjustable from 380mm to 430mm (15" - 17") with stepless adjustment on armrests, legrests and backrest) (NC 2138 + 2142)	<input type="radio"/>
ARJ2122	<input type="checkbox"/>	Set at 430mm (17") (adjustable from 430mm to 480mm (17" - 19) with stepless adjustment on armrests, legrests and backrest) (NC 2138 + 2142)	<input type="radio"/>
ARJ2123	<input type="checkbox"/>	Set at 480mm (19") (adjustable from 480mm to 530mm (19" - 21") with stepless adjustment on armrests, legrests and backrest) (NC 2138 + 2142)	<input type="radio"/>
ARJ2124	<input type="checkbox"/>	Set at 530mm (21") (adjustable from 530mm to 580mm (21" - 23") with stepless adjustment on armrests, legrests and backrest)	<input type="radio"/>
SEAT DEPTH RAIL LENGTH*			
ARJ2129	<input type="checkbox"/>	Standard Depth Package (16" - 20", 410mm - 510mm simple stepless adjustment from the side) (NC 2138, 2196)	<input type="radio"/>
ARJ2130	<input type="checkbox"/>	Extended Depth Package (20" - 22", 510mm - 560mm simple stepless adjustment from the side) (NC 2132, 2134)	<input type="radio"/>
SEAT DEPTH SETTING*			
ARJ2132	<input type="checkbox"/>	410mm (16") (410mm - 510mm simple stepless adjustment from the side) (NC 2196)	<input type="radio"/>
ARJ2134	<input type="checkbox"/>	460mm (18") (410mm - 510mm simple stepless adjustment from the side)	<input type="radio"/>
ARJ2136	<input type="checkbox"/>	510mm (20") (410mm - 510mm simple stepless adjustment from the side)	<input type="radio"/>
ARJ2138	<input type="checkbox"/>	560mm (22") (510mm - 560mm simple stepless adjustment from the side) (Cushion (option 2142) is only available with Seat Width (option 2124) 530mm)	<input type="radio"/>
SEAT CUSHION*			
ARJ2142	<input type="checkbox"/>	Matrx Flo-tech Lite	<input checked="" type="checkbox"/>
BACKREST SUPPORT*			
ARJ2162	<input type="checkbox"/>	Backposts 540mm (Tall) must be selected for Flex3 adjustable backrest plates	<input checked="" type="checkbox"/>
PUSH HANDLES FOR BACKPOSTS*			
ARJ2107	<input type="checkbox"/>	Backposts WITH push handles	<input checked="" type="checkbox"/>
BACKREST UPHOLSTERY / CUSHION*			
ARJ2173	<input type="checkbox"/>	Flex3 adjustable Backrest Plates with Laguna Backrest cushion	<input checked="" type="checkbox"/>
UPHOLSTERY COLOUR*			
ARJ2152	<input type="checkbox"/>	Black Black multi stretch PU coated fabric Upholstery	<input type="radio"/>
ARJ2155	<input type="checkbox"/>	Black Spacer fabric	<input type="radio"/>
ARMREST*			
ARJ2191	<input type="checkbox"/>	T-Armrest Standard, Range of adjustment from seat plate 245mm to 315mm	<input checked="" type="checkbox"/>
ARMPAD*			
ARJ2196	<input type="checkbox"/>	Long Armpad (375mm) (NC 2132)	<input checked="" type="checkbox"/>
HEADREST / NECKREST - Adequate headrest option must be selected*			
ARJ2324	<input type="checkbox"/>	Elan Standard Pad 355mm (10") wide with standard cover complete with Elan Headrest Hardware	<input checked="" type="checkbox"/>
SEAT LIFTER / SEAT TILT*			
ARJ1094	<input type="checkbox"/>	Seat lifter with Electrical Tilt (300mm, 30°) with CoG shift Seat Height 450mm - 480mm (Also choose option 1265)	<input checked="" type="checkbox"/>
SEAT HEIGHT SETTING FOR 20 DEGREES SEAT TILT*			
ARZ1003	<input type="checkbox"/>	Seat height set at 460mm (18") (from floor to front of seat plate at 0° tilt)	<input checked="" type="checkbox"/>
BACKREST RECLINER*			
ARJ0975	<input type="checkbox"/>	Manual Settable Recliner (0° - 30°, in 3.75° steps) Foldable	<input checked="" type="checkbox"/>

<b>FOOTREST HANGER TYPE*</b>			
ARJ0556	<input type="checkbox"/>	Centre Mounted: Manual : Hanger, centre mounted, angle settable (0° - 90°) : knee to heel length 390 - 510mm, (including calf pad and angle adjustable footplate)	◆
<b>FOOTPLATE TYPE*</b>			
ARJ0559	<input type="checkbox"/>	Angle Adjustable Footplate for Centre Mounted Hanger	◆
<b>CONTROLS*</b>			
ARJ2204	<input type="checkbox"/>	REM216 Drive, Actuators & Lights	◆
<b>CONTROL BRACKETRY*</b>			
ARJ1707	<input type="checkbox"/>	Swing Away RH Bracket	○
ARJ1708	<input type="checkbox"/>	Swing Away LH Bracket	○
<b>SHROUD COLOUR*</b>			
ARJ0228	<input type="checkbox"/>	Wet Black	○
ARJ2601	<input type="checkbox"/>	SpaceStation Silver	○
ARJ2602	<input type="checkbox"/>	Island Blue	○
ARJ2603	<input type="checkbox"/>	Grape JellyBean (Purple)	○
ARJ2607	<input type="checkbox"/>	Rockstar Red	○
ARJ2608	<input type="checkbox"/>	Lights Out Black - Matt	○
ARJ2609	<input type="checkbox"/>	Invacare Blue	○
<b>RIM COLOUR*</b>			
ARJ2611	<input type="checkbox"/>	SpaceStation Silver	○
ARJ2612	<input type="checkbox"/>	Island Blue	○
ARJ2613	<input type="checkbox"/>	Grape JellyBean (Purple)	○
ARJ2617	<input type="checkbox"/>	Rockstar Red	○
ARJ2618	<input type="checkbox"/>	Lights Out Black - Matt	○
ARJ2619	<input type="checkbox"/>	Invacare Blue	○
ARJ2620	<input type="checkbox"/>	Wet Black	○
ARJ2621	<input type="checkbox"/>	Without Rim insert	○
<b>FORKS*</b>			
ARJ2632	<input type="checkbox"/>	Double sided fork	◆
<b>CENTRE WHEELS (DRIVING WHEELS)*</b>			
ARJ0733	<input type="checkbox"/>	355mm x 76mm (14" x 3") Black Tyres with Gel Foam Inserts	◆
<b>MOTOR TYPE*</b>			
ARJ1908	<input type="checkbox"/>	10km/h (6mph) (4 pole) (includes LED light system and indicators)	◆
<b>BATTERY*</b>			
ARJ1260	<input type="checkbox"/>	73.5 Ah (Gel)	◆
<b>CHARGER*</b>			
ARJ1343	<input type="checkbox"/>	OptiCharge Charger 8A (200 - 250V)	◆
<b>TRANSPORTATION KIT*</b>			
ARJ1494	<input type="checkbox"/>	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory, see options below)	◆
<b>FIXATION BELT*</b>			
ARJ2113	<input type="checkbox"/>	Lapbelt (adjustable both sides) adjustable mounting on seat rail extrusion	◆

**REAR VIEW MIRROR\***

ARJ1415	<input type="checkbox"/>	Rear view Mirror RH - Must be selected with 10 km/h Motor	<input type="radio"/>
ARJ1339	<input type="checkbox"/>	Rear view Mirror LH - Must be selected with 10 km/h Motor	<input type="radio"/>

**TOTAL PRICE**

**Non Conforming Options:**

ARJ2132	410mm (16") (410mm - 510mm simple stepless adjustment from the side)
ARJ2134	460mm (18") (410mm - 510mm simple stepless adjustment from the side)
ARJ2138	560mm (22") (510mm - 560mm simple stepless adjustment from the side) (Cushion (option 2142) is only available with Seat Width (option 2124) 530mm)
ARJ2196	Long Armpad (375mm)



For a more complex customer requirement please refer to our Invacare Unique special order form. email [specialsuk@invacare.com](mailto:specialsuk@invacare.com) or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Description of special request (please include relevant option code or part number):