



<b>SEAT DEPTH SETTING SETTABLE IN 3 STEPS*</b>			
ARH0396	<input type="checkbox"/>	Seat Depth 560mm	<input type="radio"/>
ARH0397	<input type="checkbox"/>	Seat Depth 595mm	<input type="radio"/>
ARH0398	<input type="checkbox"/>	Seat Depth 630mm	<input type="radio"/>
<b>SEAT CUSHION*</b>			
ARH1484	<input type="checkbox"/>	MAX Seat Cushion	<input type="radio"/>
<b>ADAPTOR TO REDUCE THE SEAT DEPTH BY 130MM</b>			
ARH0399	<input type="checkbox"/>	Adaptor to shorten the seat depth by 130mm (Required if seat width 490mm is selected in combination with seat depth 595mm or 630mm)	<input type="radio"/>
<b>BACKREST UPHOLSTERY*</b>			
ARH1106	<input type="checkbox"/>	Tension Adjustable	<input checked="" type="checkbox"/>
<b>UPHOLSTERY COLOUR*</b>			
ARH0278	<input type="checkbox"/>	Black multi stretch PU coated fabric	<input checked="" type="checkbox"/>
<b>ARMREST INCLUDING REINFORCED DEPTH ADJUSTABLE BRACKET*</b>			
ARH0635	<input type="checkbox"/>	Fixed Armrest, 250mm - 340mm tool adjustable in height, detachable	<input checked="" type="checkbox"/>
<b>ARMPADS*</b>			
ARH0649	<input type="checkbox"/>	Long armpad (375mm)	<input checked="" type="checkbox"/>
<b>PUSH BAR*</b>			
ARH1307	<input type="checkbox"/>	Push bar	<input checked="" type="checkbox"/>
<b>HEADREST - Adequate headrest option must be selected*</b>			
ARH1432	<input type="checkbox"/>	Headrest: height, depth and sideways adjustable (Azalea type)	<input checked="" type="checkbox"/>
<b>SEAT TILT*</b>			
ARH1095	<input type="checkbox"/>	Electric heavy duty basic Tilt (0° - 20°) select either 470mm or 495mm seat to ground height. (Max user weight 180kg)	<input checked="" type="checkbox"/>
<b>SEAT HEIGHT SETTING*</b>			
ARH1024	<input type="checkbox"/>	Fixed seat height 470mm (18.5")	<input type="radio"/>
ARH1046	<input type="checkbox"/>	Fixed seat height 495mm (19.5")	<input type="radio"/>
<b>BACKREST RECLINER*</b>			
ARH0958	<input type="checkbox"/>	Manually settable recliner by spindle (0° - 25°)	<input checked="" type="checkbox"/>
<b>FOOTREST HANGER TYPE*</b>			
ARH0503	<input type="checkbox"/>	Side Mounted : Vari F : Angle adjustable with tools 0° - 70°, length 290mm - 460mm	<input checked="" type="checkbox"/>
<b>FOOTPLATE TYPE*</b>			
ARH0512	<input type="checkbox"/>	Fixed Footplate	<input checked="" type="checkbox"/>
<b>CONTROLS*</b>			
ARH2203	<input type="checkbox"/>	REM211 Drive & Actuators	<input checked="" type="checkbox"/>
<b>CONTROL BRACKETRY*</b>			
ARH1707	<input type="checkbox"/>	Swing Away RH Bracket	<input type="radio"/>
ARH1708	<input type="checkbox"/>	Swing Away LH Bracket	<input type="radio"/>
<b>SHROUD COLOUR*</b>			
ARH0228	<input type="checkbox"/>	Wet Black	<input type="radio"/>
ARH2601	<input type="checkbox"/>	SpaceStation Silver	<input type="radio"/>
ARH2602	<input type="checkbox"/>	Island Blue	<input type="radio"/>
ARH2603	<input type="checkbox"/>	Grape JellyBean (Purple)	<input type="radio"/>
ARH2607	<input type="checkbox"/>	Rockstar Red	<input type="radio"/>
ARH2608	<input type="checkbox"/>	Lights Out Black - Matt	<input type="radio"/>
ARH2609	<input type="checkbox"/>	Invacare Blue	<input type="radio"/>

<b>RIM COLOUR*</b>		
ARH2611	<input type="checkbox"/>	SpaceStation Silver <span style="float: right;">○</span>
ARH2612	<input type="checkbox"/>	Island Blue <span style="float: right;">○</span>
ARH2613	<input type="checkbox"/>	Grape JellyBean (Purple) <span style="float: right;">○</span>
ARH2617	<input type="checkbox"/>	Rockstar Red <span style="float: right;">○</span>
ARH2618	<input type="checkbox"/>	Lights Out Black - Matt <span style="float: right;">○</span>
ARH2619	<input type="checkbox"/>	Invacare Blue <span style="float: right;">○</span>
ARH2620	<input type="checkbox"/>	Wet Black <span style="float: right;">○</span>
ARH2621	<input type="checkbox"/>	Without Rim insert <span style="float: right;">○</span>
<b>FORKS*</b>		
ARH2632	<input type="checkbox"/>	Double sided fork <span style="float: right;">◆</span>
<b>CENTRE WHEELS (DRIVING WHEELS)*</b>		
ARH0733	<input type="checkbox"/>	355mm x 76mm (14" x 3") Black Tyres with Gel Foam Inserts <span style="float: right;">◆</span>
<b>MOTOR TYPE*</b>		
ARH1282	<input type="checkbox"/>	6 km/h (4mph) (4 pole) <span style="float: right;">◆</span>
<b>BATTERY*</b>		
ARH1260	<input type="checkbox"/>	73.5 Ah <span style="float: right;">◆</span>
<b>CHARGER*</b>		
ARH1343	<input type="checkbox"/>	OptiCharge Charger 8A (200 - 250V) <span style="float: right;">◆</span>
<b>LIGHTS*</b>		
ARH1254	<input type="checkbox"/>	Without light system <span style="float: right;">◆</span>
<b>TRANSPORTATION KIT</b>		
ARH1494	<input type="checkbox"/>	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory, see options below) <span style="float: right;">◆</span>
<b>FIXATION BELT*</b>		
ARH1329	<input type="checkbox"/>	MAX Lapbelt (adjustable both sides) fitted to base of backrest <span style="float: right;">◆</span>

<b>TOTAL PRICE</b>	<input style="width: 100px; height: 20px;" type="text"/>
--------------------	--

**Non Conforming Options:**



For a more complex customer requirement please refer to our Invacare Unique special order form. email [specialsuk@invacare.com](mailto:specialsuk@invacare.com) or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Description of special request (please include relevant option code or part number):