



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Motability Retail Prescription Form

MAY 2024
Customer Ref:

Customer Ref\*, Contact No\*, Delivery Address\*, Account No\*, Order No\*, Order Qty\*, Order date\*

TDX SP2 NB M2 724002 - Weekly Price £50.05

Standard Features

- Max User Weight 150kg (23.5 Stone)
Crash test ISO7176 Part-19
Modulite Seating
Flo tech lite cushion
Flex3 adjustable backrest plates with Laguna backrest cushion
Flip-up armrests, (height, and depth adjustable) with standard armpad
Elan Headrest Standard Pad
Powered Tilt 0° - 20°
Manual centre mounted hanger, with angle adjustable footplates
LiNX control REM211 with swing away bracket
solid Black tyres
6km/h (4mph) 4 pole motors
4 Point tie-down system (Q'Straint QLK-150 compatible)
Lapbelt with autolock (adjustable both sides)

Standard Base price: £6,826



Please make your selection with a cross in the box [ ] Free of charge [O] standard feature [◆] POA special request

Customer weight = KG

Table with columns: Code, Description. Rows include BASE PRODUCT\* (ARM2100), USER WEIGHT\* (ARM1920, ARM1921, ARM1922), and BASE MODEL\* (ARM0000).

<b>SEATING TYPE*</b>			
ARM2100	<input type="checkbox"/>	Modulite Seating Unit	◆
<b>SEAT FRAME*</b>			
ARM2111	<input type="checkbox"/>	Telescopic Seat Frame	◆
<b>SEAT WIDTH*</b>			
ARM2121	<input type="checkbox"/>	380mm (15")	○
ARM2122	<input type="checkbox"/>	430mm (17")	○
ARM2123	<input type="checkbox"/>	480mm (19")	○
ARM2124	<input type="checkbox"/>	530mm (21")	○
<b>SEAT DEPTH RAIL LENGTH*</b>			
ARM2129	<input type="checkbox"/>	Standard 410mm - 510mm	(NC 2138) ○
ARM2190	<input type="checkbox"/>	Extended 510mm - 560mm	(NC 2132, 2134) ○
<b>SEAT DEPTH*</b>			
ARM2132	<input type="checkbox"/>	410mm (16")	(NC 2196) ○
ARM2134	<input type="checkbox"/>	460mm (18")	○
ARM2136	<input type="checkbox"/>	510mm (20")	○
ARM2138	<input type="checkbox"/>	560mm (22")	○
<b>SEAT CUSHION*</b>			
ARM2142	<input type="checkbox"/>	Matrix Flo-Tech Lite	◆
<b>BACKREST SUPPORT*</b>			
ARM2162	<input type="checkbox"/>	Tall 540mm	◆
<b>BACKREST UPHOLSTERY / CUSHION*</b>			
ARM2173	<input type="checkbox"/>	Flex3 Backrest plates with Laguna cushion	◆
<b>UPHOLSTERY TYPE*</b>			
ARM2152	<input type="checkbox"/>	Black Black multi stretch PU coated fabric Upholstery	○
ARM2155	<input type="checkbox"/>	Black Spacer fabric	○
<b>ARMREST*</b>			
ARM2191	<input type="checkbox"/>	T-Armrest standard	◆
<b>ARMPAD*</b>			
ARM2196	<input type="checkbox"/>	Long 375mm Armpad	◆
<b>HEADREST / NECKREST - Adequate headrest option must be selected*</b>			
ARM1350	<input type="checkbox"/>	Neckrest height and depth adjustable	◆
<b>SEAT LIFTER / SEAT TILT*</b>			
ARM1092	<input type="checkbox"/>	Power Tilt (max user weight 136kg)	◆
<b>SEAT HEIGHT SETTING FOR 20° SEAT TILT*</b>			
ARM0354	<input type="checkbox"/>	Seat height 450mm (17.75")	○
ARM0364	<input type="checkbox"/>	Seat height 475mm (18.75")	○
<b>BACKREST RECLINER*</b>			
ARM0975	<input type="checkbox"/>	Manual settable Recliner	◆
<b>FOOTREST HANGER TYPE*</b>			
ARM0556	<input type="checkbox"/>	Centre mounted : Manual	◆
<b>FOOTPLATE TYPE*</b>			
ARM0559	<input type="checkbox"/>	Angle Adjustable Footplate Medium 140mm x 240mm (5.5" x 9.5") for Centre Mounted Hanger only	◆
<b>CONTROLS*</b>			
ARM2203	<input type="checkbox"/>	REM211	◆

<b>CONTROL BRACKETRY*</b>			
ARM1707	<input type="checkbox"/>	Swing away RH Bracket	<input type="radio"/>
ARM1708	<input type="checkbox"/>	Swing away LH Bracket	<input type="radio"/>
<b>SHROUD COLOUR*</b>			
ARM0228	<input type="checkbox"/>	Wet Black	<input type="radio"/>
ARM2601	<input type="checkbox"/>	SpaceStation Silver	<input type="radio"/>
ARM2602	<input type="checkbox"/>	Island Blue	<input type="radio"/>
ARM2603	<input type="checkbox"/>	Grape JellyBean (Purple)	<input type="radio"/>
ARM2607	<input type="checkbox"/>	Rockstar Red	<input type="radio"/>
ARM2608	<input type="checkbox"/>	Lights Out Black - Matt	<input type="radio"/>
ARM2609	<input type="checkbox"/>	Invacare Blue	<input type="radio"/>
<b>RIM COLOUR*</b>			
ARM2611	<input type="checkbox"/>	SpaceStation Silver	<input type="radio"/>
ARM2612	<input type="checkbox"/>	Island Blue	<input type="radio"/>
ARM2613	<input type="checkbox"/>	Grape JellyBean (Purple)	<input type="radio"/>
ARM2617	<input type="checkbox"/>	Rockstar Red	<input type="radio"/>
ARM2618	<input type="checkbox"/>	Lights Out Black - Matt	<input type="radio"/>
ARM2619	<input type="checkbox"/>	Invacare Blue	<input type="radio"/>
ARM2620	<input type="checkbox"/>	Wet Black	<input type="radio"/>
ARM2621	<input type="checkbox"/>	Without rim insert	<input type="radio"/>
<b>FORKS*</b>			
ARM2632	<input type="checkbox"/>	Double sided	<input type="checkbox"/>
<b>CENTRE WHEELS (DRIVING WHEELS)*</b>			
ARM0733	<input type="checkbox"/>	Black tyres with gel foam inserts (solid)	<input type="checkbox"/>
<b>MOTOR TYPE*</b>			
ARM1282	<input type="checkbox"/>	6 km/h (4mph)	<input type="checkbox"/>
<b>BATTERY*</b>			
ARM1257	<input type="checkbox"/>	50 Ah Gel	<input type="checkbox"/>
<b>CHARGER*</b>			
ARM1343	<input type="checkbox"/>	OptiCharge 8A (200 - 250V)	<input type="checkbox"/>
<b>TRANSPORTATION KIT*</b>			
ARM1494	<input type="checkbox"/>	4-Point tie-down	<input type="checkbox"/>
<b>FIXATION BELT*</b>			
ARM1329	<input type="checkbox"/>	Lap Belt fitted to base of backrest	<input type="radio"/>
ARM2113	<input type="checkbox"/>	Lap Belt adjustable mounting on seat rail extrusion	<input type="checkbox"/>

<b>TOTAL PRICE</b>	<input type="text"/>
--------------------	----------------------

**Non Conforming Options:**

ARM2132	410mm (16")
ARM2134	460mm (18")
ARM2138	560mm (22")
ARM2196	Long 375mm Armpad



For a more complex customer requirement please refer to our Invacare Unique special order form. email [specialsuk@invacare.com](mailto:specialsuk@invacare.com) or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Description of special request (please include relevant option code or part number):