



Motability Retail Prescription Form

Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

MAY 2024
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

**TDX SP2 NB M4 1001336 - Weekly Price
£46.97**

Standard Features

- Max User Weight 150kg (23.5 Stone)
- Crash test ISO7176 Part-19
- Modulite Seating
- Matrix Flo-tech cushion
- Flex3 adjustable backrest plates with Laguna backrest cushion
- Standard T-armrests, detachable, Long armpad (325mm)
- Manual settable seat Tilt (0° - 20°)
- Manual settable Recliner (0° - 30° in 3.75° steps)
- 80° hangers with fixed footplate
- LiNX control REM211 with swing away bracket
- Solid Black tyres
- 6km/h (4mph) motors
- 4 Point tie-down brackets (Q'Straint QLK-150 compatible)
- Lap belt with auto lock (flexible mount in side rails)

Standard Base price: £6,121



Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

Code	Description	
BASE PRODUCT*		
ARM2100	<input type="checkbox"/> TDX SP2 NARROW BASE	£6,121
USER WEIGHT*		
ARM1920	<input type="checkbox"/> User weight up to 125kg	<input type="radio"/>
ARM1921	<input type="checkbox"/> User weight 126kg to 136kg	<input type="radio"/>
ARM1922	<input type="checkbox"/> User weight 137kg to 150kg	<input type="radio"/>
BASE MODEL*		
ARM0000	<input type="checkbox"/> TDX SP2 Narrow Base	<input checked="" type="checkbox"/>

SEATING TYPE*			
ARM2100	<input type="checkbox"/>	Modulite Seating Unit	◆
SEAT FRAME*			
ARM2111	<input type="checkbox"/>	Telescopic Seat Frame	◆
SEAT WIDTH*			
ARM2121	<input type="checkbox"/>	380mm (15")	○
ARM2122	<input type="checkbox"/>	430mm (17")	○
ARM2123	<input type="checkbox"/>	480mm (19")	○
ARM2124	<input type="checkbox"/>	530mm (21")	○
SEAT DEPTH RAIL LENGTH*			
ARM2129	<input type="checkbox"/>	Standard 410mm - 510mm	(NC 2138) ○
ARM2190	<input type="checkbox"/>	Extended 510mm - 560mm	(NC 2132, 2134) ○
SEAT DEPTH*			
ARM2132	<input type="checkbox"/>	410mm (16")	(NC 2196) ○
ARM2134	<input type="checkbox"/>	460mm (18")	○
ARM2136	<input type="checkbox"/>	510mm (20")	○
ARM2138	<input type="checkbox"/>	560mm (22")	○
SEAT CUSHION*			
ARM2142	<input type="checkbox"/>	Matrix Flo-Tech Lite	◆
BACKREST SUPPORT*			
ARM2162	<input type="checkbox"/>	Tall 540mm	◆
BACKREST UPHOLSTERY / CUSHION*			
ARM2171	<input type="checkbox"/>	Tension Adjustable Backrest	◆
UPHOLSTERY TYPE*			
ARM2151	<input type="checkbox"/>	Black fabric	◆
ARMREST*			
ARM2191	<input type="checkbox"/>	T-Armrest standard	◆
ARMPAD*			
ARM2196	<input type="checkbox"/>	Long 375mm Armpad	◆
HEADREST / NECKREST - Adequate headrest option must be selected*			
ARM1495	<input type="checkbox"/>	Push bar 22mm	◆
SEAT LIFTER / SEAT TILT*			
ARM1091	<input type="checkbox"/>	Manually settable	◆
SEAT HEIGHT SETTING FOR 20° SEAT TILT*			
ARM0354	<input type="checkbox"/>	Seat height 450mm (17.75")	○
ARM0364	<input type="checkbox"/>	Seat height 475mm (18.75")	○
BACKREST RECLINER*			
ARM0975	<input type="checkbox"/>	Manual settable Recliner	◆
FOOTREST HANGER TYPE*			
ARM0526	<input type="checkbox"/>	Side mounted : Standard 80° hanger	◆
FOOTPLATE TYPE*			
ARM0512	<input type="checkbox"/>	Fixed Footplate	◆
CONTROLS*			
ARM2204	<input type="checkbox"/>	REM216	◆
CONTROL BRACKETRY*			
ARM1707	<input type="checkbox"/>	Swing away RH Bracket	○
ARM1708	<input type="checkbox"/>	Swing away LH Bracket	○

SHROUD COLOUR*			
ARM0228	<input type="checkbox"/>	Wet Black	<input type="radio"/>
ARM2601	<input type="checkbox"/>	SpaceStation Silver	<input type="radio"/>
ARM2602	<input type="checkbox"/>	Island Blue	<input type="radio"/>
ARM2603	<input type="checkbox"/>	Grape JellyBean (Purple)	<input type="radio"/>
ARM2607	<input type="checkbox"/>	Rockstar Red	<input type="radio"/>
ARM2608	<input type="checkbox"/>	Lights Out Black - Matt	<input type="radio"/>
ARM2609	<input type="checkbox"/>	Invacare Blue	<input type="radio"/>
RIM COLOUR*			
ARM2611	<input type="checkbox"/>	SpaceStation Silver	<input type="radio"/>
ARM2612	<input type="checkbox"/>	Island Blue	<input type="radio"/>
ARM2613	<input type="checkbox"/>	Grape JellyBean (Purple)	<input type="radio"/>
ARM2617	<input type="checkbox"/>	Rockstar Red	<input type="radio"/>
ARM2618	<input type="checkbox"/>	Lights Out Black - Matt	<input type="radio"/>
ARM2619	<input type="checkbox"/>	Invacare Blue	<input type="radio"/>
ARM2620	<input type="checkbox"/>	Wet Black	<input type="radio"/>
ARM2621	<input type="checkbox"/>	Without rim insert	<input type="radio"/>
FORKS*			
ARM2632	<input type="checkbox"/>	Double sided	<input type="checkbox"/>
CENTRE WHEELS (DRIVING WHEELS)*			
ARM0733	<input type="checkbox"/>	Black tyres with gel foam inserts (solid)	<input type="checkbox"/>
MOTOR TYPE*			
ARM1908	<input type="checkbox"/>	10 km/h (6mph)	<input type="checkbox"/>
BATTERY*			
ARM1257	<input type="checkbox"/>	50 Ah Gel	<input type="checkbox"/>
CHARGER*			
ARM1343	<input type="checkbox"/>	OptiCharge 8A (200 - 250V)	<input type="checkbox"/>
TRANSPORTATION KIT*			
ARM1494	<input type="checkbox"/>	4-Point tie-down	<input type="checkbox"/>
FIXATION BELT*			
ARM1329	<input type="checkbox"/>	Lap Belt fitted to base of backrest	<input type="radio"/>
ARM2113	<input type="checkbox"/>	Lap Belt adjustable mounting on seat rail extrusion	<input type="checkbox"/>
REAR VIEW MIRROR			
ARM1415	<input type="checkbox"/>	RH Mirror	<input type="radio"/>
ARM1339	<input type="checkbox"/>	LH Mirror	<input type="radio"/>

TOTAL PRICE	<input type="text"/>
--------------------	----------------------

Non Conforming Options:

ARM2132	410mm (16")
ARM2134	460mm (18")
ARM2138	560mm (22")
ARM2196	Long 375mm Armpad



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):