



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Motability Retail Prescription Form

JULY 2024
Customer Ref:

Customer Ref* Account No*
Contact No* Order No*
Delivery Address* Order Qty*
Order date*

Storm4 Modulite M32 1001076 - Weekly
Price £62.35

Standard Features

- Max user weight 150Kg (23.5 stone)
Crash tested ISO 7176 -Part 19
Modulite seating system
Flex3 Laguna Backrest
Standard T-Armrests, detachable with long armpad
Manual settable seat Tilt 0° - 20° Electrical Seat Riser
Manual settable Recliner (0° - 30° in 3.75° steps)
VARI F with Fixed Footplate
LinX Control REM216 with swing away bracket
Pneumatic Castors
Pneumatic Drive Wheels
73.5 Ah Batteries
12 km/h (4mph) Motors
LED Light system with indicators
4 Point tie down system
Lapbelt with autolock (adjustable both sides)
Anti-Tippers



Standard Base price: £8,950

Please make your selection with a cross in the box
Free of charge
standard feature
POA special request

Customer weight = KG

Table with columns: Code, Description, Price. Rows include CHAIR TYPE* (Storm 4 Modulite, Special Storm 4 Modulite), SEATING* (Modulite Seating Unit), and SEATING OPTION* (Telescopic Seat Frame).

SEAT WIDTH*			
ATBXXXX	<input type="checkbox"/>	Seat width 330mm (13") seat width achieved on a 380mm seat frame (2121) with inrigged armrests (2183)	<input type="radio"/>
ATB2121	<input type="checkbox"/>	Set at 380mm (adjustable from 380mm to 430mm (15" - 17") with stepless adjustment on armrests, legrests and backrest)	<input type="radio"/>
ATB2122	<input type="checkbox"/>	Set at 430mm (adjustable from 430mm to 480mm (17" - 19) with stepless adjustment on armrests, legrests and backrest)	<input type="radio"/>
ATB2123	<input type="checkbox"/>	Set at 480mm (adjustable from 480mm to 530mm (19" - 21") with stepless adjustment on armrests, legrests and backrest)	<input type="radio"/>
ATB2124	<input type="checkbox"/>	Set at 530mm (adjustable from 530mm to 580mm (21" - 23") with stepless adjustment on armrests, legrests and backrest)	<input type="radio"/>
SEAT DEPTH*			
ATB2132	<input type="checkbox"/>	410mm (16") (adjustable 410mm - 510mm stepless adjustment) (NC 2196)	<input type="radio"/>
ATB2134	<input type="checkbox"/>	460mm (18") (adjustable 410mm - 510mm stepless adjustment)	<input type="radio"/>
ATB2136	<input type="checkbox"/>	510mm (20") (adjustable 410mm - 510mm stepless adjustment)	<input type="radio"/>
SEAT CUSHION*			
ATB2142	<input type="checkbox"/>	Matrx Flo-tech Lite	<input checked="" type="checkbox"/>
BACKREST SUPPORT*			
ATB2162	<input type="checkbox"/>	Backposts 540mm (Tall) must be selected for Flex3 Laguna backrest	<input checked="" type="checkbox"/>
PUSH HANDLES FOR BACKPOSTS*			
ATB2107	<input type="checkbox"/>	Backposts WITH push handles	<input checked="" type="checkbox"/>
BACKREST UPHOLSTERY / CUSHION*			
ATB2173	<input type="checkbox"/>	Flex3 adjustable Backrest Plates with Laguna Backrest cushion	<input checked="" type="checkbox"/>
UPHOLSTERY TYPE*			
ATB2152	<input type="checkbox"/>	Black multi stretch PU coated fabric	<input type="radio"/>
ATB2155	<input type="checkbox"/>	Black Spacer fabric	<input type="radio"/>
ARMREST (PAIR)*			
ATB2191	<input type="checkbox"/>	T-Armrest Standard, Range of adjustment from seat plate 245mm to 315mm	<input checked="" type="checkbox"/>
ARMPAD (PAIR)*			
ATB2196	<input type="checkbox"/>	Long Armpad (375mm) (NC 2132)	<input checked="" type="checkbox"/>
HEADREST / NECKREST - Adequate headrest option must be selected*			
ATB1350	<input type="checkbox"/>	neckrest: height and depth adjustable	<input checked="" type="checkbox"/>
SEAT HEIGHT AND SEAT RISER*			
ATB1023	<input type="checkbox"/>	Fixed Seat Height 450mm (17.75")	<input type="radio"/>
ATB1025	<input type="checkbox"/>	Fixed Seat Height 480mm (19")	<input type="radio"/>
ATB1030	<input type="checkbox"/>	Electrical Seat Riser (400mm - 650mm) (16" - 25.5")	<input type="radio"/>
ATB1040	<input type="checkbox"/>	Electrical Seat Riser (440mm - 690mm) incl. adapter (17.25" - 27")	<input type="radio"/>
SEAT TILT*			
ATB1002	<input type="checkbox"/>	Manual Seat Tilt by spindle	<input checked="" type="checkbox"/>
BACKREST RECLINER*			
ATB0975	<input type="checkbox"/>	Manual Settable Recliner (0° - 30°, in 3.75° steps) Foldable	<input checked="" type="checkbox"/>
FOOTREST HANGER TYPE (PAIR)*			
ATB0503	<input type="checkbox"/>	Vari F : Angle adjustable with tools 0° - 70°, length 290mm - 460mm	<input checked="" type="checkbox"/>

FOOTPLATE TYPE*			
ATB0512	<input type="checkbox"/>	Fixed Footplate	○
CONTROLS*			
ATB2204	<input type="checkbox"/>	REM216 Drive, Actuators & Lights	◆
CONTROL BRACKETRY*			
ATB1707	<input type="checkbox"/>	Swing Away RH Bracket	○
ATB1708	<input type="checkbox"/>	Swing Away LH Bracket	○
SHROUD COLOUR*			
ATB0202	<input type="checkbox"/>	Blue	○
ATB0228	<input type="checkbox"/>	Black	○
ATB0220	<input type="checkbox"/>	Champagne	○
ATB0204	<input type="checkbox"/>	Transparent Cherry Red	○
FRONT WHEEL*			
ATB0851	<input type="checkbox"/>	3.00 - 6 Grey Pneumatic	◆
REAR WHEEL*			
ATB0771	<input type="checkbox"/>	3.00 - 8 Grey Pneumatic	◆
MOTOR TYPE*			
ATB1909	<input type="checkbox"/>	12 km/h (8mph) (includes LED light system with indicators)	◆
BATTERY*			
ATB1364	<input type="checkbox"/>	70 Ah Gel	◆
SUSPENSION*			
ATB1240	<input type="checkbox"/>	Suspension low (up to 80kg user weight)	○
ATB1241	<input type="checkbox"/>	Suspension middle (from 70kg up to 100kg user weight)	○
ATB1242	<input type="checkbox"/>	Suspension hard (from 90kg up to 150kg user weight)	○
CHARGER*			
ATB1343	<input type="checkbox"/>	OptiCharge Charger 8A (200 - 250V)	◆
ANTI TIPPER*			
ATB1205	<input type="checkbox"/>	Anti-tipper (standard version)	○
TRANSPORTATION KIT*			
ATB1494	<input type="checkbox"/>	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory, see options below)	◆
FIXATION BELT*			
ATB2113	<input type="checkbox"/>	Lap belt with auto lock (flexible mount in side rails)	◆
REAR VIEW MIRROR			
ATB1415	<input type="checkbox"/>	Rear view Mirror RH - Must be selected with 10 km/h Motor	○
ATB1339	<input type="checkbox"/>	Rear view Mirror LH - Must be selected with 10 km/h Motor	○

TOTAL PRICE	
--------------------	--

Non Conforming Options:

ATB2132 410mm (16") (adjustable 410mm - 510mm stepless adjustment)

ATB2196 Long Armpad (375mm)



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):