



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

JULY 2024
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

TDX SP2 Max M15 1001086 - Weekly Price £48.21

Standard Features

- Max user weight 180kg (28 stone)
- Crash tested ISO 7176 -Part 19
- Tension Adjustable backrest
- Armrest with reinforced depth adjustable bracket with Short armpad
- Powered 20 degree seat Tilt
- Manual settable backrest Recline
- VARI F hangers with fixed footplate
- LiNX Control REM211 with swing away bracket
- Solid Black Tyres
- 6 km/h (4mph) 4 Pole Motors
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)

Standard Base price: £5,882



Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

Code	Description	
CHAIR TYPE*		
ARH0000	TDX SP2 Max	£5,882
ARZ0000	Special TDX SP2 Max (Please provide details in special section)	POA
SEATING TYPE*		
ARH1145	MAX Seating system	◆
SEAT WIDTH SETTING*		
ARH0378	Seat Width of Seat Frame 550mm (+100mm adjustment in arms)	○
ARH0380	Seat Width of Seat Frame 610mm (+100mm adjustment in arms)	○

SEAT DEPTH SETTING SETTABLE IN 3 STEPS*			
ARH0396	<input type="checkbox"/>	Seat Depth 560mm	<input type="radio"/>
ARH0397	<input type="checkbox"/>	Seat Depth 595mm	<input type="radio"/>
ARH0398	<input type="checkbox"/>	Seat Depth 630mm	<input type="radio"/>
SEAT CUSHION*			
ARH1484	<input type="checkbox"/>	MAX Seat Cushion	<input type="radio"/>
ADAPTOR TO REDUCE THE SEAT DEPTH BY 130MM			
ARH0399	<input type="checkbox"/>	Adaptor to shorten the seat depth by 130mm (Required if seat width 490mm is selected in combination with seat depth 595mm or 630mm)	<input type="radio"/>
BACKREST UPHOLSTERY*			
ARH1106	<input type="checkbox"/>	Tension Adjustable	<input checked="" type="checkbox"/>
UPHOLSTERY COLOUR*			
ARH0278	<input type="checkbox"/>	Black multi stretch PU coated fabric	<input checked="" type="checkbox"/>
ARMREST INCLUDING REINFORCED DEPTH ADJUSTABLE BRACKET*			
ARH0635	<input type="checkbox"/>	Fixed Armrest, 250mm - 340mm tool adjustable in height, detachable	<input checked="" type="checkbox"/>
ARMPADS*			
ARH0649	<input type="checkbox"/>	Long armpad (375mm)	<input checked="" type="checkbox"/>
PUSH BAR*			
ARH1307	<input type="checkbox"/>	Push bar	<input checked="" type="checkbox"/>
HEADREST - Adequate headrest option must be selected*			
ARH1432	<input type="checkbox"/>	Headrest: height, depth and sideways adjustable (Azalea type)	<input checked="" type="checkbox"/>
SEAT TILT*			
ARH1095	<input type="checkbox"/>	Electric heavy duty basic Tilt (0° - 20°) select either 470mm or 495mm seat to ground height. (Max user weight 180kg)	<input checked="" type="checkbox"/>
SEAT HEIGHT SETTING*			
ARH1024	<input type="checkbox"/>	Fixed seat height 470mm (18.5")	<input type="radio"/>
ARH1046	<input type="checkbox"/>	Fixed seat height 495mm (19.5")	<input type="radio"/>
BACKREST RECLINER*			
ARH0958	<input type="checkbox"/>	Manually settable recliner by spindle (0° - 25°)	<input checked="" type="checkbox"/>
FOOTREST HANGER TYPE*			
ARH0503	<input type="checkbox"/>	Side Mounted : Vari F : Angle adjustable with tools 0° - 70°, length 290mm - 460mm	<input checked="" type="checkbox"/>
FOOTPLATE TYPE*			
ARH0512	<input type="checkbox"/>	Fixed Footplate	<input checked="" type="checkbox"/>
CONTROLS*			
ARH2203	<input type="checkbox"/>	REM211 Drive & Actuators	<input checked="" type="checkbox"/>
CONTROL BRACKETRY*			
ARH1707	<input type="checkbox"/>	Swing Away RH Bracket	<input type="radio"/>
ARH1708	<input type="checkbox"/>	Swing Away LH Bracket	<input type="radio"/>
SHROUD COLOUR*			
ARH0228	<input type="checkbox"/>	Wet Black	<input type="radio"/>
ARH2601	<input type="checkbox"/>	SpaceStation Silver	<input type="radio"/>
ARH2602	<input type="checkbox"/>	Island Blue	<input type="radio"/>
ARH2603	<input type="checkbox"/>	Grape JellyBean (Purple)	<input type="radio"/>
ARH2607	<input type="checkbox"/>	Rockstar Red	<input type="radio"/>
ARH2608	<input type="checkbox"/>	Lights Out Black - Matt	<input type="radio"/>
ARH2609	<input type="checkbox"/>	Invacare Blue	<input type="radio"/>

RIM COLOUR*			
ARH2611	<input type="checkbox"/>	SpaceStation Silver	<input type="radio"/>
ARH2612	<input type="checkbox"/>	Island Blue	<input type="radio"/>
ARH2613	<input type="checkbox"/>	Grape JellyBean (Purple)	<input type="radio"/>
ARH2617	<input type="checkbox"/>	Rockstar Red	<input type="radio"/>
ARH2618	<input type="checkbox"/>	Lights Out Black - Matt	<input type="radio"/>
ARH2619	<input type="checkbox"/>	Invacare Blue	<input type="radio"/>
ARH2620	<input type="checkbox"/>	Wet Black	<input type="radio"/>
ARH2621	<input type="checkbox"/>	Without Rim insert	<input type="radio"/>
FORKS*			
ARH2632	<input type="checkbox"/>	Double sided fork	<input type="checkbox"/>
CENTRE WHEELS (DRIVING WHEELS)*			
ARH0733	<input type="checkbox"/>	355mm x 76mm (14" x 3") Black Tyres with Gel Foam Inserts	<input type="checkbox"/>
MOTOR TYPE*			
ARH1282	<input type="checkbox"/>	6 km/h (4mph) (4 pole)	<input type="checkbox"/>
BATTERY*			
ARH1260	<input type="checkbox"/>	73.5 Ah	<input type="checkbox"/>
CHARGER*			
ARH1343	<input type="checkbox"/>	OptiCharge Charger 8A (200 - 250V)	<input type="checkbox"/>
LIGHTS*			
ARH1254	<input type="checkbox"/>	Without light system	<input type="checkbox"/>
TRANSPORTATION KIT			
ARH1494	<input type="checkbox"/>	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory, see options below)	<input type="checkbox"/>
FIXATION BELT*			
ARH1329	<input type="checkbox"/>	MAX Lapbelt (adjustable both sides) fitted to base of backrest	<input type="checkbox"/>

TOTAL PRICE	<input type="text"/>
--------------------	----------------------

Non Conforming Options:



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):