



Motability Retail Prescription Form

Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

JULY 2024
Customer Ref:

| | |
|--------------------------|--------------------|
| Customer Ref* | Account No* |
| Contact No* | Order No* |
| Delivery Address* | Order Qty* |
| | Order date* |

**TDX SP2 HD M16 1001587 - Weekly Price
£50.90**

Standard Features

- Max user weight 180kg (28 stone)
- Crash tested ISO 7176 -Part 19
- Modulite Seating System
- Matrix Flo-tech Cushion
- Flex 3 adjustable Backrest plates with Laguna backrest cushion
- Standard T-armrests, detachable with long armpad
- Neckrest
- Powereed tilt 20°
- Manual settable backrest Recline 0° - 30° in 3.75° steps
- 80° hangers with fixed footplate
- LiNX control REM211 swing away bracket
- Solid Black tyres
- 6km/h (4mph) 4 pole motors
- 4-Point tie-down system (Q'Straint QLK-150 compatible)
- Lapbelt with autolock (adjustable both sides)



Standard Base price: £6,517

Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

| Code | Description | |
|---------------------|-------------------------------------|-------------------------------------|
| CHAIR TYPE* | | |
| ARJ0000 | TDX SP2 Heavy Duty | <input checked="" type="checkbox"/> |
| User Weight* | | |
| ARJ1918 | User weight less 100 kg | <input type="checkbox"/> |
| ARJ1919 | User weight 101 to 125 kg | <input type="checkbox"/> |
| ARJ1921 | User weight 126 to 136 kg | <input type="checkbox"/> |
| ARJ1922 | User weight 137 to 150 kg (NC 1092) | <input type="checkbox"/> |
| ARJ1923 | User weight 151 to 160 kg (NC 1092) | <input type="checkbox"/> |
| ARJ1924 | User weight 161 to 180 kg (NC 1092) | <input type="checkbox"/> |

| BASE TYPE* | | | |
|---|--------------------------|--|-------------------------------------|
| ARJ2100 | <input type="checkbox"/> | Modulite Seating unit | £6,517 |
| SEAT WIDTH* | | | |
| ARJ2121 | <input type="checkbox"/> | Set at 380mm (15") (adjustable from 380mm to 430mm (15" - 17") with stepless adjustment on armrests, legrests and backrest) (NC 2138 + 2142) | <input type="radio"/> |
| ARJ2122 | <input type="checkbox"/> | Set at 430mm (17") (adjustable from 430mm to 480mm (17" - 19) with stepless adjustment on armrests, legrests and backrest) (NC 2138 + 2142) | <input type="radio"/> |
| ARJ2123 | <input type="checkbox"/> | Set at 480mm (19") (adjustable from 480mm to 530mm (19" - 21") with stepless adjustment on armrests, legrests and backrest) (NC 2138 + 2142) | <input type="radio"/> |
| ARJ2124 | <input type="checkbox"/> | Set at 530mm (21") (adjustable from 530mm to 580mm (21" - 23") with stepless adjustment on armrests, legrests and backrest) | <input type="radio"/> |
| SEAT DEPTH RAIL LENGTH* | | | |
| ARJ2129 | <input type="checkbox"/> | Standard Depth Package (16" - 20", 410mm - 510mm simple stepless adjustment from the side) (NC 2138, 2196) | <input type="radio"/> |
| ARJ2130 | <input type="checkbox"/> | Extended Depth Package (20" - 22", 510mm - 560mm simple stepless adjustment from the side) (NC 2132, 2134) | <input type="radio"/> |
| SEAT DEPTH SETTING* | | | |
| ARJ2132 | <input type="checkbox"/> | 410mm (16") (410mm - 510mm simple stepless adjustment from the side) (NC 2196) | <input type="radio"/> |
| ARJ2134 | <input type="checkbox"/> | 460mm (18") (410mm - 510mm simple stepless adjustment from the side) | <input type="radio"/> |
| ARJ2136 | <input type="checkbox"/> | 510mm (20") (410mm - 510mm simple stepless adjustment from the side) | <input type="radio"/> |
| ARJ2138 | <input type="checkbox"/> | 560mm (22") (510mm - 560mm simple stepless adjustment from the side) (Cushion (option 2142) is only available with Seat Width (option 2124) 530mm) | <input type="radio"/> |
| SEAT CUSHION* | | | |
| ARJ2142 | <input type="checkbox"/> | Matrx Flo-tech Lite | <input checked="" type="checkbox"/> |
| BACKREST SUPPORT* | | | |
| ARJ2162 | <input type="checkbox"/> | Backposts 540mm (Tall) must be selected for Flex3 adjustable backrest plates | <input checked="" type="checkbox"/> |
| PUSH HANDLES FOR BACKPOSTS* | | | |
| ARJ2107 | <input type="checkbox"/> | Backposts WITH push handles | <input checked="" type="checkbox"/> |
| BACKREST UPHOLSTERY / CUSHION* | | | |
| ARJ2173 | <input type="checkbox"/> | Flex3 adjustable Backrest Plates with Laguna Backrest cushion | <input checked="" type="checkbox"/> |
| UPHOLSTERY COLOUR* | | | |
| ARJ2152 | <input type="checkbox"/> | Black Black multi stretch PU coated fabric Upholstery | <input type="radio"/> |
| ARJ2155 | <input type="checkbox"/> | Black Spacer fabric | <input type="radio"/> |
| ARMREST* | | | |
| ARJ2191 | <input type="checkbox"/> | T-Armrest Standard, Range of adjustment from seat plate 245mm to 315mm | <input checked="" type="checkbox"/> |
| ARMPAD* | | | |
| ARJ2196 | <input type="checkbox"/> | Long Armpad (375mm) (NC 2132) | <input checked="" type="checkbox"/> |
| HEADREST / NECKREST - Adequate headrest option must be selected* | | | |
| ARJ1350 | <input type="checkbox"/> | Neckrest height and depth adjustable | <input checked="" type="checkbox"/> |
| SEAT LIFTER / SEAT TILT* | | | |
| ARJ1092 | <input type="checkbox"/> | Electric Tilt (0° - 20°) Seat Height 450mm | <input checked="" type="checkbox"/> |
| SEAT HEIGHT SETTING FOR 20 DEGREES SEAT TILT | | | |
| ARJ0354 | <input type="checkbox"/> | Seat height set at 450mm (17.75") (from floor to front of seat plate at 0° tilt) | <input type="radio"/> |
| ARJ0364 | <input type="checkbox"/> | Seat height set at 475mm (18.75") (from floor to front of seat plate at 0° tilt) | <input type="radio"/> |

| | | | |
|--|--------------------------|---|---|
| BACKREST RECLINER* | | | |
| ARJ0975 | <input type="checkbox"/> | Manual Settable Recliner (0° - 30°, in 3.75° steps) Foldable | ◆ |
| FOOTREST HANGER TYPE* | | | |
| ARJ0526 | <input type="checkbox"/> | Side Mounted : Standard 80° hangers | ◆ |
| FOOTPLATE TYPE* | | | |
| ARJ0512 | <input type="checkbox"/> | Fixed Footplate | ◆ |
| CONTROLS* | | | |
| ARJ2203 | <input type="checkbox"/> | REM211 Drive & Actuators | ◆ |
| CONTROL BRACKETRY* | | | |
| ARJ1707 | <input type="checkbox"/> | Swing Away RH Bracket | ○ |
| ARJ1708 | <input type="checkbox"/> | Swing Away LH Bracket | ○ |
| SHROUD COLOUR* | | | |
| ARJ0228 | <input type="checkbox"/> | Wet Black | ○ |
| ARJ2601 | <input type="checkbox"/> | SpaceStation Silver | ○ |
| ARJ2602 | <input type="checkbox"/> | Island Blue | ○ |
| ARJ2603 | <input type="checkbox"/> | Grape JellyBean (Purple) | ○ |
| ARJ2607 | <input type="checkbox"/> | Rockstar Red | ○ |
| ARJ2608 | <input type="checkbox"/> | Lights Out Black - Matt | ○ |
| ARJ2609 | <input type="checkbox"/> | Invacare Blue | ○ |
| RIM COLOUR* | | | |
| ARJ2611 | <input type="checkbox"/> | SpaceStation Silver | ○ |
| ARJ2612 | <input type="checkbox"/> | Island Blue | ○ |
| ARJ2613 | <input type="checkbox"/> | Grape JellyBean (Purple) | ○ |
| ARJ2617 | <input type="checkbox"/> | Rockstar Red | ○ |
| ARJ2618 | <input type="checkbox"/> | Lights Out Black - Matt | ○ |
| ARJ2619 | <input type="checkbox"/> | Invacare Blue | ○ |
| ARJ2620 | <input type="checkbox"/> | Wet Black | ○ |
| ARJ2621 | <input type="checkbox"/> | Without Rim insert | ○ |
| FORKS* | | | |
| ARJ2632 | <input type="checkbox"/> | Double sided fork | ◆ |
| CENTRE WHEELS (DRIVING WHEELS)* | | | |
| ARJ0733 | <input type="checkbox"/> | 355mm x 76mm (14" x 3") Black Tyres with Gel Foam Inserts | ◆ |
| MOTOR TYPE* | | | |
| ARJ1282 | <input type="checkbox"/> | 6 km/h (4mph) (4 pole) | ◆ |
| BATTERY* | | | |
| ARJ1260 | <input type="checkbox"/> | 73.5 Ah (Gel) | ◆ |
| CHARGER* | | | |
| ARJ1343 | <input type="checkbox"/> | OptiCharge Charger 8A (200 - 250V) | ◆ |
| TRANSPORTATION KIT* | | | |
| ARJ1494 | <input type="checkbox"/> | 4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory, see options below) | ◆ |
| FIXATION BELT* | | | |
| ARJ2113 | <input type="checkbox"/> | Lapbelt (adjustable both sides) adjustable mounting on seat rail extrusion | ◆ |

| | |
|--------------------|--|
| TOTAL PRICE | |
|--------------------|--|

Non Conforming Options:

| | |
|---------|--|
| ARJ2132 | 410mm (16") (410mm - 510mm simple stepless adjustment from the side) |
| ARJ2134 | 460mm (18") (410mm - 510mm simple stepless adjustment from the side) |
| ARJ2138 | 560mm (22") (510mm - 560mm simple stepless adjustment from the side) (Cushion (option 2142) is only available with Seat Width (option 2124) 530mm) |
| ARJ2196 | Long Armpad (375mm) |
| ARJ1092 | Electric Tilt (0°- 20°) Seat Height 450mm |



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):