



# Motability Retail Prescription Form

**Yes, you can:**

Tel: 01656 776222 Fax: 01656 776220  
Email: ordersuk@invacare.com  
Online Spares available at: www.invacare.co.uk

JULY 2024  
Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## Esprit Action Junior M120 706812 - Weekly Price £38.59

### Standard Features

- Max User Weight : 80kg (12.6 Stone)
- Single crossbar (Width Adjustable)
- Seat widths 305mm to 380mm 12" - 15"
- Seat depth 325mm (13") to 375mm (15") adjustable (set at 350mm)
- Seat height 460mm (18")
- Nylon seat upholstery
- FIXED 7° BACKREST 430mm backrest
- Standard nylon sling type backrest upholstery
- 80° swing in / swing out legrest hangers
- Flip up angle and depth adjustable footplates
- Height adjustable, removable arm armrests
- Choice of 3 standard colours
- Anti tippers With jack up function
- Buckle belt



Please make your selection with a cross in the box     Free of charge     standard feature    POA special request

Customer weight =                      KG

Code	Description	
<b>BASE PRODUCT*</b>		
SPD0050	<input type="checkbox"/> ESPRIT ACTION JUNIOR (SINGLE CROSSBAR) (NC 2300)	£3,982
<b>CHAIR TYPE*</b>		
SPD0100	<input type="checkbox"/> BASE MODEL ESPRIT ACTION	<input checked="" type="checkbox"/>
<b>CHASSIS*</b>		
SPD0200	<input type="checkbox"/> SINGLE CROSSBAR (SW 305 MM TO 380 MM / SD 325 MM TO 375 MM) USER MAX WEIGHT = 80Kg	<input checked="" type="checkbox"/>
<b>BRACKET TYPE*</b>		
SPD0300	<input type="checkbox"/> STANDARD BRACKET	<input type="radio"/>
SPD0310	<input type="checkbox"/> ALBER E35 SWIVEL ARM	£360

<b>USER REMOTE*</b>			
SPD0400		USER REMOTE CONTROL RH	○
SPD0410		USER REMOTE CONTROL LH	○
<b>SPEED*</b>			
SPD0500		6 KM/H	◆
<b>BATTERY*</b>			
SPD0600		BATTERY 7,65AH/36V (NON SPILLABLE LEAD ACID BATTERY)	◆
<b>POWER OPTIONS</b>			
SPD0610		ADDITIONAL BATTERY 7,65AH/36V	£271
<b>SEAT WIDTH*</b>			
SPD1000		305 MM SEAT WIDTH	○
SPD1010		330 MM SEAT WIDTH	○
SPD1020		355 MM SEAT WIDTH	○
SPD1030		380 MM SEAT WIDTH	○
<b>SEAT DEPTH*</b>			
SPD1200		325 TO 375 MM ADJUSTABLE SEAT DEPTH SET AT 350 MM	◆
<b>SEAT HEIGHT*</b>			
SPD1300		460 MM SEAT HEIGHT	○
SPD1310		485 MM SEAT HEIGHT	○
<b>SEAT UPHOLSTERY TYPE*</b>			
SPD1400		NYLON WITH LOOP BANDS STRIPS SEAT UPHOLSTERY, WIDTH & DEPTH ADJUSTABLE (SINGLE CROSSBAR ONLY)	◆
<b>BACKREST TYPE*</b>			
SPD2000		FIXED 7° BACKREST 2 HEIGHTS 400/430 MM (DELIVERED 430 MM)	◆
SPD2020		FOLDING BACKREST 430 MM HEIGHT	£32
SPD2040		ANGLE ADJUSTABLE AND FOLD DOWN (-15° TO + 15°) FIXED 7° HEIGHT 425-450 MM / 17"-18" SET AT 425 MM / 17" (2 POSITIONS; 25 MM / 1" STEPS THROUGH UPHOLSTERY)	£97
SPD2050		RECLINING/FOLDING (GAS SPRING) 0/30° BACKREST 3 HEIGHTS 390/410/430 MM (DELIVERED 430 MM)	£64
SPD2060		RECLINING (RATCHET) 0/30° BACKREST 3 HEIGHTS 390/410/430 MM (DELIVERED 430 MM)	£124
SPD2080		RECLINING (GAS SPRING) 0/30° BACKREST 3 HEIGHTS 390/410/430 MM (DELIVERED 430 MM)	£51
<b>BACKREST UPHOLSTERY TYPE*</b>			
SPD2300		LUMBAR SUPPORT (+1 STRAP) BI MATERIAL PADDED COVER	◆
SPD2320		CONTOUR ADJUSTABLE (FULL STRAPS) BI MATERIAL COVER	○
<b>HANGERS RH*</b>			
SPD3000		80° SHORT, UL FROM 200 TO 330MM	◆
<b>HANGERS LH*</b>			
SPD3100		80° SHORT, UL FROM 200 TO 330MM	◆
<b>FOOTPLATES RH*</b>			
SPD3200		FLIP UP RH FOOTPLATE INCLUDING HEEL STRAP	○
SPD3210		FLIP UP ANGLE AND DEPTH ADJUSTABLE RH FOOTPLATE, INCLUDING HEEL STRAP	£42
SPD3220		FLIP UP ONE PIECE FOOTPLATE, DEPTH & ANGLE ADJUSTABLE (NC 3310)	£76

FOOTPLATES LH*			
SPD3300	<input type="checkbox"/>	FLIP UP LH FOOTPLATE, INCLUDING HEEL STRAP	○
SPD3310	<input type="checkbox"/>	FLIP UP ANGLE AND DEPTH ADJUSTABLE LH FOOTPLATE, INCLUDING HEEL STRAP	£42
ARMRESTS RH*			
SPD4000	<input type="checkbox"/>	REMOVABLE RH DESK, ADJ IN HEIGHT (235-280 for FSTF 460 & 210-280 for FSTF 485MM) DEPTH (+35 MM FRONTWARD) & WIDTH (+25MM OUTWARDS)	○
SPD4010	<input type="checkbox"/>	REMOVABLE RH FULL LENGTH, ADJ IN HEIGHT (210 TO 280 MM) DEPTH (+35 MM BACKWARD) & WIDTH (+25MM OUTWARDS)	◆
ARMRESTS LH*			
SPD4100	<input type="checkbox"/>	REMOVABLE LH DESK, ADJ IN HEIGHT (235-280 for FSTF 460 & 210-280 for FSTF 485MM) DEPTH (+35 MM FRONTWARD) & WIDTH (+25MM OUTWARDS)	○
SPD4110	<input type="checkbox"/>	REMOVABLE LH FULL LENGTH, ADJ IN HEIGHT (210 TO 280 MM) DEPTH (35 MM BACKWARD) & WIDTH (+25MM OUTWARDS)	◆
SEAT CUSHION*			
SPD5000	<input type="checkbox"/>	ESSENTIAL CUSHION 375 MM DEPTH	○
SEATING OPTIONS			
SPD5110	<input type="checkbox"/>	HEADREST BLACK PU TR26	£153
SPD5120	<input type="checkbox"/>	HEADREST WITH CHEEK SUPPORT BLACK PU TR26 (INCLUDING BACKBRACE EXCEPT FOR LAGUNA)	£177
SPD5140	<input type="checkbox"/>	LATERAL POSITIONING PADS (PAIR)	£59
SPD5150	<input type="checkbox"/>	ADJUSTABLE CALF STRAP	○
BACK BRACE			
SPD5160	<input type="checkbox"/>	BACK BRACE WIDTH ADJUSTABLE	○
COLOURS*			
SPD6000	<input type="checkbox"/>	ANTHRACITE BLACK	○
SPD6010	<input type="checkbox"/>	DARK PURPLE	○
SPD6020	<input type="checkbox"/>	WATER BLUE	○
ACCENT COLOUR*			
SPD6100	<input type="checkbox"/>	BLACK	○
SPD6120	<input type="checkbox"/>	LIME	○
FRONT WHEELS*			
SPD7010	<input type="checkbox"/>	200 X 50 MM (8" X 2") BLACK PNEUMATIC TYRE	◆
REAR WHEELS*			
SPD7110	<input type="checkbox"/>	355 X 70 MM (14" X 27") BLACK PNEUMATIC TYRE	◆
ANTI TIPPER*			
SPD8010	<input type="checkbox"/>	REAR ANTI-TIPPERS REMOVABLE WITH CRUTCHES (E25), PAIR	◆
LAP BELT*			
SPD8020	<input type="checkbox"/>	POSTURE BELT	◆
TRANSPORTATION BRACKET*			
SPD8030	<input type="checkbox"/>	TRANSPORTATION BRACKET	◆
<b>TOTAL PRICE</b>			<input type="text"/>

**Non Conforming Options:**

SPD2300 LUMBAR SUPPORT (+1 STRAP) BI MATERIAL PADDED COVER

SPD3310 FLIP UP ANGLE AND DEPTH ADJUSTABLE LH FOOTPLATE, INCLUDING HEEL STRAP



For a more complex customer requirement please refer to our Invacare Unique special order form. email [specialsuk@invacare.com](mailto:specialsuk@invacare.com) or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Description of special request (please include relevant option code or part number):