



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Motability Retail Prescription Form

JULY 2024
Customer Ref:

Customer Ref\* Account No\*
Contact No\* Order No\*
Delivery Address\* Order Qty\*
Order date\*

AVIVA RX20 Standard Seat M20 707043 - Weekly Price £39.41

Standard Features

- Max user weight 136kg
Crash tested ISO 7176 -Part 19
One piece seat plate
Flo-tech lite cushion
Tension adjustable backrest
Neckrest
Standard T-Armrests, detachable with Long armpad
Powered "fixed pivot" Tilt 20°
Manual settable backrest Recline (0°- 30°, in 3.75°steps)
80° hangers with fixed footplate
Remote LiNX DLX-REM211 with swing away Bracket
8" Castor, puncture proof solid tyre
Solid black tyres
6 km/h 4mph motors
60Ah batteries
Anti tippers
4 Point tie down system
Lapbelt with autolock

Standard Base price: £4,280



Please make your selection with a cross in the box [ ] Free of charge [O] standard feature [◆] POA special request

Customer weight = KG

Table with columns: Code, Description, Price. Rows include BASE MODEL\* (JBB0000), CHASSIS WIDTH\* (JBL2810), SEATING TYPE\* (JBB2105), SEAT FRAME\* (JBB2110).

<b>SEAT WIDTH*</b>			
JBB2121	<input type="checkbox"/>	380 mm Seat Plate	<input type="radio"/>
JBB2122	<input type="checkbox"/>	430 mm Seat Plate	<input type="radio"/>
<b>SEAT DEPTH*</b>			
JBB2132	<input type="checkbox"/>	Seat Depth 410 mm	<input type="radio"/>
JBB2134	<input type="checkbox"/>	Seat Depth 460 mm	<input type="radio"/>
JBB2136	<input type="checkbox"/>	Seat Depth 510 mm	<input type="radio"/>
<b>SEAT CUSHION*</b>			
JBB2142	<input type="checkbox"/>	Matrx Flo-tech Cushion Flo-tech Lite	<input checked="" type="radio"/>
<b>BACKREST SUPPORT*</b>			
JBB2161	<input type="checkbox"/>	Backrest frame, detachable, height 480 mm For Tension Adjustable Upholstery/Matrx Backs	<input checked="" type="radio"/>
<b>BACKREST UPHOLSTERY / CUSHION*</b>			
JBB2171	<input type="checkbox"/>	Tension Adjustable Upholstery nylon overlay short = 480 mm / 7 straps measured from seat plate to top of upholstery 480 mm	<input checked="" type="radio"/>
<b>UPHOLSTERY TYPE*</b>			
JBB2151	<input type="checkbox"/>	Black Fabric Upholstery (Nylon w 2141; 2-way-stretch material w 2142, 2143)	<input checked="" type="radio"/>
<b>ARMREST (Pair)*</b>			
JBB2191	<input type="checkbox"/>	T-Arm detachable, Standard Adjustment range from seat plate: 245 - 315 mm	<input checked="" type="radio"/>
<b>ARMPAD (Pair)*</b>			
JBB2196	<input type="checkbox"/>	Long Armpad (375 mm)	<input checked="" type="radio"/>
<b>HEADREST / NECKREST - Adequate headrest option must be selected*</b>			
JBB1350	<input type="checkbox"/>	Neckrest depth adjustable	<input checked="" type="radio"/>
<b>SEAT LIFTER / SEAT TILT*</b>			
JBB1092	<input type="checkbox"/>	Powered "fixed pivot" Tilt 20°	<input checked="" type="radio"/>
<b>SEAT T GROUND HEIGHT*</b>			
JBB2857	<input type="checkbox"/>	Seat-to-Floor Height 435 mm	<input type="radio"/>
JBB2858	<input type="checkbox"/>	Seat-to-Floor Height 460 mm	<input type="radio"/>
JBB2859	<input type="checkbox"/>	Seat-to-Floor Height 485 mm	<input type="radio"/>
<b>BACKREST RECLINER*</b>			
JBB0975	<input type="checkbox"/>	Manual Settable Recliner (in 3.75° steps) 90° - 120°	<input checked="" type="radio"/>
<b>FOOTREST HANGER TYPE*</b>			
JBB0526	<input type="checkbox"/>	Standard hanger 80° pair Lower leg length: 350 - 450 mm	<input checked="" type="radio"/>
<b>FOOTPLATE TYPE*</b>			
JBB0512	<input type="checkbox"/>	Fixed Footplate pair	<input checked="" type="radio"/>
<b>CONTROLS*</b>			
JBB2203	<input type="checkbox"/>	Remote LiNX DLX-REM211 = Drive & Actuators 3 Drive functions and power module incl. 2 Actuators	<input checked="" type="radio"/>
<b>CONTROL BRACKETRY*</b>			
JBB1707	<input type="checkbox"/>	Swing Away Bracket RH	<input type="radio"/>
JBB1708	<input type="checkbox"/>	Swing Away Bracket LH	<input type="radio"/>

SHROUD COLOUR*			
JBB2871	<input type="checkbox"/>	Mid range - Black/Anthracite Shroud	<input type="radio"/>
JBB2872	<input type="checkbox"/>	Mid range - Silver Shroud	<input type="radio"/>
JBB2873	<input type="checkbox"/>	Mid range - Bengal Red Shroud	<input type="radio"/>
JBB2874	<input type="checkbox"/>	Mid range - Mercury Blue Shroud	<input type="radio"/>
JBB2875	<input type="checkbox"/>	Mid range - Midnight Purple Shroud	<input type="radio"/>
RIM COLOUR*			
JBB2882	<input type="checkbox"/>	Mid range - Black Rim Inserts	<input type="radio"/>
JBB2883	<input type="checkbox"/>	Mid range - Silver Rim Inserts	<input type="radio"/>
JBB2884	<input type="checkbox"/>	Mid range - Bengal Red Rim Inserts	<input type="radio"/>
JBB2885	<input type="checkbox"/>	Mid range - Mercury Blue Rim Inserts	<input type="radio"/>
JBB2886	<input type="checkbox"/>	Mid range - Midnight Purple Rim Inserts	<input type="radio"/>
FORKS*			
JBB2632	<input type="checkbox"/>	Double sided fork Without castor rim insert (enlarges wide base to 656 mm)	<input checked="" type="checkbox"/>
FRONT WHEEL *			
JBB0828	<input type="checkbox"/>	8" (200 mm) x 2.50-4 (65 mm) Castor, puncture proof solid tyre	<input checked="" type="checkbox"/>
DRIVE WHEEL *			
JBB0705	<input type="checkbox"/>	DW 14" (355 mm) x 3" (75 mm), Solid inlay Dual density foam insert	<input checked="" type="checkbox"/>
C.T.C SUSPENSION*			
JBB2870	<input type="checkbox"/>	CTC Suspension Control - Traction - Comfort	<input checked="" type="checkbox"/>
SUSPENSION SETTING*			
JBB1910	<input type="checkbox"/>	Soft (suggested for up to 70 kg User Weight)	<input checked="" type="checkbox"/>
JBB1912	<input type="checkbox"/>	Medium (suggested for up to 90 kg User Weight)	<input type="radio"/>
JBB1911	<input type="checkbox"/>	Firm (suggested for up to 110 kg User Weight)	<input type="radio"/>
JBB1915	<input type="checkbox"/>	Very firm (suggested for up to 136 kg User Weight)	<input type="radio"/>
MOTOR TYPE*			
JBB1282	<input type="checkbox"/>	6 km/h	<input checked="" type="checkbox"/>
BATTERY TYPE*			
JBB1259	<input type="checkbox"/>	60 Ah Gel	<input checked="" type="checkbox"/>
CHARGER*			
JBB1343	<input type="checkbox"/>	Opti Charge Charger 8A (IP44)	<input checked="" type="checkbox"/>
ANTI-TIPPERS*			
JBB1205	<input type="checkbox"/>	Anti-tipper	<input checked="" type="checkbox"/>
TRANSPORTATION KIT*			
JBB1494	<input type="checkbox"/>	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19) Adequate headrest is compulsory	<input checked="" type="checkbox"/>
LAP BELT*			
JBB1329	<input type="checkbox"/>	Lap belt with Auto Lock Standard mount on backrest	<input checked="" type="checkbox"/>

<b>TOTAL PRICE</b>	
--------------------	--

**Non Conforming Options:**



For a more complex customer requirement please refer to our Invacare Unique special order form. email [specialsuk@invacare.com](mailto:specialsuk@invacare.com) or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Description of special request (please include relevant option code or part number):