



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Motability Retail Prescription Form

JULY 2024
Customer Ref:

Table with 2 columns: Customer Ref*, Contact No*, Delivery Address* and Account No*, Order No*, Order Qty*, Order date*

E-Fix & Action 3NG M82 707362 - Weekly
Price £46.40

Standard Features

- Action 3NG + E-Fix
Max user weight 125kg (19.7 stone)
Crash tested ISO- 7176-Part 19
Seat widths 380mm to 505mm (15" -20")
Seat depths 400mm or 450mm (16" or 18")
Seat heights 460mm to 510mm (18" -20")
Nylon with hook and loop strips seat upholstery
Fixed 7°Backrest 2 heights 400mm/430mm
Lumbar support padded backrest with shaped cover
80°swing in/swing out hangers
Flip up footplates
Removable full arms
Basic 50mm foam cushion
Adjustable calf strap
Choice of 8 standard colours
8" Black solid front wheels, with 24" Black solid tyre spoked rim
Aluminium Pre-anodized handrim
Quick release axle
Push to lock brakes
Transportation kit (including buckle belt)
Alber Efix E35 1566293 (24" wheels)
Brackets
Alber anti-tippers
Adapter Kit
Mobility Plus package
Fitting at Invacare



Please make your selection with a cross in the box [] Free of charge [O] standard feature [◆] POA special request

Customer weight = KG

Table with 2 columns: Code, Description. Includes BASE PRODUCT* and PROPULSION* sections.

SEAT WIDTH*			
CEE0180	<input type="checkbox"/>	15" / 380 MM SEAT WIDTH (lapbelt adjustment 420mm - 1520mm)	○
CEE0190	<input type="checkbox"/>	16" / 405 MM SEAT WIDTH (lapbelt adjustment 420mm - 1520mm)	○
CEE0200	<input type="checkbox"/>	17" / 430 MM SEAT WIDTH (lapbelt adjustment 420mm - 1520mm)	◆
CEE0210	<input type="checkbox"/>	18" / 455 MM SEAT WIDTH (lapbelt adjustment 420mm - 1880mm)	○
CEE0220	<input type="checkbox"/>	19" / 480 MM SEAT WIDTH (lapbelt adjustment 420mm - 1880mm)	○
CEE0230	<input type="checkbox"/>	20" / 505 MM SEAT WIDTH (lapbelt adjustment 420mm - 1880mm)	○
SEAT DEPTH*			
CEE0240	<input type="checkbox"/>	16" / 400 MM SEAT DEPTH (NC 0879)	○
CEE0250	<input type="checkbox"/>	18" / 450 MM SEAT DEPTH (NC 0878)	◆
SEAT HEIGHT*			
CEE0260	<input type="checkbox"/>	20" / 510 MM SEAT HEIGHT	○
CEE0270	<input type="checkbox"/>	19" / 485 MM SEAT HEIGHT	○
CEE0280	<input type="checkbox"/>	18" / 460 MM SEAT HEIGHT	◆
SEAT UPHOLSTERY*			
CEE0330	<input type="checkbox"/>	NYLON W/O HOOK AND LOOP STRIPS	◆
BACKREST*			
CEE0410	<input type="checkbox"/>	FIXED 7° BACKREST 2 HEIGHTS 16" / 400 or 17" / 430 MM FIX PUSH HANDLES	◆
BACKREST UPHOLSTERY*			
CEE0474	<input type="checkbox"/>	LUMBAR SUPPORT BI MATERIAL PADDED COVER	◆
LEGREST HANGER RH*			
CEE0510	<input type="checkbox"/>	80° RH SWING IN/SWING OUT	◆
LEGREST HANGER LH*			
CEE0560	<input type="checkbox"/>	80° LH SWING IN/SWING OUT	◆
FOOTREST RH*			
CEE0600	<input type="checkbox"/>	FLIP UP RH FOOTPLATE	◆
FOOTREST LH*			
CEE0620	<input type="checkbox"/>	FLIP UP LH FOOTPLATE	◆
ARMREST RH*			
CEE0640	<input type="checkbox"/>	REMOVABLE RH (N°1) DESK ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED (NC 0760, 0839)	○
CEE0650	<input type="checkbox"/>	REMOVABLE RH (N°1) FULL LENGTH ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED (NC 0760, 0839)	○
CEE0660	<input type="checkbox"/>	SWING AWAY DETACHABLE RH (N°2) DESK 240 MM HEIGHT (NC 0740, 0750, 0839)	◆
ARMREST LH*			
CEE0740	<input type="checkbox"/>	REMOVABLE LH (N°1) DESK ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED	○
CEE0750	<input type="checkbox"/>	REMOVABLE LH (N°1) FULL LENGTH ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED	○
CEE0760	<input type="checkbox"/>	SWING AWAY DETACHABLE LH (N°2) DESK 240 MM HEIGHT	◆
CEE0839	<input type="checkbox"/>	TUBULAR ARMREST LH, SWING AWAY PADDED TUBE WITH REMOVABLE MUDGUARD LH (N°6) (CUSHION IS RECOMMENDED FOR SEAT HEIGHT 410 AND 435 MM) (NC 1200, 1345)	○

SEAT CUSHION			
CEE0849	<input type="checkbox"/>	BASIC 50 MM FOAM CUSHION NYLON AND ANTISLIP COVER	◆
SEAT CUSHION DEPTH			
CEE0878	<input type="checkbox"/>	FOR 400 MM SEAT DEPTH	○
CEE0879	<input type="checkbox"/>	FOR 450 MM SEAT DEPTH	○
COLOUR*			
CEE0961	<input type="checkbox"/>	SLATE GREY (RAL7024 FINE TEXTURE)	○
CEE0964	<input type="checkbox"/>	PEARL WHITE	○
CEE0965	<input type="checkbox"/>	ELECTRIC GREEN	○
CEE0966	<input type="checkbox"/>	MATT BLACK	○
CEE0974	<input type="checkbox"/>	BLUE	○
CEE0976	<input type="checkbox"/>	DARK PURPLE	○
CEE0969	<input type="checkbox"/>	RED MATT	○
CEE0972	<input type="checkbox"/>	WATER BLUE	○
CEE0980	<input type="checkbox"/>	ANTHRACITE BLACK	○
CASTOR HOUSING*			
CEE0995	<input type="checkbox"/>	STANDARD	◆
FRONT CASTORS*			
CEE1010	<input type="checkbox"/>	8" X 1.1/4" BLACK SOLID TYRE	◆
CEE1030	<input type="checkbox"/>	8" X 2" BLACK SOLID TYRE	○
REAR WHEEL SIZE*			
CEE1068	<input type="checkbox"/>	24"	◆
REAR WHEELS*			
CEE1071	<input type="checkbox"/>	STANDARD (BASIC) (BLACK HUB, COLORLESS RIM AND ZINC PLATED SPOKES)	◆
TYRES*			
CEE1082	<input type="checkbox"/>	SOLID BLACK 1.3/8" / 35MM	◆
CEE1088	<input type="checkbox"/>	PNEUMATIC LOW PROFILE BLACK 1" / 25MM	○
HANDRIM*			
CEE1140	<input type="checkbox"/>	ALUMINIUM PREANODIZED HANDRIM	◆
REAR WHEEL AXLE*			
CEE1180	<input type="checkbox"/>	Q/R AXLE	◆
REAR WHEEL POSITION*			
CEE1190	<input type="checkbox"/>	STANDARD REAR WHEELS POSITION	◆
CEE1195	<input type="checkbox"/>	MID WHEEL PROPULSION REAR WHEEL POSITION	○
CEE1200	<input type="checkbox"/>	PASSIVE REAR WHEELS POSITION	○
BRAKES*			
CEE1270	<input type="checkbox"/>	PUSH TO LOCK BRAKES	◆
TRANSPORTATION KIT*			
CEE1345	<input type="checkbox"/>	TRANSPORTATION BRACKETS	◆
SAFETY OPTIONS*			
CEE1360	<input type="checkbox"/>	POSTURE BELT	◆

TOTAL PRICE	<input type="text"/>
--------------------	----------------------

Non Conforming Options:

CEE0740	REMOVABLE LH (N°1) DESK ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED
CEE0750	REMOVABLE LH (N°1) FULL LENGTH ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED
CEE0760	SWING AWAY DETACHABLE LH (N°2) DESK 240 MM HEIGHT
CEE0839	TUBULAR ARMREST LH, SWING AWAY PADDED TUBE WITH REMOVABLE MUDGUARD LH (N°6) (CUSHION IS RECOMMENDED FOR SEAT HEIGHT 410 AND 435 MM)
CEE0878	FOR 400 MM SEAT DEPTH
CEE0879	FOR 450 MM SEAT DEPTH
CEE1195	MID WHEEL PROPULSION REAR WHEEL POSITION
CEE1200	PASSIVE REAR WHEELS POSITION
CEE1345	TRANSPORTATION BRACKETS



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):