



# Motability Retail Prescription Form

**Yes, you can:**

Tel: 01656 776222 Fax: 01656 776220  
Email: ordersuk@invacare.com  
Online Spares available at: www.invacare.co.uk

JULY 2024  
Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## Storm4 X-Plore Modulite M34 707458 - Weekly Price £58.10

### Standard Features

- Max user weight 150Kg (23.5 stone)
- Crash tested ISO7176-Part 19
- Modulite seating system
- Flex3 Laguna backrest
- Matrix Flo-tech Lite cushion
- Backposts 540mm (Tall)
- Standard T-armrests, detachable with Long armpads
- Manual settable seat Tilt
- Manual settable backrest Recline
- Angle adjustable legrest hangers (VARI F) with Fixed footplates
- LINX REM216 control with swing away bracket
- Puncture proof (sealed) drive wheels and front castors
- 12 km/h (8mph) includes LED light system with indicators, Mirror
- Lapbelt (adjustable both sides)

**Standard Base price: £7,885**



Please make your selection with a cross in the box     Free of charge     standard feature    POA special request

Customer weight =                      KG

Code	Description	
<b>CHAIR TYPE*</b>		
ATE0017	Storm 4 X-plore	£7,885
ATZ0017	Special Storm 4 X-plore	POA
<b>SEATING TYPE*</b>		
ATE2100	Invacare Modulite 'Flex3'	◆
<b>SEAT FRAME*</b>		
ATE2111	Telescopic Seat Frame with metal seat plate	◆

<b>SEAT WIDTH*</b>			
ATEXXXX	<input type="checkbox"/>	Seat width 330mm (13")	<input type="radio"/>
ATE2121	<input type="checkbox"/>	Set at 380mm	<input type="radio"/>
ATE2122	<input type="checkbox"/>	Set at 430mm	<input type="radio"/>
ATE2123	<input type="checkbox"/>	Set at 480mm	<input type="radio"/>
ATE2124	<input type="checkbox"/>	Set at 530mm	<input type="radio"/>
<b>SEAT DEPTH*</b>			
ATE2132	<input type="checkbox"/>	410mm (16")	(NC 2196) <input type="radio"/>
ATE2134	<input type="checkbox"/>	460mm (18")	<input type="radio"/>
ATE2136	<input type="checkbox"/>	510mm (20")	<input type="radio"/>
<b>SEAT CUSHION*</b>			
ATE2142	<input type="checkbox"/>	Matrx Flo-tech Lite	<input checked="" type="radio"/>
<b>BACKREST SUPPORT*</b>			
ATE2162	<input type="checkbox"/>	Backposts 540mm (Tall)	<input checked="" type="radio"/>
<b>PUSH HANDLES FOR BACKPOSTS*</b>			
ATE2107	<input type="checkbox"/>	Backposts WITH push handles	<input checked="" type="radio"/>
<b>BACKREST UPHOLSTREY / CUSHION*</b>			
ATE2173	<input type="checkbox"/>	Flex3 adjustable Backrest Plates with Laguna Backrest cushion	<input checked="" type="radio"/>
<b>UPHOLSTERY COLOUR*</b>			
ATE2152	<input type="checkbox"/>	Black multi stretch PU coated fabric	<input type="radio"/>
ATE2155	<input type="checkbox"/>	Black Spacer fabric	<input type="radio"/>
<b>ARMREST*</b>			
ATE2191	<input type="checkbox"/>	T-Armrest Standard	<input checked="" type="radio"/>
<b>ARMPAD*</b>			
ATE2196	<input type="checkbox"/>	Long Armpad (375mm)	<input checked="" type="radio"/>
<b>HEADREST / NECKREST - Adequate headrest option must be selected*</b>			
ATE1495	<input type="checkbox"/>	Push bar	<input checked="" type="radio"/>
<b>SEAT HEIGHT SETTING &amp; SEAT RISER*</b>			
ATE1023	<input type="checkbox"/>	Fixed Seat Height 450mm	<input type="radio"/>
ATE1025	<input type="checkbox"/>	Fixed Seat Height 480mm	<input type="radio"/>
<b>SEAT TILT*</b>			
ATE1002	<input type="checkbox"/>	Manual Seat Tilt by spindle	<input checked="" type="radio"/>
<b>BACKREST RECLINER*</b>			
ATE0975	<input type="checkbox"/>	Manual Settable Recliner (0° - 30°, in 3.75° steps) Foldable	<input checked="" type="radio"/>
<b>FOOTREST HANGER TYPE (PAIR)*</b>			
ATE0503	<input type="checkbox"/>	Vari F : Angle adjustable with tools 0° - 70°, length 290mm - 460mm	<input checked="" type="radio"/>
<b>FOOTPLATE TYPE*</b>			
ATE0512	<input type="checkbox"/>	Fixed Footplate	<input checked="" type="radio"/>
<b>CONTROLS*</b>			
ATE2204	<input type="checkbox"/>	REM216 Drive, Actuators & Lights	<input checked="" type="radio"/>
<b>CONTROL BRACKETRY*</b>			
ATE1707	<input type="checkbox"/>	Swing Away RH Bracket	<input type="radio"/>
ATE1708	<input type="checkbox"/>	Swing Away LH Bracket	<input type="radio"/>

<b>SHROUD COLOUR*</b>			
ATE0202	<input type="checkbox"/>	Blue	<input type="radio"/>
ATE0228	<input type="checkbox"/>	Black	<input type="radio"/>
ATE0220	<input type="checkbox"/>	Champagne	<input type="radio"/>
ATE0204	<input type="checkbox"/>	Transparent Cherry Red	<input type="radio"/>
<b>CASTOR WHEEL*</b>			
ATE0853	<input type="checkbox"/>	3.00 - 6 Grey Puncture Proof (sealed)	<input checked="" type="radio"/>
<b>DRIVE WHEEL*</b>			
ATE0773	<input type="checkbox"/>	3.00 - 8 Grey Puncture Proof (sealed)	<input checked="" type="radio"/>
<b>MOTOR TYPE*</b>			
ATE1909	<input type="checkbox"/>	12 km/h (8mph) includes LED light system with indicators	<input checked="" type="radio"/>
<b>BATTERY*</b>			
ATE1364	<input type="checkbox"/>	73.5 Ah GEL	<input checked="" type="radio"/>
<b>REAR SUSPENSION*</b>			
ATE1233	<input type="checkbox"/>	Hitec spring X-plore	<input type="radio"/>
<b>FRONT SUSPENSION*</b>			
ATE1238	<input type="checkbox"/>	Front suspension (suggested for up to 100kg)	<input type="radio"/>
ATE1239	<input type="checkbox"/>	Front suspension (suggested for above 100kg)	<input type="radio"/>
<b>CHARGER*</b>			
ATE1343	<input type="checkbox"/>	OptiCharge Charger 8A (200 - 250V)	<input checked="" type="radio"/>
<b>ANTI-TIPPERS*</b>			
ATE1205	<input type="checkbox"/>	Anti-tipper (standard version)	<input checked="" type="radio"/>
<b>TRANSPORTATION*</b>			
ATE1494	<input type="checkbox"/>	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory).	<input checked="" type="radio"/>
<b>FIXATION BELT*</b>			
ATE1329	<input type="checkbox"/>	Lapbelt (adjustable both sides) fitted to base of backrest	<input checked="" type="radio"/>
<b>REAR VIEW MIRROR*</b>			
ATE1415	<input type="checkbox"/>	Rear view Mirror RH - Must be selected with 10 km/h Motor	<input type="radio"/>
ATE1339	<input type="checkbox"/>	Rear view Mirror LH - Must be selected with 10 km/h Motor	<input type="radio"/>

<b>TOTAL PRICE</b>	
--------------------	--

**Non Conforming Options:**

ATE2196          Long Armpad (375mm)

---



For a more complex customer requirement please refer to our Invacare Unique special order form. email [specialsuk@invacare.com](mailto:specialsuk@invacare.com) or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Description of special request (please include relevant option code or part number):