



# Motability Retail Prescription Form

**Yes, you can:**

Tel: 01656 776222 Fax: 01656 776220  
Email: ordersuk@invacare.com  
Online Spares available at: www.invacare.co.uk

JULY 2024  
Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

**TDX SP2 NB M1 744349 - Weekly Price  
£48.98**

### Standard Features

- Max User Weight 150kg (23.5 Stone)
- Crash test ISO7176 Part-19
- Modulite Seating
- Matrix Flo-tech cushion
- Flex3 adjustable backrest plates with Laguna backrest cushion
- Standard T-armrests, detachable, Long armpad (325mm)
- Manual settable seat Tilt (0° - 20° )
- Manual settable Recliner (0° - 30° in 3.75° steps)
- 80° hangers with fixed footplate
- LiNX control REM211 with swing away bracket
- Solid Black tyres
- 6km/h (4mph) motors
- 4 Point tie-down brackets (Q'Straint QLK-150 compatible)
- Lap belt with auto lock (flexible mount in side rails)

**Standard Base price: £5,964**



Please make your selection with a cross in the box     Free of charge     standard feature    POA    special request

Customer weight =                      KG

Code	Description	
<b>BASE PRODUCT*</b>		
ARM2100	<input type="checkbox"/> TDX SP2 NARROW BASE	£5,964
<b>USER WEIGHT*</b>		
ARM1920	<input type="checkbox"/> User weight up to 125kg	<input type="radio"/>
ARM1921	<input type="checkbox"/> User weight 126kg to 136kg	<input type="radio"/>
ARM1922	<input type="checkbox"/> User weight 137kg to 150kg	<input type="radio"/>
<b>BASE MODEL*</b>		
ARM0000	<input type="checkbox"/> TDX SP2 Narrow Base	<input checked="" type="checkbox"/>

<b>SEATING TYPE*</b>			
ARM2100	<input type="checkbox"/>	Modulite Seating Unit	◆
<b>SEAT FRAME*</b>			
ARM2111	<input type="checkbox"/>	Telescopic Seat Frame	◆
<b>SEAT WIDTH*</b>			
ARM2121	<input type="checkbox"/>	380mm (15")	○
ARM2122	<input type="checkbox"/>	430mm (17")	○
ARM2123	<input type="checkbox"/>	480mm (19")	○
ARM2124	<input type="checkbox"/>	530mm (21")	○
<b>SEAT DEPTH RAIL LENGTH*</b>			
ARM2129	<input type="checkbox"/>	Standard 410mm - 510mm	(NC 2138) ○
ARM2190	<input type="checkbox"/>	Extended 510mm - 560mm	(NC 2132, 2134) ○
<b>SEAT DEPTH*</b>			
ARM2132	<input type="checkbox"/>	410mm (16")	(NC 2196) ○
ARM2134	<input type="checkbox"/>	460mm (18")	○
ARM2136	<input type="checkbox"/>	510mm (20")	○
ARM2138	<input type="checkbox"/>	560mm (22")	○
<b>SEAT CUSHION*</b>			
ARM2142	<input type="checkbox"/>	Matrix Flo-Tech Lite	◆
<b>BACKREST SUPPORT*</b>			
ARM2162	<input type="checkbox"/>	Tall 540mm	◆
<b>BACKREST UPHOLSTERY / CUSHION*</b>			
ARM2173	<input type="checkbox"/>	Flex3 Backrest plates with Laguna cushion	◆
<b>UPHOLSTERY TYPE*</b>			
ARM2152	<input type="checkbox"/>	Black Black multi stretch PU coated fabric Upholstery	○
ARM2155	<input type="checkbox"/>	Black Spacer fabric	○
<b>ARMREST*</b>			
ARM2191	<input type="checkbox"/>	T-Armrest standard	◆
<b>ARMPAD*</b>			
ARM2196	<input type="checkbox"/>	Long 375mm Armpad	◆
<b>HEADREST / NECKREST - Adequate headrest option must be selected*</b>			
ARM1495	<input type="checkbox"/>	Push bar 22mm	◆
<b>SEAT LIFTER / SEAT TILT*</b>			
ARM1091	<input type="checkbox"/>	Manually settable	◆
<b>SEAT HEIGHT SETTING FOR 20° SEAT TILT*</b>			
ARM0354	<input type="checkbox"/>	Seat height 450mm (17.75")	○
ARM0364	<input type="checkbox"/>	Seat height 475mm (18.75")	○
<b>BACKREST RECLINER*</b>			
ARM0975	<input type="checkbox"/>	Manual settable Recliner	◆
<b>FOOTREST HANGER TYPE*</b>			
ARM0526	<input type="checkbox"/>	Side mounted : Standard 80° hanger	◆
<b>FOOTPLATE TYPE*</b>			
ARM0512	<input type="checkbox"/>	Fixed Footplate	◆
<b>CONTROLS*</b>			
ARM2203	<input type="checkbox"/>	REM211	◆

<b>CONTROL BRACKETRY*</b>		
ARM1707	<input type="checkbox"/>	Swing away RH Bracket <span style="float: right;">○</span>
ARM1708	<input type="checkbox"/>	Swing away LH Bracket <span style="float: right;">○</span>
<b>SHROUD COLOUR*</b>		
ARM0228	<input type="checkbox"/>	Wet Black <span style="float: right;">○</span>
ARM2601	<input type="checkbox"/>	SpaceStation Silver <span style="float: right;">○</span>
ARM2602	<input type="checkbox"/>	Island Blue <span style="float: right;">○</span>
ARM2603	<input type="checkbox"/>	Grape JellyBean (Purple) <span style="float: right;">○</span>
ARM2607	<input type="checkbox"/>	Rockstar Red <span style="float: right;">○</span>
ARM2608	<input type="checkbox"/>	Lights Out Black - Matt <span style="float: right;">○</span>
ARM2609	<input type="checkbox"/>	Invacare Blue <span style="float: right;">○</span>
<b>RIM COLOUR*</b>		
ARM2611	<input type="checkbox"/>	SpaceStation Silver <span style="float: right;">○</span>
ARM2612	<input type="checkbox"/>	Island Blue <span style="float: right;">○</span>
ARM2613	<input type="checkbox"/>	Grape JellyBean (Purple) <span style="float: right;">○</span>
ARM2617	<input type="checkbox"/>	Rockstar Red <span style="float: right;">○</span>
ARM2618	<input type="checkbox"/>	Lights Out Black - Matt <span style="float: right;">○</span>
ARM2619	<input type="checkbox"/>	Invacare Blue <span style="float: right;">○</span>
ARM2620	<input type="checkbox"/>	Wet Black <span style="float: right;">○</span>
ARM2621	<input type="checkbox"/>	Without rim insert <span style="float: right;">○</span>
<b>FORKS*</b>		
ARM2632	<input type="checkbox"/>	Double sided <span style="float: right;">◆</span>
<b>CENTRE WHEELS (DRIVING WHEELS)*</b>		
ARM0733	<input type="checkbox"/>	Black tyres with gel foam inserts (solid) <span style="float: right;">◆</span>
<b>MOTOR TYPE*</b>		
ARM1282	<input type="checkbox"/>	6 km/h (4mph) <span style="float: right;">◆</span>
<b>BATTERY*</b>		
ARM1257	<input type="checkbox"/>	50 Ah Gel <span style="float: right;">◆</span>
<b>CHARGER*</b>		
ARM1343	<input type="checkbox"/>	OptiCharge 8A (200 - 250V) <span style="float: right;">◆</span>
<b>TRANSPORTATION KIT*</b>		
ARM1494	<input type="checkbox"/>	4-Point tie-down <span style="float: right;">◆</span>
<b>FIXATION BELT*</b>		
ARM1329	<input type="checkbox"/>	Lap Belt fitted to base of backrest <span style="float: right;">○</span>
ARM2113	<input type="checkbox"/>	Lap Belt adjustable mounting on seat rail extrusion <span style="float: right;">◆</span>

<b>TOTAL PRICE</b>	<input style="width: 100px; height: 20px;" type="text"/>
--------------------	--

**Non Conforming Options:**

ARM2132	410mm (16")
ARM2134	460mm (18")
ARM2138	560mm (22")
ARM2196	Long 375mm Armpad



For a more complex customer requirement please refer to our Invacare Unique special order form. email [specialsuk@invacare.com](mailto:specialsuk@invacare.com) or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Description of special request (please include relevant option code or part number):