



Motability Retail Prescription Form

Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

JULY 2024
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

**TDX SP2 HD M10 900092 - Weekly Price
£47.38**

Standard Features

- Max user weight 180kg (23.6 stone)
- Crash tested ISO 7176 -Part 19
- Modulite Seating System
- Matrix Flo-tech Cushion
- Flex 3 adjustable Backrest plates with Laguna backrest cushion
- Standard T-armrests, detachable with Long armpad
- Manual settable seat Tilt 0° - 20°
- Manual settable backrest Recline 0° - 30° in 3.75° steps
- 80° hangers with fixed footplate
- LiNX control REM216 swing away bracket
- Solid Black tyres
- 6km/h (4mph) 4 pole motors
- 4-Point tie-down system (Q'Straint QLK-150 compatible)
- Lapbelt with autolock (adjustable both sides)

Standard Base price: £5,976



Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

Code	Description	
CHAIR TYPE*		
ARJ0000	TDX SP2 Heavy Duty	<input checked="" type="checkbox"/>
User Weight*		
ARJ1918	User weight less 100 kg	<input type="checkbox"/>
ARJ1919	User weight 101 to 125 kg	<input type="checkbox"/>
ARJ1921	User weight 126 to 136 kg	<input type="checkbox"/>
ARJ1922	User weight 137 to 150 kg	<input type="checkbox"/>
BASE TYPE*		
ARJ2100	Modulite Seating unit	£5,976

SEAT WIDTH*			
ARJ2121	<input type="checkbox"/>	Set at 380mm (15") (adjustable from 380mm to 430mm (15" - 17") with stepless adjustment on armrests, legrests and backrest) (NC 2138 + 2142)	<input type="radio"/>
ARJ2122	<input type="checkbox"/>	Set at 430mm (17") (adjustable from 430mm to 480mm (17" - 19") with stepless adjustment on armrests, legrests and backrest) (NC 2138 + 2142)	<input type="radio"/>
ARJ2123	<input type="checkbox"/>	Set at 480mm (19") (adjustable from 480mm to 530mm (19" - 21") with stepless adjustment on armrests, legrests and backrest) (NC 2138 + 2142)	<input type="radio"/>
ARJ2124	<input type="checkbox"/>	Set at 530mm (21") (adjustable from 530mm to 580mm (21" - 23") with stepless adjustment on armrests, legrests and backrest)	<input type="radio"/>
SEAT DEPTH RAIL LENGTH*			
ARJ2129	<input type="checkbox"/>	Standard Depth Package (16" - 20", 410mm - 510mm simple stepless adjustment from the side) (NC 2138, 2196)	<input type="radio"/>
ARJ2130	<input type="checkbox"/>	Extended Depth Package (20" - 22", 510mm - 560mm simple stepless adjustment from the side) (NC 2132, 2134)	<input type="radio"/>
SEAT DEPTH SETTING*			
ARJ2132	<input type="checkbox"/>	410mm (16") (410mm - 510mm simple stepless adjustment from the side) (NC 2196)	<input type="radio"/>
ARJ2134	<input type="checkbox"/>	460mm (18") (410mm - 510mm simple stepless adjustment from the side)	<input type="radio"/>
ARJ2136	<input type="checkbox"/>	510mm (20") (410mm - 510mm simple stepless adjustment from the side)	<input type="radio"/>
ARJ2138	<input type="checkbox"/>	560mm (22") (510mm - 560mm simple stepless adjustment from the side) (Cushion (option 2142) is only available with Seat Width (option 2124) 530mm)	<input type="radio"/>
SEAT CUSHION*			
ARJ2142	<input type="checkbox"/>	Matrx Flo-tech Lite	<input checked="" type="checkbox"/>
BACKREST SUPPORT*			
ARJ2162	<input type="checkbox"/>	Backposts 540mm (Tall) must be selected for Flex3 adjustable backrest plates	<input checked="" type="checkbox"/>
PUSH HANDLES FOR BACKPOSTS*			
ARJ2107	<input type="checkbox"/>	Backposts WITH push handles	<input checked="" type="checkbox"/>
BACKREST UPHOLSTERY / CUSHION*			
ARJ2173	<input type="checkbox"/>	Flex3 adjustable Backrest Plates with Laguna Backrest cushion	<input checked="" type="checkbox"/>
UPHOLSTERY COLOUR*			
ARJ2152	<input type="checkbox"/>	Black Black multi stretch PU coated fabric Upholstery	<input type="radio"/>
ARJ2155	<input type="checkbox"/>	Black Spacer fabric	<input type="radio"/>
ARMREST*			
ARJ2191	<input type="checkbox"/>	T-Armrest Standard, Range of adjustment from seat plate 245mm to 315mm	<input checked="" type="checkbox"/>
ARMPAD*			
ARJ2196	<input type="checkbox"/>	Long Armpad (375mm) (NC 2132)	<input checked="" type="checkbox"/>
HEADREST / NECKREST - Adequate headrest option must be selected*			
ARJ1495	<input type="checkbox"/>	Push bar dia. 22mm in conjunction with telescopic seat frame/backposts, No Headrest - no mounting hardware.	<input checked="" type="checkbox"/>
SEAT LIFTER / SEAT TILT*			
ARJ1091	<input type="checkbox"/>	Manually settable Seat Tilt (0° - 20°) Seat Height 450mm	<input checked="" type="checkbox"/>
SEAT HEIGHT SETTING FOR 20 DEGREES SEAT TILT*			
ARJ0354	<input type="checkbox"/>	Seat height set at 450mm (17.75") (from floor to front of seat plate at 0° tilt)	<input checked="" type="checkbox"/>
BACKREST RECLINER*			
ARJ0975	<input type="checkbox"/>	Manual Settable Recliner (0° - 30°, in 3.75° steps) Foldable	<input checked="" type="checkbox"/>

FOOTREST HANGER TYPE*			
ARJ0526	<input type="checkbox"/>	Side Mounted : Standard 80° hangers	◆
FOOTPLATE TYPE*			
ARJ0512	<input type="checkbox"/>	Fixed Footplate	◆
CONTROLS*			
ARJ2201	<input type="checkbox"/>	REM110 Drive only	◆
CONTROL BRACKETRY*			
ARJ1707	<input type="checkbox"/>	Swing Away RH Bracket	○
ARJ1708	<input type="checkbox"/>	Swing Away LH Bracket	○
SHROUD COLOUR*			
ARJ0228	<input type="checkbox"/>	Wet Black	○
ARJ2601	<input type="checkbox"/>	SpaceStation Silver	○
ARJ2602	<input type="checkbox"/>	Island Blue	○
ARJ2603	<input type="checkbox"/>	Grape JellyBean (Purple)	○
ARJ2607	<input type="checkbox"/>	Rockstar Red	○
ARJ2608	<input type="checkbox"/>	Lights Out Black - Matt	○
ARJ2609	<input type="checkbox"/>	Invacare Blue	○
RIM COLOUR*			
ARJ2611	<input type="checkbox"/>	SpaceStation Silver	○
ARJ2612	<input type="checkbox"/>	Island Blue	○
ARJ2613	<input type="checkbox"/>	Grape JellyBean (Purple)	○
ARJ2617	<input type="checkbox"/>	Rockstar Red	○
ARJ2618	<input type="checkbox"/>	Lights Out Black - Matt	○
ARJ2619	<input type="checkbox"/>	Invacare Blue	○
ARJ2620	<input type="checkbox"/>	Wet Black	○
ARJ2621	<input type="checkbox"/>	Without Rim insert	○
FORKS*			
ARJ2632	<input type="checkbox"/>	Double sided fork	○
CENTRE WHEELS (DRIVING WHEELS)*			
ARJ0733	<input type="checkbox"/>	355mm x 76mm (14" x 3") Black Tyres with Gel Foam Inserts	○
MOTOR TYPE*			
ARJ1282	<input type="checkbox"/>	6 km/h (4mph) (4 pole)	◆
BATTERY*			
ARJ1260	<input type="checkbox"/>	73.5 Ah (Gel)	◆
CHARGER*			
ARJ1343	<input type="checkbox"/>	OptiCharge Charger 8A (200 - 250V)	◆
TRANSPORTATION KIT*			
ARJ1494	<input type="checkbox"/>	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory, see options below)	◆
FIXATION BELT*			
ARJ2113	<input type="checkbox"/>	Lapbelt (adjustable both sides) adjustable mounting on seat rail extrusion	○

TOTAL PRICE	
--------------------	--

Non Conforming Options:

ARJ2132	410mm (16") (410mm - 510mm simple stepless adjustment from the side)
ARJ2134	460mm (18") (410mm - 510mm simple stepless adjustment from the side)
ARJ2138	560mm (22") (510mm - 560mm simple stepless adjustment from the side) (Cushion (option 2142) is only available with Seat Width (option 2124) 530mm)
ARJ2196	Long Armpad (375mm)



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):