



**Yes, you can:**

Tel: 01656 776222  
 Email: ordersuk@invacare.com  
 Online Spares available at: www.invacare.co.uk

OCTOBER 2024  
 Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

**E-Motion / DuoDrive with Action 3NG M123 & M124 706257 / 706256**

**Standard Features**

- WEEKLY PRICE:
- E-Motion M25 & Action 3NG £51.60
- E-Motion DuoDrive & Action 3NG £54.00
- 
- Max user weight 125kg (19.7 stone)
- Crash tested ISO- 7176-Part 19
- Self Propulsion with 24" Quick release rear wheels
- Choice of 6 seat widths 380mm to 505mm (15" to 20")
- Choice of 2 seat depths 400mm or 450mm (16" or 18")
- Choice of 6 seat heights 385mm to 510mm (15" to 20")
- Black nylon upholstery
- Fixed 7° 430mm backrest
- Lumbar support padded backrest with shaped cover
- 80° Swing in / Swing out legrest hangers
- Flip up footplates
- Swing away and detachable desk armrests
- Choice of 8 colours
- Push or pull to lock brakes (configuration dependant)
- Transportation kit and buckle belt (adjustable both sides)
- Adjustable heel straps
- User tool kit
- 
- Alber E-Motion M25/ DuoDrive
- 24" Stainless steel handrims
- Schwalbe Righrun Black
- Anti-tippers with jack-up function (pair)
- Fitted at Invacare Included



Please make your selection with a cross in the box     Free of charge     standard feature    POA special request

Customer weight =                      KG

**Code**                      **Description**

**BASE PRODUCT\***

CEE0003        ACTION 3NG SELF PROPEL

PROPULSION*			
CEE0100		SELF PROPULSION	◆
CEE9250		E-Motion M25 Weekly Price £47.44 Including Axle adaption kit, Anti tippers with Jack up function M124	£6,478
CEE9251		E-Motion DuoDrive Weekly Price £51.43 Including Axle adaption kit, Anti tippers with Jack up function M125	£6,937
SEAT WIDTH*			
CEE0180		15" / 380 MM SEAT WIDTH (lapbelt adjustment 420mm - 1520mm)	○
CEE0190		16" / 405 MM SEAT WIDTH (lapbelt adjustment 420mm - 1520mm)	○
CEE0200		17" / 430 MM SEAT WIDTH (lapbelt adjustment 420mm - 1520mm)	◆
CEE0210		18" / 455 MM SEAT WIDTH (lapbelt adjustment 420mm - 1880mm)	○
CEE0220		19" / 480 MM SEAT WIDTH (lapbelt adjustment 420mm - 1880mm)	○
CEE0230		20" / 505 MM SEAT WIDTH (lapbelt adjustment 420mm - 1880mm)	○
SEAT DEPTH*			
CEE0240		16" / 400 MM SEAT DEPTH (NC 0879)	○
CEE0250		18" / 450 MM SEAT DEPTH (NC 0878)	◆
SEAT HEIGHT*			
CEE0260		20" / 510 MM SEAT HEIGHT (NC 0511, 0561)	○
CEE0270		19" / 485 MM SEAT HEIGHT (NC 0511, 0561)	○
CEE0280		18" / 460 MM SEAT HEIGHT	◆
CEE0290		17" / 435 MM SEAT HEIGHT (NC 1010)	○
CEE0300		16" / 410 MM SEAT HEIGHT (NC 1010)	○
CEE0310		15" / 385 MM SEAT HEIGHT , (NC 0511, 0561, 0660, 0670, 0760, 0770, 0839, 1010, 1068)	○
SEAT UPHOLSTERY*			
CEE0330		NYLON W/O HOOK AND LOOP STRIPS	◆
CEE0340		NYLON WITH HOOK AND LOOP STRIPS	○
CEE0345		TENSION ADJUSTABLE SEAT UPHOLSTERY & DEPTH ADJUSTABLE -50MM	○
BACKREST*			
CEE0380		FIXED 0° BACKREST 2 HEIGHTS 16"/ 400 MM or 17" / 430 MM	○
CEE0390		FIXED 0° BACKREST 20" / 510 MM HEIGHT	○
CEE0410		FIXED 7° BACKREST 2 HEIGHTS 16" / 400 or 17" / 430 MM FIX PUSH HANDLES	○
CEE0420		FIXED 7° BACKREST 20" / 510 MM HEIGHT FIX PUSH HANDLES	○
CEE0428		FIXED 7° BACKREST ACTIVE 14" / 350 MM HT ADJ TO 16" / 400 MM FIX PUSH HANDLES	○
CEE0455		FOLDING 7° BACKREST 17"	○
CEE0460		FOLDING 7° BACKREST 20"	○
BACKREST UPHOLSTERY*			
CEE0465		WITHOUT BACK UPHOLSTERY	○
CEE0474		LUMBAR SUPPORT BI MATERIAL PADDED COVER	○
LEGREST HANGER RH*			
CEE0510		80° RH SWING IN/SWING OUT (NC 66666)	◆
CEE0511		90° RH SWING IN/SWING OUT (NC 0600, 1010, 66666)	○

LEGREST HANGER LH*			
CEE0560	<input type="checkbox"/>	80° LH SWING IN/SWING OUT	◆
CEE0561	<input type="checkbox"/>	90° LH SWING IN/SWING OUT	○ (NC 0620, 1010)
FOOTREST RH*			
CEE0600	<input type="checkbox"/>	FLIP UP RH FOOTPLATE	◆
CEE66666	<input type="checkbox"/>	NO FOOTREST	○
FOOTREST LH*			
CEE0620	<input type="checkbox"/>	FLIP UP LH FOOTPLATE	◆
CEE55555	<input type="checkbox"/>	NO FOOTREST	○
ARMREST RH*			
CEE0640	<input type="checkbox"/>	REMOVABLE RH (N°1) DESK ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED	○ (NC 0760, 0770, 0839)
CEE0650	<input type="checkbox"/>	REMOVABLE RH (N°1) FULL LENGTH ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED	○ (NC 0760, 0770, 0839)
CEE0660	<input type="checkbox"/>	SWING AWAY DETACHABLE RH (N°2) DESK 240 MM HEIGHT	◆ (NC 0740, 0750, 0839, 0955)
CEE0670	<input type="checkbox"/>	SWING AWAY DETACHABLE RH (N°2) FULL LENGTH 240 MM HEIGHT	○ (NC 0740, 0750, 0839, 0955)
ARMREST LH*			
CEE0740	<input type="checkbox"/>	REMOVABLE LH (N°1) DESK ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED	○
CEE0750	<input type="checkbox"/>	REMOVABLE LH (N°1) FULL LENGTH ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED	○
CEE0760	<input type="checkbox"/>	SWING AWAY DETACHABLE LH (N°2) DESK 240 MM HEIGHT	◆ (NC 0955)
CEE0770	<input type="checkbox"/>	SWING AWAY DETACHABLE LH (N°2) FULL LENGTH 240 MM HEIGHT	○ (NC 0955)
CEE0839	<input type="checkbox"/>	TUBULAR ARMREST LH, SWING AWAY PADDED TUBE WITH REMOVABLE MUDGUARD LH (N°6) (CUSHION IS RECOMMENDED FOR SEAT HEIGHT 410 AND 435 MM)	○ (NC 0955, 1200, 1345)
SEAT CUSHION DEPTH			
CEE0878	<input type="checkbox"/>	FOR 400 MM SEAT DEPTH	○
CEE0879	<input type="checkbox"/>	FOR 450 MM SEAT DEPTH	○
OTHER SEATING OPTIONS			
CEE0930	<input type="checkbox"/>	CALF STRAP	○
CEE0955	<input type="checkbox"/>	W/O AUTOLOCK SAFETY DEVICE FOR ARMREST N°1 ONLY	○
COLOUR*			
CEE0961	<input type="checkbox"/>	SLATE GREY (RAL7024 FINE TEXTURE)	○
CEE0964	<input type="checkbox"/>	PEARL WHITE	○
CEE0965	<input type="checkbox"/>	ELECTRIC GREEN	○
CEE0966	<input type="checkbox"/>	MATT BLACK	○
CEE0974	<input type="checkbox"/>	BLUE	○
CEE0976	<input type="checkbox"/>	DARK PURPLE	○
CEE0969	<input type="checkbox"/>	RED MATT	○
CEE0972	<input type="checkbox"/>	WATER BLUE	○
CEE0980	<input type="checkbox"/>	ANTHRACITE BLACK	○
CASTOR HOUSING*			
CEE0995	<input type="checkbox"/>	STANDARD	◆
FRONT CASTORS*			
CEE1010	<input type="checkbox"/>	8" X 1.1/4" BLACK SOLID TYRE	◆

REAR WHEEL SIZE*		
CEE1068	<input type="checkbox"/>	24" <span style="float: right;">◆</span>
REAR WHEELS*		
CEE1071	<input type="checkbox"/>	STANDARD (BASIC) (BLACK HUB, COLORLESS RIM AND ZINC PLATTED SPOKES ) <span style="float: right;">◆</span>
TYRES*		
CEE1082	<input type="checkbox"/>	SOLID BLACK 1.3/8" / 35MM <span style="float: right;">◆</span>
CEE1088	<input type="checkbox"/>	PNEUMATIC LOW PROFILE BLACK 1" / 25MM <span style="float: right;">○</span>
HANDRIM*		
CEE1140	<input type="checkbox"/>	ALUMINIUM PREANODIZED HANDRIM <span style="float: right;">◆</span>
REAR WHEEL AXLE*		
CEE1180	<input type="checkbox"/>	Q/R AXLE <span style="float: right;">◆</span>
REAR WHEEL POSITION*		
CEE1190	<input type="checkbox"/>	STANDARD REAR WHEELS POSITION <span style="float: right;">◆</span>
CEE1200	<input type="checkbox"/>	PASSIVE REAR WHEELS POSITION <span style="float: right;">○</span>
BRAKES*		
CEE1270	<input type="checkbox"/>	PUSH TO LOCK BRAKES <span style="float: right;">◆</span>
CEE1275	<input type="checkbox"/>	PULL TO LOCK BRAKES <span style="float: right;">○</span>
TRANSPORTATION KIT*		
CEE1345	<input type="checkbox"/>	TRANSPORTATION BRACKETS <span style="float: right;">◆</span>
SAFETY OPTIONS*		
CEE1360	<input type="checkbox"/>	POSTURE BELT <span style="float: right;">◆</span>

<b>TOTAL PRICE</b>	<input style="width: 100px; height: 20px;" type="text"/>
--------------------	--

**Non Conforming Options:**

CEE0511	90° RH SWING IN/SWING OUT
CEE0561	90° LH SWING IN/SWING OUT
CEE0600	FLIP UP RH FOOTPLATE
CEE66666	NO FOOTREST
CEE0620	FLIP UP LH FOOTPLATE
CEE0660	SWING AWAY DETACHABLE RH (N°2) DESK 240 MM HEIGHT
CEE0670	SWING AWAY DETACHABLE RH (N°2) FULL LENGTH 240 MM HEIGHT
CEE0740	REMOVABLE LH (N°1) DESK ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED
CEE0750	REMOVABLE LH (N°1) FULL LENGTH ADJ. 200 TO 280 MM HEIGHT AUTOLOCK SAFETY DEVICE INCLUDED
CEE0760	SWING AWAY DETACHABLE LH (N°2) DESK 240 MM HEIGHT
CEE0770	SWING AWAY DETACHABLE LH (N°2) FULL LENGTH 240 MM HEIGHT
CEE0839	TUBULAR ARMREST LH, SWING AWAY PADDED TUBE WITH REMOVABLE MUDGUARD LH (N°6) (CUSHION IS RECOMMENDED FOR SEAT HEIGHT 410 AND 435 MM)
CEE0878	FOR 400 MM SEAT DEPTH
CEE0879	FOR 450 MM SEAT DEPTH
CEE0955	W/O AUTOLOCK SAFETY DEVICE FOR ARMREST N°1 ONLY
CEE1010	8" X 1.1/4" BLACK SOLID TYRE
CEE1068	24"
CEE1200	PASSIVE REAR WHEELS POSITION
CEE1345	TRANSPORTATION BRACKETS



For a more complex customer requirement  
please refer to our Invacare Unique special order form.  
email [specialsuk@invacare.com](mailto:specialsuk@invacare.com)

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):