



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Motability Retail Prescription Form

OCTOBER 2024
Customer Ref:

| | |
|--------------------------|--------------------|
| Customer Ref* | Account No* |
| Contact No* | Order No* |
| Delivery Address* | Order Qty* |
| | Order date* |

Emotion M25 or Duo Drive with Kuschall K series 2.0 M121 & M125 1001887 / 706258

Standard Features

- Max User Weight 130kg (20.5 stone)
- Crash Tested ISO 7176-Part 19
- Active Rigid Wheelchair with Height Adjustable Rigid Footrest
- Multiple adjustments, Front and Rear Seat Heights and Centre of Gravity
- Full Hydroformed Aluminium Front Frame 75° Front Angle
- 24" Self Propel Q/R Wheels with Grey Aluminium Handrims and Pneumatic Tyres
- 5 position rear wheel extension
- Choice of 10 Seat Widths 320mm to 500mm
- Choice of 8 Seat Depths 350mm to 525mm
- Range of Front Seat Heights from 450mm to 520mm
- Range of Rear Seat Heights from 390mm to 490mm
- Lockdown Folding, Angle and Height Adjustable Backrest
- Black Nylon with Reflective Trim, Tension Adjustable Upholstery
- Without Push Handles, Carbon Clothesguards
- Black Aluminium Axle Tube with 0° Camber
- Black Comfort Front Castor
- Standard Brake Lever Bent
- Tool Kit



Standard Base price: £8,802

Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

| Code | Description | |
|----------------------|---|--------|
| BASE PRODUCT* | | |
| DKA0067 | KÜSCHALL K-SERIES | ◆ |
| DKA9250 | E-Motion M25 Weekly Price £61.34 Including Anti tippers with Jack up function M121 | £8,343 |
| DKA9251 | E-Motion DuoDrive Weekly Price £63.48 Including Anti tippers with Jack up function M125 | £8,802 |

| FRAME* | | | | |
|---------------|--|---|---|---|
| DKA1040 | | V Style Standard Frame Angle 75° Aluminium | (NC 1475, 1490, 1495) | ◆ |
| DKA1050 | | Tapered Style Standard Frame Angle 75° Aluminium | (NC 1475, 1490, 1495, 1740, 1745, 1750) | ○ |
| DKA1060 | | V Style +5cm Frame Angle 75° Aluminium. Frame to accommodate a longer seat depth and seat height for taller clients | (NC 1410, 1420, 1435, 1740, 1745, 1750, 1755, 1760) | ○ |
| DKA1070 | | V Style Standard Frame Angle 90° Aluminium | (NC 1475, 1490, 1495) | ○ |
| DKA1080 | | Tapered Style Standard Frame Angle 90° Aluminium. Taper before the frame bend creating a parallel style front frame with a wider footrest space | (NC 1475, 1490, 1495, 1740, 1745, 1750, 5300, 5350) | ○ |
| DKA1090 | | V Style +5cm Frame Angle 90° Aluminium. Frame to accommodate a longer seat depth and seat height for taller clients | (NC 1410, 1420, 1435, 1740, 1745, 1750, 1755, 1760) | ○ |

| SEAT WIDTH* | | | | |
|--------------------|--|--------|-----------|---|
| DKA1310 | | 320 mm | (NC 3420) | ○ |
| DKA1315 | | 340 mm | (NC 3420) | ○ |
| DKA1320 | | 360 mm | | ○ |
| DKA1325 | | 380 mm | | ○ |
| DKA1330 | | 400 mm | | ◆ |
| DKA1335 | | 420 mm | | ○ |
| DKA1340 | | 440 mm | | ○ |
| DKA1345 | | 460 mm | | ○ |
| DKA1350 | | 480 mm | | ○ |
| DKA1355 | | 500 mm | | ○ |

| SEAT DEPTH* | | | | |
|--------------------|--|--------|-----------------|---|
| DKA1410 | | 350 mm | (NC 1500, 1510) | ○ |
| DKA1420 | | 375 mm | | ○ |
| DKA1435 | | 400 mm | | ◆ |
| DKA1450 | | 425 mm | | ○ |
| DKA1465 | | 450 mm | | ○ |
| DKA1475 | | 475 mm | | ○ |
| DKA1490 | | 500 mm | | ○ |
| DKA1495 | | 525 mm | | ○ |

| SEAT POSITION. Seat module can be moved forward and backward along the frame to achieve a more active or stable rear wheel position* | | | | |
|---|--|--|--|---|
| DKA1500 | | Position 1 = 140 mm MOST ACTIVE POSITION | | ○ |
| DKA1510 | | Position 2 = 125 mm | | ○ |
| DKA1520 | | Position 3 = 110 mm | | ◆ |
| DKA1530 | | Position 4 = 95 mm | | ○ |
| DKA1540 | | Position 5 = 80 mm LEAST ACTIVE POSITION | | ○ |

BACK HEIGHT. Transportation kit is only possible with backrest heights of 405mm and above*

| | | | | |
|---------|--|---|-----------------|---|
| DKA1600 | | Back Height for 3rd party backrest. Final height depends on the height of the positioning backrest, select if MatrX/Fixed backrest is required. | (NC 2100, 2220) | ○ |
| DKA1605 | | 270 mm | (NC 2150, 2240) | ○ |
| DKA1610 | | 285 mm | (NC 2150, 2240) | ○ |
| DKA1615 | | 300 mm | (NC 2150, 2240) | ○ |
| DKA1620 | | 315 mm | (NC 2150, 2240) | ○ |
| DKA1625 | | 330 mm | (NC 2150, 2240) | ○ |
| DKA1630 | | 345 mm | (NC 2150, 2240) | ○ |
| DKA1635 | | 360 mm | (NC 2150, 2240) | ○ |
| DKA1640 | | 375 mm | (NC 2150, 2240) | ○ |
| DKA1645 | | 390 mm | (NC 2150, 2240) | ○ |
| DKA1650 | | 405 mm | (NC 2150, 2240) | ○ |
| DKA1655 | | 420 mm | (NC 2150, 2240) | ○ |
| DKA1660 | | 435 mm | (NC 2150, 2240) | ○ |
| DKA1665 | | 450 mm | (NC 2150, 2240) | ○ |
| DKA1670 | | 465 mm | (NC 2150, 2240) | ○ |
| DKA1675 | | 480 mm | (NC 2150, 2240) | ○ |

FRONT SEAT HEIGHT. For 75° frames, the maximum lower leg must be front seat height minus 30mm and minus 40mm for 90° frames*

| | | | | |
|---------|--|--------|-----------------|---|
| DKA1740 | | 450 mm | (NC 5300, 5350) | ○ |
| DKA1745 | | 460 mm | (NC 5300, 5350) | ○ |
| DKA1750 | | 470 mm | | ○ |
| DKA1755 | | 480 mm | | ○ |
| DKA1760 | | 490 mm | | ○ |
| DKA1765 | | 500 mm | | ◆ |
| DKA1770 | | 510 mm | | ○ |
| DKA1775 | | 520 mm | | ○ |

REAR SEAT HEIGHT*

| | | | | |
|---------|--|--------|-----------|---|
| DKA1815 | | 390 mm | (NC 6040) | ○ |
| DKA1820 | | 400 mm | (NC 6040) | ○ |
| DKA1825 | | 410 mm | | ○ |
| DKA1830 | | 420 mm | | ○ |
| DKA1835 | | 430 mm | | ○ |
| DKA1840 | | 440 mm | | ◆ |
| DKA1845 | | 450 mm | | ○ |
| DKA1850 | | 460 mm | | ○ |
| DKA1855 | | 470 mm | | ○ |
| DKA1860 | | 480 mm | | ○ |
| DKA1865 | | 490 mm | (NC 6410) | ○ |

| LOWER LEG LENGTH. Cushion height must be deducted* | | | |
|---|--|--|---|
| DKA1918 | | 290 mm | ○ |
| DKA1920 | | 300 mm | ○ |
| DKA1922 | | 310 mm | ○ |
| DKA1924 | | 320 mm | ○ |
| DKA1926 | | 330 mm | ○ |
| DKA1928 | | 340 mm | ○ |
| DKA1930 | | 350 mm | ○ |
| DKA1932 | | 360 mm | ○ |
| DKA1934 | | 370 mm | ○ |
| DKA1936 | | 380 mm | ○ |
| DKA1938 | | 390 mm | ○ |
| DKA1940 | | 400 mm | ○ |
| DKA1942 | | 410 mm | ○ |
| DKA1944 | | 420 mm | ○ |
| DKA1946 | | 430 mm | ◆ |
| DKA1948 | | 440 mm | ○ |
| DKA1950 | | 450 mm | ○ |
| DKA1952 | | 460 mm | ○ |
| DKA1954 | | 470 mm | ○ |
| DKA1956 | | 480 mm | ○ |
| BACKREST MATERIAL* | | | |
| DKA2000 | | Aluminium | ◆ |
| BACKREST TYPE* | | | |
| DKA2100 | | Standard Backrest with Backbrace, lockdown folding for easy transportation (NC 2240) | ◆ |
| DKA2150 | | Backrest with Lower Backbrace for use with Matrix fixed back. Backrest Tube Height is 285 mm, lockdown folding for easy transportation (NC 2220) | ○ |
| BACKREST UPHOLSTERY* | | | |
| DKA2220 | | Tension Adjustable Backrest, with Premium Black Nylon and Reflective Trim Upholstery Cover, Back Pocket/hook and loop fastener | ◆ |
| DKA2240 | | Without Backrest Upholstery Select if fitting a Matrix/Fixed backrest, no upholstery required | ○ |
| PUSH HANDLES* | | | |
| DKA2400 | | Without Push Handles | ◆ |
| BACKREST ANGLE* | | | |
| DKA2520 | | 78° angle from the seat | ○ |
| DKA2540 | | 82° angle from the seat | ○ |
| DKA2560 | | 86° angle from the seat | ○ |
| DKA2580 | | 90° angle from the seat | ◆ |
| DKA2590 | | 94° angle from the seat | ○ |
| SEAT UPHOLSTERY* | | | |
| DKA2600 | | Standard Nylon Seat Upholstery with one 5cm central band (Loop) | ◆ |
| FOOTREST* | | | |
| DKA3200 | | Footrest Tube, Rigid - Aluminium Hydroformed Profile | ◆ |

| CALF STRAP* | | | |
|--------------------------|--------------------------|---|---|
| DKA3300 | <input type="checkbox"/> | Calf strap mounted on the frame with hook and loop | ◆ |
| DKA3305 | <input type="checkbox"/> | Calf strap mounted with bracket support | ○ |
| CLOTHESGUARDS/MUDGUARDS* | | | |
| DKA3420 | <input type="checkbox"/> | Carbon Clothesguards, Flip up, Fold in | ◆ |
| FRAME COLOUR* | | | |
| DKA4000 | <input type="checkbox"/> | Space Black (RAL 9005) | ○ |
| DKA4010 | <input type="checkbox"/> | Anthracite Grey (RAL 7016) | ○ |
| DKA4020 | <input type="checkbox"/> | Pure White (RAL 9016) (NC 4220) | ○ |
| DKA4030 | <input type="checkbox"/> | Carmine Red (RAL 3002) | ○ |
| DKA4040 | <input type="checkbox"/> | Azure Blue (RAL 5017) | ○ |
| FRAME DECAL TYPE* | | | |
| DKA4195 | <input type="checkbox"/> | Tag (Small Kuschall brand transfer on front of frame) (NC 4200, 4210, 4220, 4230, 4240, 4250, 2222) | ◆ |
| DKA4197 | <input type="checkbox"/> | Without Tag (Kuschall branding only) (NC 4200, 4210, 4220, 4230, 4240, 4250, 4300, 4310) | ○ |
| BOLD DECAL COLOUR* | | | |
| DKA4200 | <input type="checkbox"/> | Black | ○ |
| DKA4210 | <input type="checkbox"/> | Dark Grey | ○ |
| DKA4220 | <input type="checkbox"/> | Light Grey | ○ |
| DKA4230 | <input type="checkbox"/> | Red | ○ |
| DKA4240 | <input type="checkbox"/> | Blue | ○ |
| DKA4250 | <input type="checkbox"/> | Gold | ○ |
| DKA1111 | <input type="checkbox"/> | No bold colour, select only with Option DKA4195 Tag or DKA4197 Without Tag | ○ |
| TAG COLOUR* | | | |
| DKA4300 | <input type="checkbox"/> | Black | ○ |
| DKA4310 | <input type="checkbox"/> | Red | ○ |
| DKA2222 | <input type="checkbox"/> | No tag colour, select only with Option DKA4190 Bold Decal or DKA4197 Without Tag (NC 4195) | ○ |
| FRONT CASTORS* | | | |
| DKA5000 | <input type="checkbox"/> | 75mm x 25mm (3" x 1") Supersport | ○ |
| DKA5050 | <input type="checkbox"/> | 100mm x 30mm (4" x 1.2") Everyday Black hub and tyre | ○ |
| DKA5100 | <input type="checkbox"/> | 100mm x 30mm (4" x 1.2") Comfort Black hub and tyre | ○ |
| DKA5300 | <input type="checkbox"/> | 125mm x 30mm (5" x 1.2") Everyday Black hub and tyre | ○ |
| DKA5350 | <input type="checkbox"/> | 125mm x 30mm (5" x 1.2") Comfort Black hub and tyre | ○ |
| REAR WHEEL POSITION* | | | |
| DKA6020 | <input type="checkbox"/> | Standard, Fixed Rear Wheel Position (NC 6100, 6110, 6120, 6130, 6140) | ◆ |
| DKA6040 | <input type="checkbox"/> | Rear Wheel Extension (5 Positions) | ◆ |
| CAMBER AXLE TYPE* | | | |
| DKA6060 | <input type="checkbox"/> | Standard Aluminium Axle | ◆ |
| POSITION OF AXLE | | | |
| DKA6100 | <input type="checkbox"/> | 1st Position from the front. If 2nd axle is selected, 1st axle must be placed in 1st position and 2nd in the 5th position | ○ |
| DKA6110 | <input type="checkbox"/> | 2nd Position from the front | ○ |
| DKA6120 | <input type="checkbox"/> | 3rd Position from the front | ○ |
| DKA6130 | <input type="checkbox"/> | 4th Position from the front | ○ |
| DKA6140 | <input type="checkbox"/> | 5th Position from the front | ○ |

| CAMBER ANGLE* | | | |
|------------------|--------------------------|--|---|
| DKA6200 | <input type="checkbox"/> | 0° Camber | ◆ |
| REAR WHEEL SIZE* | | | |
| DKA6410 | <input type="checkbox"/> | 24" | ◆ |
| REAR WHEEL TYPE* | | | |
| DKA6510 | <input type="checkbox"/> | 24", 26" Kuschall Starec Wheel, 24" Black Cross Spokes, Rim & Hub. 26" Silver Cross-Spokes, Colourless hub and Black rim | ◆ |
| REAR WHEEL TYRE* | | | |
| DKA6600 | <input type="checkbox"/> | 24" Maxxis Lite guard Pneumatic tyre. Herringbone tread, Silvershield puncture protection | ◆ |
| HANDRIM* | | | |
| DKA6780 | <input type="checkbox"/> | 24", 25", 26" Light Aluminium Grey, long tab with 2 holes for 3 cm or 4cm mountings, mounted at 3cm distance | ◆ |
| BRAKES* | | | |
| DKA7210 | <input type="checkbox"/> | Standard Push to Lock Brakes, Lever Bent | ◆ |
| OTHER OPTIONS | | | |
| DKA8550 | <input type="checkbox"/> | Tool Kit | ◆ |

| | |
|--------------------|--|
| TOTAL PRICE | |
|--------------------|--|

Non Conforming Options:

| | |
|---------|--|
| DKA1410 | 350 mm |
| DKA1420 | 375 mm |
| DKA1435 | 400 mm |
| DKA1475 | 475 mm |
| DKA1490 | 500 mm |
| DKA1495 | 525 mm |
| DKA1500 | Position 1 = 140 mm MOST ACTIVE POSITION |
| DKA1510 | Position 2 = 125 mm |
| DKA1740 | 450 mm |
| DKA1745 | 460 mm |
| DKA1750 | 470 mm |
| DKA1755 | 480 mm |
| DKA1760 | 490 mm |
| DKA2100 | Standard Backrest with Backbrace, lockdown folding for easy transportation |
| DKA2150 | Backrest with Lower Backbrace for use with Matrix fixed back. Backrest Tube Height is 285 mm, lockdown folding for easy transportation |
| DKA2220 | Tension Adjustable Backrest, with Premium Black Nylon and Reflective Trim Upholstery Cover, Back Pocket/hook and loop fastener |
| DKA2240 | Without Backrest Upholstery Select if fitting a Matrix/Fixed backrest, no upholstery required |
| DKA3420 | Carbon Clothesguards, Flip up, Fold in |
| DKA4195 | Tag (Small Kuschall brand transfer on front of frame) |
| DKA4200 | Black |
| DKA4210 | Dark Grey |
| DKA4220 | Light Grey |
| DKA4230 | Red |
| DKA4240 | Blue |
| DKA4250 | Gold |
| DKA4300 | Black |
| DKA4310 | Red |
| DKA2222 | No tag colour, select only with Option DKA4190 Bold Decal or DKA4197 Without Tag |
| DKA5300 | 125mm x 30mm (5" x 1.2") Everyday Black hub and tyre |
| DKA5350 | 125mm x 30mm (5" x 1.2") Comfort Black hub and tyre |
| DKA6040 | Rear Wheel Extension (5 Positions) |
| DKA6100 | 1st Position from the front. If 2nd axle is selected, 1st axle must be placed in 1st position and 2nd in the 5th position |
| DKA6110 | 2nd Position from the front |
| DKA6120 | 3rd Position from the front |
| DKA6130 | 4th Position from the front |
| DKA6140 | 5th Position from the front |
| DKA6410 | 24" |



For a more complex customer requirement
please refer to our Invacare Unique special order form.
email specialsuk@invacare.com

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____

Description of special request (please include relevant option code or part number):

