



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Motability Retail Prescription Form

OCTOBER 2024
Customer Ref:

Customer Ref*, Contact No*, Delivery Address*, Account No*, Order No*, Order Qty*, Order date*

Esprit Action Junior M120 706812 - Weekly
Price £38.59

Standard Features

- Max User Weight : 80kg (12.6 Stone)
Single crossbar (Width Adjustable)
Seat widths 305mm to 380mm 12" - 15"
Seat depth 325mm (13") to 375mm (15") adjustable (set at 350mm)
Seat height 460mm (18")
Nylon seat upholstery
FIXED 7° BACKREST 430mm backrest
Standard nylon sling type backrest upholstery
80° swing in / swing out legrest hangers
Flip up angle and depth adjustable footplates
Height adjustable, removable arm armrests
Choice of 3 standard colours
Anti tippers With jack up function
Buckle belt



Please make your selection with a cross in the box
Free of charge
standard feature
POA
special request

Customer weight = KG

Table with columns: Code, Description, Price. Rows include BASE PRODUCT*, CHAIR TYPE*, CHASSIS*, and BRACKET TYPE* with various options and prices.

USER REMOTE*			
SPD0400		USER REMOTE CONTROL RH	○
SPD0410		USER REMOTE CONTROL LH	○
SPEED*			
SPD0500		6 KM/H	◆
BATTERY*			
SPD0600		BATTERY 7,65AH/36V (NON SPILLABLE LEAD ACID BATTERY)	◆
POWER OPTIONS			
SPD0610		ADDITIONAL BATTERY 7,65AH/36V	£271
SEAT WIDTH*			
SPD1000		305 MM SEAT WIDTH	○
SPD1010		330 MM SEAT WIDTH	○
SPD1020		355 MM SEAT WIDTH	○
SPD1030		380 MM SEAT WIDTH	○
SEAT DEPTH*			
SPD1200		325 TO 375 MM ADJUSTABLE SEAT DEPTH SET AT 350 MM	◆
SEAT HEIGHT*			
SPD1300		460 MM SEAT HEIGHT	○
SPD1310		485 MM SEAT HEIGHT	○
SEAT UPHOLSTERY TYPE*			
SPD1400		NYLON WITH LOOP BANDS STRIPS SEAT UPHOLSTERY, WIDTH & DEPTH ADJUSTABLE (SINGLE CROSSBAR ONLY)	◆
BACKREST TYPE*			
SPD2000		FIXED 7° BACKREST 2 HEIGHTS 400/430 MM (DELIVERED 430 MM)	◆
SPD2020		FOLDING BACKREST 430 MM HEIGHT	£32
SPD2040		ANGLE ADJUSTABLE AND FOLD DOWN (-15° TO + 15°) FIXED 7° HEIGHT 425-450 MM / 17"-18" SET AT 425 MM / 17" (2 POSITIONS; 25 MM / 1" STEPS THROUGH UPHOLSTERY)	£97
SPD2050		RECLINING/FOLDING (GAS SPRING) 0/30° BACKREST 3 HEIGHTS 390/410/430 MM (DELIVERED 430 MM)	£64
SPD2060		RECLINING (RATCHET) 0/30° BACKREST 3 HEIGHTS 390/410/430 MM (DELIVERED 430 MM)	£124
SPD2080		RECLINING (GAS SPRING) 0/30° BACKREST 3 HEIGHTS 390/410/430 MM (DELIVERED 430 MM)	£51
BACKREST UPHOLSTERY TYPE*			
SPD2300		LUMBAR SUPPORT (+1 STRAP) BI MATERIAL PADDED COVER	◆
SPD2320		CONTOUR ADJUSTABLE (FULL STRAPS) BI MATERIAL COVER	○
HANGERS RH*			
SPD3000		80° SHORT, UL FROM 200 TO 330MM	◆
HANGERS LH*			
SPD3100		80° SHORT, UL FROM 200 TO 330MM	◆
FOOTPLATES RH*			
SPD3200		FLIP UP RH FOOTPLATE INCLUDING HEEL STRAP	○
SPD3210		FLIP UP ANGLE AND DEPTH ADJUSTABLE RH FOOTPLATE, INCLUDING HEEL STRAP	£42
SPD3220		FLIP UP ONE PIECE FOOTPLATE, DEPTH & ANGLE ADJUSTABLE (NC 3310)	£76

FOOTPLATES LH*			
SPD3300	<input type="checkbox"/>	FLIP UP LH FOOTPLATE, INCLUDING HEEL STRAP	○
SPD3310	<input type="checkbox"/>	FLIP UP ANGLE AND DEPTH ADJUSTABLE LH FOOTPLATE, INCLUDING HEEL STRAP	£42
ARMRESTS RH*			
SPD4000	<input type="checkbox"/>	REMOVABLE RH DESK, ADJ IN HEIGHT (235-280 for FSTF 460 & 210-280 for FSTF 485MM) DEPTH (+35 MM FRONTWARD) & WIDTH (+25MM OUTWARDS)	○
SPD4010	<input type="checkbox"/>	REMOVABLE RH FULL LENGTH, ADJ IN HEIGHT (210 TO 280 MM) DEPTH (+35 MM BACKWARD) & WIDTH (+25MM OUTWARDS)	◆
ARMRESTS LH*			
SPD4100	<input type="checkbox"/>	REMOVABLE LH DESK, ADJ IN HEIGHT (235-280 for FSTF 460 & 210-280 for FSTF 485MM) DEPTH (+35 MM FRONTWARD) & WIDTH (+25MM OUTWARDS)	○
SPD4110	<input type="checkbox"/>	REMOVABLE LH FULL LENGTH, ADJ IN HEIGHT (210 TO 280 MM) DEPTH (35 MM BACKWARD) & WIDTH (+25MM OUTWARDS)	◆
SEAT CUSHION*			
SPD5000	<input type="checkbox"/>	ESSENTIAL CUSHION 375 MM DEPTH	○
SEATING OPTIONS			
SPD5110	<input type="checkbox"/>	HEADREST BLACK PU TR26	£153
SPD5120	<input type="checkbox"/>	HEADREST WITH CHEEK SUPPORT BLACK PU TR26 (INCLUDING BACKBRACE EXCEPT FOR LAGUNA)	£177
SPD5140	<input type="checkbox"/>	LATERAL POSITIONING PADS (PAIR)	£59
SPD5150	<input type="checkbox"/>	ADJUSTABLE CALF STRAP	○
BACK BRACE			
SPD5160	<input type="checkbox"/>	BACK BRACE WIDTH ADJUSTABLE	○
COLOURS*			
SPD6000	<input type="checkbox"/>	ANTHRACITE BLACK	○
SPD6010	<input type="checkbox"/>	DARK PURPLE	○
SPD6020	<input type="checkbox"/>	WATER BLUE	○
ACCENT COLOUR*			
SPD6100	<input type="checkbox"/>	BLACK	○
SPD6120	<input type="checkbox"/>	LIME	○
FRONT WHEELS*			
SPD7010	<input type="checkbox"/>	200 X 50 MM (8" X 2") BLACK PNEUMATIC TYRE	◆
REAR WHEELS*			
SPD7110	<input type="checkbox"/>	355 X 70 MM (14" X 27") BLACK PNEUMATIC TYRE	◆
ANTI TIPPER*			
SPD8010	<input type="checkbox"/>	REAR ANTI-TIPPERS REMOVABLE WITH CRUTCHES (E25), PAIR	◆
LAP BELT*			
SPD8020	<input type="checkbox"/>	POSTURE BELT	◆
TRANSPORTATION BRACKET*			
SPD8030	<input type="checkbox"/>	TRANSPORTATION BRACKET	◆

TOTAL PRICE	
--------------------	--

Non Conforming Options:

SPD2300 LUMBAR SUPPORT (+1 STRAP) BI MATERIAL PADDED COVER

SPD3310 FLIP UP ANGLE AND DEPTH ADJUSTABLE LH FOOTPLATE, INCLUDING HEEL STRAP



For a more complex customer requirement
please refer to our Invacare Unique special order form.
email specialsuk@invacare.com

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):