



Yes, you can:

# Motability Retail Prescription Form

Tel: 01656 776222  
Email: ordersuk@invacare.com  
Online Spares available at: www.invacare.co.uk

OCTOBER 2024  
Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

**Storm4 Modulite M31 1001116 - Weekly  
Price £59.24**

### Standard Features

- Max user weight 150Kg (23.5 stone)
- Crash tested ISO 7176 -Part 19
- Modulite seating system
- Flex 3 Laguna Backrest
- Standard T-Armrests, detachable with Long armpad
- Electrical Tilt 0° to 25°
- Manual settable Recliner (0° - 30° in 3.75° steps)
- VARI F with Fixed Footplate
- LinX Control REM216 with swing away bracket
- Pneumatic Castors
- Pneumatic Drive Wheels
- 12 km/h (8mph) Motors, LED light system with indicators
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)
- Anti-Tippers



**Standard Base price: £8,334**

Please make your selection with a cross in the box     Free of charge     standard feature     POA     special request

Customer weight =                      KG

Code	Description	
<b>CHAIR TYPE*</b>		
ATB0000	<input type="checkbox"/> Storm 4 Modulite	£8,334
<b>SEATING*</b>		
ATB2100	<input type="checkbox"/> Modulite Seating Unit (Flex3 seating unit)	◆
<b>SEATING OPTION*</b>		
ATB2111	<input type="checkbox"/> Telescopic Seat Frame with metal seat plate (with stepless adjustment on armrests, legrests and backrest)	◆

<b>SEAT WIDTH*</b>			
ATBXXXX	<input type="checkbox"/>	Seat width 330mm (13") seat width achieved on a 380mm seat frame (2121) with inrigged armrests (2183)	<input type="radio"/>
ATB2121	<input type="checkbox"/>	Set at 380mm (adjustable from 380mm to 430mm (15" - 17") with stepless adjustment on armrests, legrests and backrest)	<input type="radio"/>
ATB2122	<input type="checkbox"/>	Set at 430mm (adjustable from 430mm to 480mm (17" - 19) with stepless adjustment on armrests, legrests and backrest)	<input type="radio"/>
ATB2123	<input type="checkbox"/>	Set at 480mm (adjustable from 480mm to 530mm (19" - 21") with stepless adjustment on armrests, legrests and backrest)	<input type="radio"/>
ATB2124	<input type="checkbox"/>	Set at 530mm (adjustable from 530mm to 580mm (21" - 23") with stepless adjustment on armrests, legrests and backrest)	<input type="radio"/>
<b>SEAT DEPTH*</b>			
ATB2132	<input type="checkbox"/>	410mm (16") (adjustable 410mm - 510mm stepless adjustment) (NC 2196)	<input type="radio"/>
ATB2134	<input type="checkbox"/>	460mm (18") (adjustable 410mm - 510mm stepless adjustment)	<input type="radio"/>
ATB2136	<input type="checkbox"/>	510mm (20") (adjustable 410mm - 510mm stepless adjustment)	<input type="radio"/>
<b>SEAT CUSHION*</b>			
ATB2142	<input type="checkbox"/>	Matrx Flo-tech Lite	<input checked="" type="radio"/>
<b>BACKREST SUPPORT*</b>			
ATB2162	<input type="checkbox"/>	Backposts 540mm (Tall) must be selected for Flex3 Laguna backrest	<input checked="" type="radio"/>
<b>PUSH HANDLES FOR BACKPOSTS*</b>			
ATB2107	<input type="checkbox"/>	Backposts WITH push handles	<input checked="" type="radio"/>
<b>BACKREST UPHOLSTERY / CUSHION*</b>			
ATB2173	<input type="checkbox"/>	Flex3 adjustable Backrest Plates with Laguna Backrest cushion	<input checked="" type="radio"/>
<b>UPHOLSTERY TYPE*</b>			
ATB2152	<input type="checkbox"/>	Black multi stretch PU coated fabric	<input type="radio"/>
ATB2155	<input type="checkbox"/>	Black Spacer fabric	<input type="radio"/>
<b>ARMREST (PAIR)*</b>			
ATB2191	<input type="checkbox"/>	T-Armrest Standard, Range of adjustment from seat plate 245mm to 315mm	<input checked="" type="radio"/>
<b>ARMPAD (PAIR)*</b>			
ATB2196	<input type="checkbox"/>	Long Armpad (375mm) (NC 2132)	<input checked="" type="radio"/>
<b>HEADREST / NECKREST - Adequate headrest option must be selected*</b>			
ATB1350	<input type="checkbox"/>	neckrest: height and depth adjustable	<input checked="" type="radio"/>
<b>SEAT HEIGHT AND SEAT RISER*</b>			
ATB1023	<input type="checkbox"/>	Fixed Seat Height 450mm (17.75")	<input type="radio"/>
ATB1025	<input type="checkbox"/>	Fixed Seat Height 480mm (19")	<input type="radio"/>
<b>SEAT TILT*</b>			
ATB1014	<input type="checkbox"/>	Electrical Tilt 0° to 25°	<input checked="" type="radio"/>
<b>BACKREST RECLINER*</b>			
ATB0975	<input type="checkbox"/>	Manual Settable Recliner (0° - 30°, in 3.75° steps) Foldable	<input checked="" type="radio"/>
<b>FOOTREST HANGER TYPE (PAIR)*</b>			
ATB0503	<input type="checkbox"/>	Vari F : Angle adjustable with tools 0° - 70°, length 290mm - 460mm	<input type="radio"/>
<b>FOOTPLATE TYPE*</b>			
ATB0512	<input type="checkbox"/>	Fixed Footplate	<input type="radio"/>

<b>CONTROLS*</b>			
ATB2204	<input type="checkbox"/>	REM216 Drive, Actuators & Lights	◆
<b>CONTROL BRACKETRY*</b>			
ATB1707	<input type="checkbox"/>	Swing Away RH Bracket	○
ATB1708	<input type="checkbox"/>	Swing Away LH Bracket	○
<b>SHROUD COLOUR*</b>			
ATB0202	<input type="checkbox"/>	Blue	○
ATB0228	<input type="checkbox"/>	Black	○
ATB0220	<input type="checkbox"/>	Champagne	○
ATB0204	<input type="checkbox"/>	Transparent Cherry Red	○
<b>FRONT WHEEL*</b>			
ATB0851	<input type="checkbox"/>	3.00 - 6 Grey Pneumatic	◆
<b>REAR WHEEL*</b>			
ATB0771	<input type="checkbox"/>	3.00 - 8 Grey Pneumatic	◆
<b>MOTOR TYPE*</b>			
ATB1909	<input type="checkbox"/>	12 km/h (8mph) (includes LED light system with indicators)	(NC 1437) ◆
<b>BATTERY*</b>			
ATB1364	<input type="checkbox"/>	70 Ah Gel	◆
<b>SUSPENSION*</b>			
ATB1240	<input type="checkbox"/>	Suspension low (up to 80kg user weight)	○
ATB1241	<input type="checkbox"/>	Suspension middle (from 70kg up to 100kg user weight)	○
ATB1242	<input type="checkbox"/>	Suspension hard (from 90kg up to 150kg user weight)	○
<b>CHARGER*</b>			
ATB1343	<input type="checkbox"/>	OptiCharge Charger 8A (200 - 250V)	◆
<b>LIGHTS*</b>			
ATB1437	<input type="checkbox"/>	LED Light system with indicators	(NC 1909) ◆
<b>ANTI TIPPER*</b>			
ATB1205	<input type="checkbox"/>	Anti-tipper (standard version)	◆
<b>TRANSPORTATION KIT*</b>			
ATB1494	<input type="checkbox"/>	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory, see options below)	◆
<b>FIXATION BELT*</b>			
ATB2113	<input type="checkbox"/>	Lap belt with auto lock (flexible mount in side rails)	◆
<b>REAR VIEW MIRROR</b>			
ATB1415	<input type="checkbox"/>	Rear view Mirror RH - Must be selected with 10 km/h Motor	○
ATB1339	<input type="checkbox"/>	Rear view Mirror LH - Must be selected with 10 km/h Motor	○

<b>TOTAL PRICE</b>	
--------------------	--

**Non Conforming Options:**

ATB2132 410mm (16") (adjustable 410mm - 510mm stepless adjustment)

ATB2196 Long Armpad (375mm)

ATB1909 12 km/h (8mph) (includes LED light system with indicators)

ATB1437 LED Light system with indicators



For a more complex customer requirement  
please refer to our Invacare Unique special order form.  
email [specialsuk@invacare.com](mailto:specialsuk@invacare.com)

**Special Instructions**

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_

Description of special request (please include relevant option code or part number):  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_