



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Motability Retail Prescription Form

JANUARY 2025
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

**TDX SP2 HD M16 1001587 - Weekly Price
£50.90**

Standard Features

- Max user weight 180kg (28 stone)
- Crash tested ISO 7176 -Part 19
- Modulite Seating System
- Matrix Flo-tech Cushion
- Flex 3 adjustable Backrest plates with Laguna backrest cushion
- Standard T-armrests, detachable with long armpad
- Neckrest
- Powereed tilt 20°
- Manual settable backrest Recline 0° - 30° in 3.75° steps
- 80° hangers with fixed footplate
- LiNX control REM211 swing away bracket
- Solid Black tyres
- 6km/h (4mph) 4 pole motors
- 4-Point tie-down system (Q'Straint QLK-150 compatible)
- Lapbelt with autolock (adjustable both sides)



Standard Base price: £6,517

Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

Code	Description	
CHAIR TYPE*		
ARJ0000	TDX SP2 Heavy Duty	◆
User Weight*		
ARJ1918	User weight less 100 kg	○
ARJ1919	User weight 101 to 125 kg	○
ARJ1921	User weight 126 to 136 kg	○
ARJ1922	User weight 137 to 150 kg (NC 1092)	○
ARJ1923	User weight 151 to 160 kg (NC 1092)	○
ARJ1924	User weight 161 to 180 kg (NC 1092)	○

BASE TYPE*			
ARJ2100	<input type="checkbox"/>	Modulite Seating unit	£6,517
SEAT WIDTH*			
ARJ2121	<input type="checkbox"/>	Set at 380mm (15") (adjustable from 380mm to 430mm (15" - 17") with stepless adjustment on armrests, legrests and backrest) (NC 2138 + 2142)	<input type="radio"/>
ARJ2122	<input type="checkbox"/>	Set at 430mm (17") (adjustable from 430mm to 480mm (17" - 19") with stepless adjustment on armrests, legrests and backrest) (NC 2138 + 2142)	<input type="radio"/>
ARJ2123	<input type="checkbox"/>	Set at 480mm (19") (adjustable from 480mm to 530mm (19" - 21") with stepless adjustment on armrests, legrests and backrest) (NC 2138 + 2142)	<input type="radio"/>
ARJ2124	<input type="checkbox"/>	Set at 530mm (21") (adjustable from 530mm to 580mm (21" - 23") with stepless adjustment on armrests, legrests and backrest)	<input type="radio"/>
SEAT DEPTH RAIL LENGTH*			
ARJ2129	<input type="checkbox"/>	Standard Depth Package (16" - 20", 410mm - 510mm simple stepless adjustment from the side) (NC 2138, 2196)	<input type="radio"/>
ARJ2130	<input type="checkbox"/>	Extended Depth Package (20" - 22", 510mm - 560mm simple stepless adjustment from the side) (NC 2132, 2134)	<input type="radio"/>
SEAT DEPTH SETTING*			
ARJ2132	<input type="checkbox"/>	410mm (16") (410mm - 510mm simple stepless adjustment from the side) (NC 2196)	<input type="radio"/>
ARJ2134	<input type="checkbox"/>	460mm (18") (410mm - 510mm simple stepless adjustment from the side)	<input type="radio"/>
ARJ2136	<input type="checkbox"/>	510mm (20") (410mm - 510mm simple stepless adjustment from the side)	<input type="radio"/>
ARJ2138	<input type="checkbox"/>	560mm (22") (510mm - 560mm simple stepless adjustment from the side) (Cushion (option 2142) is only available with Seat Width (option 2124) 530mm)	<input type="radio"/>
SEAT CUSHION*			
ARJ2142	<input type="checkbox"/>	Matrx Flo-tech Lite	<input checked="" type="checkbox"/>
BACKREST SUPPORT*			
ARJ2162	<input type="checkbox"/>	Backposts 540mm (Tall) must be selected for Flex3 adjustable backrest plates	<input checked="" type="checkbox"/>
PUSH HANDLES FOR BACKPOSTS*			
ARJ2107	<input type="checkbox"/>	Backposts WITH push handles	<input checked="" type="checkbox"/>
BACKREST UPHOLSTERY / CUSHION*			
ARJ2173	<input type="checkbox"/>	Flex3 adjustable Backrest Plates with Laguna Backrest cushion	<input checked="" type="checkbox"/>
UPHOLSTERY COLOUR*			
ARJ2152	<input type="checkbox"/>	Black Black multi stretch PU coated fabric Upholstery	<input type="radio"/>
ARJ2155	<input type="checkbox"/>	Black Spacer fabric	<input type="radio"/>
ARMREST*			
ARJ2191	<input type="checkbox"/>	T-Armrest Standard, Range of adjustment from seat plate 245mm to 315mm	<input checked="" type="checkbox"/>
ARMPAD*			
ARJ2196	<input type="checkbox"/>	Long Armpad (375mm) (NC 2132)	<input checked="" type="checkbox"/>
HEADREST / NECKREST - Adequate headrest option must be selected*			
ARJ1350	<input type="checkbox"/>	Neckrest height and depth adjustable	<input checked="" type="checkbox"/>
SEAT LIFTER / SEAT TILT*			
ARJ1092	<input type="checkbox"/>	Electric Tilt (0° - 20°) Seat Height 450mm	<input checked="" type="checkbox"/>
ARJ1095	<input type="checkbox"/>	Electric HD Tilt (0° - 20°) Seat Height 450mm (upto 180kg Max User Weight)	<input type="radio"/>

SEAT HEIGHT SETTING FOR 20 DEGREES SEAT TILT			
ARJ0354	<input type="checkbox"/>	Seat height set at 450mm (17.75") (from floor to front of seat plate at 0° tilt)	<input type="radio"/>
ARJ0364	<input type="checkbox"/>	Seat height set at 475mm (18.75") (from floor to front of seat plate at 0° tilt)	<input type="radio"/>
BACKREST RECLINER*			
ARJ0975	<input type="checkbox"/>	Manual Settable Recliner (0° - 30°, in 3.75° steps) Foldable	<input checked="" type="checkbox"/>
FOOTREST HANGER TYPE*			
ARJ0526	<input type="checkbox"/>	Side Mounted : Standard 80° hangers	<input checked="" type="checkbox"/>
FOOTPLATE TYPE*			
ARJ0512	<input type="checkbox"/>	Fixed Footplate	<input checked="" type="checkbox"/>
CONTROLS*			
ARJ2203	<input type="checkbox"/>	REM211 Drive & Actuators	<input checked="" type="checkbox"/>
CONTROL BRACKETRY*			
ARJ1707	<input type="checkbox"/>	Swing Away RH Bracket	<input type="radio"/>
ARJ1708	<input type="checkbox"/>	Swing Away LH Bracket	<input type="radio"/>
SHROUD COLOUR*			
ARJ2640	<input type="checkbox"/>	Obsidian Black Frame Shrouds	<input type="radio"/>
ARJ2641	<input type="checkbox"/>	Ruby Red Frame Shrouds	<input type="radio"/>
ARJ2642	<input type="checkbox"/>	Purple Passion Frame Shrouds	<input type="radio"/>
ARJ2643	<input type="checkbox"/>	Peacock Blue Frame Shrouds	<input type="radio"/>
ARJ2644	<input type="checkbox"/>	Arctic Blue Frame Shrouds	<input type="radio"/>
RIM COLOUR*			
ARJ2650	<input type="checkbox"/>	Obsidian Black Frame Shrouds	<input type="radio"/>
ARJ2651	<input type="checkbox"/>	Ruby Red Frame Shrouds	<input type="radio"/>
ARJ2652	<input type="checkbox"/>	Purple Passion Frame Shrouds	<input type="radio"/>
ARJ2653	<input type="checkbox"/>	Peacock Blue Frame Shrouds	<input type="radio"/>
ARJ2654	<input type="checkbox"/>	Arctic Blue Frame Shrouds	<input type="radio"/>
FORKS*			
ARJ2632	<input type="checkbox"/>	Double sided fork	<input checked="" type="checkbox"/>
CENTRE WHEELS (DRIVING WHEELS)*			
ARJ0705	<input type="checkbox"/>	DW 14" (355 mm) x 3" (75 mm), Solid inlay, dual density foam insert	<input checked="" type="checkbox"/>
MOTOR TYPE*			
ARJ1282	<input type="checkbox"/>	6 km/h (4mph) (4 pole)	<input checked="" type="checkbox"/>
BATTERY*			
ARJ1260	<input type="checkbox"/>	73.5 Ah (Gel)	<input checked="" type="checkbox"/>
CHARGER*			
ARJ1343	<input type="checkbox"/>	OptiCharge Charger 8A (200 - 250V)	<input checked="" type="checkbox"/>
TRANSPORTATION KIT*			
ARJ1494	<input type="checkbox"/>	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory, see options below)	<input checked="" type="checkbox"/>
FIXATION BELT*			
ARJ2113	<input type="checkbox"/>	Lapbelt (adjustable both sides) adjustable mounting on seat rail extrusion	<input checked="" type="checkbox"/>

TOTAL PRICE	
--------------------	--

Non Conforming Options:

ARJ2132	410mm (16") (410mm - 510mm simple stepless adjustment from the side)
ARJ2134	460mm (18") (410mm - 510mm simple stepless adjustment from the side)
ARJ2138	560mm (22") (510mm - 560mm simple stepless adjustment from the side) (Cushion (option 2142) is only available with Seat Width (option 2124) 530mm)
ARJ2196	Long Armpad (375mm)
ARJ1092	Electric Tilt (0°- 20°) Seat Height 450mm



For a more complex customer requirement
please refer to our Invacare Unique special order form.
email specialsuk@invacare.com

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____

Description of special request (please include relevant option code or part number):

