



Yes, you can:

Tel: 01656 776222  
Email: ordersuk@invacare.com  
Online Spares available at: www.invacare.co.uk

# Motability Retail Prescription Form

JANUARY 2025  
Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## AVIVA RX40 HD Modulite M28 707010 - Weekly Price £60.48

### Standard Features

- Max user weight 160kg (25 Stone)
- Crash tested ISO 7176 -Part 19
- Modulite adjustable seating system
- Flex3 Laguna backrest
- Neckrest depth adjustable
- Standard T-Armrests, detachable with Long armpad
- Powered "fixed pivot" Tilt 20°
- Manual settable backrest Recline (0°- 30° in 3.75° steps)
- 80° hangers with fixed footplate
- LiNX control REM216 with swing away bracket
- 9" Castor, puncture proof solid tyre
- Solid inlay black drive wheels
- 12 km/h Ultra High Speed Package combined with LED Light system with Indicators, 30 Lux headlights
- 73Ah batteries
- Anti tippers
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)

**Standard Base price: £8,315**



Please make your selection with a cross in the box     Free of charge     standard feature    POA    special request

Customer weight =                      KG

Code	Description	
<b>BASE MODEL*</b>		
JBH0002	AVIVA RX 40 HD Base Model Max user weight up to 160 kg HD Version, includes RX 40 Performance Package Wide Base only	£8,315
<b>CHASSIS WIDTH*</b>		
JBH2830	Chassis Wide (635 mm/25")	◆
<b>SEATING TYPE*</b>		
JBH2100	Invacare® Modulite Seat Unit	◆

<b>SEAT FRAME*</b>			
JBH2111	<input type="checkbox"/>	Telescopic Seat Frame and Plate with Side Rails width 380 to 530 mm Additional 50 mm by arms	◆
<b>SEAT WIDTH*</b>			
JBH2121	<input type="checkbox"/>	380 mm Seat Plate Can be reduced to 330 mm in combination with inrigged armrests	○
JBH2122	<input type="checkbox"/>	430 mm Seat Plate Can be reduced to 380 mm in combination with inrigged armrests	○
JBH2123	<input type="checkbox"/>	480 mm Seat Plate	○
JBH2124	<input type="checkbox"/>	530 mm Seat Plate	○
<b>SEAT DEPTH*</b>			
JBH2132	<input type="checkbox"/>	Seat Depth 410 mm (NC 2196)	○
JBH2134	<input type="checkbox"/>	Seat Depth 460 mm	○
JBH2136	<input type="checkbox"/>	Seat Depth 510 mm	○
<b>SEAT CUSHION*</b>			
JBH2142	<input type="checkbox"/>	Matrx Flo-tech Cushion Flo-tech Lite	◆
<b>BACKREST SUPPORT*</b>			
JBH2162	<input type="checkbox"/>	Backposts tall (540 mm) For Flex3 backrests/Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery Width adjustable cast bracket	◆
<b>BACKREST UPHOLSTERY / CUSHION*</b>			
JBH2173	<input type="checkbox"/>	Flex3 adjustable Backrest Plates with Laguna Backrest cushion Flex3 Backrest plate width and height adjustable 560 mm - 790 mm	◆
<b>UPHOLSTERY TYPE*</b>			
JBH2151	<input type="checkbox"/>	Black Fabric Upholstery (Nylon with option 2141; 2-way-stretch material with options 2142, 2143)	○
JBH2152	<input type="checkbox"/>	Black multi stretch PU coated fabric Upholstery	○
JBH2155	<input type="checkbox"/>	Black Spacer Fabric Upholstery	○
<b>ARMREST (Pair)*</b>			
JBH2191	<input type="checkbox"/>	T-Arm detachable, Standard Adjustment range from seat plate: 245 - 315 mm	◆
<b>ARMPAD (Pair)*</b>			
JBH2196	<input type="checkbox"/>	Long Armpad (375 mm) (NC 2132)	◆
<b>HEADREST / NECKREST - Adequate headrest option must be selected*</b>			
JBH1350	<input type="checkbox"/>	Neckrest depth adjustable	◆
<b>SEAT LIFTER / SEAT TILT*</b>			
JBH1092	<input type="checkbox"/>	Powered "fixed pivot" Tilt 20°	◆
<b>SEAT TO GROUND HEIGHT*</b>			
JBH2857	<input type="checkbox"/>	Seat-to-Floor Height 435 mm	○
JBH2858	<input type="checkbox"/>	Seat-to-Floor Height 460 mm	○
JBH2859	<input type="checkbox"/>	Seat-to-Floor Height 485 mm	○
<b>BACKREST RECLINER*</b>			
JBH0975	<input type="checkbox"/>	Manual Settable Recliner (in 3.75° steps) 90° - 120°	◆
<b>FOOTREST HANGER TYPE*</b>			
JBH0526	<input type="checkbox"/>	Standard hanger 80° pair Lower leg length: 350 - 450 mm	◆
<b>FOOTPLATE TYPE*</b>			
JBH0512	<input type="checkbox"/>	Fixed Footplate pair	◆
<b>CONTROLS*</b>			
JBH2204	<input type="checkbox"/>	Remote LiNX DLX-REM216 = Drive, Actuators & Lights 3 Drive functions and power module incl. 2 Actuators & Lights	◆

<b>CONTROL BRACKETRY*</b>			
JBH1707	<input type="checkbox"/>	Swing Away Bracket RH	<input type="radio"/>
JBH1708	<input type="checkbox"/>	Swing Away Bracket LH	<input type="radio"/>
<b>SHROUD COLOUR*</b>			
JBH2871	<input type="checkbox"/>	Mid range - Black/Anthracite Shroud	<input type="radio"/>
JBH2872	<input type="checkbox"/>	Mid range - Silver Shroud	<input type="radio"/>
JBH2873	<input type="checkbox"/>	Mid range - Bengal Red Shroud	<input type="radio"/>
JBH2874	<input type="checkbox"/>	Mid range - Mercury Blue Shroud	<input type="radio"/>
JBH2875	<input type="checkbox"/>	Mid range - Midnight Purple Shroud	<input type="radio"/>
<b>RIM COLOUR*</b>			
JBH2882	<input type="checkbox"/>	Mid range - Black Rim Inserts	<input type="radio"/>
JBH2883	<input type="checkbox"/>	Mid range - Silver Rim Inserts	<input type="radio"/>
JBH2884	<input type="checkbox"/>	Mid range - Bengal Red Rim Inserts	<input type="radio"/>
JBH2885	<input type="checkbox"/>	Mid range - Mercury Blue Rim Inserts	<input type="radio"/>
JBH2886	<input type="checkbox"/>	Mid range - Midnight Purple Rim Inserts	<input type="radio"/>
JBH2888	<input type="checkbox"/>	Premium - Racing Green Rim Inserts	<input type="radio"/>
JBH2889	<input type="checkbox"/>	Premium - Deep Ocean Blue Satin Rim Inserts	<input type="radio"/>
JBH2890	<input type="checkbox"/>	Premium - Frosty Blue Rim Inserts	<input type="radio"/>
JBH2891	<input type="checkbox"/>	Premium - Ultimate Red Rim Inserts	<input type="radio"/>
JBH2892	<input type="checkbox"/>	Premium - Falcon White Rim Inserts	<input type="radio"/>
<b>FORKS*</b>			
JBH2632	<input type="checkbox"/>	Double sided fork Without castor rim insert (enlarges wide base to 656 mm)	<input checked="" type="radio"/>
<b>FRONT WHEEL *</b>			
JBH0836	<input type="checkbox"/>	9" (225 mm) x 2.80 (70 mm) Castor, puncture proof solid tyre	<input checked="" type="radio"/>
<b>DRIVE WHEEL *</b>			
JBH0705	<input type="checkbox"/>	DW 14" (355 mm) x 3" (75 mm), Solid inlay Dual density foam insert	<input checked="" type="radio"/>
<b>C.T.C SUSPENSION*</b>			
JBH2870	<input type="checkbox"/>	CTC Suspension Control - Traction - Comfort	<input checked="" type="radio"/>
<b>SUSPENSION SETTING*</b>			
JBH1916	<input type="checkbox"/>	Suspension up to 160 kg User Weight For RX 40 HD only	<input checked="" type="radio"/>
<b>MOTOR TYPE*</b>			
JBH1909	<input type="checkbox"/>	12 km/h Ultra High Speed Package Always combined with LED Light system with Indicators, 30 Lux headlights	<input checked="" type="radio"/>
<b>BATTERY TYPE*</b>			
JBH1260	<input type="checkbox"/>	73,5 Ah Gel	<input checked="" type="radio"/>
<b>CHARGER*</b>			
JBH1343	<input type="checkbox"/>	Opti Charge Charger 8A (IP44)	<input checked="" type="radio"/>
<b>ANTI-TIPPERS*</b>			
JBH1205	<input type="checkbox"/>	Anti-tipper	<input checked="" type="radio"/>
<b>TRANSPORTATION KIT*</b>			
JBH1494	<input type="checkbox"/>	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19) Adequate headrest is compulsory, Not compatible with Seat lifter	<input checked="" type="radio"/>
<b>LAP BELT*</b>			
JBH2113	<input type="checkbox"/>	Lap belt with Auto Lock Flexible mount in side rails	<input checked="" type="radio"/>

**REAR VIEW MIRROR\***

JBH1415		Rear view mirror RH	○
JBH1339		Rear view mirror LH	○

<b>TOTAL PRICE</b>			
--------------------	--	--	--

**Non Conforming Options:**

JBH2132      Seat Depth 410 mm

JBH2196      Long Armpad (375 mm)



For a more complex customer requirement  
please refer to our Invacare Unique special order form.  
email [specialsuk@invacare.com](mailto:specialsuk@invacare.com)

**Special Instructions**

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_

Description of special request (please include relevant option code or part number):  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_