



# Motability Retail Prescription Form

**Yes, you can:**

Tel: 01656 776222  
Email: ordersuk@invacare.com  
Online Spares available at: www.invacare.co.uk

JANUARY 2025  
Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## AVIVA RX40 Modulite M24 707026 - weekly Price £53.95

### Standard Features

- Max user weight 136kg
- Crash tested ISO 7176 -Part 19
- Modulite seating system
- Tension adjustable backrest
- Standard T-Armrests, detachable with long armpad
- Manually settable Seat Tilt 0 - 15°
- Manual settable backrest Recline (0°- 30°, in 3.75°steps)
- 80° hangers with fixed footplate
- LiNX Control REM216 with swing away Bracket
- 9" Castor, puncture proof solid tyre
- Solid inlay Dual density foam insert
- 12 km/h Ultra High Speed Package
- LED Light system with Indicators, 30 Lux headlights
- 73,5 Ah Gel Batteries
- Anti tippers
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)



**Standard Base price: £7,424**

Please make your selection with a cross in the box     Free of charge     standard feature    POA special request

Customer weight =                      KG

**Code**                      **Description**

#### BASE MODEL\*

JBF0001	<input type="checkbox"/>	AVIVA RX 40 Base Model Max user weight up to 136 kg includes Performance Package: 120A Power Module + minimum 60 Ah Batteries	£7,424
---------	--------------------------	--	--------

#### CHASSIS WIDTH\*

JBL2810	<input type="checkbox"/>	Chassis Narrow (560 mm/22")	<input checked="" type="checkbox"/>
---------	--------------------------	-----------------------------	-------------------------------------

#### SEATING TYPE\*

JBF2100	<input type="checkbox"/>	Invacare® Modulite Seat Unit	<input checked="" type="checkbox"/>
---------	--------------------------	------------------------------	-------------------------------------

SEAT FRAME*			
JBF2111	<input type="checkbox"/>	Telescopic Seat Frame and Plate with Side Rails width 380 to 530 mm Additional 50 mm by arms	○
JBF2112	<input type="checkbox"/>	Telescopic Seat Frame + Sling Seat with Side Rails width 380 to 530 mm; + 30 mm seat height Additional 50 mm by arms	○
SEAT WIDTH*			
JBF2121	<input type="checkbox"/>	380 mm Seat Plate Can be reduced to 330 mm in combination with inrigged armrests	○
JBF2122	<input type="checkbox"/>	430 mm Seat Plate Can be reduced to 380 mm in combination with inrigged armrests	○
JBF2123	<input type="checkbox"/>	480 mm Seat Plate	○
JBF2124	<input type="checkbox"/>	530 mm Seat Plate	○
SEAT DEPTH*			
JBF2132	<input type="checkbox"/>	Seat Depth 410 mm	○
JBF2134	<input type="checkbox"/>	Seat Depth 460 mm	○
JBF2136	<input type="checkbox"/>	Seat Depth 510 mm	○
SEAT CUSHION*			
JBF2142	<input type="checkbox"/>	Matrx Flo-tech Cushion Flo-tech Lite	◆
BACKREST SUPPORT*			
JBF2163	<input type="checkbox"/>	Backposts short (480 mm) For Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery Width adjustable cast bracket	◆
JBF2162	<input type="checkbox"/>	Backposts tall (540 mm) For Flex3 backrests/Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery Width adjustable cast bracket	○
BACKREST UPHOLSTERY / CUSHION*			
JBF2171	<input type="checkbox"/>	Tension Adjustable Upholstery nylon overlay short = 480 mm / 7 straps measured from seat plate to top of upholstery 480 mm	◆
UPHOLSTERY TYPE*			
JBF2151	<input type="checkbox"/>	Black Fabric Upholstery (Nylon w 2141; 2-way-stretch material w 2142, 2143)	◆
ARMREST (Pair)*			
JBF2191	<input type="checkbox"/>	T-Arm detachable, Standard Adjustment range from seat plate: 245 - 315 mm	◆
ARMPAD (Pair)*			
JBF2196	<input type="checkbox"/>	Long Armpad (375 mm)	◆
SEAT LIFTER / SEAT TILT*			
JBF1091	<input type="checkbox"/>	Manually settable Seat Tilt 0 - 15°	◆
SEAT TO GROUND HEIGHT*			
JBF2857	<input type="checkbox"/>	Seat-to-Floor Height 435 mm	○
JBF2858	<input type="checkbox"/>	Seat-to-Floor Height 460 mm	○
JBF2859	<input type="checkbox"/>	Seat-to-Floor Height 485 mm	○
BACKREST RECLINER*			
JBF0975	<input type="checkbox"/>	Manual Settable Recliner (in 3.75° steps) 90° - 120°	◆
FOOTREST HANGER TYPE*			
JBF0526	<input type="checkbox"/>	Standard hanger 80° pair Lower leg length: 350 - 450 mm	◆
FOOTPLATE TYPE*			
JBF0512	<input type="checkbox"/>	Fixed Footplate pair	◆

<b>CONTROLS*</b>			
JBF2204	<input type="checkbox"/>	Remote LiNX DLX-REM216 = Drive, Actuators & Lights 3 Drive functions and power module incl. 2 Actuators & Lights	◆
<b>CONTROL BRACKETRY*</b>			
JBF1707	<input type="checkbox"/>	Swing Away Bracket RH	○
JBF1708	<input type="checkbox"/>	Swing Away Bracket LH	○
<b>SHROUD COLOUR*</b>			
JBF2871	<input type="checkbox"/>	Mid range - Black/Anthracite Shroud	○
JBF2872	<input type="checkbox"/>	Mid range - Silver Shroud	○
JBF2873	<input type="checkbox"/>	Mid range - Bengal Red Shroud	○
JBF2874	<input type="checkbox"/>	Mid range - Mercury Blue Shroud	○
JBF2875	<input type="checkbox"/>	Mid range - Midnight Purple Shroud	○
<b>RIM COLOUR*</b>			
JBF2882	<input type="checkbox"/>	Mid range - Black Rim Inserts	○
JBF2883	<input type="checkbox"/>	Mid range - Silver Rim Inserts	○
JBF2884	<input type="checkbox"/>	Mid range - Bengal Red Rim Inserts	○
JBF2885	<input type="checkbox"/>	Mid range - Mercury Blue Rim Inserts	○
JBF2886	<input type="checkbox"/>	Mid range - Midnight Purple Rim Inserts	○
<b>FORKS*</b>			
JBF2632	<input type="checkbox"/>	Double sided fork Without castor rim insert (enlarges wide base to 656 mm)	◆
<b>FRONT WHEEL*</b>			
JBF0836	<input type="checkbox"/>	9" (225 mm) x 2.80 (70 mm) Castor, puncture proof solid tyre	◆
<b>DRIVE WHEEL*</b>			
JBF0705	<input type="checkbox"/>	DW 14" (355 mm) x 3" (75 mm), Solid inlay Dual density foam insert	◆
<b>C.T.C SUSPENSION*</b>			
JBF2870	<input type="checkbox"/>	CTC Suspension Control - Traction - Comfort	◆
<b>SUSPENSION SETTING*</b>			
JBF1910	<input type="checkbox"/>	Soft (suggested for up to 70 kg User Weight)	○
JBF1912	<input type="checkbox"/>	Medium (suggested for up to 90 kg User Weight)	○
JBF1911	<input type="checkbox"/>	Firm (suggested for up to 110 kg User Weight)	○
JBF1915	<input type="checkbox"/>	Very firm (suggested for up to 136 kg User Weight)	○
<b>MOTOR TYPE*</b>			
JBF1909	<input type="checkbox"/>	12 km/h Ultra High Speed Package Always combined with LED Light system with Indicators, 30 Lux headlights	◆
<b>BATTERY TYPE*</b>			
JBF1260	<input type="checkbox"/>	73,5 Ah Gel	◆
<b>CHARGER*</b>			
JBF1343	<input type="checkbox"/>	Opti Charge Charger 8A (IP44)	◆
<b>ANTI-TIPPERS*</b>			
JBF1205	<input type="checkbox"/>	Anti-tipper	◆
<b>TRANSPORTATION KIT*</b>			
JBF1494	<input type="checkbox"/>	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19) Adequate headrest is compulsory	◆
<b>LAP BELT*</b>			
JBF1329	<input type="checkbox"/>	Lap belt with Auto Lock Standard mount on backrest	◆

**REAR VIEW MIRROR\***

JBF1415		Rear view mirror RH	○
JBF1339		Rear view mirror LH	○

**TOTAL PRICE**

**Non Conforming Options:**



For a more complex customer requirement  
please refer to our Invacare Unique special order form.  
email [specialsuk@invacare.com](mailto:specialsuk@invacare.com)

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):