



**Yes, you can:**

Tel: 01656 776222  
Email: ordersuk@invacare.com  
Online Spares available at: www.invacare.co.uk

JANUARY 2025  
Customer Ref:

|                          |                    |
|--------------------------|--------------------|
| <b>Customer Ref*</b>     | <b>Account No*</b> |
| <b>Contact No*</b>       | <b>Order No*</b>   |
| <b>Delivery Address*</b> | <b>Order Qty*</b>  |
|                          | <b>Order date*</b> |

## AVIVA RX40 Modulite M26 707017 - Weekly Price £60.49

### Standard Features

- Max user weight 136kg (21 Stone)
- Crash tested ISO 7176 -Part 19
- Modulite seating system
- Flex3 Tension Adjustable Upholstery
- Standard T-Armrests, detachable with long armpad
- Neckrest depth adjustable
- Powered "fixed pivot" Tilt 20°
- Manual Recliner
- Standard hanger 80° pair with Fixed Footplates
- LiNX control REM216 with swing away bracket
- 8" Castor, puncture proof solid tyre
- Solid inlay Dual density foam insert
- 12 km/h Ultra High Speed Package with LED Light system with Indicators, 30 Lux headlights & Mirror
- 73Ah batteries
- Anti tippers
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)



**Standard Base price: £8,062**

Please make your selection with a cross in the box     Free of charge     standard feature    POA    special request

Customer weight =                      KG

| Code                  | Description  |        |
|-----------------------|--|--------|
| <b>BASE MODEL*</b>    |  |        |
| JBG0001               | AVIVA RX 40 Base Model Max user weight up to 136 kg<br>includes Performance Package: 120A Power Module + minimum 60 Ah Batteries | £8,062 |
| <b>CHASSIS WIDTH*</b> |  |        |
| JBL2810               | Chassis Narrow (560 mm/22")  | ◆      |
| <b>SEATING TYPE*</b>  |  |        |
| JBG2100               | Invacare® Modulite Seat Unit   | ◆      |

|   |                          |   |   |
|---|--------------------------|---|---|
| <b>SEAT FRAME*</b>  |                          |   |   |
| JBG2111   | <input type="checkbox"/> | Telescopic Seat Frame and Plate with Side Rails width 380 to 530 mm Additional 50 mm by arms  | ◆ |
| <b>SEAT WIDTH*</b>  |                          |   |   |
| JBG2121   | <input type="checkbox"/> | 380 mm Seat Plate Can be reduced to 330 mm in combination with inrigged armrests  | ○ |
| JBG2122   | <input type="checkbox"/> | 430 mm Seat Plate Can be reduced to 380 mm in combination with inrigged armrests  | ○ |
| JBG2123   | <input type="checkbox"/> | 480 mm Seat Plate   | ○ |
| JBG2124   | <input type="checkbox"/> | 530 mm Seat Plate   | ○ |
| <b>SEAT DEPTH*</b>  |                          |   |   |
| JBG2132   | <input type="checkbox"/> | Seat Depth 410 mm   | ○ |
| JBG2134   | <input type="checkbox"/> | Seat Depth 460 mm   | ○ |
| JBG2136   | <input type="checkbox"/> | Seat Depth 510 mm   | ○ |
| <b>SEAT CUSHION*</b>  |                          |   |   |
| JBG2142   | <input type="checkbox"/> | Matrx Flo-tech Cushion Flo-tech Lite  | ◆ |
| <b>BACKREST SUPPORT*</b>  |                          |   |   |
| JBG2162   | <input type="checkbox"/> | Backposts tall (540 mm) For Flex3 backrests/Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery<br>Width adjustable cast bracket | ◆ |
| <b>PUSH HANDLES FOR BACKPOSTS</b>                                       |                          |   |   |
| JBG2108   | <input type="checkbox"/> | Backposts without push handles  | ◆ |
| <b>BACKREST UPHOLSTERY / CUSHION*</b>                                   |                          |   |   |
| JBG2173   | <input type="checkbox"/> | Flex3 adjustable Backrest Plates with Laguna Backrest cushion   | ◆ |
| <b>UPHOLSTERY TYPE*</b>   |                          |   |   |
| JBG2152   | <input type="checkbox"/> | Black multi stretch PU coated fabric Upholstery   | ○ |
| JBG2155   | <input type="checkbox"/> | Black Spacer Fabric Upholstery  | ○ |
| <b>ARMREST (Pair)*</b>  |                          |   |   |
| JBG2191   | <input type="checkbox"/> | T-Arm detachable, Standard Adjustment range from seat plate:<br>245 - 315 mm  | ◆ |
| <b>ARMPAD (Pair)*</b>   |                          |   |   |
| JBG2196   | <input type="checkbox"/> | Long Armpad (375 mm)  | ◆ |
| <b>HEADREST / NECKREST - Adequate headrest option must be selected*</b> |                          |   |   |
| JBG1350   | <input type="checkbox"/> | Neckrest depth adjustable   | ◆ |
| <b>SEAT LIFTER / SEAT TILT*</b>   |                          |   |   |
| JBG1092   | <input type="checkbox"/> | Powered "fixed pivot" Tilt 20°  | ◆ |
| <b>SEAT TO GROUND HEIGHT*</b>   |                          |   |   |
| JBG2857   | <input type="checkbox"/> | Seat-to-Floor Height 435 mm   | ○ |
| JBG2858   | <input type="checkbox"/> | Seat-to-Floor Height 460 mm   | ○ |
| JBG2859   | <input type="checkbox"/> | Seat-to-Floor Height 485 mm   | ○ |
| <b>BACKREST RECLINER*</b>   |                          |   |   |
| JBG0975   | <input type="checkbox"/> | Manual Settable Recliner (in 3.75° steps) 90° - 120°  | ◆ |
| <b>FOOTREST HANGER TYPE*</b>  |                          |   |   |
| JBG0526   | <input type="checkbox"/> | Standard hanger 80° pair Lower leg length: 350 - 450 mm   | ◆ |
| <b>FOOTPLATE TYPE*</b>  |                          |   |   |
| JBG0512   | <input type="checkbox"/> | Fixed Footplate pair  | ◆ |
| <b>CONTROLS*</b>  |                          |   |   |
| JBG2204   | <input type="checkbox"/> | DLX-REM216 = Drive, Actuators & Lights 3 Drive functions and power module incl. 2 Actuators & Lights  | ◆ |

| <b>CONTROL BRACKETRY*</b>  |                          |  |                                  |
|----------------------------|--------------------------|--|----------------------------------|
| JBG1707                    | <input type="checkbox"/> | Swing Away Bracket RH  | <input type="radio"/>            |
| JBG1708                    | <input type="checkbox"/> | Swing Away Bracket LH  | <input type="radio"/>            |
| <b>SHROUD COLOUR*</b>      |                          |  |                                  |
| JBG2871                    | <input type="checkbox"/> | Mid range - Black/Anthracite Shroud  | <input type="radio"/>            |
| JBG2872                    | <input type="checkbox"/> | Mid range - Silver Shroud  | <input type="radio"/>            |
| JBG2873                    | <input type="checkbox"/> | Mid range - Bengal Red Shroud  | <input type="radio"/>            |
| JBG2874                    | <input type="checkbox"/> | Mid range - Mercury Blue Shroud  | <input type="radio"/>            |
| JBG2875                    | <input type="checkbox"/> | Mid range - Midnight Purple Shroud   | <input type="radio"/>            |
| <b>RIM COLOUR*</b>         |                          |  |                                  |
| JBG2882                    | <input type="checkbox"/> | Mid range - Black Rim Inserts  | <input type="radio"/>            |
| JBG2883                    | <input type="checkbox"/> | Mid range - Silver Rim Inserts   | <input type="radio"/>            |
| JBG2884                    | <input type="checkbox"/> | Mid range - Bengal Red Rim Inserts   | <input type="radio"/>            |
| JBG2885                    | <input type="checkbox"/> | Mid range - Mercury Blue Rim Inserts   | <input type="radio"/>            |
| JBG2886                    | <input type="checkbox"/> | Mid range - Midnight Purple Rim Inserts  | <input type="radio"/>            |
| <b>FORKS*</b>              |                          |  |                                  |
| JBG2632                    | <input type="checkbox"/> | Double sided fork Without castor rim insert<br>(enlarges wide base to 656 mm)  | <input checked="" type="radio"/> |
| <b>FRONT WHEEL*</b>        |                          |  |                                  |
| JBG0836                    | <input type="checkbox"/> | 9" (225 mm) x 2.80 (70 mm) Castor, puncture proof solid tyre   | <input checked="" type="radio"/> |
| <b>DRIVE WHEEL*</b>        |                          |  |                                  |
| JBG0705                    | <input type="checkbox"/> | DW 14" (355 mm) x 3" (75 mm), Solid inlay Dual density foam insert   | <input checked="" type="radio"/> |
| <b>C.T.C SUSPENSION*</b>   |                          |  |                                  |
| JBG2870                    | <input type="checkbox"/> | CTC Suspension Control - Traction - Comfort  | <input checked="" type="radio"/> |
| <b>SUSPENSION SETTING*</b> |                          |  |                                  |
| JBG1910                    | <input type="checkbox"/> | Soft (suggested for up to 70 kg User Weight)   | <input type="radio"/>            |
| JBG1912                    | <input type="checkbox"/> | Medium (suggested for up to 90 kg User Weight)   | <input type="radio"/>            |
| JBG1911                    | <input type="checkbox"/> | Firm (suggested for up to 110 kg User Weight)  | <input type="radio"/>            |
| JBG1915                    | <input type="checkbox"/> | Very firm (suggested for up to 136 kg User Weight)   | <input type="radio"/>            |
| <b>MOTOR TYPE*</b>         |                          |  |                                  |
| JBG1909                    | <input type="checkbox"/> | 12 km/h Ultra High Speed Package<br>Always combined with LED Light system with Indicators, 30 Lux headlights                         | <input checked="" type="radio"/> |
| <b>BATTERY TYPE*</b>       |                          |  |                                  |
| JBG1260                    | <input type="checkbox"/> | 73,5 Ah Gel  | <input checked="" type="radio"/> |
| <b>CHARGER*</b>            |                          |  |                                  |
| JBG1343                    | <input type="checkbox"/> | Opti Charge Charger 8A (IP44)  | <input checked="" type="radio"/> |
| <b>ANTI-TIPPERS*</b>       |                          |  |                                  |
| JBG1205                    | <input type="checkbox"/> | Anti-tipper  | <input checked="" type="radio"/> |
| <b>TRANSPORTATION KIT*</b> |                          |  |                                  |
| JBG1494                    | <input type="checkbox"/> | 4-Point tie-down<br>(Tested version for using wheelchair as a vehicle seat according to ISO 7176-19) Adequate headrest is compulsory | <input checked="" type="radio"/> |
| <b>LAP BELT*</b>           |                          |  |                                  |
| JBG2113                    | <input type="checkbox"/> | Lap belt with Auto Lock Flexible mount in side rails   | <input checked="" type="radio"/> |
| <b>REAR VIEW MIRROR*</b>   |                          |  |                                  |
| JBG1415                    | <input type="checkbox"/> | Rear view mirror RH  | <input type="radio"/>            |
| JBG1339                    | <input type="checkbox"/> | Rear view mirror LH  | <input type="radio"/>            |

---

|                    |  |
|--------------------|--|
| <b>TOTAL PRICE</b> |  |
|--------------------|--|

**Non Conforming Options:**



For a more complex customer requirement  
please refer to our Invacare Unique special order form.  
email [specialsuk@invacare.com](mailto:specialsuk@invacare.com)

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):