



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Motability Retail Prescription Form

JANUARY 2025
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Storm4 Modulite M30 708645 - Weekly Price £54.88

Standard Features

- Max user weight 150Kg (23.5 stone)
- Crash tested ISO 7176 -Part 19
- Modulite seating system
- Matrix Flo-tech Lite Cushion
- Tension Adjustable Backrest
- Standard T-Armrests, detachable with long armpad
- Manual settable seat Tilt 0° - 20°
- Manual settable Recliner (0° - 30° in 3.75° steps)
- VARI F hangers with Fixed Footplate
- LiNX Control REM216 withswing away bracket
- Pneumatic Castors
- Pneumatic Drive Wheels
- 73.5 Ah Batteries
- 12km/h (8mph) Motors
- LED light system with indicators and Mirror
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)
- Anti-Tippers



Standard Base price: £7,361

Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

Code **Description**

CHAIR TYPE*

ATB0000	<input type="checkbox"/>	Storm 4 Modulite	£7,361
---------	--------------------------	------------------	--------

SEATING TYPE*

ATB2100	<input type="checkbox"/>	Modulite Seating Unit (Flex3 seating unit)	◆
---------	--------------------------	--	---

SEATING OPTIONS*

ATB2111	<input type="checkbox"/>	Telescopic Seat Frame with metal seat plate (with stepless adjustment on armrests, legrests and backrest)	◆
---------	--------------------------	---	---

SEAT WIDTH*			
ATBXXXX	<input type="checkbox"/>	Seat width 330mm (13") seat width achieved on a 380mm seat frame (2121) with inrigged armrests (2183)	<input type="radio"/>
ATB2121	<input type="checkbox"/>	Set at 380mm (adjustable from 380mm to 430mm (15" - 17") with stepless adjustment on armrests, legrests and backrest)	<input type="radio"/>
ATB2122	<input type="checkbox"/>	Set at 430mm (adjustable from 430mm to 480mm (17" - 19) with stepless adjustment on armrests, legrests and backrest)	<input type="radio"/>
ATB2123	<input type="checkbox"/>	Set at 480mm (adjustable from 480mm to 530mm (19" - 21") with stepless adjustment on armrests, legrests and backrest)	<input type="radio"/>
ATB2124	<input type="checkbox"/>	Set at 530mm (adjustable from 530mm to 580mm (21" - 23") with stepless adjustment on armrests, legrests and backrest)	<input type="radio"/>
SEAT DEPTH*			
ATB2132	<input type="checkbox"/>	410mm (16") (adjustable 410mm - 510mm stepless adjustment) (NC 2196)	<input type="radio"/>
ATB2134	<input type="checkbox"/>	460mm (18") (adjustable 410mm - 510mm stepless adjustment)	<input type="radio"/>
ATB2136	<input type="checkbox"/>	510mm (20") (adjustable 410mm - 510mm stepless adjustment)	<input type="radio"/>
SEAT CUSHION*			
ATB2142	<input type="checkbox"/>	Matrx Flo-tech Lite	<input checked="" type="checkbox"/>
BACKREST SUPPORT*			
ATB2162	<input type="checkbox"/>	Backposts 540mm (Tall) must be selected for Flex3 Laguna backrest	<input type="radio"/>
ATB2163	<input type="checkbox"/>	Backposts (Short)	<input type="radio"/>
PUSH HANDLES FOR BACKPOSTS*			
ATB2107	<input type="checkbox"/>	Backposts WITH push handles	<input checked="" type="checkbox"/>
BACKREST UPHOLSTERY / CUSHION*			
ATB2171	<input type="checkbox"/>	Basic Tension Adjustable Backrest (540mm or 480mm upholstery height) Straps easily accessible and adjustable from the rear, Black fabric only	<input checked="" type="checkbox"/>
UPHOLSTREY TYPE*			
ATB2151	<input type="checkbox"/>	Black fabric	<input checked="" type="checkbox"/>
ARMREST (PAIR)*			
ATB2191	<input type="checkbox"/>	T-Armrest Standard, Range of adjustment from seat plate 245mm to 315mm	<input checked="" type="checkbox"/>
ARMPAD (PAIR)*			
ATB2196	<input type="checkbox"/>	Long Armpad (375mm) (NC 2132)	<input checked="" type="checkbox"/>
SEAT HEIGHT AND SEAT RISER*			
ATB1023	<input type="checkbox"/>	Fixed Seat Height 450mm (17.75")	<input type="radio"/>
ATB1025	<input type="checkbox"/>	Fixed Seat Height 480mm (19")	<input type="radio"/>
SEAT TILT*			
ATB1002	<input type="checkbox"/>	Manual Seat Tilt by spindle	<input checked="" type="checkbox"/>
BACKREST RECLINER*			
ATB0975	<input type="checkbox"/>	Manual Settable Recliner (0° - 30°, in 3.75° steps) Foldable	<input checked="" type="checkbox"/>
FOOTREST HANGER TYPE (PAIR)*			
ATB0503	<input type="checkbox"/>	Vari F : Angle adjustable with tools 0° - 70°, length 290mm - 460mm	<input checked="" type="checkbox"/>
FOOTPLATE TYPE*			
ATB0512	<input type="checkbox"/>	Fixed Footplate	<input checked="" type="checkbox"/>
CONTROLS*			
ATB2204	<input type="checkbox"/>	REM216 Drive, Actuators & Lights	<input checked="" type="checkbox"/>

CONTROL BRACKETRY*			
ATB1707	<input type="checkbox"/>	Swing Away RH Bracket	<input type="radio"/>
ATB1708	<input type="checkbox"/>	Swing Away LH Bracket	<input type="radio"/>
SHROUD COLOUR*			
ATB0202	<input type="checkbox"/>	Blue	<input type="radio"/>
ATB0228	<input type="checkbox"/>	Black	<input type="radio"/>
ATB0204	<input type="checkbox"/>	Transparent Cherry Red	<input type="radio"/>
FRONT WHEEL*			
ATB0851	<input type="checkbox"/>	3.00 - 6 Grey Pneumatic	<input checked="" type="radio"/>
REAR WHEEL*			
ATB0771	<input type="checkbox"/>	3.00 - 8 Grey Pneumatic	<input checked="" type="radio"/>
MOTOR TYPE*			
ATB1909	<input type="checkbox"/>	12 km/h (8mph) (includes LED light system with indicators)	<input checked="" type="radio"/>
BATTERY*			
ATB1364	<input type="checkbox"/>	70 Ah Gel	<input checked="" type="radio"/>
SUSPENSION*			
ATB1240	<input type="checkbox"/>	Suspension low (up to 80kg user weight)	<input type="radio"/>
ATB1241	<input type="checkbox"/>	Suspension middle (from 70kg up to 100kg user weight)	<input type="radio"/>
ATB1242	<input type="checkbox"/>	Suspension hard (from 90kg up to 150kg user weight)	<input type="radio"/>
CHARGER*			
ATB1343	<input type="checkbox"/>	OptiCharge Charger 8A (200 - 250V)	<input checked="" type="radio"/>
ANTI TIPPER*			
ATB1205	<input type="checkbox"/>	Anti-tipper (standard version)	<input checked="" type="radio"/>
TRANSPORTATION KIT*			
ATB1494	<input type="checkbox"/>	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory, see options below)	<input checked="" type="radio"/>
FIXATION BELT*			
ATB1329	<input type="checkbox"/>	Lap Belt with Auto Lock	<input checked="" type="radio"/>
REAR VIEW MIRROR*			
ATB1415	<input type="checkbox"/>	Rear view Mirror RH - Must be selected with 10 km/h Motor	<input type="radio"/>
ATB1339	<input type="checkbox"/>	Rear view Mirror LH - Must be selected with 10 km/h Motor	<input type="radio"/>

TOTAL PRICE	
--------------------	--

Non Conforming Options:

ATB2132 410mm (16") (adjustable 410mm - 510mm stepless adjustment)

ATB2196 Long Armpad (375mm)



For a more complex customer requirement
please refer to our Invacare Unique special order form.
email specialsuk@invacare.com

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):