



Yes, you can:

Motability Retail Prescription Form

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

JANUARY 2025
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Storm4 Modulite M32 1001076 - Weekly Price £62.35

Standard Features

- Max user weight 150Kg (23.5 stone)
- Crash tested ISO 7176 -Part 19
- Modulite seating system
- Flex3 Laguna Backrest
- Standard T-Armrests, detachable with long armpad
- Manual settable seat Tilt 0° - 20° Electrical Seat Riser
- Manual settable Recliner (0° - 30° in 3.75° steps)
- VARI F with Fixed Footplate
- LinX Control REM216 with swing away bracket
- Pneumatic Castors
- Pneumatic Drive Wheels
- 73.5 Ah Batteries
- 12 km/h (4mph) Motors
- LED Light system with indicators
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)
- Anti-Tippers



Standard Base price: £8,950

Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

Code	Description	
CHAIR TYPE*		
ATB0000	Storm 4 Modulite	£8,950
ATZ0000	Special Storm 4 Modulite (Please provide details in special section)	POA
SEATING*		
ATB2100	Modulite Seating Unit (Flex3 seating unit)	◆
SEATING OPTION*		
ATB2111	Telescopic Seat Frame with metal seat plate (with stepless adjustment on armrests, legrests and backrest)	◆

SEAT WIDTH*			
ATBXXXX	<input type="checkbox"/>	Seat width 330mm (13") seat width achieved on a 380mm seat frame (2121) with inrigged armrests (2183)	<input type="radio"/>
ATB2121	<input type="checkbox"/>	Set at 380mm (adjustable from 380mm to 430mm (15" - 17") with stepless adjustment on armrests, legrests and backrest)	<input type="radio"/>
ATB2122	<input type="checkbox"/>	Set at 430mm (adjustable from 430mm to 480mm (17" - 19) with stepless adjustment on armrests, legrests and backrest)	<input type="radio"/>
ATB2123	<input type="checkbox"/>	Set at 480mm (adjustable from 480mm to 530mm (19" - 21") with stepless adjustment on armrests, legrests and backrest)	<input type="radio"/>
ATB2124	<input type="checkbox"/>	Set at 530mm (adjustable from 530mm to 580mm (21" - 23") with stepless adjustment on armrests, legrests and backrest)	<input type="radio"/>
SEAT DEPTH*			
ATB2132	<input type="checkbox"/>	410mm (16") (adjustable 410mm - 510mm stepless adjustment) (NC 2196)	<input type="radio"/>
ATB2134	<input type="checkbox"/>	460mm (18") (adjustable 410mm - 510mm stepless adjustment)	<input type="radio"/>
ATB2136	<input type="checkbox"/>	510mm (20") (adjustable 410mm - 510mm stepless adjustment)	<input type="radio"/>
SEAT CUSHION*			
ATB2142	<input type="checkbox"/>	Matrx Flo-tech Lite	<input checked="" type="checkbox"/>
BACKREST SUPPORT*			
ATB2162	<input type="checkbox"/>	Backposts 540mm (Tall) must be selected for Flex3 Laguna backrest	<input checked="" type="checkbox"/>
PUSH HANDLES FOR BACKPOSTS*			
ATB2107	<input type="checkbox"/>	Backposts WITH push handles	<input checked="" type="checkbox"/>
BACKREST UPHOLSTERY / CUSHION*			
ATB2173	<input type="checkbox"/>	Flex3 adjustable Backrest Plates with Laguna Backrest cushion	<input checked="" type="checkbox"/>
UPHOLSTERY TYPE*			
ATB2152	<input type="checkbox"/>	Black multi stretch PU coated fabric	<input type="radio"/>
ATB2155	<input type="checkbox"/>	Black Spacer fabric	<input type="radio"/>
ARMREST (PAIR)*			
ATB2191	<input type="checkbox"/>	T-Armrest Standard, Range of adjustment from seat plate 245mm to 315mm	<input checked="" type="checkbox"/>
ARMPAD (PAIR)*			
ATB2196	<input type="checkbox"/>	Long Armpad (375mm) (NC 2132)	<input checked="" type="checkbox"/>
HEADREST / NECKREST - Adequate headrest option must be selected*			
ATB1350	<input type="checkbox"/>	neckrest: height and depth adjustable	<input checked="" type="checkbox"/>
SEAT HEIGHT AND SEAT RISER*			
ATB1023	<input type="checkbox"/>	Fixed Seat Height 450mm (17.75")	<input type="radio"/>
ATB1025	<input type="checkbox"/>	Fixed Seat Height 480mm (19")	<input type="radio"/>
ATB1030	<input type="checkbox"/>	Electrical Seat Riser (400mm - 650mm) (16" - 25.5")	<input type="radio"/>
ATB1040	<input type="checkbox"/>	Electrical Seat Riser (440mm - 690mm) incl. adapter (17.25" - 27")	<input type="radio"/>
SEAT TILT*			
ATB1002	<input type="checkbox"/>	Manual Seat Tilt by spindle	<input checked="" type="checkbox"/>
BACKREST RECLINER*			
ATB0975	<input type="checkbox"/>	Manual Settable Recliner (0° - 30°, in 3.75° steps) Foldable	<input checked="" type="checkbox"/>
FOOTREST HANGER TYPE (PAIR)*			
ATB0503	<input type="checkbox"/>	Vari F : Angle adjustable with tools 0° - 70°, length 290mm - 460mm	<input checked="" type="checkbox"/>

FOOTPLATE TYPE*		
ATB0512	<input type="checkbox"/>	Fixed Footplate ○
CONTROLS*		
ATB2204	<input type="checkbox"/>	REM216 Drive, Actuators & Lights ◆
CONTROL BRACKETRY*		
ATB1707	<input type="checkbox"/>	Swing Away RH Bracket ○
ATB1708	<input type="checkbox"/>	Swing Away LH Bracket ○
SHROUD COLOUR*		
ATB0202	<input type="checkbox"/>	Blue ○
ATB0228	<input type="checkbox"/>	Black ○
ATB0204	<input type="checkbox"/>	Transparent Cherry Red ○
FRONT WHEEL*		
ATB0851	<input type="checkbox"/>	3.00 - 6 Grey Pneumatic ◆
REAR WHEEL*		
ATB0771	<input type="checkbox"/>	3.00 - 8 Grey Pneumatic ◆
MOTOR TYPE*		
ATB1909	<input type="checkbox"/>	12 km/h (8mph) (includes LED light system with indicators) ◆
BATTERY*		
ATB1364	<input type="checkbox"/>	70 Ah Gel ◆
SUSPENSION*		
ATB1240	<input type="checkbox"/>	Suspension low (up to 80kg user weight) ○
ATB1241	<input type="checkbox"/>	Suspension middle (from 70kg up to 100kg user weight) ○
ATB1242	<input type="checkbox"/>	Suspension hard (from 90kg up to 150kg user weight) ○
CHARGER*		
ATB1343	<input type="checkbox"/>	OptiCharge Charger 8A (200 - 250V) ◆
ANTI TIPPER*		
ATB1205	<input type="checkbox"/>	Anti-tipper (standard version) ○
TRANSPORTATION KIT*		
ATB1494	<input type="checkbox"/>	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory, see options below) ◆
FIXATION BELT*		
ATB2113	<input type="checkbox"/>	Lap belt with auto lock (flexible mount in side rails) ◆
REAR VIEW MIRROR		
ATB1415	<input type="checkbox"/>	Rear view Mirror RH - Must be selected with 10 km/h Motor ○
ATB1339	<input type="checkbox"/>	Rear view Mirror LH - Must be selected with 10 km/h Motor ○

TOTAL PRICE	
--------------------	--

Non Conforming Options:

ATB2132 410mm (16") (adjustable 410mm - 510mm stepless adjustment)

ATB2196 Long Armpad (375mm)



For a more complex customer requirement
please refer to our Invacare Unique special order form.
email specialsuk@invacare.com

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):