



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Motability Retail Prescription Form

JANUARY 2025
Customer Ref:

Customer Ref*, Account No*, Contact No*, Order No*, Delivery Address*, Order Qty*, Order date*

Storm4 X-Plore Modulite M34 707458 - Weekly Price £58.10

Standard Features

- Max user weight 150Kg (23.5 stone)
Crash tested ISO7176-Part 19
Modulite seating system
Flex3 Laguna backrest
Matrix Flo-tech Lite cushion
Backposts 540mm (Tall)
Standard T-armrests, detachable with Long armpads
Manual settable seat Tilt
Manual settable backrest Recline
Angle adjustable legrest hangers (VARI F) with Fixed footplates
LINX REM216 control with swing away bracket
Puncture proof (sealed) drive wheels and front castors
12 km/h (8mph) includes LED light system with indicators, Mirror
Lapbelt (adjustable both sides)

Standard Base price: £7,885



Please make your selection with a cross in the box [] Free of charge [O] standard feature [◆] POA special request

Customer weight = KG

Table with columns: Code, Description, Price. Rows include CHAIR TYPE* (ATE0017, ATZ0017), SEATING TYPE* (ATE2100), and SEAT FRAME* (ATE2111).

| SEAT WIDTH* | | | |
|---|--------------------------|--|----------------------------------|
| ATEXXXX | <input type="checkbox"/> | Seat width 330mm (13") | <input type="radio"/> |
| ATE2121 | <input type="checkbox"/> | Set at 380mm | <input type="radio"/> |
| ATE2122 | <input type="checkbox"/> | Set at 430mm | <input type="radio"/> |
| ATE2123 | <input type="checkbox"/> | Set at 480mm | <input type="radio"/> |
| ATE2124 | <input type="checkbox"/> | Set at 530mm | <input type="radio"/> |
| SEAT DEPTH* | | | |
| ATE2132 | <input type="checkbox"/> | 410mm (16") | <input type="radio"/> |
| | | (NC 2196) | |
| ATE2134 | <input type="checkbox"/> | 460mm (18") | <input type="radio"/> |
| ATE2136 | <input type="checkbox"/> | 510mm (20") | <input type="radio"/> |
| SEAT CUSHION* | | | |
| ATE2142 | <input type="checkbox"/> | Matrx Flo-tech Lite | <input checked="" type="radio"/> |
| BACKREST SUPPORT* | | | |
| ATE2162 | <input type="checkbox"/> | Backposts 540mm (Tall) | <input checked="" type="radio"/> |
| PUSH HANDLES FOR BACKPOSTS* | | | |
| ATE2107 | <input type="checkbox"/> | Backposts WITH push handles | <input checked="" type="radio"/> |
| BACKREST UPHOLSTREY / CUSHION* | | | |
| ATE2173 | <input type="checkbox"/> | Flex3 adjustable Backrest Plates with Laguna Backrest cushion | <input checked="" type="radio"/> |
| UPHOLSTERY COLOUR* | | | |
| ATE2152 | <input type="checkbox"/> | Black multi stretch PU coated fabric | <input type="radio"/> |
| ATE2155 | <input type="checkbox"/> | Black Spacer fabric | <input type="radio"/> |
| ARMREST* | | | |
| ATE2191 | <input type="checkbox"/> | T-Armrest Standard | <input checked="" type="radio"/> |
| ARMPAD* | | | |
| ATE2196 | <input type="checkbox"/> | Long Armpad (375mm) | <input checked="" type="radio"/> |
| HEADREST / NECKREST - Adequate headrest option must be selected* | | | |
| ATE1495 | <input type="checkbox"/> | Push bar | <input checked="" type="radio"/> |
| SEAT HEIGHT SETTING & SEAT RISER* | | | |
| ATE1023 | <input type="checkbox"/> | Fixed Seat Height 450mm | <input type="radio"/> |
| ATE1025 | <input type="checkbox"/> | Fixed Seat Height 480mm | <input type="radio"/> |
| SEAT TILT* | | | |
| ATE1002 | <input type="checkbox"/> | Manual Seat Tilt by spindle | <input checked="" type="radio"/> |
| BACKREST RECLINER* | | | |
| ATE0975 | <input type="checkbox"/> | Manual Settable Recliner (0° - 30°, in 3.75° steps) Foldable | <input checked="" type="radio"/> |
| FOOTREST HANGER TYPE (PAIR)* | | | |
| ATE0503 | <input type="checkbox"/> | Vari F : Angle adjustable with tools 0° - 70°, length 290mm - 460mm | <input checked="" type="radio"/> |
| FOOTPLATE TYPE* | | | |
| ATE0512 | <input type="checkbox"/> | Fixed Footplate | <input checked="" type="radio"/> |
| CONTROLS* | | | |
| ATE2204 | <input type="checkbox"/> | REM216 Drive, Actuators & Lights | <input checked="" type="radio"/> |
| CONTROL BRACKETRY* | | | |
| ATE1707 | <input type="checkbox"/> | Swing Away RH Bracket | <input type="radio"/> |
| ATE1708 | <input type="checkbox"/> | Swing Away LH Bracket | <input type="radio"/> |

| SHROUD COLOUR* | | | |
|--------------------------|--------------------------|---|----------------------------------|
| ATE0202 | <input type="checkbox"/> | Blue | <input type="radio"/> |
| ATE0228 | <input type="checkbox"/> | Black | <input type="radio"/> |
| ATE0204 | <input type="checkbox"/> | Transparent Cherry Red | <input type="radio"/> |
| CASTOR WHEEL* | | | |
| ATE0853 | <input type="checkbox"/> | 3.00 - 6 Grey Puncture Proof (sealed) | <input checked="" type="radio"/> |
| DRIVE WHEEL* | | | |
| ATE0773 | <input type="checkbox"/> | 3.00 - 8 Grey Puncture Proof (sealed) | <input checked="" type="radio"/> |
| MOTOR TYPE* | | | |
| ATE1909 | <input type="checkbox"/> | 12 km/h (8mph) includes LED light system with indicators | <input checked="" type="radio"/> |
| BATTERY* | | | |
| ATE1364 | <input type="checkbox"/> | 73.5 Ah GEL | <input checked="" type="radio"/> |
| REAR SUSPENSION* | | | |
| ATE1233 | <input type="checkbox"/> | Hitec spring X-plore | <input type="radio"/> |
| FRONT SUSPENSION* | | | |
| ATE1238 | <input type="checkbox"/> | Front suspension (suggested for up to 100kg) | <input type="radio"/> |
| ATE1239 | <input type="checkbox"/> | Front suspension (suggested for above 100kg) | <input type="radio"/> |
| CHARGER* | | | |
| ATE1343 | <input type="checkbox"/> | OptiCharge Charger 8A (200 - 250V) | <input checked="" type="radio"/> |
| ANTI-TIPPERS* | | | |
| ATE1205 | <input type="checkbox"/> | Anti-tipper (standard version) | <input checked="" type="radio"/> |
| TRANSPORTATION* | | | |
| ATE1494 | <input type="checkbox"/> | 4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory). | <input checked="" type="radio"/> |
| FIXATION BELT* | | | |
| ATE1329 | <input type="checkbox"/> | Lapbelt (adjustable both sides) fitted to base of backrest | <input checked="" type="radio"/> |
| REAR VIEW MIRROR* | | | |
| ATE1415 | <input type="checkbox"/> | Rear view Mirror RH - Must be selected with 10 km/h Motor | <input type="radio"/> |
| ATE1339 | <input type="checkbox"/> | Rear view Mirror LH - Must be selected with 10 km/h Motor | <input type="radio"/> |

| | |
|--------------------|----------------------|
| TOTAL PRICE | <input type="text"/> |
|--------------------|----------------------|

Non Conforming Options:

ATE2196 Long Armpad (375mm)



For a more complex customer requirement
please refer to our Invacare Unique special order form.
email specialsuk@invacare.com

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):