



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Motability Retail Prescription Form

JANUARY 2025
Customer Ref:

| | |
|--------------------------|--------------------|
| Customer Ref* | Account No* |
| Contact No* | Order No* |
| Delivery Address* | Order Qty* |
| | Order date* |

Storm4 X-Plore Modulite M37 1001120 - Weekly Price £62.26

Standard Features

- Max user weight 150Kg (23.5 stone)
- Crash tested ISO7176-Part 19
- Modulite seating system
- Flex3 Laguna Backrest
- Standard T-Armrests, detachable with Long Armpads
- Neckrest
- Electrical Seat Riser, Electrical Tilt from 0° to 25°
- Manual settable Recliner (0° - 30° in 3.75° steps)
- VARI F hangers with Fixed Footplate
- LiNX control REM216 with swing away bracket
- Puncture proof castors
- Puncture proof Drive Wheels
- 12 km/h (8mph) Motors, LED light system with indicators
- 4 Point-tie down system
- Lapbelt with autolock (adjustable both sides)
- Anti-Tippers
- X-plore suspension 50mm all round

Standard Base price: £9,653



Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

| Code | Description | |
|----------------------|--|--------|
| CHAIR TYPE* | | |
| ATE0017 | <input type="checkbox"/> Storm 4 X-plore | £9,653 |
| SEATING TYPE* | | |
| ATE2100 | <input type="checkbox"/> Invacare Modulite 'Flex3' | ◆ |
| SEAT FRAME* | | |
| ATE2111 | <input type="checkbox"/> Telescopic Seat Frame with metal seat plate | ◆ |

| SEAT WIDTH* | | | |
|--|--------------------------|--|----------------------------------|
| ATEXXXX | <input type="checkbox"/> | Seat width 330mm (13") | <input type="radio"/> |
| ATE2121 | <input type="checkbox"/> | Set at 380mm | <input type="radio"/> |
| ATE2122 | <input type="checkbox"/> | Set at 430mm | <input type="radio"/> |
| ATE2123 | <input type="checkbox"/> | Set at 480mm | <input type="radio"/> |
| ATE2124 | <input type="checkbox"/> | Set at 530mm | <input type="radio"/> |
| SEAT DEPTH* | | | |
| ATE2132 | <input type="checkbox"/> | 410mm (16") (NC 2196) | <input type="radio"/> |
| ATE2134 | <input type="checkbox"/> | 460mm (18") | <input type="radio"/> |
| ATE2136 | <input type="checkbox"/> | 510mm (20") | <input type="radio"/> |
| SEAT CUSHION* | | | |
| ATE2142 | <input type="checkbox"/> | Matrx Flo-tech Lite | <input checked="" type="radio"/> |
| BACKREST SUPPORT* | | | |
| ATE2162 | <input type="checkbox"/> | Backposts 540mm (Tall) | <input checked="" type="radio"/> |
| PUSH HANDLES FOR BACKPOSTS* | | | |
| ATE2107 | <input type="checkbox"/> | Backposts WITH push handles | <input checked="" type="radio"/> |
| BACKREST UPHOLSTREY / CUSHION* | | | |
| ATE2173 | <input type="checkbox"/> | Flex3 adjustable Backrest Plates with Laguna Backrest cushion | <input checked="" type="radio"/> |
| UPHOLSTERY COLOUR* | | | |
| ATE2152 | <input type="checkbox"/> | Black multi stretch PU coated fabric | <input type="radio"/> |
| ATE2155 | <input type="checkbox"/> | Black Spacer fabric | <input type="radio"/> |
| ARMREST* | | | |
| ATE2191 | <input type="checkbox"/> | T-Armrest Standard | <input checked="" type="radio"/> |
| ARMPAD* | | | |
| ATE2196 | <input type="checkbox"/> | Long Armpad (375mm) | <input checked="" type="radio"/> |
| HEADREST / NECKREST - Adequate headrest option must be selected* | | | |
| ATE1350 | <input type="checkbox"/> | Neckrest: height and depth Adjustable | <input checked="" type="radio"/> |
| SEAT HEIGHT SETTING & SEAT RISER* | | | |
| ATE1040 | <input type="checkbox"/> | Electrical Seat Riser (440mm - 690mm) | <input checked="" type="radio"/> |
| SEAT TILT* | | | |
| ATE1014 | <input type="checkbox"/> | Electrical Tilt from 0° to 25° | <input checked="" type="radio"/> |
| BACKREST RECLINER* | | | |
| ATE0975 | <input type="checkbox"/> | Manual Settable Recliner (0° - 30°, in 3.75° steps) Foldable | <input checked="" type="radio"/> |
| FOOTREST HANGER TYPE (PAIR)* | | | |
| ATE0503 | <input type="checkbox"/> | Vari F : Angle adjustable with tools 0° - 70°, length 290mm - 460mm | <input checked="" type="radio"/> |
| FOOTPLATE TYPE* | | | |
| ATE0512 | <input type="checkbox"/> | Fixed Footplate | <input checked="" type="radio"/> |
| CONTROLS* | | | |
| ATE2204 | <input type="checkbox"/> | REM216 Drive, Actuators & Lights | <input checked="" type="radio"/> |
| CONTROL BRACKETRY* | | | |
| ATE1707 | <input type="checkbox"/> | Swing Away RH Bracket | <input type="radio"/> |
| ATE1708 | <input type="checkbox"/> | Swing Away LH Bracket | <input type="radio"/> |
| SHROUD COLOUR* | | | |
| ATE0202 | <input type="checkbox"/> | Blue | <input type="radio"/> |
| ATE0228 | <input type="checkbox"/> | Black | <input type="radio"/> |
| ATE0204 | <input type="checkbox"/> | Transparent Cherry Red | <input type="radio"/> |

| | | |
|--------------------------|--------------------------|---|
| CASTOR WHEEL* | | |
| ATE0853 | <input type="checkbox"/> | 3.00 - 6 Grey Puncture Proof (sealed) ◆ |
| DRIVE WHEEL* | | |
| ATE0773 | <input type="checkbox"/> | 3.00 - 8 Grey Puncture Proof (sealed) ◆ |
| MOTOR TYPE* | | |
| ATE1909 | <input type="checkbox"/> | 12 km/h (8mph) includes LED light system with indicators ◆ |
| BATTERY* | | |
| ATE1364 | <input type="checkbox"/> | 73.5 Ah GEL ◆ |
| REAR SUSPENSION* | | |
| ATE1233 | <input type="checkbox"/> | Hitec spring X-plore ○ |
| FRONT SUSPENSION* | | |
| ATE1238 | <input type="checkbox"/> | Front suspension (suggested for up to 100kg) ○ |
| ATE1239 | <input type="checkbox"/> | Front suspension (suggested for above 100kg) ○ |
| CHARGER* | | |
| ATE1343 | <input type="checkbox"/> | OptiCharge Charger 8A (200 - 250V) ◆ |
| LIGHTS* | | |
| ATZ1006 | <input type="checkbox"/> | Lights Included with Motor Package 10 & 12kmh ◆ |
| ANTI-TIPPERS* | | |
| ATE1205 | <input type="checkbox"/> | Anti-tipper (standard version) ◆ |
| TRANSPORTATION* | | |
| ATE1494 | <input type="checkbox"/> | 4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory). ◆ |
| FIXATION BELT* | | |
| ATE2113 | <input type="checkbox"/> | Lapbelt (adjustable both sides) adjustable mounting on seat rail extrusion ◆ |
| REAR VIEW MIRROR | | |
| ATE1415 | <input type="checkbox"/> | Rear view Mirror RH - Must be selected with 10 km/h Motor ○ |
| ATE1339 | <input type="checkbox"/> | Rear view Mirror LH - Must be selected with 10 km/h Motor ○ |

| | |
|--------------------|--|
| TOTAL PRICE | |
|--------------------|--|

Non Conforming Options:

ATE2196 Long Armpad (375mm)



For a more complex customer requirement
please refer to our Invacare Unique special order form.
email specialsuk@invacare.com

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):