



Yes, you can:

Motability Retail Prescription Form

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

APRIL 2025
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

AVIVA RX40 Modulite M24 707026 - weekly
Price £53.95

Standard Features

- Max user weight 136kg
- Crash tested ISO 7176 -Part 19
- Modulite seating system
- Tension adjustable backrest
- Standard T-Armrests, detachable with long armpad
- Manually settable Seat Tilt 0 - 15°
- Manual settable backrest Recline (0°- 30°, in 3.75°steps)
- 80° hangers with fixed footplate
- LiNX Control REM216 with swing away Bracket
- 9" Castor, puncture proof solid tyre
- Solid inlay Dual density foam insert
- 12 km/h Ultra High Speed Package
- LED Light system with Indicators, 30 Lux headlights
- 73,5 Ah Gel Batteries
- Anti tippers
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)



Standard Base price: £7,424

Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

Code	Description	
BASE MODEL*		
JBF0001	AVIVA RX 40 Base Model Max user weight up to 136 kg includes Performance Package: 120A Power Module + minimum 60 Ah Batteries	£7,424
CHASSIS WIDTH*		
JBL2810	Chassis Narrow (560 mm/22")	◆
SEATING TYPE*		
JBF2100	Invacare® Modulite Seat Unit	◆

SEAT FRAME*			
JBF2111	<input type="checkbox"/>	Telescopic Seat Frame and Plate with Side Rails width 380 to 530 mm Additional 50 mm by arms	<input type="radio"/>
JBF2112	<input type="checkbox"/>	Telescopic Seat Frame + Sling Seat with Side Rails width 380 to 530 mm; + 30 mm seat height Additional 50 mm by arms	<input type="radio"/>
SEAT WIDTH*			
JBF2121	<input type="checkbox"/>	380 mm Seat Plate Can be reduced to 330 mm in combination with inrigged armrests	<input type="radio"/>
JBF2122	<input type="checkbox"/>	430 mm Seat Plate Can be reduced to 380 mm in combination with inrigged armrests	<input type="radio"/>
JBF2123	<input type="checkbox"/>	480 mm Seat Plate	<input type="radio"/>
JBF2124	<input type="checkbox"/>	530 mm Seat Plate	<input type="radio"/>
SEAT DEPTH*			
JBF2132	<input type="checkbox"/>	Seat Depth 410 mm	<input type="radio"/>
JBF2134	<input type="checkbox"/>	Seat Depth 460 mm	<input type="radio"/>
JBF2136	<input type="checkbox"/>	Seat Depth 510 mm	<input type="radio"/>
SEAT CUSHION*			
JBF2142	<input type="checkbox"/>	Matrx Flo-tech Cushion Flo-tech Lite	<input checked="" type="radio"/>
BACKREST SUPPORT*			
JBF2163	<input type="checkbox"/>	Backposts short (480 mm) For Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery Width adjustable cast bracket	<input checked="" type="radio"/>
JBF2162	<input type="checkbox"/>	Backposts tall (540 mm) For Flex3 backrests/Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery Width adjustable cast bracket	<input type="radio"/>
BACKREST UPHOLSTERY / CUSHION*			
JBF2171	<input type="checkbox"/>	Tension Adjustable Upholstery nylon overlay short = 480 mm / 7 straps measured from seat plate to top of upholstery 480 mm	<input checked="" type="radio"/>
UPHOLSTERY TYPE*			
JBF2151	<input type="checkbox"/>	Black Fabric Upholstery (Nylon w 2141; 2-way-stretch material w 2142, 2143)	<input checked="" type="radio"/>
ARMREST (Pair)*			
JBF2191	<input type="checkbox"/>	T-Arm detachable, Standard Adjustment range from seat plate: 245 - 315 mm	<input checked="" type="radio"/>
ARMPAD (Pair)*			
JBF2196	<input type="checkbox"/>	Long Armpad (375 mm)	<input checked="" type="radio"/>
MANUAL / POWERED SEATING OPTIONS*			
JBF1091	<input type="checkbox"/>	Manually settable Seat Tilt 0 - 15°	<input checked="" type="radio"/>
SEAT TO GROUND HEIGHT*			
JBF2857	<input type="checkbox"/>	Seat-to-Floor Height 435 mm	<input type="radio"/>
JBF2858	<input type="checkbox"/>	Seat-to-Floor Height 460 mm	<input type="radio"/>
JBF2859	<input type="checkbox"/>	Seat-to-Floor Height 485 mm	<input type="radio"/>
BACKREST RECLINER*			
JBF0975	<input type="checkbox"/>	Manual Settable Recliner (in 3.75° steps) 90° - 120°	<input checked="" type="radio"/>
LEG REST TYPE*			
JBF0526	<input type="checkbox"/>	Standard hanger 80° pair Lower leg length: 350 - 450 mm	<input checked="" type="radio"/>
FOOTPLATE TYPE*			
JBF0512	<input type="checkbox"/>	Fixed Footplate pair	<input checked="" type="radio"/>

LINX CONTROLS*			
JBF2204	<input type="checkbox"/>	Remote LiNX DLX-REM216 = Drive, Actuators & Lights 3 Drive functions and power module incl. 2 Actuators & Lights	◆
CONTROL BRACKETRY*			
JBF1707	<input type="checkbox"/>	Swing Away Bracket RH	○
JBF1708	<input type="checkbox"/>	Swing Away Bracket LH	○
SHROUD COLOUR*			
JBF2871	<input type="checkbox"/>	Mid range - Black/Anthracite Shroud	○
JBF2872	<input type="checkbox"/>	Mid range - Silver Shroud	○
JBF2873	<input type="checkbox"/>	Mid range - Bengal Red Shroud	○
JBF2874	<input type="checkbox"/>	Mid range - Mercury Blue Shroud	○
JBF2875	<input type="checkbox"/>	Mid range - Midnight Purple Shroud	○
RIM COLOUR*			
JBF2882	<input type="checkbox"/>	Mid range - Black Rim Inserts	○
JBF2883	<input type="checkbox"/>	Mid range - Silver Rim Inserts	○
JBF2884	<input type="checkbox"/>	Mid range - Bengal Red Rim Inserts	○
JBF2885	<input type="checkbox"/>	Mid range - Mercury Blue Rim Inserts	○
JBF2886	<input type="checkbox"/>	Mid range - Midnight Purple Rim Inserts	○
FORKS*			
JBF2632	<input type="checkbox"/>	Double sided fork Without castor rim insert (enlarges wide base to 656 mm)	◆
FRONT WHEEL*			
JBF0836	<input type="checkbox"/>	9" (225 mm) x 2.80 (70 mm) Castor, puncture proof solid tyre	◆
DRIVE WHEEL*			
JBF0705	<input type="checkbox"/>	DW 14" (355 mm) x 3" (75 mm), Solid inlay Dual density foam insert	◆
C.T.C SUSPENSION*			
JBF2870	<input type="checkbox"/>	CTC Suspension Control - Traction - Comfort	◆
SUSPENSION SETTING*			
JBF1910	<input type="checkbox"/>	Soft (suggested for up to 70 kg User Weight)	○
JBF1912	<input type="checkbox"/>	Medium (suggested for up to 90 kg User Weight)	○
JBF1911	<input type="checkbox"/>	Firm (suggested for up to 110 kg User Weight)	○
JBF1915	<input type="checkbox"/>	Very firm (suggested for up to 136 kg User Weight)	○
MOTOR TYPE*			
JBF1909	<input type="checkbox"/>	12 km/h Ultra High Speed Package Always combined with LED Light system with Indicators, 30 Lux headlights	◆
BATTERY TYPE*			
JBF1260	<input type="checkbox"/>	73,5 Ah Gel	◆
CHARGER*			
JBF1343	<input type="checkbox"/>	Opti Charge Charger 8A (IP44)	◆
ANTI-TIPPERS*			
JBF1205	<input type="checkbox"/>	Anti-tipper	◆
TRANSPORTATION KIT*			
JBF1494	<input type="checkbox"/>	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19) Adequate headrest is compulsory	◆
LAP BELT*			
JBF1329	<input type="checkbox"/>	Lap belt with Auto Lock Standard mount on backrest	◆

REAR VIEW MIRROR*

JBF1415	<input type="checkbox"/>	Rear view mirror RH	<input type="radio"/>
JBF1339	<input type="checkbox"/>	Rear view mirror LH	<input type="radio"/>

TOTAL PRICE

Non Conforming Options:



For a more complex customer requirement
please refer to our Invacare Unique special order form.
email specialsuk@invacare.com

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):