



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

APRIL 2025
Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

**Erivo R10 M38 1001072 - Weekly Price
£33.78**

Standard Features

- ERIVO R10
- 11Ah Battery
- 17KM range
- Max user weight 150kg (23.5 Stone)
- 6kmh
- High quality foam seat cushion
- Solid seat plate
- Padded backrest with large built in pocket
- Fold-up armrests
- Fold up and height adjustable footplate (3 positions)
- LiNX REM060 remote with simple swivel arm
- 8" PU front castor
- Rear wheel suspension with 12.5" drive wheel, solid tyre



Standard Base price: £2,900

Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

Code	Description	
ERIVO*		
1595187	<input type="checkbox"/> R10, 11Ah Battery, 17KM range	£2,900

TOTAL PRICE	
--------------------	--