



**Yes, you can:**

Tel: 01656 776222  
Email: [ordersuk@invacare.com](mailto:ordersuk@invacare.com)  
Online Spares available at: [www.invacare.co.uk](http://www.invacare.co.uk)

# Motability Retail Prescription Form

APRIL 2025  
Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## Storm4 Modulite M32 1001076 - Weekly Price £62.35

### Standard Features

- Max user weight 150Kg (23.5 stone)
- Crash tested ISO 7176 -Part 19
- Modulite seating system
- Flex3 Laguna Backrest
- Standard T-Armrests, detachable with long armpad
- Manual settable seat Tilt 0° - 20° Electrical Seat Riser
- Manual settable Recliner (0° - 30° in 3.75° steps)
- VARI F with Fixed Footplate
- LiNX Control REM216 with swing away bracket
- Pneumatic Castors
- Pneumatic Drive Wheels
- 73.5 Ah Batteries
- 12 km/h (4mph) Motors
- LED Light system with indicators
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)
- Anti-Tippers



**Standard Base price: £8,950**

Please make your selection with a cross in the box	<input type="radio"/> Free of charge	<input checked="" type="radio"/> standard feature	<input type="radio"/> POA	<input type="radio"/> special request
--	--------------------------------------	---	---------------------------	---------------------------------------

Customer weight =                      KG

Code	Description	
<b>CHAIR TYPE*</b>		
ATB0000	Storm 4 Modulite	£8,950
ATZ0000	Special Storm 4 Modulite (Please provide details in special section)	POA
<b>SEATING*</b>		
ATB2100	Modulite Seating Unit (Flex3 seating unit)	◆
<b>SEATING OPTION*</b>		
ATB2111	Telescopic Seat Frame with metal seat plate (with stepless adjustment on armrests, legrests and backrest)	◆

SEAT WIDTH*			
ATBXXXX	<input type="checkbox"/>	Seat width 330mm (13") seat width achieved on a 380mm seat frame (2121) with inrigged armrests (2183)	<input type="radio"/>
ATB2121	<input type="checkbox"/>	Set at 380mm (adjustable from 380mm to 430mm (15" - 17") with stepless adjustment on armrests, legrests and backrest)	<input type="radio"/>
ATB2122	<input type="checkbox"/>	Set at 430mm (adjustable from 430mm to 480mm (17" - 19) with stepless adjustment on armrests, legrests and backrest)	<input type="radio"/>
ATB2123	<input type="checkbox"/>	Set at 480mm (adjustable from 480mm to 530mm (19" - 21") with stepless adjustment on armrests, legrests and backrest)	<input type="radio"/>
ATB2124	<input type="checkbox"/>	Set at 530mm (adjustable from 530mm to 580mm (21" - 23") with stepless adjustment on armrests, legrests and backrest)	<input type="radio"/>
SEAT DEPTH*			
ATB2132	<input type="checkbox"/>	410mm (16") (adjustable 410mm - 510mm stepless adjustment) (NC 2196)	<input type="radio"/>
ATB2134	<input type="checkbox"/>	460mm (18") (adjustable 410mm - 510mm stepless adjustment)	<input type="radio"/>
ATB2136	<input type="checkbox"/>	510mm (20") (adjustable 410mm - 510mm stepless adjustment)	<input type="radio"/>
SEAT CUSHION*			
ATB2142	<input type="checkbox"/>	Matrx Flo-tech Lite	<input checked="" type="checkbox"/>
BACKREST SUPPORT*			
ATB2162	<input type="checkbox"/>	Backposts 540mm (Tall) must be selected for Flex3 Laguna backrest	<input checked="" type="checkbox"/>
PUSH HANDLES FOR BACKPOSTS*			
ATB2107	<input type="checkbox"/>	Backposts WITH push handles	<input checked="" type="checkbox"/>
BACKREST UPHOLSTERY / CUSHION*			
ATB2173	<input type="checkbox"/>	Flex3 adjustable Backrest Plates with Laguna Backrest cushion	<input checked="" type="checkbox"/>
UPHOLSTERY TYPE*			
ATB2152	<input type="checkbox"/>	Black multi stretch PU coated fabric	<input type="radio"/>
ATB2155	<input type="checkbox"/>	Black Spacer fabric	<input type="radio"/>
ARMREST (PAIR)*			
ATB2191	<input type="checkbox"/>	T-Armrest Standard, Range of adjustment from seat plate 245mm to 315mm	<input checked="" type="checkbox"/>
ARMPAD (PAIR)*			
ATB2196	<input type="checkbox"/>	Long Armpad (375mm) (NC 2132)	<input checked="" type="checkbox"/>
HEADREST / NECKREST (Adequate headrest option must be selected)*			
ATB1350	<input type="checkbox"/>	neckrest: height and depth adjustable	<input checked="" type="checkbox"/>
SEAT HEIGHT AND SEAT RISER*			
ATB1023	<input type="checkbox"/>	Fixed Seat Height 450mm (17.75")	<input type="radio"/>
ATB1025	<input type="checkbox"/>	Fixed Seat Height 480mm (19")	<input type="radio"/>
ATB1030	<input type="checkbox"/>	Electrical Seat Riser (400mm - 650mm) (16" - 25.5")	<input type="radio"/>
ATB1040	<input type="checkbox"/>	Electrical Seat Riser (440mm - 690mm) incl. adapter (17.25" - 27")	<input type="radio"/>
SEAT TILT*			
ATB1002	<input type="checkbox"/>	Manual Seat Tilt by spindle	<input checked="" type="checkbox"/>
BACKREST RECLINER*			
ATB0975	<input type="checkbox"/>	Manual Settable Recliner (0° - 30°, in 3.75° steps) Foldable	<input checked="" type="checkbox"/>
LEG REST TYPE*			
ATB0503	<input type="checkbox"/>	Vari F : Angle adjustable with tools 0° - 70°, length 290mm - 460mm	<input checked="" type="checkbox"/>

<b>FOOTPLATE TYPE*</b>			
ATB0512	<input type="checkbox"/>	Fixed Footplate	<input type="radio"/>
<b>LINX CONTROLS*</b>			
ATB2204	<input type="checkbox"/>	REM216 Drive, Actuators & Lights	<input checked="" type="radio"/>
<b>CONTROL BRACKETRY*</b>			
ATB1707	<input type="checkbox"/>	Swing Away RH Bracket	<input type="radio"/>
ATB1708	<input type="checkbox"/>	Swing Away LH Bracket	<input type="radio"/>
<b>SHROUD COLOUR*</b>			
ATB0202	<input type="checkbox"/>	Blue	<input type="radio"/>
ATB0228	<input type="checkbox"/>	Black	<input type="radio"/>
ATB0204	<input type="checkbox"/>	Transparent Cherry Red	<input type="radio"/>
<b>FRONT WHEEL*</b>			
ATB0851	<input type="checkbox"/>	3.00 - 6 Grey Pneumatic	<input checked="" type="radio"/>
<b>REAR WHEEL*</b>			
ATB0771	<input type="checkbox"/>	3.00 - 8 Grey Pneumatic	<input checked="" type="radio"/>
<b>MOTOR TYPE*</b>			
ATB1909	<input type="checkbox"/>	12 km/h (8mph) (includes LED light system with indicators)	<input checked="" type="radio"/>
<b>BATTERY*</b>			
ATB1364	<input type="checkbox"/>	70 Ah Gel	<input checked="" type="radio"/>
<b>SUSPENSION*</b>			
ATB1240	<input type="checkbox"/>	Suspension low (up to 80kg user weight)	<input type="radio"/>
ATB1241	<input type="checkbox"/>	Suspension middle (from 70kg up to 100kg user weight)	<input type="radio"/>
ATB1242	<input type="checkbox"/>	Suspension hard (from 90kg up to 150kg user weight)	<input type="radio"/>
<b>CHARGER*</b>			
ATB1343	<input type="checkbox"/>	OptiCharge Charger 8A (200 - 250V)	<input checked="" type="radio"/>
<b>ANTI TIPPER*</b>			
ATB1205	<input type="checkbox"/>	Anti-tipper (standard version)	<input type="radio"/>
<b>TRANSPORTATION KIT*</b>			
ATB1494	<input type="checkbox"/>	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory, see options below)	<input checked="" type="radio"/>
<b>FIXATION BELT*</b>			
ATB2113	<input type="checkbox"/>	Lap belt with auto lock (flexible mount in side rails)	<input checked="" type="radio"/>
<b>REAR VIEW MIRROR</b>			
ATB1415	<input type="checkbox"/>	Rear view Mirror RH - Must be selected with 10 km/h Motor	<input type="radio"/>
ATB1339	<input type="checkbox"/>	Rear view Mirror LH - Must be selected with 10 km/h Motor	<input type="radio"/>

<b>TOTAL PRICE</b>			
--------------------	--	--	--

**Non Conforming Options:**

ATB2132	410mm (16") (adjustable 410mm - 510mm stepless adjustment)
ATB2196	Long Armpad (375mm)



For a more complex customer requirement  
please refer to our Invacare Unique special order form.  
email [specialsuk@invacare.com](mailto:specialsuk@invacare.com)

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):