



Yes, you can:

Tel: 01656 776222
Email: ordersuk@invacare.com
Online Spares available at: www.invacare.co.uk

Motability Retail Prescription Form

APRIL 2025
Customer Ref:

Customer Ref*, Contact No*, Delivery Address*, Account No*, Order No*, Order Qty*, Order date*

TDX SP2 Low Rider M14 712116 - Weekly
Price £60.00

Standard Features

- Max user weight 160kg (25.2 Stone)
Crash tested ISO7176-Part 19
Modulite seating
Matrix Flo-tech Lite cushion
Without push handles
Flex 3 Adjustable backrest plates with laguna backrest cushion
Elan Headrest
Standard T-armrests, detachable with Long armpad
Seat Lifter with powered Tilt (300mm, 30°) CoG Shift
Seat height either 403mm
Manual settable Recliner 0° - 30° in 3.75° steps
Standard 80° Hangers with Angle and Depth Adjustable footplates(adapt)
LiNX control REM211 with swing away bracket
Solid Black tyres
6km/h (4mph) 4 Pole motors
4-Point tie-down system
Lapbelt with autolock (adjustable both sides)

Standard Base price: £8,358



Please make your selection with a cross in the box
Free of charge
standard feature
POA
special request

Customer weight = KG

Table with columns: Code, Description, Selection. Rows include BASE PRODUCT* (ARQ0072) and USER WEIGHT* (ARQ1918-ARQ1923).

SEATING TYPE*		
ARQ2100	<input type="checkbox"/>	Modulite Seating Unit £8,358
SEAT FRAME*		
ARQ2111	<input type="checkbox"/>	Telescopic Seat Frame with metal seat plate ◆
SEAT WIDTH*		
ARQ2122	<input type="checkbox"/>	Set at 430mm (17") (adjustable from 430mm to 480mm (17" - 19)) ○
ARQ2123	<input type="checkbox"/>	Set at 480mm(19") (adjustable from 480mm to 530mm (19" - 21")) ○
ARQ2124	<input type="checkbox"/>	Set at 530mm (21") (adjustable from 530mm to 580mm (21" - 23")) ○
SEAT DEPTH RAIL LENGTH*		
ARQ2129	<input type="checkbox"/>	Standard Depth Package (16" - 20", 410mm - 510mm) ○
ARQ2130	<input type="checkbox"/>	Extended Depth Package (20" - 22", 510mm - 560mm) (NC 2134) ○
SEAT DEPTH*		
ARQ2132	<input type="checkbox"/>	410mm (16") (410mm - 510mm) (NC 2196) ○
ARQ2134	<input type="checkbox"/>	460mm (18") (410mm - 510mm) ○
ARQ2136	<input type="checkbox"/>	510mm (20") (410mm - 510mm) ○
ARQ2138	<input type="checkbox"/>	560mm (22") (510mm - 560mm) ○
SEAT CUSHION*		
ARQ2142	<input type="checkbox"/>	Matrx Flo-tech Lite ◆
BACKREST SUPPORT*		
ARQ2162	<input type="checkbox"/>	Backposts 540mm (Tall) must be selected for Flex3 adjustable backrest plates ◆
BACKREST UPHOLSTERY / CUSHION*		
ARQ2173	<input type="checkbox"/>	Flex3 Adjustable Backrest Plates with Laguna Backrest cushion ◆
UPHOLSTERY TYPE*		
ARQ2152	<input type="checkbox"/>	Black multi stretch PU coated fabric ○
ARQ2155	<input type="checkbox"/>	Black Spacer fabric ○
ARMREST*		
ARQ2191	<input type="checkbox"/>	T-Armrest Standard, Range of adjustment ◆
ARMPAD*		
ARQ2196	<input type="checkbox"/>	Long Armpad (375mm) ◆
HEADREST / NECKREST*		
ARQ2324	<input type="checkbox"/>	Elan Standard Pad 355mm (10") wide ◆
MANUAL / POWERED SEATING OPTIONS*		
ARQ1094	<input type="checkbox"/>	Seat lifter with Electrical Tilt ◆
BACKREST RECLINER*		
ARQ0975	<input type="checkbox"/>	Manual Settable Recliner ◆
LEG REST TYPE*		
ARQ0526	<input type="checkbox"/>	Side Mounted : Standard 80° hangers ◆
FOOTPLATE TYPE*		
ARQ0512	<input type="checkbox"/>	Fixed Footplate ◆
LINX CONTROLS*		
ARQ2203	<input type="checkbox"/>	REM211 Drive & Actuators ◆
CONTROL BRACKETRY*		
ARQ1707	<input type="checkbox"/>	Swing Away RH Bracket ○
ARQ1708	<input type="checkbox"/>	Swing Away LH Bracket ○

SHROUD COLOUR*			
ARQ2640	<input type="checkbox"/>	Obsidian Black Frame Shrouds	<input type="radio"/>
ARQ2641	<input type="checkbox"/>	Ruby Red Frame Shrouds	<input type="radio"/>
ARQ2642	<input type="checkbox"/>	Purple Passion Frame Shrouds	<input type="radio"/>
ARQ2643	<input type="checkbox"/>	Peacock Blue Frame Shrouds	<input type="radio"/>
ARQ2644	<input type="checkbox"/>	Arctic Blue Frame Shrouds	<input type="radio"/>
RIM COLOUR*			
ARQ2650	<input type="checkbox"/>	Obsidian Black Rim Inserts	<input type="radio"/>
ARQ2651	<input type="checkbox"/>	Ruby Red Rim Inserts	<input type="radio"/>
ARQ2652	<input type="checkbox"/>	Purple Passion Rim Inserts	<input type="radio"/>
ARQ2653	<input type="checkbox"/>	Peacock Blue Rim Inserts	<input type="radio"/>
ARQ2654	<input type="checkbox"/>	Arctic Blue Rim Inserts	<input type="radio"/>
FORKS*			
ARQ2632	<input type="checkbox"/>	Double sided fork	<input type="checkbox"/>
CENTRE WHEELS (DRIVE WHEELS)*			
ARQ0705	<input type="checkbox"/>	DW 14" (355 mm) x 3" (75 mm), Solid inlay, dual density foam insert	<input type="checkbox"/>
MOTOR TYPE*			
ARQ1282	<input type="checkbox"/>	6 km/h (4mph) (4 pole)	<input type="checkbox"/>
CHARGER*			
ARQ1343	<input type="checkbox"/>	OptiCharge Charger 8A (200 - 250V)	<input type="checkbox"/>
LIGHTS*			
ARQ1254	<input type="checkbox"/>	Without Light system	<input type="checkbox"/>
TRANSPORTATION KIT*			
ARQ1494	<input type="checkbox"/>	4-Point tie-down	<input type="checkbox"/>
FIXATION BELT*			
ARQ2113	<input type="checkbox"/>	Lapbelt (adjustable both sides)	<input type="checkbox"/>

TOTAL PRICE	
--------------------	--

Non Conforming Options:

ARQ2134 460mm (18") (410mm - 510mm)

ARQ2196 Long Armpad (375mm)



For a more complex customer requirement
please refer to our Invacare Unique special order form.
email specialsuk@invacare.com

Special Instructions
Please supply the following information:

Name: _____ Email: _____

Phone: _____

Description of special request (please include relevant option code or part number):

