



**Yes, you can:**

Tel: 01656 776222  
Email: [ordersuk@invacare.com](mailto:ordersuk@invacare.com)  
Online Spares available at: [www.invacare.co.uk](http://www.invacare.co.uk)

# Motability Retail Prescription Form

APRIL 2025  
Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

**TDX SP2 Max M15 1001086 - Weekly Price  
£48.21**

## Standard Features

- Max user weight 180kg (28 stone)
- Crash tested ISO 7176 -Part 19
- Tension Adjustable backrest
- Armrest with reinforced depth adjustable bracket with Short armpad
- Powered 20 degree seat Tilt
- Manual settable backrest Recline
- VARI F hangers with fixed footplate
- LiNX Control REM211 with swing away bracket
- Solid Black Tyres
- 6 km/h (4mph) 4 Pole Motors
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)

**Standard Base price: £5,882**



Please make your selection with a cross in the box    ☐ Free of charge    ☒ standard feature    POA special request

Customer weight =                      KG

Code	Description	
<b>CHAIR TYPE*</b>		
ARH0000	TDX SP2 Max	£5,882
ARZ0000	Special TDX SP2 Max (Please provide details in special section)	POA
<b>SEATING TYPE*</b>		
ARH1145	MAX Seating system	<input checked="" type="radio"/>
<b>SEAT WIDTH SETTING*</b>		
ARH0378	Seat Width of Seat Frame 550mm (+100mm adjustment in arms)	<input type="radio"/>
ARH0380	Seat Width of Seat Frame 610mm (+100mm adjustment in arms)	<input type="radio"/>

SEAT DEPTH SETTING SETTABLE IN 3 STEPS*			
ARH0396	<input type="checkbox"/>	Seat Depth 560mm	<input type="radio"/>
ARH0397	<input type="checkbox"/>	Seat Depth 595mm	<input type="radio"/>
ARH0398	<input type="checkbox"/>	Seat Depth 630mm	<input type="radio"/>
SEAT CUSHION*			
ARH1484	<input type="checkbox"/>	MAX Seat Cushion	<input type="radio"/>
ADAPTOR TO REDUCE THE SEAT DEPTH BY 130MM			
ARH0399	<input type="checkbox"/>	Adaptor to shorten the seat depth by 130mm (Required if seat width 490mm is selected in combination with seat depth 595mm or 630mm)	<input type="radio"/>
BACKREST UPHOLSTERY*			
ARH1106	<input type="checkbox"/>	Tension Adjustable	<input checked="" type="radio"/>
UPHOLSTERY COLOUR*			
ARH0278	<input type="checkbox"/>	Black multi stretch PU coated fabric	<input checked="" type="radio"/>
ARMREST INCLUDING REINFORCED DEPTH ADJUSTABLE BRACKET*			
ARH0635	<input type="checkbox"/>	Fixed Armrest, 250mm - 340mm tool adjustable in height, detachable	<input checked="" type="radio"/>
ARMPADS*			
ARH0649	<input type="checkbox"/>	Long armpad (375mm)	<input checked="" type="radio"/>
PUSH BAR*			
ARH1307	<input type="checkbox"/>	Push bar	<input checked="" type="radio"/>
HEADREST (Adequate headrest option must be selected)*			
ARH1432	<input type="checkbox"/>	Headrest: height, depth and sideways adjustable (Azalea type)	<input checked="" type="radio"/>
MANUAL / POWERED SEATING OPTIONS*			
ARH1095	<input type="checkbox"/>	Electric heavy duty basic Tilt (0° - 20°) select either 470mm or 495mm seat to ground height. (Max user weight 180kg)	<input checked="" type="radio"/>
SEAT HEIGHT SETTING*			
ARH1024	<input type="checkbox"/>	Fixed seat height 470mm (18.5")	<input type="radio"/>
ARH1046	<input type="checkbox"/>	Fixed seat height 495mm (19.5")	<input type="radio"/>
BACKREST RECLINER*			
ARH0958	<input type="checkbox"/>	Manually settable recliner by spindle (0° - 25°)	<input checked="" type="radio"/>
LEG REST TYPE*			
ARH0503	<input type="checkbox"/>	Side Mounted : Vari F : Angle adjustable with tools 0° - 70°, length 290mm - 460mm	<input checked="" type="radio"/>
FOOTPLATE TYPE*			
ARH0512	<input type="checkbox"/>	Fixed Footplate	<input checked="" type="radio"/>
LINX CONTROLS*			
ARH2203	<input type="checkbox"/>	REM211 Drive & Actuators	<input checked="" type="radio"/>
CONTROL BRACKETRY*			
ARH1707	<input type="checkbox"/>	Swing Away RH Bracket	<input type="radio"/>
ARH1708	<input type="checkbox"/>	Swing Away LH Bracket	<input type="radio"/>
SHROUD COLOUR*			
ARH2640	<input type="checkbox"/>	Obsidian Black Frame Shrouds	<input type="radio"/>
ARH2641	<input type="checkbox"/>	Ruby Red Frame Shrouds	<input type="radio"/>
ARH2642	<input type="checkbox"/>	Purple Passion Frame Shrouds	<input type="radio"/>
ARH2643	<input type="checkbox"/>	Peacock Blue Frame Shrouds	<input type="radio"/>
ARH2644	<input type="checkbox"/>	Arctic Blue Frame Shrouds	<input type="radio"/>

RIM COLOUR*			
ARH2650	<input type="checkbox"/>	Obsidian Black Rim Inserts	<input type="radio"/>
ARH2651	<input type="checkbox"/>	Ruby Red Rim Inserts	<input type="radio"/>
ARH2652	<input type="checkbox"/>	Purple Passion Rim Inserts	<input type="radio"/>
ARH2653	<input type="checkbox"/>	Peacock Blue Rim Inserts	<input type="radio"/>
ARH2654	<input type="checkbox"/>	Arctic Blue Rim Inserts	<input type="radio"/>
FORKS*			
ARH2632	<input type="checkbox"/>	Double sided fork	<input checked="" type="checkbox"/>
CENTRE WHEELS (DRIVE WHEELS)*			
ARH0705	<input type="checkbox"/>	DW 14" (355 mm) x 3" (75 mm), Solid inlay, dual density foam	<input checked="" type="checkbox"/>
MOTOR TYPE*			
ARH1282	<input type="checkbox"/>	6 km/h (4mph) (4 pole)	<input checked="" type="checkbox"/>
BATTERY*			
ARH1260	<input type="checkbox"/>	73.5 Ah	<input checked="" type="checkbox"/>
CHARGER*			
ARH1343	<input type="checkbox"/>	OptiCharge Charger 8A (200 - 250V)	<input checked="" type="checkbox"/>
LIGHTS*			
ARH1254	<input type="checkbox"/>	Without light system	<input checked="" type="checkbox"/>
TRANSPORTATION KIT			
ARH1494	<input type="checkbox"/>	4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory, see options below)	<input checked="" type="checkbox"/>
FIXATION BELT*			
ARH1329	<input type="checkbox"/>	MAX Lapbelt (adjustable both sides) fitted to base of backrest	<input checked="" type="checkbox"/>

TOTAL PRICE			
-------------	--	--	--

Non Conforming Options:



For a more complex customer requirement  
please refer to our Invacare Unique special order form.  
email [specialsuk@invacare.com](mailto:specialsuk@invacare.com)

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):