



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
 Email: ordersuk@invacare.com
 Online Spares available at: www.invacare.co.uk

Power Retail Prescription Form

SEPTEMBER 2024
 Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

AVIVA RX40 HD

Standard Features

- Max user weight 160kg (25 Stone)
- Crash tested ISO 7176 -Part 19
- Modulite adjustable seating system
- Telescopic adjustable seat widths from 330mm - 530mm (13" - 21") with an additional 50mm, (2") by arms
- Seat depths 410mm - 510mm stepless adjustment (16" - 20")
- Tension adjustable backrest
- Standard T-Armrests, detachable with short armpad
- Manual settable seat Tilt (0° - 20°)
- Manual settable backrest Recline (0° - 30° in 3.75° steps)
- 80° hangers with fixed footplate
- LiNX control REM110 with standard bracket
- Choice of 5 Shroud colours, 5 Rim insert colours
- 9" (225 mm) x 2.80 (70 mm) Castor, puncture proof solid tyre
- 14" (355 mm) x 3" (75 mm), Solid inlay black drive wheels
- CTC Suspension Control - Traction - Comfort
- 6 km/h (4mph) 4 pole and 60Ah batteries
- 8A (200 - 250V) charger
- Anti tippers
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)

Standard Base price: £6,234



Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

Code	Description	
BASE MODEL*		
JBH0002	AVIVA RX 40 HD Base Model Max user weight up to 160 kg HD Version, includes RX 40 Performance Package Wide Base only	(NC 2193) £6,234
JBX1300	Special Extended lead time: + 10 days	POA
CHASSIS WIDTH*		
JBH2830	Chassis Wide (635 mm/25")	◆

SEATING TYPE*			
JBH2100	<input type="checkbox"/>	Invacare® Modulite Seat Unit	○
SEAT FRAME*			
JBH2111	<input type="checkbox"/>	Telescopic Seat Frame and Plate with Side Rails width 380 to 530 mm Additional 50 mm by arms	○
JBH2112	<input type="checkbox"/>	Telescopic Seat Frame + Sling Seat with Side Rails width 380 to 530 mm; + 30 mm seat height Additional 50 mm by arms	○
SEAT WIDTH*			
JBH2121	<input type="checkbox"/>	380 mm Seat Plate Can be reduced to 330 mm in combination with inrigged armrests	(NC 2174, 2175, 2176) ○
JBH2122	<input type="checkbox"/>	430 mm Seat Plate Can be reduced to 380 mm in combination with inrigged armrests	○
JBH2123	<input type="checkbox"/>	480 mm Seat Plate	(NC 2702, 2704, 2706) ○
JBH2124	<input type="checkbox"/>	530 mm Seat Plate	(NC 2174, 2175, 2176, 2177, 2178, 2702, 2704, 2706, 2106) ○
Body Fit Kit			
JBH2106	<input type="checkbox"/>	Body fit Kit Outrigging Backrest Brackets and Pivot Point (outrig: 1 inch/2,54 cm per side)	○
SEAT DEPTH*			
JBH2132	<input type="checkbox"/>	Seat Depth 410 mm	(NC 2196, 2197, 2198) ○
JBH2134	<input type="checkbox"/>	Seat Depth 460 mm	○
JBH2136	<input type="checkbox"/>	Seat Depth 510 mm	○
SEAT CUSHION*			
JBH2142	<input type="checkbox"/>	Matrx Flo-tech Cushion Flo-tech Lite	£117
JBH2147	<input type="checkbox"/>	Matrx Libra Cushion, 90 mm thick Extended lead time: + 4 days	£388
JBH2140	<input type="checkbox"/>	Without seat cushion	○
BACKREST SUPPORT*			
JBH2163	<input type="checkbox"/>	Backposts short (480 mm) For Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery Width adjustable cast bracket	(NC 2172, 2173) ○
JBH2162	<input type="checkbox"/>	Backposts tall (540 mm) For Flex3 backrests/Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery Width adjustable cast bracket	◆
PUSH HANDLES FOR BACKPOSTS			
JBH2108	<input type="checkbox"/>	Backposts without push handles	(NC 2219, 2220, 2227, 2228, 2245, 2291, 2292) ○

BACKREST UPHOLSTERY / CUSHION*			
JBH2171		Tension Adjustable Upholstery nylon overlay short = 480 mm / 7 straps measured from seat plate to top of upholstery 480 mm	(NC 2155, 2230, 2247) ◆
JBH2179		Flex3 Tension Adjustable Upholstery height 540 mm (seat plate to top of upholstery)	(NC 2151, 2192, 2193, 2194, 2199, 2155, 2230, 2247) £154
JBH2172		Flex3 adjustable Backrest Plates (without Backrest Cushion) Flex3 Backrest plate width and height adjustable 560 mm - 790 mm	(NC 2163, 2128) £250
JBH2173		Flex3 adjustable Backrest Plates with Laguna Backrest cushion Flex3 Backrest plate width and height adjustable 560 mm - 790 mm	(NC 2163, 2128) £406
JBH2174		Matrx Posture Back, height 400 mm (Mesh fabric) Corresponds to SW setting chosen (up to 136 kg max user weight)	(NC 2121) £557
JBH2175		Matrx Posture Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen (up to 136 kg max user weight)	(NC 2121) £557
JBH2176		Matrx Elite E2 Back , height 500 mm (Mesh fabric) Corresponds to SW setting chosen (up to 136 kg max user weight)	(NC 2121) £557
JBH2177		Matrx Elite TR Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen	£557
JBH2178		Matrx Elite E2 Deep Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen (up to 136 kg max user weight)	£647
JBH2170		Without Backrest Upholstery (Back Posts/Back Frame only)	(NC 2324, 2325, 2327, 2230, 2247) ○
UPHOLSTERY TYPE*			
JBH2151		Black Fabric Upholstery (Nylon with option 2141; 2- way-stretch material with options 2142, 2143)	◆
JBH2152		Black multi stretch PU coated fabric Upholstery	(NC 2171) ○
JBH2155		Black Spacer Fabric Upholstery	(NC 2171) ○
ARMREST (Pair)*			
JBH2191		T-Arm detachable, Standard Adjustment range from seat plate: 245 - 315 mm	◆
JBH2181		T-Arm detachable, Extended height Adjustment range from seat plate: 295 - 365 mm	£41
JBH2198		Padded Armrest with integrated cloth guard	(NC 2132, 2195, 2196, 2197, 2199) £110
JBH2193		Flip-up Armrest Up to 150 kg max user weight	(NC 1445, 2179) £65
JBH2192		Following Armrest Height adjustable from 230 - 300 mm	(NC 1445, 2179) £184
JBH2194		Following Armrest Height adjustable from 300 - 360 mm	(NC 1445, 2179) £184
ARMPAD (Pair)*			
JBH2195		Short Armpad (325 mm)	(NC 2198) ◆
JBH2196		Long Armpad (375 mm)	(NC 2132, 2198) ○
JBH2197		Long Ergonomic Armpad (397 mm)	(NC 2132, 2198, 2230, 2247) ○
JBH2199		Soft Armpad In combination with 10-way switch, the switch 2251 has to be mounted on the opposite side of the remote	(NC 2198) £95
HIP SUPPORTS LH			
JBH2184		Hip Support LH, pad size 140 x 100 mm	(NC 2198, 2701, 2702, 2703, 2704, 2705, 2706) £252
JBH2186		Hip Support LH, pad size 230 x 100 mm	(NC 2198, 2701, 2702, 2703, 2704, 2705, 2706) £252
JBH2188		Hip Support LH, pad size 330 x 100 mm	(NC 2198, 2701, 2702, 2703, 2704, 2705, 2706) £252
HIP SUPPORTS RH			
JBH2185		Hip Support RH, pad size 140 x 100 mm	(NC 2198, 2701, 2702, 2703, 2704, 2705, 2706) £252
JBH2187		Hip Support RH, pad size 230 x 100 mm	(NC 2198, 2701, 2702, 2703, 2704, 2705, 2706) £252
JBH2189		Hip Support RH, pad size 330 x 100 mm	(NC 2198, 2701, 2702, 2703, 2704, 2705, 2706) £252

MOUNTING HARDWARE FOR HIP SUPPORTS			
JBH2164		Mounting hardware for hip supports (6,1" / 155 mm high)	○
JBH2165		Mounting hardware for hip supports (8,1" / 206 mm high)	○
HEADREST / NECKREST - Adequate headrest option must be selected*			
JBH1350		Neckrest depth adjustable (NC 2230, 2247)	£81
JBH1432		Headrest sideways adjustable (NC 2230, 2247)	£226
JBH1427		Headrest with cheek support (NC 2230, 2247)	£226
JBH2324		Elan Standard Pad 255 mm / 10" wide w/standard cover (ESP26) Complete with Elan Headrest Hardware (NC 2230, 2247, 2170)	£361
JBH2325		Elan Occipital Pad 230x180 mm / 9"x7" w/standard cover (EOP23) Complete with Elan Headrest Hardware (NC 2230, 2247, 2170)	£361
JBH2327		Elan 4 Point Headrest Pad 280x255 mm / 11"x10" w/standard Complete with Elan Headrest Hardware (NC 2230, 2247, 2170)	£489
JBH1496		Non Invacare Headrest but profile mount fitting Flex3 hard back (NC 2170, 2171, 2174, 2175, 2176, 2177, 2178, 1288, 1289)	○
JBH1495		Non Invacare Headrest - no mounting hardware With Backbrace/Pushbar dia. 22mm in conjunction with telescopic seat frame/backposts	○
LATERAL SUPPORTS			
JBH1445		Lateral Support with extended arm, multi-adjustable pad (NC 2170, 2171, 2174, 2175, 2176, 2177, 2178, 2192, 2193, 2194)	£168
SEAT LIFTER / SEAT TILT*			
JBH1091		Manually settable Seat Tilt 0 - 15° (NC 1401)	○
JBH1092		Powered "fixed pivot" Tilt 20° (NC 2201)	£404
JBH1093		Powered Tilt (0°- 30°) with CoG shift (NC 2201)	£1,073
JBH1094		Seat lifter with powered tilt (300 mm, 30°) with CoG shift (NC 2201, 2220)	£1,893
SEAT TO GROUND HEIGHT*			
JBH2857		Seat-to-Floor Height 435 mm	○
JBH2858		Seat-to-Floor Height 460 mm	○
JBH2859		Seat-to-Floor Height 485 mm	○
BACKREST RECLINER*			
JBH0975		Manual Settable Recliner (in 3.75° steps) 90° - 120°	◆
JBH0977		Powered Recliner 30°	£505
CENTRE MOUNTED LEGREST*			
JBH0556		Centre mount legrest, angle settable Length: 280 - 385 mm (NC 0517, 1244, 1416)	£373
JBH0589		LNX Power centre mount legrest, length compensating Length: 340 - 410 mm (NC 2201, 1244, 1416)	£1,640
JBH2855		Without Legrests	○
CENTRE MOUNT FOOTPLATES / FOOT PLATFORM*			
JBH2365		Angle & width adjustable footplates S (127 W x 190 mm D / 5"x7.5")	£103
JBH0559		Angle & width adjustable footplates M (140 W x 240 mm D / 5.5"x9.5")	○
JBH0516		Angle & width adjustable footplates L (150 W x 290 mm D / 6"x11.5")	£135
JBH0517		Angle adjustable foot platform M (305 W x 275 mm D / 12"x10.75")	£270
JBX3001		Please select If Centre mount footplate is not required. (NC 0589)	○

FOOTREST HANGER TYPE*			
JBH0503		Hanger Vari F pair Lower leg length: 290 - 460 mm Angle adjustment: 70 - 0° with tools	◆
JBH0533		Hanger Vari A pair Lower leg length: 290 - 460 mm Angle adjustment: 70 - 0° with handwheel and calf pad	£233
JBH0531		Hanger ADE pair Lower leg length: 290 - 460 mm Powered Angle adjustment: 80 - 0° with calf pad	(NC 2201) £669
JBH0501		Without Hanger pair	(NC 0504, 0513) ○
FOOTPLATE TYPE*			
JBH0512		Fixed Footplate pair	(NC 0533, 0531) ◆
JBH0504		Angle Adjustable Footplate pair	£32
JBH0513		Angle and depth adjustable Footplate pair	£118
JBX3003		No Side Mounted Footplate required	(NC 0503, 0533, 0531) ○
CONTROLS*			
JBH2201		Remote LiNX DLX-REM110 = Drive only 1 Drive function with Drive only power module	(NC 2837, 1092, 1093, 1094, 0977, 0589, 2231, 2232, 2230, 2247, 2257, 2233, 2240, 2258, 2219, 2245, 2270, 1908) ◆
JBH2203		Remote LiNX DLX-REM211 = Drive & Actuators 3 Drive functions and power module incl. 2 Actuators	(NC 1437, 2231, 2232, 2233, 2230, 2247, 2257, 2240, 2270, 2258, 1908, 1909) ○
JBH2204		Remote LiNX DLX-REM216 = Drive, Actuators & Lights 3 Drive functions and power module incl. 2 Actuators & Lights	(NC 2232, 2233, 2230, 2247, 2257, 2240, 2270, 2258) £151
JBH2205		Rehab Remote LiNX, touch-screen user interface DLX-REM400 Configurable jack sockets	(NC 2220) £341
JBH2206		Rehab Remote LiNX without joystick for alternative access users DLX-REM500 Requires a user input/alternative joystick: 2231, 2232, 2229, 2230, 2247, 2233	(NC 2220, 2270, 1709, 1710) £776
REM400 TOGGLE OPTION			
JBH2270		Assignable Toggle Switch Module f. REM400 (DLX- SW400TGL)	(NC 2203, 2204, 2206) £316
CONTROL BRACKETRY*			
JBH1207		Standard Bracket RH	(NC 2285, 2283, 2287, 2207, 2265, 2275, 1404) ○
JBH1208		Standard Bracket LH	(NC 2286, 2284, 2288, 2208, 2266, 2276, 1405) ○
JBH1707		Swing Away Bracket RH	(NC 2285, 2283, 2287, 2207, 2265, 2275, 1404) £133
JBH1708		Swing Away Bracket LH	(NC 2286, 2284, 2288, 2208, 2266, 2276, 1405) £133
JBH1709		Matrx Resolve Swing in/out bracket RH	(NC 2283, 2287, 2207, 2265, 2275, 1404) £219
JBH1710		Matrx Resolve Swing in/out bracket LH	(NC 2284, 2288, 2208, 2266, 2276, 1405) £219
JBH2213		Centre Mount, swing-away to RH, Nucleus	(NC 2283, 2287, 2207, 2265, 2275, 1415, 2217, 1404, 1405) £498
JBH2214		Centre Mount, swing-away to LH, Nucleus	(NC 2284, 2288, 2208, 2266, 2276, 1339, 2218, 1404, 1405) £498
JBH2227		Bracket RH, mounted for Attendant	(NC 1287, 2220, 2128, 2108, 2219, 1221, 1312) £64
JBH2228		Bracket LH, mounted for Attendant	(NC 1286, 2220, 2128, 2108, 2219, 1222, 1363) £64
JBX2239		Individual mounting option per request (Special) Requires additional information/specification	POA
DUAL CONTROL			
JBH2220		Dual Control LiNX (DLX-REM050 as secondary remote)	(NC 2273, 2108, 2205, 2206, 1221, 1222, 2227, 2228) £395
JBH2219		Compact Dual Control LiNX (DLX-ACU200)	(NC 2108, 1221, 1222, 2227, 2228) £395
INVACARE LINX IDC DUAL CONTROL			
JBH2245		Invacare LiNX IDC Dual control	(NC 2108, 2259, 2260) £750
DUAL CONTROL BRACKETRY			
JBH2259		Bracket RH, Attendant	(NC 2291, 2245, 1312) ○
JBH2260		Bracket LH, Attendant	(NC 2292, 2245, 1363) ○
INVACARE LINX IDC DUAL CONTROL BRACKETRY			
JBH1221		Bracket RH, detachable	(NC 2291, 2220, 2219, 2227, 1312) ○
JBH1222		Bracket LH, detachable	(NC 2292, 2220, 2219, 2228, 1363) ○

LINX COMPACT REMOTE				
JBH2231		Proportional Compact Remote DLX-CR400	(NC 2203, 2204)	£719
JBH2232		Proportional Compact Remote Low Force DLX-CR-400LF	(NC 2203, 2204)	£960
LINX COMPACT REMOTE BRACKET				
JBH2207		RH Mount	(NC 1707, 2213, 2217, 1404)	○
JBH2208		LH Mount	(NC 1708, 2214, 2218, 1405)	○
JBH2265		Centre Mount, swing-away to RH (Nucleus) If several joysticks (or Joystick + REM500) should be mounted on one Nucleus, please select special mounting option 2215!	(NC 1707, 2213, 2217, 1404, 1405, 1415)	£498
JBH2266		Centre Mount, swing-away to LH (Nucleus) If several joysticks (or Joystick + REM500) should be mounted on one Nucleus, please select special mounting option 2215!	(NC 1708, 2214, 2218, 1404, 1405, 1339)	£498
JBX2215		Individual mounting option per request (Special) Requires additional information/specification		POA
CHIN CONTROL				
UNI9250		Mo-vis Micro joystick		POA
UNI9251		Mo-vis Multi joystick		POA
CHIN CONTROL BRACKETRY MANUAL				
JBH1289		Swing Away Bracket RH, side manual (ASL628) Requires Elan Headrest	(NC 1350, 1432, 1427, 1496)	○
JBH1288		Swing Away Bracket LH, side manual (ASL628) Requires Elan Headrest	(NC 1350, 1432, 1427, 1496)	○
CHIN CONTROL BRACKETRY POWERED				
JBH1287		Swing Away Bracket RH, powered	(NC 2259, 1221, 2291, 1312, 1415)	£2,183
JBH1286		Swing Away Bracket LH, powered	(NC 2260, 1222, 2292, 1363, 1339)	£2,183
LINX HEAD CONTROL				
JBH2230		Head Control with ultrasonic proximity sensors (ATOM, ASL104) with standard wings M3	(NC 2203, 2204, 2233, 2324, 2325, 2327, 2257, 2240, 2175, 2171, 2170, 2197, 1350, 1432, 1427)	£3,719
JBH2247		Head Control with ultrasonic proximity sensors (PROTON, ASL104P) with PROTON wings	(NC 2203, 2204, 2233, 2324, 2325, 2327, 2257, 2240, 2175, 2171, 2170, 2197, 1350, 1432, 1427)	£4,356
LINX FOUR -SWITCH PROXIMITY ARRAY				
JBH2257		Four-Switch Proximity Array (ASL106)	(NC 2203, 2204, 2230, 2247, 2233, 2240, 2128)	£5,760
JBH2275		Centre Mount, swing-away to RH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm	(NC 1707, 2213, 2217, 1404, 1405, 1415)	£498
JBH2276		Centre Mount, swing-away to LH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm	(NC 1708, 2214, 2218, 1404, 1405, 1339)	£498
JBX2269		Individual mounting option per request (Special) Requires additional information/specification		POA
LINX WIRELESS ACCESSORIES FOR HEAD CONTROL AND PROXIMITY CONTROL				
JBH2249		Wireless Triple Switch Receiver (ASL557-3)		£341
JBH2250		Wireless Mouse Emulator (ASL558)		£817
LINX MICRO EXTREMITY CONTROL BRACKET				
JBH2233		MEC Proportional Joystick (Micro Extremity Control ASL130)	(NC 2203, 2204, 2230, 2247, 2257, 2128)	£3,612
LINX HEAVY DUTY REMOTE				
JBH2240		Proportional Heavy Duty Remote Joystick w mode switch (ASL133)	(NC 2203, 2204, 2230, 2247, 2257, 2128)	£2,013

LINX MICRO EXTREMITY AND HEAVY DUTY CONTROL BRACKET				
JBH2285		Tray Mount RH (100x150 mm)	(NC 1707, 2213, 2217, 1404, 1405, 1415)	○
JBH2286		Tray Mount LH (100x150 mm)	(NC 1708, 2214, 2218, 1404, 1405, 1339)	○
JBH2283		Centre Mount, swing-away to RH, Nucleus	(NC 1707, 2213, 2217, 1404, 1405, 1415)	£498
JBH2284		Centre Mount, swing-away to LH, Nucleus	(NC 1708, 2214, 2218, 1404, 1405, 1339)	£498
JBH2287		Centre Mount, swing-away to RH, Nucleus Incl. Tray ASL603, 100x150 mm	(NC 1707, 2213, 2217, 1404, 1405, 1415)	£498
JBH2288		Centre Mount, swing-away to LH, Nucleus Incl. Tray ASL603, 100x150 mm	(NC 1708, 2214, 2218, 1404, 1405, 1339)	£498
JBX2289		Individual mounting option per request (Special) Requires additional information/specification		POA
LINX 10-WAY SWITCH				
JBH2251		10-Way Switch DLX-SW10-A for direct Powered Seating Access		£505
LINX 10-WAY SWITCH BRACKET				
JBH2217		Bracket RH, mounted for user	(NC 2207, 2287, 2283, 2275, 2213, 2265, 1404)	○
JBH2218		Bracket LH, mounted for user	(NC 2208, 2288, 2284, 2286, 2276, 2214, 2266, 1405)	○
JBH2291		Bracket RH, mounted for attendant	(NC 2259, 1221, 1287, 2108)	£64
JBH2292		Bracket LH, mounted for attendant	(NC 2260, 1222, 1288, 2108)	£64
LINX USB CHARGER PORT				
JBH2253		USB Charger Port DLX-USB02-A	(NC 2702, 2703, 2704, 2705, 2706)	£165
LINX USB CHARGER PORT BRACKET				
JBH2281		RH Mount		○
JBH2282		LH Mount		○
REMOTE STOP SWITCH				
JBH2258		Remote Stop Switch (ASL504)	(NC 2203, 2204)	£1,050
LINX IN200 MODULE				
JBH2273		IN200 Module (Jack port for external power button)		£241
LINX INPUT MODULE				
UNI9252		IN500 Module (for DB9 Drive input)	(NC 2201, 2203, 2203)	£785
FURTHER LINX COMPONENTS / ACCESSORIES				
JBH2300		To be enclosed only		○
JBX2301		To be connected and programmed (Special)		POA
BUDDY BUTTONS				
JBH1685		Red Buddy Button		£88
JBH1686		Yellow Buddy Button		£88
JBH1687		Green Buddy Button		£88
JBH1690		Blue Buddy Button		£88
EGG SWITCH				
JBH2254		Egg Switch, black		£88
PIKO BUTTONS				
JBH1711		Red Piko Button		£158
JBH1712		Yellow Piko Button		£158
JBH1713		Green Piko Button		£158
JBH1714		Orange Piko Button		£158
JBH1715		Blue Piko Button		£158
JBH1716		Black Piko Button		£158

BODYPOINT JOYSTICK KNOBS			
JBH2721		3 inch U-Shaped, Joystick Handle, Bodypoint PC101A	£63
JBH2722		4 inch U-Shaped, Joystick Handle, Bodypoint PC102A	£63
JBH2723		Domed Rubber, Joystick Handle, Bodypoint PC107A	£63
JBH2724		Mushroom Grooved, Joystick Handle, Bodypoint PC110A	£63
5:1 SWITCH ADAPTOR			
JBH1681		5:1 Switch Adapter (9-pin sub D to 5 jack ports)	£121
SHROUD COLOUR*			
JBH2871		Mid range - Black/Anthracite Shroud	○
JBH2872		Mid range - Silver Shroud	○
JBH2873		Mid range - Bengal Red Shroud	○
JBH2874		Mid range - Mercury Blue Shroud	○
JBH2875		Mid range - Midnight Purple Shroud	○
JBX0299		Premium - Individual Shroud Colour (includes coloured spreader bars on Modulite back) RAL number required Extended lead lime: + 15 days	£214
JBH2865		Premium - Racing Green Shroud	£104
JBH2866		Premium - Deep Ocean Blue Satin Shroud	£104
JBH2867		Premium - Frosty Blue Shroud	£104
JBH2868		Premium - Ultimate Red Shroud	£104
JBH2869		Premium - Falcon White Shroud	£104
RIM COLOUR*			
JBH2882		Mid range - Black Rim Inserts	○
JBH2883		Mid range - Silver Rim Inserts	○
JBH2884		Mid range - Bengal Red Rim Inserts	○
JBH2885		Mid range - Mercury Blue Rim Inserts	○
JBH2886		Mid range - Midnight Purple Rim Inserts	○
JBX2887		Premium - Individual Rim Insert Colour RAL number required Extended lead lime: + 15 days	○
JBH2888		Premium - Racing Green Rim Inserts	○
JBH2889		Premium - Deep Ocean Blue Satin Rim Inserts	○
JBH2890		Premium - Frosty Blue Rim Inserts	○
JBH2891		Premium - Ultimate Red Rim Inserts	○
JBH2892		Premium - Falcon White Rim Inserts	○
FORKS*			
JBH2631		Single sided fork	£226
JBH2632		Double sided fork Without castor rim insert (enlarges wide base to 656 mm)	◆
FRONT WHEEL*			
JBH0835		9" (225 mm) x 2.80 (70 mm) Castor, pneumatic	○
JBH0836		9" (225 mm) x 2.80 (70 mm) Castor, puncture proof solid tyre	◆
CASTOR LOCKS			
JBH1490		Castor Lock	£197

DRIVE WHEEL*			
JBH0704		DW 14" (355 mm) x 3" (75 mm), Pneumatic	○
JBH0705		DW 14" (355 mm) x 3" (75 mm), Solid inlay Dual density foam insert	◆
C.T.C SUSPENSION*			
JBH2870		CTC Suspension Control - Traction - Comfort	◆
SUSPENSION SETTING*			
JBH1916		Suspension up to 160 kg User Weight For RX 40 HD only	◆
MOTOR TYPE*			
JBH1294		Indoor setting, reduced speed approx. 3 km/h	○
JBH1282		6 km/h	◆
JBH1908		10 km/h Package Always combined with LED Light system with Indicators, 30 Lux headlights	(NC 2203, 1437, 1254, 2837) £543
JBH1909		12 km/h Ultra High Speed Package Always combined with LED Light system with Indicators, 30 Lux headlights	(NC 2201, 2203, 1437, 1254, 2837) £653
BATTERY TYPE*			
JBH1259		60 Ah Gel	◆
JBH1260		73,5 Ah Gel	(NC 1450) £92
G TRAC			
JBH1215		G-Trac Module	£324
CHARGER*			
JBH1343		Opti Charge Charger 8A (IP44)	◆
JBH1313		Opti Charge Charger 10A (IP44)	£151
ANTI-TIPPERS*			
JBH1205		Anti-tipper	◆
LIGHTS			
JBH1437		LED Light system with Indicators 30 Lux headlights	(NC 2203, 1908, 1909) £479
JBH2837		LED Light system extra bright with indicators 60 Lux headlights	(NC 2203) £143
JBH1254		Without Light system	(NC 1909) ○
KERB RIDER			
JBH1244		Kerb Rider With park position feature	(NC 0556, 0589, 1294) £212
TRANSPORTATION KIT*			
JBH1450		Preparation for DAHL-Docking station (Adequate headrest is compulsory) Includes a mounting bracket in the battery box Docking Station for car implementation has to be ordered separately with DAHL Engineering	£586
JBH1494		4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19) Adequate headrest is compulsory, Not compatible with Seat lifter	◆
JBH1401		6-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19)) required for variants with combined Lift & Tilt module Adequate headrest is compulsory	(NC 1091, 1092, 1093) ○
EASY ACCESS FRONT TIE DOWNS			
JBH2894		Easily accessible front tie-down points, additionally mounted for 4- or 6- point-tie-down For transport worthy versions (1494/1401) only	£49
LAP BELT*			
JBH1329		Lap belt with Auto Lock Standard mount on backrest	○
JBH2113		Lap belt with Auto Lock Flexible mount in side rails	○

POSITIONING BELT				
JBH2701		BP Evoflex Pelvic Stabilizer, Rehab Latch Size L Padded length 620 mm	(NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253)	£77
JBH2702		BP Evoflex Pelvic Stabilizer, Rehab Latch Size M Padded length 460 mm	(NC 2123, 2124, 2184, 2185, 2186, 2187, 2188, 2189, 2253)	£77
JBH2703		BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size L Padded length 620 mm	(NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253)	£79
JBH2704		BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size M Padded length 460 mm	(NC 2123, 2124, 2184, 2185, 2186, 2187, 2188, 2189, 2253)	£79
JBH2705		BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size L Padded length 620 mm	(NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253)	£79
JBH2706		BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size M Padded length 460 mm	(NC 2123, 2124, 2184, 2185, 2186, 2187, 2188, 2189, 2253)	£79
REAR VIEW MIRROR				
JBH1415		Rear view mirror RH	(NC 2217, 2285, 2213, 2265, 2275, 2283, 2287)	£42
JBH1339		Rear view mirror LH	(NC 2218, 2286, 2214, 2266, 2276, 2284, 2288)	£42
OTHER OPTIONS				
JBH2128		Bag Hooks on Modulite Seat Max load up to 6 kg = impact on max user weight	(NC 2172, 2173, 2227, 2228, 2233, 2240, 2257, 2291, 2292)	£77
JBH1312		Cane holder RH	(NC 2227, 2259, 1221)	£42
JBH1363		Cane holder LH	(NC 2228, 2260, 1222)	£42
JBH1416		Fender for Castors		£88
JBH1417		Fender for Drive wheels		£169
JBH1404		Tray, Swing Away RH	(NC 1707, 2275, 2276, 2287, 2288, 2283, 2284, 2285, 2286, 2265, 2266, 2207, 2213, 2214, 2217)	£167
JBH1405		Tray, Swing Away LH	(NC 1708, 2275, 2276, 2287, 2288, 2283, 2284, 2285, 2286, 2265, 2266, 2208, 2213, 2214, 2218)	£167

TOTAL PRICE	
--------------------	--

Non Conforming Options:

JBH2121	380 mm Seat Plate Can be reduced to 330 mm in combination with inrigged armrests
JBH2122	430 mm Seat Plate Can be reduced to 380 mm in combination with inrigged armrests
JBH2123	480 mm Seat Plate
JBH2124	530 mm Seat Plate
JBH2106	Body fit Kit Outrigging Backrest Brackets and Pivot Point (outrig: 1 inch/2,54 cm per side)
JBH2132	Seat Depth 410 mm
JBH2163	Backposts short (480 mm) For Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery Width adjustable cast bracket
JBH2108	Backposts without push handles
JBH2171	Tension Adjustable Upholstery nylon overlay short = 480 mm / 7 straps measured from seat plate to top of upholstery 480 mm
JBH2179	Flex3 Tension Adjustable Upholstery height 540 mm (seat plate to top of upholstery)
JBH2172	Flex3 adjustable Backrest Plates (without Backrest Cushion) Flex3 Backrest plate width and height adjustable 560 mm - 790 mm
JBH2173	Flex3 adjustable Backrest Plates with Laguna Backrest cushion Flex3 Backrest plate width and height adjustable 560 mm - 790 mm
JBH2174	Matrx Posture Back, height 400 mm (Mesh fabric) Corresponds to SW setting chosen (up to 136 kg max user weight)
JBH2175	Matrx Posture Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen (up to 136 kg max user weight)
JBH2176	Matrx Elite E2 Back , height 500 mm (Mesh fabric) Corresponds to SW setting chosen (up to 136 kg max user weight)
JBH2177	Matrx Elite TR Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen
JBH2178	Matrx Elite E2 Deep Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen (up to 136 kg max user weight)
JBH2170	Without Backrest Upholstery (Back Posts/Back Frame only)
JBH2151	Black Fabric Upholstery (Nylon with option 2141; 2-way-stretch material with options 2142, 2143)
JBH2155	Black Spacer Fabric Upholstery
JBH2198	Padded Armrest with integrated cloth guard
JBH2193	Flip-up Armrest Up to 150 kg max user weight
JBH2192	Following Armrest Height adjustable from 230 - 300 mm
JBH2194	Following Armrest Height adjustable from 300 - 360 mm
JBH2195	Short Armpad (325 mm)
JBH2196	Long Armpad (375 mm)
JBH2197	Long Ergonomic Armpad (397 mm)
JBH2199	Soft Armpad In combination with 10-way switch, the switch 2251 has to be mounted on the opposite side of the remote
JBH2184	Hip Support LH, pad size 140 x 100 mm
JBH2186	Hip Support LH, pad size 230 x 100 mm
JBH2188	Hip Support LH, pad size 330 x 100 mm
JBH2185	Hip Support RH, pad size 140 x 100 mm
JBH2187	Hip Support RH, pad size 230 x 100 mm
JBH2189	Hip Support RH, pad size 330 x 100 mm

JBH1350	Neckrest depth adjustable
JBH1432	Headrest sideways adjustable
JBH1427	Headrest with cheek support
JBH2324	Elan Standard Pad 255 mm / 10" wide w/standard cover (ESP26) Complete with Elan Headrest Hardware
JBH2325	Elan Occipital Pad 230x180 mm / 9"x7" w/standard cover (EOP23) Complete with Elan Headrest Hardware
JBH2327	Elan 4 Point Headrest Pad 280x255 mm / 11"x10" w/standard Complete with Elan Headrest Hardware
JBH1496	Non Invacare Headrest but profile mount fitting Flex3 hard back
JBH1445	Lateral Support with extended arm, multi-adjustable pad
JBH1091	Manually settable Seat Tilt 0 - 15°
JBH1092	Powered "fixed pivot" Tilt 20°
JBH1093	Powered Tilt (0° - 30°) with CoG shift
JBH1094	Seat lifter with powered tilt (300 mm, 30°) with CoG shift
JBH0977	Powered Recliner 30°
JBH0556	Centre mount legrest, angle settable Length: 280 - 385 mm
JBH0589	LNX Power centre mount legrest, length compensating Length: 340 - 410 mm
JBH0517	Angle adjustable foot platform M (305 W x 275 mm D / 12"x10.75")
JBH0503	Hanger Vari F pair Lower leg length: 290 - 460 mm Angle adjustment: 70 - 0° with tools
JBH0533	Hanger Vari A pair Lower leg length: 290 - 460 mm Angle adjustment: 70 - 0° with handwheel and calf pad
JBH0531	Hanger ADE pair Lower leg length: 290 - 460 mm Powered Angle adjustment: 80 - 0° with calf pad
JBH0504	Angle Adjustable Footplate pair
JBH0513	Angle and depth adjustable Footplate pair
JBH2201	Remote LiNX DLX-REM110 = Drive only 1 Drive function with Drive only power module
JBH2203	Remote LiNX DLX-REM211 = Drive & Actuators 3 Drive functions and power module incl. 2 Actuators
JBH2204	Remote LiNX DLX-REM216 = Drive, Actuators & Lights 3 Drive functions and power module incl. 2 Actuators & Lights
JBH2205	Rehab Remote LiNX, touch-screen user interface DLX-REM400 Configurable jack sockets
JBH2206	Rehab Remote LiNX without joystick for alternative access users DLX-REM500 Requires a user input/alternative joystick: 2231, 2232, 2229, 2230, 2247, 2233
JBH2270	Assignable Toggle Switch Module f. REM400 (DLX-SW400TGL)
JBH1707	Swing Away Bracket RH
JBH1708	Swing Away Bracket LH
JBH1709	Matrx Resolve Swing in/out bracket RH
JBH1710	Matrx Resolve Swing in/out bracket LH
JBH2213	Centre Mount, swing-away to RH, Nucleus
JBH2214	Centre Mount, swing-away to LH, Nucleus
JBH2227	Bracket RH, mounted for Attendant
JBH2228	Bracket LH, mounted for Attendant
JBH2220	Dual Control LiNX (DLX-REM050 as secondary remote)
JBH2219	Compact Dual Control LiNX (DLX-ACU200)
JBH2245	Invacare LiNX IDC Dual control
JBH2259	Bracket RH, Attendant
JBH2260	Bracket LH, Attendant

JBH1221	Bracket RH, detachable
JBH1222	Bracket LH, detachable
JBH2231	Proportional Compact Remote DLX-CR400
JBH2232	Proportional Compact Remote Low Force DLX-CR-400LF
JBH2207	RH Mount
JBH2208	LH Mount
JBH2265	Centre Mount, swing-away to RH (Nucleus) If several joysticks (or Joystick + REM500) should be mounted on one Nucleus, please select special mounting option 2215!
JBH2266	Centre Mount, swing-away to LH (Nucleus) If several joysticks (or Joystick + REM500) should be mounted on one Nucleus, please select special mounting option 2215!
JBH1289	Swing Away Bracket RH, side manual (ASL628) Requires Elan Headrest
JBH1288	Swing Away Bracket LH, side manual (ASL628) Requires Elan Headrest
JBH1287	Swing Away Bracket RH, powered
JBH1286	Swing Away Bracket LH, powered
JBH2230	Head Control with ultrasonic proximity sensors (ATOM, ASL104) with standard wings M3
JBH2247	Head Control with ultrasonic proximity sensors (PROTON, ASL104P) with PROTON wings
JBH2257	Four-Switch Proximity Array (ASL106)
JBH2275	Centre Mount, swing-away to RH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)
JBH2276	Centre Mount, swing-away to LH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)
JBH2233	MEC Proportional Joystick (Micro Extremity Control ASL130)
JBH2240	Proportional Heavy Duty Remote Joystick w mode switch (ASL133)
JBH2285	Tray Mount RH (100x150 mm)
JBH2286	Tray Mount LH (100x150 mm)
JBH2283	Centre Mount, swing-away to RH, Nucleus
JBH2284	Centre Mount, swing-away to LH, Nucleus
JBH2287	Centre Mount, swing-away to RH, Nucleus Incl. Tray ASL603, 100x150 mm
JBH2288	Centre Mount, swing-away to LH, Nucleus Incl. Tray ASL603, 100x150 mm
JBH2217	Bracket RH, mounted for user
JBH2218	Bracket LH, mounted for user
JBH2291	Bracket RH, mounted for attendant
JBH2292	Bracket LH, mounted for attendant
JBH2253	USB Charger Port DLX-USB02-A
JBH2258	Remote Stop Switch (ASL504)
JBH2273	IN200 Module (Jack port for external power button)
JBH1294	Indoor setting, reduced speed approx. 3 km/h
JBH1908	10 km/h Package Always combined with LED Light system with Indicators, 30 Lux headlights
JBH1909	12 km/h Ultra High Speed Package Always combined with LED Light system with Indicators, 30 Lux headlights
JBH1437	LED Light system with Indicators 30 Lux headlights
JBH2837	LED Light system extra bright with indicators 60 Lux headlights
JBH1254	Without Light system
JBH1244	Kerb Rider With park position feature
JBH1450	Preparation for DAHL-Docking station (Adequate headrest is compulsory) Includes a mounting bracket in the battery box Docking Station for car implementation has to be ordered separately with DAHL Engineering

JBH1401	6-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19)) required for variants with combined Lift & Tilt module Adequate headrest is compulsory
JBH2701	BP Evoflex Pelvic Stabilizer, Rehab Latch Size L Padded length 620 mm
JBH2702	BP Evoflex Pelvic Stabilizer, Rehab Latch Size M Padded length 460 mm
JBH2703	BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size L Padded length 620 mm
JBH2704	BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size M Padded length 460 mm
JBH2705	BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size L Padded length 620 mm
JBH2706	BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size M Padded length 460 mm
JBH1415	Rear view mirror RH
JBH1339	Rear view mirror LH
JBH2128	Bag Hooks on Modulite Seat Max load up to 6 kg = impact on max user weight
JBH1312	Cane holder RH
JBH1363	Cane holder LH
JBH1416	Fender for Castors
JBH1404	Tray, Swing Away RH
JBH1405	Tray, Swing Away LH



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):