



Yes, you can:

# Power Retail Prescription Form

Tel: 01656 776222 Fax: 01656 776220  
Email: ordersuk@invacare.com  
Online Spares available at: www.invacare.co.uk

SEPTEMBER 2024  
Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## AVIVA RX40 Modulite

### Standard Features

- Max user weight 136kg (21 Stone)
- Crash tested ISO 7176 -Part 19
- Modulite highly adjustable seating system
- Telescopic adjustable seat widths from 330mm - 530mm (13" - 21") with an additional 50mm, (2") by arms
- Seat depths 410mm - 510mm stepless adjustment (16" - 20")
- Tension adjustable backrest
- Standard T-Armrests, detachable with short armpad
- Manual settable seat Tilt (0° - 15°)
- Manual settable backrest Recline (0° - 30° in 3.75° steps)
- 80° hangers with fixed footplate
- LiNX control REM110 with standard bracket
- Choice of 5 Shroud colours, 5 Rim insert colours
- 9" (225 mm) x 2.80 (70 mm) Castor, pneumatic
- 14" (355 mm) x 3" (75 mm), Pneumatic drive wheels
- CTC Suspension Control - Traction - Comfort
- 6 km/h (4mph) 4 pole and 52Ah batteries
- 8A (200 - 250V) charger
- Anti tippers
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)

**Standard Base price: £6,234**



Please make your selection with a cross in the box     Free of charge     standard feature    POA special request

Customer weight =                      KG

Code	Description	
<b>BASE MODEL*</b>		
JBG001	AVIVA RX 40 Base Model Max user weight up to 136 kg includes Performance Package: 120A Power Module + minimum 60 Ah Batteries	£6,234
JBX1300	Special Extended lead time: + 10 days	POA
<b>CHASSIS WIDTH*</b>		
JBL2810	Chassis Narrow (560 mm/22")	◆
<b>SEATING TYPE*</b>		
JBG2100	Invacare® Modulite Seat Unit	◆

<b>SEAT FRAME*</b>			
JBG2111		Telescopic Seat Frame and Plate with Side Rails width 380 to 530 mm Additional 50 mm by arms	○
JBG2112		Telescopic Seat Frame + Sling Seat with Side Rails width 380 to 530 mm; + 30 mm seat height Additional 50 mm by arms	○
<b>SEAT WIDTH*</b>			
JBG2121		380 mm Seat Plate Can be reduced to 330 mm in combination with inrigged armrests (NC 2174, 2175, 2176)	○
JBG2122		430 mm Seat Plate Can be reduced to 380 mm in combination with inrigged armrests	○
JBG2123		480 mm Seat Plate (NC 2704, 2706)	○
JBG2124		530 mm Seat Plate (NC 2175, 2176, 2177, 2178, 2143, 2702, 2704, 2706, 2106)	○
<b>Body Fit Kit</b>			
JBG2106		Body fit Kit Outrigging Backrest Brackets and Pivot Point (outrig: 1 inch/2,54 cm per side) (NC 2192, 2194)	○
<b>SEAT DEPTH*</b>			
JBG2132		Seat Depth 410 mm (NC 2197, 2198)	○
JBG2134		Seat Depth 460 mm	○
JBG2136		Seat Depth 510 mm	○
<b>SEAT CUSHION*</b>			
JBG2142		Matrx Flo-tech Cushion Flo-tech Lite	£118
JBG2143		Matrx Flo-tech Contour Visco; (2-way stretch material)	£165
JBG2147		Matrx Libra Cushion, 90 mm thick Extended lead time: + 4 days	£389
JBG2144		Matrx Vi Cushion (Mesh fabric)	£360
JBG2140		Without seat cushion	○
<b>BACKREST SUPPORT*</b>			
JBG2163		Backposts short (480 mm) For Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery Width adjustable cast bracket (NC 2173)	○
JBG2162		Backposts tall (540 mm) For Flex3 backrests/Tension Adjustable Upholstery/Matrx Backs, measured from seat plate to top of upholstery Width adjustable cast bracket	○
<b>PUSH HANDLES FOR BACKPOSTS</b>			
JBG2108		Backposts without push handles (NC 2219, 2220, 2227, 2228, 2245, 2291, 2292)	○

BACKREST UPHOLSTERY / CUSHION*			
JBG2171		Tension Adjustable Upholstery nylon overlay short = 480 mm / 7 straps measured from seat plate to top of upholstery 480 mm	(NC 2230, 2247) ◆
JBG2179		Flex3 Tension Adjustable Upholstery height 540 mm (seat plate to top of upholstery)	(NC 2192, 2193, 2194, 2199, 2155, 2230, 2247) £153
JBG2172		Flex3 adjustable Backrest Plates without Backrest Cushion	(NC 2128) £250
JBG2173		Flex3 adjustable Backrest Plates with Laguna Backrest cushion	(NC 2128) £406
JBG2174		Matrx Posture Back, height 400 mm (Mesh fabric) Corresponds to SW setting chosen (up to 136 kg max user weight)	£558
JBG2175		Matrx Posture Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen (up to 136 kg max user weight)	£558
JBG2176		Matrx Elite Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen (up to 136 kg max user weight)	£558
JBG2177		Matrx Elite TR Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen	£558
JBG2178		Matrx Elite Deep Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen (up to 136 kg max user weight)	£648
JBG2170		Without Backrest Upholstery (Back Posts/Back Frame only)	(NC 2324, 2325, 2327, 2230, 2247) ○
UPHOLSTERY TYPE*			
JBG2151		Black Fabric Upholstery (Nylon w 2141; 2-way-stretch material w 2142, 2143)	◆
JBG2152		Black multi stretch PU coated fabric Upholstery	○
JBX3002		No Upholstry Colour required - No backrest upholstery	(NC 2179, 2173, 2174, 2175, 2176, 2177, 2178) ○
JBG2155		Black Spacer Fabric Upholstery	(NC 2171) ○
ARMREST (Pair)*			
JBG2191		T-Arm detachable, Standard Adjustment range from seat plate: 245 - 315 mm	◆
JBG2181		T-Arm detachable, Extended height Adjustment range from seat plate: 295 - 365 mm	£41
JBG2198		Padded Armrest with padded side guard (WATER FALL TYPE)	(NC 2195, 2196, 2197, 2199) £109
JBG2183		Inrigged Armrests Reduces SW by 25 mm each side	○
JBG2193		Flip-up Armrest Up to 150 kg max user weight	(NC 1445, 2179) £65
JBG2192		Following Armrest Height adjustable from 230 - 300 mm	(NC 2179) £185
JBG2194		Following Armrest Height adjustable from 300 - 360 mm	(NC 2179) £185
ARMPAD (Pair)*			
JBG2195		Short Armpad (325 mm)	◆
JBG2196		Long Armpad (375 mm)	(NC 2198) ○
JBG2197		Long Ergonomic Armpad (397 mm)	(NC 2198, 2230, 2247) ○
JBG2199		Soft Armpad In combination with 10-way switch, the switch 2251 has to be mounted on the opposite side of the remote	£96
JBX3000		Please selct with Option 2198 - Padded Armrest with padded side guard (WATER FALL TYPE)	(NC 2181, 2183, 2193, 2192) ○
HIP SUPPORTS LH			
JBG2184		Hip Support LH, pad size 140 x 100 mm	(NC 2183, 2701, 2702, 2703, 2704, 2705, 2706) £252
JBG2186		Hip Support LH, pad size 230 x 100 mm	(NC 2183, 2701, 2702, 2703, 2704, 2705, 2706) £252
JBG2188		Hip Support LH, pad size 330 x 100 mm	(NC 2183, 2701, 2702, 2703, 2704, 2705, 2706) £252

<b>HIP SUPPORTS RH</b>				
JBG2185		Hip Support RH, pad size 140 x 100 mm	(NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	£252
JBG2187		Hip Support RH, pad size 230 x 100 mm	(NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	£252
JBG2189		Hip Support RH, pad size 330 x 100 mm	(NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	£252
<b>MOUNTING HARDWARE FOR HIP SUPPORTS</b>				
JBG2164		Mounting hardware for hip supports (6,1" / 155 mm high)		○
JBG2165		Mounting hardware for hip supports (8,1" / 206 mm high)		○
<b>HEADREST / NECKREST - Adequate headrest option must be selected*</b>				
JBG1350		Neckrest depth adjustable	(NC 2247)	£81
JBG1432		Headrest sideways adjustable	(NC 2247)	£226
JBG1427		Headrest with cheek support	(NC 2247)	£226
JBG2324		Elan Standard Pad 255 mm / 10" wide w/standard cover (ESP26) Complete with Elan Headrest Hardware	(NC 2247, 2170)	£361
JBG2325		Elan Occipital Pad 230x180 mm / 9"x7" w/standard cover (EOP23) Complete with Elan Headrest Hardware	(NC 2247, 2170)	£361
JBG2327		Elan 4 Point Headrest Pad 280x255 mm / 11"x10" w/standard Complete with Elan Headrest Hardware	(NC 2247, 2170)	£489
JBG1496		Non Invacare Headrest but profile mount fitting Flex3 hard back	(NC 2171, 2174, 2175, 2176, 2177, 2178, 1288, 1289)	○
JBG1495		Non Invacare Headrest - no mounting hardware With Backbrace/Pushbar dia. 22mm in conjunction with telescopic seat frame/backposts		○
<b>LATERAL SUPPORTS</b>				
JBG1445		Lateral Support with extended arm, multi-adjustable pad	(NC 2171, 2174, 2175, 2176, 2177, 2178, 2192, 2193, 2194)	£168
<b>SEAT LIFTER / SEAT TILT*</b>				
JBG1091		Manually settable Seat Tilt 0 - 15°		○
JBG1092		Powered "fixed pivot" Tilt 20°	(NC 2201)	£404
JBG1093		Powered Tilt (0° - 30°) with CoG shift	(NC 2201)	£1,073
JBG2052		Electric Tilt (0°-50°) with COG shift	(NC 2201, 2171)	£1,322
JBG1094		Seat lifter with powered tilt (300 mm, 30°) with CoG shift	(NC 2201, 2220)	£1,893
JBG2053		Seat lifter with Electrical Tilt (300mm, 45°) with CoG shift	(NC 2201, 2220, 1910, 1912, 2171)	£2,359
<b>SEAT TO GROUND HEIGHT*</b>				
JBG2857		Seat-to-Floor Height 435 mm		○
JBG2858		Seat-to-Floor Height 460 mm		○
JBG2859		Seat-to-Floor Height 485 mm		○
<b>BACKREST RECLINER*</b>				
JBG0975		Manual Settable Recliner (in 3.75° steps) 90° - 120°		◆
JBG0977		Powered Recliner 30°		£505
<b>CENTRE MOUNTED LEGREST*</b>				
JBG0556		Centre mount legrest, angle settable Length: 280 - 385 mm	(NC 0517, 1244, 1416)	£373
JBG0589		LNx Power centre mount legrest, length compensating Length: 340 - 410 mm	(NC 2201, 1244, 1416)	£1,640
JBG2855		Without Legrests		○

<b>CENTRE MOUNT FOOTPLATES / FOOT PLATFORM*</b>			
JBG2365		Angle & width adjustable footplates S (127 W x 190 mm D / 5"x7.5")	£103
JBG0559		Angle & width adjustable footplates M (140 W x 240 mm D / 5.5"x9.5")	○
JBG0516		Angle & width adjustable footplates L (150 W x 290 mm D / 6"x11.5)	£135
JBG0517		Angle adjustable foot platform M (305 W x 275 mm D / 12"x10.75")	£270
JBX3001		Please select if Centre mount footplate is not required. (NC 0589)	○
<b>SIDE MOUNTED HANGERS LH</b>			
JBG0555		Amputee Support LH	£194
<b>SIDE MOUNTED HANGERS RH</b>			
JBG0554		Amputee Support RH	£194
<b>FOOTREST HANGER TYPE*</b>			
JBG0526		Standard hanger 80° pair Lower leg length: 350 - 450 mm	◆
JBG0503		Hanger Vari F pair Lower leg length: 290 - 460 mm Angle adjustment: 70 - 0° with tools	£120
JBG0533		Hanger Vari A pair Lower leg length: 290 - 460 mm Angle adjustment: 70 - 0° with handwheel and calf pad	£234
JBG0531		Hanger ADE pair Lower leg length: 290 - 460 mm (NC 2201) Powered Angle adjustment: 80 - 0° with calf pad	£669
JBG0501		Without Hanger pair (NC 0504, 0513)	○
<b>FOOTPLATE TYPE*</b>			
JBG0512		Fixed Footplate pair (NC 0533, 0531)	◆
JBG0504		Angle Adjustable Footplate pair (NC 0526)	£32
JBG0513		Angle and depth adjustable Footplate pair	£119
JBX3003		No Side Mounted Footplate required (NC 0503, 0533, 0531)	○
<b>CONTROLS*</b>			
JBG2201		DLX-REM110 = Drive only 1 Drive function with Drive only power module (NC 2837, 1092, 1093, 1094, 0977, 0589, 2231, 2232, 2230, 2247, 2257, 2233, 2240, 2258, 2219, 2245, 2270, 1908)	◆
JBG2203		DLX-REM211 = Drive & Actuators 3 Drive functions and power module incl. 2 Actuators (NC 1437, 2231, 2232, 2229, 2233, 2230, 2247, 2257, 2240, 2270, 2258, 1908, 1909)	○
JBG2204		DLX-REM216 = Drive, Actuators & Lights 3 Drive functions and power module incl. 2 Actuators & Lights (NC 2232, 2229, 2233, 2230, 2247, 2257, 2240, 2270, 2258)	£151
JBG2205		DLX-REM400 Configurable jack sockets (NC 2220)	£341
JBG2206		DLX-REM500 Requires a user input/alternative joystick (NC 2220, 2270, 1709, 1710)	£776
<b>REM400 TOGGLE OPTION</b>			
JBG2270		Assignable Toggle Switch Module for REM400 (NC 2203, 2204, 2206)	£316
<b>CONTROL BRACKETRY*</b>			
JBG1207		Standard Bracket RH (NC 2285, 2283, 2287, 2207, 2265, 2275, 1404)	○
JBG1208		Standard Bracket LH (NC 2286, 2284, 2288, 2208, 2266, 2276, 1405)	○
JBG1707		Swing Away Bracket RH (NC 2285, 2283, 2287, 2207, 2265, 2275, 1404)	£132
JBG1708		Swing Away Bracket LH (NC 2286, 2284, 2288, 2208, 2266, 2276, 1405)	£132
JBG1709		Matrx Resolve Swing in/out bracket RH (NC 2283, 2287, 2207, 2265, 2275, 1404)	£219
JBG1710		Matrx Resolve Swing in/out bracket LH (NC 2284, 2288, 2208, 2266, 2276, 1405)	£219
JBG2213		Centre Mount, swing-away to RH, Nucleus (NC 2283, 2287, 2207, 2265, 2275, 1415, 2217, 1404, 1405)	£498
JBG2214		Centre Mount, swing-away to LH, Nucleus (NC 2284, 2288, 2208, 2266, 2276, 1339, 2218, 1404, 1405)	£498
JBG2227		Bracket RH, mounted for Attendant (NC 1287, 2220, 2128, 2108, 2219, 1221, 1312)	£64
JBG2228		Bracket LH, mounted for Attendant (NC 1286, 2220, 2128, 2108, 2219, 1222, 1363)	£64
JBX2239		Individual mounting option per request (Special) Requires additional information/specification	POA

DUAL CONTROL			
JBG2220		Dual Control LiNX (DLX-REM050 as secondary remote)	(NC 2273, 2108, 2205, 2206, 1221, 1222, 2227, 2228) £395
JBG2219		Compact Dual Control LiNX (DLX-ACU200)	(NC 2108, 1221, 1222, 2227, 2228) £395
INVACARE LINX IDC DUAL CONTROL			
JBG2245		Invacare LiNX IDC Dual control	(NC 2108, 2259, 2260) £751
DUAL CONTROL BRACKERTY			
JBG2259		Bracket RH, Attendant	(NC 2291, 2245, 1312) ○
JBG2260		Bracket LH, Attendant	(NC 2292, 2245, 1363) ○
INVACARE LINX IDC DUAL CONTROL BRACKETRY			
JBG1221		Bracket RH, detachable	(NC 2291, 2220, 2219, 2227, 1312) ○
JBG1222		Bracket LH, detachable	(NC 2292, 2220, 2219, 2228, 1363) ○
LINX COMPACT REMOTE			
JBG2231		Proportional Compact Remote DLX-CR400	(NC 2203, 2204) £718
JBG2232		Proportional Compact Remote Low Force DLX-CR-400LF	(NC 2203, 2204) £960
LINX COMPACT REMOTE BRACKET			
JBG2207		RH Mount	(NC 1707, 2213, 2217, 1404) ○
JBG2208		LH Mount	(NC 1708, 2214, 2218, 1405) ○
JBG2265		Centre Mount, swing-away to RH (Nucleus) If several joysticks (or Joystick + REM500) should be mounted on one Nucleus, please select special mounting option 2215!	(NC 1707, 2213, 2217, 1404, 1405, 1415) £498
JBG2266		Centre Mount, swing-away to LH (Nucleus) If several joysticks (or Joystick + REM500) should be mounted on one Nucleus, please select special mounting option 2215!	(NC 1708, 2214, 2218, 1404, 1405, 1339) £498
JBX2215		Individual mounting option per request (Special) Requires additional information/specification	POA
CHIN CONTROL			
JBG2229		Proportional Chin Joystick (ASL138) SELECTED HEADREST ALWAYS WITH STD MOUNTING HARDWARE	(NC 2201, 2203, 2204, 2230, 2247, 2233, 2257, 2240, 2128) £1,734
UNI9250		Mo-vis Micro joystick	POA
UNI9251		Mo-vis Multi joystick	POA
CHIN CONTROL BRACKETRY MANUAL			
JBG1289		Swing Away Bracket RH, side manual (ASL628) Requires Elan Headrest	(NC 1350, 1432, 1427, 1496) ○
JBG1288		Swing Away Bracket LH, side manual (ASL628) Requires Elan Headrest	(NC 1350, 1432, 1427, 1496) ○
CHIN CONTROL BRACKETRY POWERED			
JBG1287		Swing Away Bracket RH, powered	(NC 2259, 1221, 2291, 1312, 1415) £2,183
JBG1286		Swing Away Bracket LH, powered	(NC 2260, 1222, 2292, 1363, 1339) £2,183
LINX HEAD CONTROL			
JBG2230		Head Control with ultrasonic proximity sensors (ATOM, ASL104) with standard wings M3	(NC 2203, 2204, 2229, 2233, 2324, 2325, 2327, 2257, 2240, 2175, 2171, 2170, 2197, 1350, 1432, 1427) £3,719
JBG2247		Head Control with ultrasonic proximity sensors (PROTON, ASL104P) with PROTON wings	(NC 2203, 2204, 2229, 2233, 2324, 2325, 2327, 2257, 2240, 2175, 2171, 2170, 2197, 1350, 1432, 1427) £4,355
LINX FOUR -SWITCH PROXIMITY ARRAY			
JBG2257		Four-Switch Proximity Array (ASL106)	(NC 2203, 2204, 2229, 2230, 2247, 2233, 2240, 2128) £5,760
JBG2275		Centre Mount, swing-away to RH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)	(NC 1707, 2213, 2217, 1404, 1405, 1415) £498
JBG2276		Centre Mount, swing-away to LH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)	(NC 1708, 2214, 2218, 1404, 1405, 1339) £498
JBX2269		Individual mounting option per request (Special) Requires additional information/specification	POA

<b>LINX WIRELESS ACCESSORIES FOR HEAD CONTROL AND PROXIMITY CONTROL</b>			
JBG2249		Wireless Triple Switch Receiver (ASL557-3)	£341
JBG2250		Wireless Mouse Emulator (ASL558)	£817
<b>LINX MICRO EXTREMITY CONTROL BRACKET</b>			
JBG2233		MEC Proportional Joystick (Micro Extremity Control ASL130)	(NC 2203, 2204, 2229, 2230, 2247, 2257, 2128) £3,612
<b>LINX HEAVY DUTY REMOTE</b>			
JBG2240		Proportional Heavy Duty Remote Joystick w mode switch (ASL133)	(NC 2203, 2204, 2229, 2230, 2247, 2257, 2128) £2,013
<b>LINX MICRO EXTREMITY AND HEAVY DUTY CONTROL BRACKET</b>			
JBG2285		Tray Mount RH (100x150 mm)	(NC 1707, 2213, 2217, 1404, 1405, 1415) ○
JBG2286		Tray Mount LH (100x150 mm)	(NC 1708, 2214, 2218, 1404, 1405, 1339) ○
JBG2283		Centre Mount, swing-away to RH, Nucleus	(NC 1707, 2213, 2217, 1404, 1405, 1415) £498
JBG2284		Centre Mount, swing-away to LH, Nucleus	(NC 1708, 2214, 2218, 1404, 1405, 1339) £498
JBG2287		Centre Mount, swing-away to RH, Nucleus Incl. Tray ASL603, 100x150 mm	(NC 1707, 2213, 2217, 1404, 1405, 1415) £498
JBG2288		Centre Mount, swing-away to LH, Nucleus Incl. Tray ASL603, 100x150 mm	(NC 1708, 2214, 2218, 1404, 1405, 1339) £498
JBX2289		Individual mounting option per request (Special) Requires additional information/specification	POA
<b>LINX 10-WAY SWITCH</b>			
JBG2251		10-Way Switch DLX-SW10-A for direct Powered Seating Access	£505
<b>LINX 10-WAY SWITCH BRACKET</b>			
JBG2217		Bracket RH, mounted for user	(NC 2207, 2287, 2283, 2275, 2213, 2265, 1404) ○
JBG2218		Bracket LH, mounted for user	(NC 2208, 2288, 2284, 2286, 2276, 2214, 2266, 1405) ○
JBG2291		Bracket RH, mounted for attendant	(NC 2259, 1221, 1287, 2108) £64
JBG2292		Bracket LH, mounted for attendant	(NC 2260, 1222, 1288, 2108) £64
<b>LINX USB CHARGER PORT</b>			
JBG2253		USB Charger Port DLX-USB02-A	(NC 2702, 2703, 2704, 2705, 2706) £165
<b>LINX USB CHARGER PORT BRACKET</b>			
JBG2281		RH Mount	○
JBG2282		LH Mount	○
<b>REMOTE STOP SWITCH</b>			
JBG2258		Remote Stop Switch (ASL504)	(NC 2203, 2204) £1,050
<b>LINX INPUT MODULE</b>			
JBG2273		IN200 Module (Jack port for external power button )	£242
<b>LINX INPUT MODULE</b>			
UNI9252		IN500 Module (for DB9 Drive input)	(NC 1094) £785
<b>FURTHER LINX COMPONENTS / ACCESSORIES</b>			
JBG2300		To be enclosed only	○
JBX2301		To be connected and programmed (Special)	POA
<b>BUDDY BUTTONS</b>			
JBG1685		Red Buddy Button	£88
JBG1686		Yellow Buddy Button	£88
JBG1687		Green Buddy Button	£88
JBG1690		Blue Buddy Button	£88
<b>EGG SWITCH</b>			
JBG2254		Egg Switch, black	£88

PIKO BUTTONS			
JBG1711		Red Piko Button	£158
JBG1712		Yellow Piko Button	£158
JBG1713		Green Piko Button	£158
JBG1714		Orange Piko Button	£158
JBG1715		Blue Piko Button	£158
JBG1716		Black Piko Button	£158
BODYPOINT JOYSTICK KNOBS			
JBG2721		3 inch U-Shaped, Joystick Handle, Bodypoint PC101A	£63
JBG2722		4 inch U-Shaped, Joystick Handle, Bodypoint PC102A	£63
JBG2723		Domed Rubber, Joystick Handle, Bodypoint PC107A	£63
JBG2724		Mushroom Grooved, Joystick Handle, Bodypoint PC110A	£63
5:1 SWITCH ADAPTOR			
JBG1681		5:1 Switch Adapter ( 9-pin sub D to 5 jack ports)	£121
SHROUD COLOUR*			
JBG2871		Mid range - Black/Anthracite Shroud	○
JBG2872		Mid range - Silver Shroud	○
JBG2873		Mid range - Bengal Red Shroud	○
JBG2874		Mid range - Mercury Blue Shroud	○
JBG2875		Mid range - Midnight Purple Shroud	○
JBX0299		Premium - Individual Shroud Colour RAL number required Extended lead lime: + 15 days	(NC 2883, 2884, 2885, 2886, 2888, 2889, 2890, 2891) £214
JBG2865		Premium - Racing Green Shroud	(NC 2883, 2884, 2885, 2886, 2887, 2889, 2890, 2891) £104
JBG2866		Premium - Deep Ocean Blue Satin Shroud	(NC 2883, 2884, 2885, 2886, 2887, 2888, 2890, 2891) £104
JBG2867		Premium - Frosty Blue Shroud	(NC 2883, 2884, 2885, 2886, 2887, 2888, 2889, 2891) £104
JBG2868		Premium - Gold Shroud	(NC 2883, 2884, 2885, 2886, 2887, 2888, 2889, 2890) £104
RIM COLOUR*			
JBG2882		Mid range - Black Rim Inserts	○
JBG2883		Mid range - Silver Rim Inserts	○
JBG2884		Mid range - Bengal Red Rim Inserts	○
JBG2885		Mid range - Mercury Blue Rim Inserts	○
JBG2886		Mid range - Midnight Purple Rim Inserts	○
JBX2887		Premium - Individual Rim Insert Colour RAL number required Extended lead lime: + 15 days	(NC 2872, 2873, 2874, 2875, 2865, 2866, 2867, 2868) ○
JBG2888		Premium - Racing Green Rim Inserts Price Included with Premium Shroud Kit	(NC 2872, 2873, 2874, 2875, 0299, 2866, 2867, 2868) ○
JBG2889		Premium - Deep Ocean Blue Satin Rim Inserts Price Included with Premium Shroud Kit	(NC 2872, 2873, 2874, 2875, 0299, 2865, 2867, 2868) ○
JBG2890		Premium - Frosty Blue Rim Inserts Price Included with Premium Shroud Kit	(NC 2872, 2873, 2874, 2875, 0299, 2865, 2866, 2868) ○
JBG2891		Premium - Gold Rim Inserts RAL1035 Price Included with Premium Shroud Kit	(NC 2872, 2873, 2874, 2875, 0299, 2865, 2866, 2867) ○
FORKS*			
JBG2631		Single sided fork	£352
JBG2632		Double sided fork Without castor rim insert (enlarges wide base to 656 mm)	◆



<b>FRONT WHEEL*</b>			
JBG0835		9" (225 mm) x 2.80 (70 mm) Castor, pneumatic	◆
JBG0836		9" (225 mm) x 2.80 (70 mm) Castor, puncture proof solid tyre	○
JBG0838		9" (225 mm) x 2.80 (70 mm) Castor, puncture protected Gel-sealed	○
JBG0828		8" (200 mm) x 2.50-4 (65 mm) Castor, puncture proof solid tyre (NC 1909)	○
<b>CASTOR LOCKS</b>			
JBG1490		Castor Lock	£197
<b>DRIVE WHEEL*</b>			
JBG0704		DW 14" (355 mm) x 3" (75 mm), Pneumatic	◆
JBG0705		DW 14" (355 mm) x 3" (75 mm), Solid inlay Dual density foam insert	○
JBG0706		DW 14" (355 mm) x 3" (75 mm), Puncture protected Gel-sealed	○
<b>C.T.C SUSPENSION*</b>			
JBG2870		CTC Suspension Control - Traction - Comfort	◆
<b>SUSPENSION SETTING*</b>			
JBG1910		Soft (suggested for up to 70 kg User Weight) (NC 1094)	◆
JBG1912		Medium (suggested for up to 90 kg User Weight) (NC 1094)	○
JBG1911		Firm (suggested for up to 110 kg User Weight)	○
JBG1915		Very firm (suggested for up to 136 kg User Weight)	○
<b>MOTOR TYPE*</b>			
JBG1294		Indoor setting, reduced speed approx. 3 km/h	○
JBG1282		6 km/h	◆
JBG1908		10 km/h Package Always combined with LED Light system with Indicators, 30 Lux headlights (NC 2203, 1437, 1254)	£543
JBG1909		12 km/h Ultra High Speed Package Always combined with LED Light system with Indicators, 30 Lux headlights (NC 0828, 2201, 2203, 1437, 1254)	£653
<b>BATTERY TYPE*</b>			
JBG1259		60 Ah Gel	◆
JBG1260		73,5 Ah Gel (NC 1450)	£92
<b>G TRAC</b>			
JBG1215		G-Trac Module	£323
<b>CHARGER*</b>			
JBG1343		Opti Charge Charger 8A (IP44)	◆
JBG1313		Opti Charge Charger 10A (IP44)	£151
<b>ANTI-TIPPERS*</b>			
JBG1205		Anti-tipper	◆
<b>LIGHTS</b>			
JBG1437		LED Light system with Indicators 30 Lux headlights (NC 2203)	£480
JBG2837		LED Light system extra bright with indicators 60 Lux headlights (NC 2203, 1282)	£143
JBG1254		Without Light system (NC 1909)	○
<b>KERB RIDER</b>			
JBG1244		Kerb Rider With park position feature (NC 0556, 0589, 1294)	£212

TRANSPORTATION KIT*			
JBG1450		Preparation for DAHL-Docking station (Adequate headrest is compulsory) Includes a mounting bracket in the battery box Docking Station for car implementation has to be ordered separately with DAHL Engineering	£586
JBG1494		4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19) Adequate headrest is compulsory	◆
JBG1401		6-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19)) required for variants with combined Lift & Tilt module Adequate headrest is compulsory	(NC 1091, 1092, 1093) ○
EASY ACCESS FRONT TIE DOWNS			
JBG2894		Easily accessible front tie-down points, additionally mounted for 4- or 6- point-tie-down For transport worthy versions (1494/1401) only	£49
LAP BELT*			
JBG1329		Lap belt with Auto Lock Standard mount on backrest	○
JBG2113		Lap belt with Auto Lock Flexible mount in side rails	◆
POSITIONING BELT			
JBG2701		BP Evoflex Pelvic Stabilizer, Rehab Latch Size L Padded length 620 mm	(NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253) £77
JBG2702		BP Evoflex Pelvic Stabilizer, Rehab Latch Size M Padded length 460 mm	(NC 2123, 2124, 2184, 2185, 2186, 2187, 2188, 2189, 2253) £77
JBG2703		BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size L Padded length 620 mm	(NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253) £79
JBG2704		BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size M Padded length 460 mm	(NC 2123, 2124, 2184, 2185, 2186, 2187, 2188, 2189, 2253) £79
JBG2705		BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size L Padded length 620 mm	(NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253) £79
JBG2706		BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size M Padded length 460 mm	(NC 2123, 2124, 2184, 2185, 2186, 2187, 2188, 2189, 2253) £79
REAR VIEW MIRROR			
JBG1415		Rear view mirror RH	(NC 2217, 2285, 2213, 2265, 2275, 2283, 2287) £44
JBG1339		Rear view mirror LH	(NC 2218, 2286, 2214, 2266, 2276, 2284, 2288) £44
OTHER OPTIONS			
JBG2128		Bag Hooks on Modulite Seat Max load up to 6 kg = impact on max user weight	(NC 2172, 2173, 2227, 2228, 2229, 2233, 2240, 2257, 2291, 2292) £77
JBG1312		Cane holder RH	(NC 2227, 2259, 1221) £42
JBG1363		Cane holder LH	(NC 2228, 2260, 1222) £42
JBG1416		Fender for Castors	£88
JBG1417		Fender for Drive wheels	£169
JBG1404		Tray, Swing Away RH	(NC 1707, 2275, 2276, 2287, 2288, 2283, 2284, 2285, 2286, 2265, 2266, 2207, 2213, 2214, 2217) £167
JBG1405		Tray, Swing Away LH	(NC 1708, 2275, 2276, 2287, 2288, 2283, 2284, 2285, 2286, 2265, 2266, 2208, 2213, 2214, 2218) £167

<b>TOTAL PRICE</b>	
--------------------	--

**Non Conforming Options:**

JBG2121	380 mm Seat Plate Can be reduced to 330 mm in combination with inrigged armrests
JBG2122	430 mm Seat Plate Can be reduced to 380 mm in combination with inrigged armrests
JBG2123	480 mm Seat Plate
JBG2124	530 mm Seat Plate
JBG2106	Body fit Kit Outrigging Backrest Brackets and Pivot Point (outrig: 1 inch/2,54 cm per side)
JBG2132	Seat Depth 410 mm
JBG2143	Matrx Flo-tech Contour Visco; (2-way stretch material)
JBG2108	Backposts without push handles
JBG2171	Tension Adjustable Upholstery nylon overlay short = 480 mm / 7 straps measured from seat plate to top of upholstery 480 mm
JBG2179	Flex3 Tension Adjustable Upholstery height 540 mm (seat plate to top of upholstery)
JBG2172	Flex3 adjustable Backrest Plates without Backrest Cushion
JBG2173	Flex3 adjustable Backrest Plates with Laguna Backrest cushion
JBG2174	Matrx Posture Back, height 400 mm (Mesh fabric) Corresponds to SW setting chosen (up to 136 kg max user weight)
JBG2175	Matrx Posture Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen (up to 136 kg max user weight)
JBG2176	Matrx Elite Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen (up to 136 kg max user weight)
JBG2177	Matrx Elite TR Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen
JBG2178	Matrx Elite Deep Back, height 500 mm (Mesh fabric) Corresponds to SW setting chosen (up to 136 kg max user weight)
JBG2170	Without Backrest Upholstery (Back Posts/Back Frame only)
JBG2155	Black Spacer Fabric Upholstery
JBG2181	T-Arm detachable, Extended height Adjustment range from seat plate: 295 - 365 mm
JBG2198	Padded Armrest with padded side guard (WATER FALL TYPE)
JBG2183	Inrigged Armrests Reduces SW by 25 mm each side
JBG2193	Flip-up Armrest Up to 150 kg max user weight
JBG2192	Following Armrest Height adjustable from 230 - 300 mm
JBG2194	Following Armrest Height adjustable from 300 - 360 mm
JBG2195	Short Armpad (325 mm)
JBG2196	Long Armpad (375 mm)
JBG2197	Long Ergonomic Armpad (397 mm)
JBG2199	Soft Armpad In combination with 10-way switch, the switch 2251 has to be mounted on the opposite side of the remote
JBG2184	Hip Support LH, pad size 140 x 100 mm
JBG2186	Hip Support LH, pad size 230 x 100 mm
JBG2188	Hip Support LH, pad size 330 x 100 mm
JBG2185	Hip Support RH, pad size 140 x 100 mm
JBG2187	Hip Support RH, pad size 230 x 100 mm
JBG2189	Hip Support RH, pad size 330 x 100 mm

JBG1350	Neckrest depth adjustable
JBG1432	Headrest sideways adjustable
JBG1427	Headrest with cheek support
JBG2324	Elan Standard Pad 255 mm / 10" wide w/standard cover (ESP26) Complete with Elan Headrest Hardware
JBG2325	Elan Occipital Pad 230x180 mm / 9"x7" w/standard cover (EOP23) Complete with Elan Headrest Hardware
JBG2327	Elan 4 Point Headrest Pad 280x255 mm / 11"x10" w/standard Complete with Elan Headrest Hardware
JBG1496	Non Invacare Headrest but profile mount fitting Flex3 hard back
JBG1445	Lateral Support with extended arm, multi-adjustable pad
JBG1091	Manually settable Seat Tilt 0 - 15°
JBG1092	Powered "fixed pivot" Tilt 20°
JBG1093	Powered Tilt (0° - 30°) with CoG shift
JBG1094	Seat lifter with powered tilt (300 mm, 30°) with CoG shift
JBG0977	Powered Recliner 30°
JBG0556	Centre mount legrest, angle settable Length: 280 - 385 mm
JBG0589	LNX Power centre mount legrest, length compensating Length: 340 - 410 mm
JBG0517	Angle adjustable foot platform M (305 W x 275 mm D / 12"x10.75")
JBG0526	Standard hanger 80° pair Lower leg length: 350 - 450 mm
JBG0503	Hanger Vari F pair Lower leg length: 290 - 460 mm Angle adjustment: 70 - 0° with tools
JBG0533	Hanger Vari A pair Lower leg length: 290 - 460 mm Angle adjustment: 70 - 0° with handwheel and calf pad
JBG0531	Hanger ADE pair Lower leg length: 290 - 460 mm Powered Angle adjustment: 80 - 0° with calf pad
JBG0504	Angle Adjustable Footplate pair
JBG0513	Angle and depth adjustable Footplate pair
JBG2201	DLX-REM110 = Drive only 1 Drive function with Drive only power module
JBG2203	DLX-REM211 = Drive & Actuators 3 Drive functions and power module incl. 2 Actuators
JBG2204	DLX-REM216 = Drive, Actuators & Lights 3 Drive functions and power module incl. 2 Actuators & Lights
JBG2205	DLX-REM400 Configurable jack sockets
JBG2206	DLX-REM500 Requires a user input/alternative joystick
JBG2270	Assignable Toggle Switch Module for REM400
JBG1707	Swing Away Bracket RH
JBG1708	Swing Away Bracket LH
JBG1709	Matrx Resolve Swing in/out bracket RH
JBG1710	Matrx Resolve Swing in/out bracket LH
JBG2213	Centre Mount, swing-away to RH, Nucleus
JBG2214	Centre Mount, swing-away to LH, Nucleus
JBG2227	Bracket RH, mounted for Attendant
JBG2228	Bracket LH, mounted for Attendant
JBG2220	Dual Control LiNX (DLX-REM050 as secondary remote)
JBG2219	Compact Dual Control LiNX (DLX-ACU200)
JBG2245	Invacare LiNX IDC Dual control
JBG2259	Bracket RH, Attendant
JBG2260	Bracket LH, Attendant
JBG1221	Bracket RH, detachable

JBG1222	Bracket LH, detachable
JBG2231	Proportional Compact Remote DLX-CR400
JBG2232	Proportional Compact Remote Low Force DLX-CR-400LF
JBG2207	RH Mount
JBG2208	LH Mount
JBG2265	Centre Mount, swing-away to RH (Nucleus) If several joysticks (or Joystick + REM500) should be mounted on one Nucleus, please select special mounting option 2215!
JBG2266	Centre Mount, swing-away to LH (Nucleus) If several joysticks (or Joystick + REM500) should be mounted on one Nucleus, please select special mounting option 2215!
JBG2229	Proportional Chin Joystick (ASL138) SELECTED HEADREST ALWAYS WITH STD MOUNTING HARDWARE
JBG1289	Swing Away Bracket RH, side manual (ASL628) Requires Elan Headrest
JBG1288	Swing Away Bracket LH, side manual (ASL628) Requires Elan Headrest
JBG1287	Swing Away Bracket RH, powered
JBG1286	Swing Away Bracket LH, powered
JBG2230	Head Control with ultrasonic proximity sensors (ATOM, ASL104) with standard wings M3
JBG2247	Head Control with ultrasonic proximity sensors (PROTON, ASL104P) with PROTON wings
JBG2257	Four-Switch Proximity Array (ASL106)
JBG2275	Centre Mount, swing-away to RH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)
JBG2276	Centre Mount, swing-away to LH, Nucleus Incl. Eclipse Tray ASL619, 360x180 mm)
JBG2233	MEC Proportional Joystick (Micro Extremity Control ASL130)
JBG2240	Proportional Heavy Duty Remote Joystick w mode switch (ASL133)
JBG2285	Tray Mount RH (100x150 mm)
JBG2286	Tray Mount LH (100x150 mm)
JBG2283	Centre Mount, swing-away to RH, Nucleus
JBG2284	Centre Mount, swing-away to LH, Nucleus
JBG2287	Centre Mount, swing-away to RH, Nucleus Incl. Tray ASL603, 100x150 mm
JBG2288	Centre Mount, swing-away to LH, Nucleus Incl. Tray ASL603, 100x150 mm
JBG2217	Bracket RH, mounted for user
JBG2218	Bracket LH, mounted for user
JBG2291	Bracket RH, mounted for attendant
JBG2292	Bracket LH, mounted for attendant
JBG2253	USB Charger Port DLX-USB02-A
JBG2258	Remote Stop Switch (ASL504)
JBG2273	IN200 Module (Jack port for external power button )
JBG2872	Mid range - Silver Shroud
JBG2873	Mid range - Bengal Red Shroud
JBG2874	Mid range - Mercury Blue Shroud
JBG2875	Mid range - Midnight Purple Shroud
JBX0299	Premium - Individual Shroud Colour RAL number required Extended lead lime: + 15 days
JBG2865	Premium - Racing Green Shroud
JBG2866	Premium - Deep Ocean Blue Satin Shroud
JBG2867	Premium - Frosty Blue Shroud
JBG2868	Premium - Gold Shroud
JBG2883	Mid range - Silver Rim Inserts

JBG2884	Mid range - Bengal Red Rim Inserts
JBG2885	Mid range - Mercury Blue Rim Inserts
JBG2886	Mid range - Midnight Purple Rim Inserts
JBX2887	Premium - Individual Rim Insert Colour RAL number required Extended lead lime: + 15 days
JBG2888	Premium - Racing Green Rim Inserts Price Included with Premium Shroud Kit
JBG2889	Premium - Deep Ocean Blue Satin Rim Inserts Price Included with Premium Shroud Kit
JBG2890	Premium - Frosty Blue Rim Inserts Price Included with Premium Shroud Kit
JBG2891	Premium - Gold Rim Inserts RAL1035 Price Included with Premium Shroud Kit
JBG0828	8" (200 mm) x 2.50-4 (65 mm) Castor, puncture proof solid tyre
JBG1910	Soft (suggested for up to 70 kg User Weight)
JBG1912	Medium (suggested for up to 90 kg User Weight)
JBG1294	Indoor setting, reduced speed approx. 3 km/h
JBG1282	6 km/h
JBG1908	10 km/h Package Always combined with LED Light system with Indicators, 30 Lux headlights
JBG1909	12 km/h Ultra High Speed Package Always combined with LED Light system with Indicators, 30 Lux headlights
JBG1437	LED Light system with Indicators 30 Lux headlights
JBG2837	LED Light system extra bright with indicators 60 Lux headlights
JBG1254	Without Light system
JBG1244	Kerb Rider With park position feature
JBG1450	Preparation for DAHL-Docking station (Adequate headrest is compulsory) Includes a mounting bracket in the battery box Docking Station for car implementation has to be ordered separately with DAHL Engineering
JBG2701	BP Evoflex Pelvic Stabilizer, Rehab Latch Size L Padded length 620 mm
JBG2702	BP Evoflex Pelvic Stabilizer, Rehab Latch Size M Padded length 460 mm
JBG2703	BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size L Padded length 620 mm
JBG2704	BP 2-Point Hip Belt, Ctr-Pull, Push Butt, Cinch-Mount Size M Padded length 460 mm
JBG2705	BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size L Padded length 620 mm
JBG2706	BP 2-Point Hip Belt, Dual-Pull, Side-Release Butt. Cinch-Mount Size M Padded length 460 mm
JBG1415	Rear view mirror RH
JBG1339	Rear view mirror LH
JBG2128	Bag Hooks on Modulite Seat Max load up to 6 kg = impact on max user weight
JBG1312	Cane holder RH
JBG1363	Cane holder LH
JBG1416	Fender for Castors
JBG1404	Tray, Swing Away RH
JBG1405	Tray, Swing Away LH



For a more complex customer requirement please refer to our Invacare Unique special order form. email [specialsuk@invacare.com](mailto:specialsuk@invacare.com) or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Description of special request (please include relevant option code or part number):