



Yes, you can:

Tel: 01656 776222  
 Email: ordersuk@invacare.com  
 Online Spares available at: www.invacare.co.uk

# Manual Retail Prescription Form

JANUARY 2025  
 Customer Ref:

|                          |                    |
|--------------------------|--------------------|
| <b>Customer Ref*</b>     | <b>Account No*</b> |
| <b>Contact No*</b>       | <b>Order No*</b>   |
| <b>Delivery Address*</b> | <b>Order Qty*</b>  |
|                          | <b>Order date*</b> |

## Kuschall Champion 2.0

### Standard Features

- Max User Weight 120kg (19 Stone)
- Lightweight Folding Fully Hydroformed Aluminium Wheelchair
- Multiple adjustments, Front, Rear Seat Heights and Centre of Gravity
- Aluminium Front Frame with 75° or 90° Front Frame Angles
- 24" Self Propel Q/R Wheels with Grey Aluminium Handrims
- and Pneumatic Tyres
- Choice of 7 Seat Widths 360mm to 480mm (14" - 19")
- Choice of 8 Seat Depths 340mm to 480mm (13" - 19")
- Choice of 5 Standard Frame Colours
- Range of Front Seat Heights 450mm to 540mm (17" - 21")
- Range of Rear Seat Heights 390mm to 480mm (15" - 19")
- Folding Backrest, Angle Adjustable and Soft Lock Feature
- Black Nylon Light Upholstery
- Angle and Height Adjustable Rigid Footrest, with Calf Strap
- Without Push Handles, Flip up Carbon Clothesguards
- Unfolding Help Cord. Black Comfort Front Castor
- Standard Brake Lever Bent
- Tool Kit



Please make your selection with a cross in the box    ○ Free of charge    ◆ standard feature    POA special request

Customer weight =                      KG

| Code                                                                                                    | Description                                                                                             |     |
|---------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------|-----|
| <b>BASE PRODUCT*</b>                                                                                    |                                                                                                         |     |
| DHA0031                                                                                                 | <input type="checkbox"/> Kuschall Champion 2.0                                                          | ○   |
| <b>UNIQUE SPECIALS (For special modifications only) Additional lead time of 4 weeks</b>                 |                                                                                                         |     |
| DHA0000                                                                                                 | <input type="checkbox"/> UNIQUE Special Request - For use when requesting any special options or builds | POA |
| <b>UNIQUE CUSTOMISATION AND PERSONALISATION OPTIONS Additional 4 weeks lead time for Unique options</b> |                                                                                                         |     |
| UHA0001                                                                                                 | <input type="checkbox"/> UNIQUE Personalisation - For use when selecting UHA customisation options      | ○   |

| <b>FRAME*</b>                         |  |                                                                                            |                                                     |        |
|---------------------------------------|--|--------------------------------------------------------------------------------------------|-----------------------------------------------------|--------|
| DHA1150                               |  | Frame Angle 75° Aluminium                                                                  | (NC 1350)                                           | £2,683 |
| DHA1160                               |  | Frame Angle 90° Aluminium                                                                  | (NC 1350)                                           | £2,683 |
| DHA1170                               |  | Frame Angle 75° Aluminium - Tapered Front Frame                                            | (NC 1320, 1325, 7010)                               | £2,683 |
| DHA1180                               |  | Frame Angle 90° Aluminium - Tapered Front Frame                                            | (NC 1320, 1325, 7010)                               | £2,683 |
| <b>CHAMPION SL HYBRID RIGID FRAME</b> |  |                                                                                            |                                                     |        |
| UHA1260                               |  | Complete frame with fully welded front seat height & centre of gravity rear wheel supports | (NC 6180)                                           | £450   |
| <b>SEAT WIDTH*</b>                    |  |                                                                                            |                                                     |        |
| DHA1320                               |  | 360 mm                                                                                     | (NC 2710, 2748, 7000)                               | ○      |
| DHA1325                               |  | 380 mm                                                                                     | (NC 7000)                                           | ○      |
| DHA1330                               |  | 400 mm                                                                                     |                                                     | ○      |
| DHA1335                               |  | 420 mm                                                                                     |                                                     | ○      |
| DHA1340                               |  | 440 mm                                                                                     |                                                     | ○      |
| DHA1345                               |  | 460 mm                                                                                     |                                                     | ○      |
| DHA1350                               |  | 480 mm                                                                                     |                                                     | ○      |
| <b>SEAT DEPTH*</b>                    |  |                                                                                            |                                                     |        |
| DHA1405                               |  | 340 mm                                                                                     | (NC 2710, 2720, 2748, 2750, 2860, 7810, 7820, 7830) | ○      |
| DHA1415                               |  | 360 mm                                                                                     | (NC 2710, 2720, 2748, 2860, 7810, 7820, 7830)       | ○      |
| DHA1425                               |  | 380 mm                                                                                     | (NC 2720)                                           | ○      |
| DHA1435                               |  | 400 mm                                                                                     |                                                     | ○      |
| DHA1445                               |  | 420 mm                                                                                     |                                                     | ○      |
| DHA1460                               |  | 440 mm                                                                                     |                                                     | ○      |
| DHA1470                               |  | 460 mm                                                                                     |                                                     | ○      |
| DHA1480                               |  | 480 mm                                                                                     | (NC 2710, 2748)                                     | ○      |
| <b>BACK HEIGHT*</b>                   |  |                                                                                            |                                                     |        |
| DHA1615                               |  | 300 mm                                                                                     | (NC 2200, 2240, 2400)                               | ○      |
| DHA1620                               |  | 315 mm                                                                                     | (NC 2200, 2240, 2400)                               | ○      |
| DHA1625                               |  | 330 mm                                                                                     | (NC 2240)                                           | ○      |
| DHA1630                               |  | 345 mm                                                                                     | (NC 2240)                                           | ○      |
| DHA1635                               |  | 360 mm                                                                                     |                                                     | ○      |
| DHA1640                               |  | 375 mm                                                                                     |                                                     | ○      |
| DHA1645                               |  | 390 mm                                                                                     |                                                     | ○      |
| DHA1650                               |  | 405 mm                                                                                     |                                                     | ○      |
| DHA1655                               |  | 420 mm                                                                                     |                                                     | ○      |
| DHA1660                               |  | 435 mm                                                                                     |                                                     | ○      |
| DHA1665                               |  | 450 mm                                                                                     |                                                     | ○      |
| DHA1670                               |  | 465 mm                                                                                     |                                                     | ○      |
| <b>FRONT SEAT HEIGHT*</b>             |  |                                                                                            |                                                     |        |
| DHA1740                               |  | 450 mm                                                                                     | (NC 1946, 1948, 1950, 1952, 1954, 1956, 1958, 1960) | ○      |
| DHA1745                               |  | 460 mm                                                                                     | (NC 1948, 1950, 1952, 1954, 1956, 1958, 1960)       | ○      |
| DHA1750                               |  | 470 mm                                                                                     | (NC 1950, 1952, 1954, 1956, 1958, 1960)             | ○      |
| DHA1755                               |  | 480 mm                                                                                     | (NC 1952, 1954, 1956, 1958, 1960)                   | ○      |
| DHA1760                               |  | 490 mm                                                                                     | (NC 1954, 1956, 1958, 1960)                         | ○      |
| DHA1765                               |  | 500 mm                                                                                     | (NC 1956, 1958, 1960)                               | ○      |
| DHA1770                               |  | 510 mm                                                                                     | (NC 1958, 1960)                                     | ○      |
| DHA1775                               |  | 520 mm                                                                                     | (NC 1960)                                           | ○      |
| DHA1780                               |  | 530 mm                                                                                     | (NC 5000)                                           | ○      |
| DHA1785                               |  | 541 mm Only possible with 5" castors                                                       | (NC 5000, 5050, 5100, 5150, 5200, 5250)             | ○      |

**REAR SEAT HEIGHT. On 75° frame, Rear Seat Height can be 0mm - 70mm lower then front seat height and 90° can be 30mm to 100mm lower\***

|         |  |                                                 |                 |   |
|---------|--|-------------------------------------------------|-----------------|---|
| DHA1815 |  | 390 mm Only possible with 22" & 24" rear wheels | (NC 6420, 6430) | ○ |
| DHA1820 |  | 400 mm Only possible with 22" & 24" rear wheels | (NC 6420, 6430) | ○ |
| DHA1825 |  | 410 mm                                          | (NC 6430)       | ○ |
| DHA1830 |  | 420 mm                                          | (NC 6430)       | ○ |
| DHA1835 |  | 430 mm                                          |                 | ○ |
| DHA1840 |  | 440 mm                                          |                 | ○ |
| DHA1845 |  | 450 mm                                          |                 | ○ |
| DHA1850 |  | 460 mm                                          |                 | ○ |
| DHA1855 |  | 470 mm Only possible with 22" & 24" rear wheels | (NC 6410)       | ○ |
| DHA1860 |  | 480 mm Only possible with 22" & 24" rear wheels | (NC 6410)       | ○ |

**LOWER LEG LENGTH. Cushion height must be deducted\***

|         |  |        |  |   |
|---------|--|--------|--|---|
| DHA1924 |  | 320 mm |  | ○ |
| DHA1926 |  | 330 mm |  | ○ |
| DHA1928 |  | 340 mm |  | ○ |
| DHA1930 |  | 350 mm |  | ○ |
| DHA1932 |  | 360 mm |  | ○ |
| DHA1934 |  | 370 mm |  | ○ |
| DHA1936 |  | 380 mm |  | ○ |
| DHA1938 |  | 390 mm |  | ○ |
| DHA1940 |  | 400 mm |  | ○ |
| DHA1942 |  | 410 mm |  | ○ |
| DHA1944 |  | 420 mm |  | ○ |
| DHA1946 |  | 430 mm |  | ○ |
| DHA1948 |  | 440 mm |  | ○ |
| DHA1950 |  | 450 mm |  | ○ |
| DHA1952 |  | 460 mm |  | ○ |
| DHA1954 |  | 470 mm |  | ○ |
| DHA1956 |  | 480 mm |  | ○ |
| DHA1958 |  | 490 mm |  | ○ |
| DHA1960 |  | 500 mm |  | ○ |

**BACKREST\***

|         |  |                                                                                              |  |   |
|---------|--|----------------------------------------------------------------------------------------------|--|---|
| DHA2140 |  | Angle Adjustable, Forward Folding with Soft Lock. Backrest remains locked for transportation |  | ◆ |
|---------|--|----------------------------------------------------------------------------------------------|--|---|

**BACKREST UPHOLSTERY\***

|         |  |                                                                                                                                |                                                           |      |
|---------|--|--------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------|------|
| DHA2200 |  | Black Standard Light Upholstery, without Tension Adjustment                                                                    | (NC 2241, 2242, 2243, 2244, 2245, 0002, 0003, 0004, 2430) | ◆    |
| DHA2220 |  | Tension Adjustable Backrest, with Premium Black Nylon and Reflective Trim Upholstery Cover, Back Pocket/hook and loop fastener | (NC 2241, 2242, 2243, 2244, 2245, 0002, 0003, 0004)       | £137 |
| DHA2230 |  | Tension Adjustable Backrest, with Black Aero, Breathable Spacer Upholstery Cover, Back Pocket/hook and loop fastener           |                                                           | £223 |
| DHA2240 |  | Without Backrest Upholstery Select if fitting a Matrx/Fixed backrest, no upholstery required                                   | (NC 2430, 2241, 2242, 2243, 2244, 2245, 0002, 0003, 0004) | ○    |

| BACK UPHOLSTERY STITCHING COLOUR |  |                                                                                                                                                |      |
|----------------------------------|--|------------------------------------------------------------------------------------------------------------------------------------------------|------|
| UHA2241                          |  | Neon Yellow                                                                                                                                    | £75  |
| UHA2242                          |  | Neon Pink                                                                                                                                      | £75  |
| UHA2243                          |  | Blue                                                                                                                                           | £75  |
| UHA2244                          |  | White                                                                                                                                          | £75  |
| UHA2245                          |  | Red                                                                                                                                            | £75  |
| UHA0002                          |  | Your personal colour stitching. Add RAL colour of your choice in special comments box, for a personalised colour stitching on back upholstery. | £75  |
| BACK UPHOLSTERY EMBROIDERY       |  |                                                                                                                                                |      |
| UHA0003                          |  | Your personal design embroidery. Supply your design with your order, subject to copyright. Design size up to 250mm x 250mm (NC 0004)           | £199 |
| BACK UPHOLSTERY FLOCKING         |  |                                                                                                                                                |      |
| UHA0004                          |  | Your personal flock design. Supply your design with your order, subject to copyright. Design size up to 280mm x 400mm                          | £99  |
| PUSH HANDLES*                    |  |                                                                                                                                                |      |
| DHA2400                          |  | Without Push Handles                                                                                                                           | ◆    |
| DHA2410                          |  | Push Handles Fixed Short/Mini Grip                                                                                                             | £133 |
| DHA2420                          |  | Push Handles Fixed Long Grip                                                                                                                   | ○    |
| DHA2430                          |  | Push Handles Height Adjustable, Rear Set (NC 3480)                                                                                             | £254 |
| DHA2440                          |  | Push Handles Height Adjustable, Integrated in backrest                                                                                         | £219 |
| DHA2450                          |  | Push Handles Foldable                                                                                                                          | £198 |
| BACKREST ANGLE*                  |  |                                                                                                                                                |      |
| DHA2510                          |  | 76° angle from the seat                                                                                                                        | ○    |
| DHA2530                          |  | 80.5° angle from the seat                                                                                                                      | ○    |
| DHA2550                          |  | 85° angle from the seat                                                                                                                        | ◆    |
| DHA2570                          |  | 89.5° angle from the seat                                                                                                                      | ○    |
| DHA2590                          |  | 94° angle from the seat                                                                                                                        | ○    |
| SEAT UPHOLSTERY*                 |  |                                                                                                                                                |      |
| DHA2600                          |  | Standard Seat Upholstery                                                                                                                       | ◆    |
| SEAT CUSHION                     |  |                                                                                                                                                |      |
| DHA2700                          |  | Standard Seat Cushion, 60mm high Double Foam with Visco Top Layer - Comfort/Low Risk (NC 2810, 2820, 2830, 2840, 2850)                         | £66  |
| DHA2710                          |  | Flo-Tech Lite Cushion 78mm high- Comfort/Low Risk (NC 2810, 2820, 2830, 2840, 2850)                                                            | £121 |
| DHA2720                          |  | Matrx Libra Cushion 83mm high, Black Smooth Cover - Very High Risk                                                                             | £410 |
| DHA2748                          |  | Contour Visco NG Cushion 100mm high - High Risk (NC 2810, 2820, 2830, 2840, 2850)                                                              | £169 |
| DCA2750                          |  | Vicair Active O2 lightweight cushion 60mm high - High risk (NC 2810, 2820, 2830, 2840, 2850, 2860)                                             | £521 |
| SEAT BAG                         |  |                                                                                                                                                |      |
| DHA2800                          |  | Seat Bag, Click in and removable, includes the fixation clips (NC 7010)                                                                        | £74  |
| MATRX LIBRA CUSHION ACCESSORIES  |  |                                                                                                                                                |      |
| DHA2810                          |  | Fluid supplement sac (Pair)                                                                                                                    | £44  |
| DHA2820                          |  | Tapered wedge (Pair)                                                                                                                           | £30  |
| DHA2830                          |  | Abductor                                                                                                                                       | £23  |
| DHA2840                          |  | Lateral side wedge                                                                                                                             | £20  |
| DHA2850                          |  | Anterior front wedge                                                                                                                           | £27  |
| SEAT RIGIDIZER                   |  |                                                                                                                                                |      |
| DHA2860                          |  | Seat rigidizer                                                                                                                                 | £39  |

| FOOTREST*                  |  |                                                                                                                                                         |      |
|----------------------------|--|---------------------------------------------------------------------------------------------------------------------------------------------------------|------|
| DHA3240                    |  | Rigid, Plastic Footplate, Set to front, Height and Angle Adjustable. Note: This does not flip up for transfer, but folds easy for transport and storage | ◆    |
| DHA3250                    |  | Rigid, Plastic Footplate, Set to rear, Height and Angle Adjustable. Note: This does not flip up for transfer, but folds easy for transport and storage  | ○    |
| CALF STRAP*                |  |                                                                                                                                                         |      |
| DHA3320                    |  | Calf Strap Mounted on Supporter                                                                                                                         | ◆    |
| DHA3330                    |  | Calf Strap Mounted on Frame                                                                                                                             | ○    |
| CLOTHESGUARDS/MUDGUARDS*   |  |                                                                                                                                                         |      |
| DHA3405                    |  | Composite Clothesguards, Flip Back - (only possible with 24" and 25" rear wheels) (NC 6430)                                                             | ○    |
| DHA3420                    |  | Carbon Clothesguards, Flip Up                                                                                                                           | ◆    |
| DHA3425                    |  | Composite Mudguards, Fixed (only possible with 24" and 25" rear wheels) (NC 3490, 6430)                                                                 | £65  |
| DHA3430                    |  | Carbon Mudguards, Fixed (NC 3490)                                                                                                                       | £247 |
| DHA3440                    |  | Carbon Mudguards, Removable (NC 3490)                                                                                                                   | £307 |
| DHA3450                    |  | Composite Mudguards, Removable (only possible with 24" and 25" rear wheels) (NC 3490, 6430)                                                             | £65  |
| ARMREST                    |  |                                                                                                                                                         |      |
| DHA3480                    |  | Tubular Armrest. Height Adjustable, Removable, Non-locking, Swing-away and Comfort Padded, no Clothesguard included                                     | £241 |
| DHA3490                    |  | Armrest, Height Adjustable, Short Armpad (T-Arm), Non-locking, no Clothesguard included                                                                 | £106 |
| FRAME COLOUR*              |  |                                                                                                                                                         |      |
| DHA4000                    |  | Space Black (Standard) (NC 4066)                                                                                                                        | ○    |
| DHA4010                    |  | Anthracite Grey (Standard) (NC 4066)                                                                                                                    | ○    |
| DHA4020                    |  | Pure White (Standard) (NC 4066, 4220)                                                                                                                   | ○    |
| DHA4030                    |  | Carmine Red (Standard) (NC 4066)                                                                                                                        | ○    |
| DHA4040                    |  | Azure Blue (Standard) (NC 4066)                                                                                                                         | ○    |
| DHA4060                    |  | Flashy Yellow (Additional 10 day lead time) (NC 4066)                                                                                                   | £31  |
| DHA4100                    |  | Black Glossy (Additional 10 day lead time) (NC 4066)                                                                                                    | £31  |
| DHA4130                    |  | Silver Grey (Additional 10 day lead time) (NC 4066)                                                                                                     | £31  |
| DHA4145                    |  | Pearl Night Blue (Additional 15 day lead time) (NC 4066)                                                                                                | £31  |
| DHA4150                    |  | Ruby Red (Additional 10 day lead time) (NC 4066)                                                                                                        | £31  |
| DHA4170                    |  | Dynamite Orange (Additional 10 day lead time) (NC 4066)                                                                                                 | £31  |
| DHA4180                    |  | Forest Green (Additional 10 day lead time) (NC 4066)                                                                                                    | £31  |
| UHA4101                    |  | Mirror Polish with standard varnish finish (Premium ceramic glaze finish available, for this option select UHA4066)                                     | £399 |
| UHA0008                    |  | Your personal frame colour. Add RAL colour of your choice in special comments box, for a personalised colour frame                                      | £299 |
| FULL FRAME COLOUR          |  |                                                                                                                                                         |      |
| SWO9998                    |  | Front/Rear frames & Crossbrace painted in frame colour choice                                                                                           | £296 |
| CERAMIC GLAZE FRAME FINISH |  |                                                                                                                                                         |      |
| UHA4066                    |  | Ceramic glaze finish on frame, for added durability                                                                                                     | £336 |
| FRAME DECAL TYPE*          |  |                                                                                                                                                         |      |
| DHA4190                    |  | Bold Decal (Full coloured transfer on frame) (NC 4300, 4310, 1111)                                                                                      | £31  |
| DHA4195                    |  | Tag (Small Kuschall brand transfer on front of frame) (NC 4200, 4210, 4220, 4230, 4240, 4250, 2222)                                                     | ◆    |
| DHA4197                    |  | Without Tag (Kuschall branding only) (NC 4200, 4210, 4220, 4230, 4240, 4250, 4300, 4310)                                                                | ○    |

| <b>BOLD DECAL COLOUR*</b>               |  |                                                                                                                                     |     |
|-----------------------------------------|--|-------------------------------------------------------------------------------------------------------------------------------------|-----|
| DHA4200                                 |  | Black                                                                                                                               | ○   |
| DHA4210                                 |  | Dark Grey                                                                                                                           | ○   |
| DHA4220                                 |  | Light Grey                                                                                                                          | ○   |
| DHA4230                                 |  | Red                                                                                                                                 | ○   |
| DHA4240                                 |  | Blue                                                                                                                                | ○   |
| DHA4250                                 |  | Gold                                                                                                                                | ○   |
| DHA1111                                 |  | No bold colour, select only with Option DHA4195 Tag or DHA4197 Without Tag (NC 4190)                                                | ○   |
| <b>TAG COLOUR*</b>                      |  |                                                                                                                                     |     |
| DHA4300                                 |  | Black                                                                                                                               | ○   |
| DHA4310                                 |  | Red                                                                                                                                 | ○   |
| DHA2222                                 |  | No tag colour, select only with Option DHA4190 Bold Decal or DHA4197 Without Tag (NC 4195)                                          | ○   |
| <b>FRONT CASTORS*</b>                   |  |                                                                                                                                     |     |
| DHA5000                                 |  | 75mm x 25mm (3" x 1") Supersport (NC 5502, 5504, 5506, 5507, 5508, 5509, 5510)                                                      | ○   |
| DHA5050                                 |  | 100mm x 30mm (4" x 1.2") Everyday Black hub and tyre (NC 5502, 5504, 5506, 5507, 5508, 5509, 5510)                                  | ○   |
| DHA5100                                 |  | 100mm x 30mm (4" x 1.2") Comfort Black hub and tyre (NC 5502, 5504, 5506, 5507, 5508, 5509, 5510)                                   | ○   |
| DHA5150                                 |  | 100mm x 25mm (4" x 1") Starec Black hub and tyre (NC 5502, 5504, 5506, 5507, 5508, 5509, 5510)                                      | £83 |
| DHA5200                                 |  | 100mm x 36mm (4" x 1.4") Aluminium Softroll Black hub and tyre (Frogs Legs) Use this option to access 4" anodised colour options    | £93 |
| DHA5250                                 |  | 100mm x 36mm (4" x 1.4") Aluminium Softroll Silver hub, Black tyre (Frog Legs) (NC 5502, 5504, 5506, 5507, 5508, 5509)              | £93 |
| DHA5300                                 |  | 125mm x 30mm (5" x 1.2") Everyday Black hub and tyre (NC 5502, 5504, 5506, 5507, 5508, 5509, 5510)                                  | ○   |
| DHA5350                                 |  | 125mm x 30mm (5" x 1.2") Comfort Black hub and tyre (NC 5502, 5504, 5506, 5507, 5508, 5509, 5510)                                   | ◆   |
| DHA5400                                 |  | 125mm x 25mm (5" x 1") Starec Black hub and tyre (NC 5502, 5504, 5506, 5507, 5508, 5509, 5510)                                      | £83 |
| DHA5450                                 |  | 125mm x 36mm (5" x 1.4") Aluminium Softroll Black hub and tyre (Frog Legs) Use this option to access 5" anodised colour options     | £93 |
| DHA5500                                 |  | 125mm x 36mm (5" x 1.4") (Frog Legs) Aluminium Softroll Silver hub, Black tyre (Frogs Legs) (NC 5502, 5504, 5506, 5507, 5508, 5509) | £93 |
| <b>ANODISED COLOUR SOFTROLL CASTORS</b> |  |                                                                                                                                     |     |
| UHA5502                                 |  | Pink                                                                                                                                | £50 |
| UHA5504                                 |  | Orange                                                                                                                              | £50 |
| UHA5506                                 |  | Blue                                                                                                                                | £50 |
| UHA5507                                 |  | Purple                                                                                                                              | £50 |
| UHA5508                                 |  | Red                                                                                                                                 | £50 |
| UHA5509                                 |  | Gold                                                                                                                                | £50 |
| <b>CERAMIC CASTOR BEARINGS</b>          |  |                                                                                                                                     |     |
| UHA5510                                 |  | Ceramic castor bearings for Aluminium soft roll castors                                                                             | £58 |

| REAR WHEEL POSITION (CENTRE OF GRAVITY)*  |  |                                                                                                                                                                                                                  |        |
|-------------------------------------------|--|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------|
| DHA6100                                   |  | 1st from front Most Active Position                                                                                                                                                                              | ○      |
| DHA6110                                   |  | 2nd from front                                                                                                                                                                                                   | ○      |
| DHA6120                                   |  | 3rd from front                                                                                                                                                                                                   | ○      |
| DHA6130                                   |  | 4th from front                                                                                                                                                                                                   | ○      |
| DHA6140                                   |  | 5th from front                                                                                                                                                                                                   | ◆      |
| DHA6150                                   |  | 6th from front                                                                                                                                                                                                   | ○      |
| DHA6160                                   |  | 7th from front                                                                                                                                                                                                   | ○      |
| DHA6170                                   |  | 8th from front                                                                                                                                                                                                   | ○      |
| DHA6180                                   |  | 9th from front                                                                                                                                                                                                   | ○      |
| QUAD QUICK RELEASE AXLE                   |  |                                                                                                                                                                                                                  |        |
| DHA6300                                   |  | Quad Quick Release/Tetra Clip (NC 6550, 6551, 6552, 6553, 6554, 0907)                                                                                                                                            | £85    |
| REAR WHEEL SIZE*                          |  |                                                                                                                                                                                                                  |        |
| DHA6410                                   |  | 24"                                                                                                                                                                                                              | ◆      |
| DHA6420                                   |  | 25" (NC 6510, 6600, 6630, 6660, 6700, 6710, 6720, 6760, 6770, 6900)                                                                                                                                              | ○      |
| DHA6430                                   |  | 26" (NC 6520, 6550, 6551, 6552, 6553, 6554, 6600, 6630, 6660, 6720, 6740, 6750, 6760, 6770)                                                                                                                      | ○      |
| REAR WHEEL TYPE*                          |  |                                                                                                                                                                                                                  |        |
| DHA6510                                   |  | 24", 26" Kuschall Starec Wheel, 24" & 26" Black Cross Spokes, Rim & Hub. (NC 6557, 6558, 6559, 6560, 6561, 6562, 6563, 6564, 6565, 6566, 6567, 6568, 6569, 6570, 6571, 6572, 6573, 6574, 6575, 6576, 6577, 6579) | ◆      |
| DHA6520                                   |  | 24", 25" High Performance Wheel Black hub and rim (NC 6557, 6558, 6559, 6560, 6561, 6562, 6563, 6564, 6565, 6566, 6567, 6568, 6569, 6570, 6571, 6572, 6573, 6574, 6575, 6576, 6577, 6579, 6660)                  | £182   |
| DHA6550                                   |  | 24", 25" Spinergy LXK, Black hub, rim and 18 spokes in Black (Unique coloured spokes available) (NC 6557, 6558, 6559, 6560, 6561, 6562, 6563, 6564, 6565, 6566, 6567, 6568, 6569, 6660)                          | £611   |
| UHA6551                                   |  | 24", 25" Spinergy Lite Extreme LX, Weight 687 grams per wheel (Unique coloured rims, hub and spokes available) (NC 6660)                                                                                         | £785   |
| UHA6552                                   |  | 24", 25" Spinergy Lite Extreme X-laced XLX, Weight 821 grams per wheel (Unique coloured rims, hub and spokes available) (NC 6660)                                                                                | £844   |
| UHA6553                                   |  | 24", 25" Spinergy X-Laced Sport Lite Extreme XSLX, Weight 757 grams per wheel (Unique coloured rims, hub and spokes available) (NC 6660)                                                                         | £893   |
| UHA6554                                   |  | 24", 25" Spinergy Carbon Fibre CLX, Weight 530 grams per wheel (Unique coloured hub and spokes available) (NC 6557, 6558, 6559, 6560, 6561, 6562, 6563, 6564, 6565, 6566)                                        | £1,084 |
| SWO0907                                   |  | 24", 25" Spinergy SPOX, Chrome hub and 18 spokes in Black. (Unique coloured spokes available) (NC 6557, 6558, 6559, 6560, 6561, 6562, 6563, 6564, 6565, 6566, 6567, 6568, 6569, 6660)                            | £486   |
| ANODISED COLOUR SPINERGY RIM (DESIGN LAB) |  |                                                                                                                                                                                                                  |        |
| UHA6557                                   |  | Blue                                                                                                                                                                                                             | £179   |
| UHA6558                                   |  | Gold                                                                                                                                                                                                             | £179   |
| UHA6559                                   |  | Green                                                                                                                                                                                                            | £179   |
| UHA6560                                   |  | Orange                                                                                                                                                                                                           | £179   |
| UHA6561                                   |  | Red                                                                                                                                                                                                              | £179   |
| UHA6562                                   |  | Purple                                                                                                                                                                                                           | £179   |
| ANODISED COLOUR SPINERGY HUB (DESIGN LAB) |  |                                                                                                                                                                                                                  |        |
| UHA6563                                   |  | Blue                                                                                                                                                                                                             | £179   |
| UHA6564                                   |  | Gold                                                                                                                                                                                                             | £179   |
| UHA6565                                   |  | Silver                                                                                                                                                                                                           | £179   |
| UHA6566                                   |  | Green                                                                                                                                                                                                            | £179   |
| UHA6567                                   |  | Orange                                                                                                                                                                                                           | £179   |
| UHA6568                                   |  | Red                                                                                                                                                                                                              | £179   |
| UHA6569                                   |  | Purple                                                                                                                                                                                                           | £179   |

| <b>SPINERGY COLOURED SPOKES</b>   |  |                                                                                                                                                                                                                                 |      |
|-----------------------------------|--|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------|
| UHA6570                           |  | Blue                                                                                                                                                                                                                            | £84  |
| UHA6571                           |  | Green                                                                                                                                                                                                                           | £84  |
| UHA6572                           |  | Orange                                                                                                                                                                                                                          | £84  |
| UHA6573                           |  | Pink                                                                                                                                                                                                                            | £84  |
| UHA6574                           |  | Red                                                                                                                                                                                                                             | £84  |
| UHA6575                           |  | Purple                                                                                                                                                                                                                          | £84  |
| UHA6576                           |  | White                                                                                                                                                                                                                           | £84  |
| UHA6577                           |  | Yellow                                                                                                                                                                                                                          | £84  |
| <b>CERAMIC BEARING REAR WHEEL</b> |  |                                                                                                                                                                                                                                 |      |
| UHA6579                           |  | Ceramic Bearing for Spinergy rear wheels                                                                                                                                                                                        | £170 |
| <b>REAR WHEEL TYRE*</b>           |  |                                                                                                                                                                                                                                 |      |
| DHA6600                           |  | 24" x 1" Maxxis Lite guard Pneumatic tyre. Herringbone tread, Silvershield puncture protection                                                                                                                                  | ◆    |
| DHA6625                           |  | 24", 25", 26" x 1" Maxxis Alpha Protect. Lightweight, Pneumatic and highly puncture resistant.                                                                                                                                  | £100 |
| DHA6630                           |  | 24" x 1" MBL Speedlite Super Lightweight (Dual layer very high puncture resistance)                                                                                                                                             | £97  |
| DHA6660                           |  | 24" x 1 3/8" Solid Profile Black Tyre (NC 7230)                                                                                                                                                                                 | £66  |
| <b>HANDRIM*</b>                   |  |                                                                                                                                                                                                                                 |      |
| DHA6700                           |  | 24", 26" Supergripp, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance                                                                                                                                      | £223 |
| DHA6710                           |  | 26" Steel, short tab with 1 hole, mounted at 3 cm distance                                                                                                                                                                      | £62  |
| DHA6720                           |  | 24" Titanium, long tab with 2 holes for 3cm or 4cm mounting, mounted at 3cm distance                                                                                                                                            | £192 |
| DHA6740                           |  | 24", 25" Carbolife Gekko, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3 cm distance. Unique shape for low-effort propulsion, with fully integrated rubber strip for extra grip                                     | £320 |
| DHA6750                           |  | 24", 25" Carbolife Gekko S, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion, with fully integrated rubber strip for extra grip, smaller profile for smaller hands | £320 |
| DHA6760                           |  | 24" Carbolife Curve L, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion                                                                                            | £284 |
| DHA6770                           |  | 24" Carbolife Quadro, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion, with full anti-slip coating over the entire handrim                                        | £356 |
| DHA6780                           |  | 24", 25", 26" Light Aluminium Grey, long tab with 2 holes for 3cm or 4cm mounting, mounted at 3cm distance                                                                                                                      | ◆    |
| DHA6790                           |  | 24", 25", 26" Light Aluminium Black, short tab with 1 hole for 3cm mounting                                                                                                                                                     | £42  |
| <b>REAR WHEEL CAMBER*</b>         |  |                                                                                                                                                                                                                                 |      |
| DHA6820                           |  | 0° Camber                                                                                                                                                                                                                       | ◆    |
| DHA6830                           |  | 1° Camber                                                                                                                                                                                                                       | ○    |
| DHA6840                           |  | 3° Camber (NC 7230)                                                                                                                                                                                                             | ○    |
| <b>SPOKEGUARD</b>                 |  |                                                                                                                                                                                                                                 |      |
| DHA6900                           |  | 24", 26" Transparent Spokeguards                                                                                                                                                                                                | £126 |



| <b>FIXATION CLIP*</b>                      |  |                                                                                                                                                                   |      |
|--------------------------------------------|--|-------------------------------------------------------------------------------------------------------------------------------------------------------------------|------|
| DHA7000                                    |  | Fixation Strap (Hook and loop strap under the seat holds the frame together when folded for easy transportation, designed for clients with limited hand function) | ○    |
| DHA7010                                    |  | Fixation Clip (Plastic clip on either side of the seat holds the frame together when folded for easy transportation)                                              | ◆    |
| DHA3333                                    |  | No fixation clip, select with Tapered frame options 1170, 1180                                                                                                    | ○    |
| <b>UNFOLDING HELP*</b>                     |  |                                                                                                                                                                   |      |
| DHA7100                                    |  | Unfolding Help (Cord to help open the folding mechanism from the front of the chair)                                                                              | ◆    |
| <b>UNIQUE CUSTOMISATION COLOUR PACKAGE</b> |  |                                                                                                                                                                   |      |
| UHA6960                                    |  | UNIQUE KIT 1 - Includes custom colour forks and bearing cover. Select colour from the 10 choices below                                                            | £111 |
| UHA6961                                    |  | UNIQUE KIT 2 - Includes custom colour forks and bearing cover, castor fork supporter, rear wheel supporter, brake clamps. Select colour from the 10 choices below | £285 |
| <b>UNIQUE CUSTOMISATION COLOURS</b>        |  |                                                                                                                                                                   |      |
| UHA6963                                    |  | Red                                                                                                                                                               | ○    |
| UHA6964                                    |  | Purple                                                                                                                                                            | ○    |
| UHA6965                                    |  | Copper                                                                                                                                                            | ○    |
| UHA6966                                    |  | Green                                                                                                                                                             | ○    |
| UHA6967                                    |  | Orange                                                                                                                                                            | ○    |
| UHA6968                                    |  | Gold                                                                                                                                                              | ○    |
| UHA6969                                    |  | Champagne                                                                                                                                                         | ○    |
| UHA6970                                    |  | Pink                                                                                                                                                              | ○    |
| UHA6971                                    |  | Blue                                                                                                                                                              | ○    |
| UHA6972                                    |  | Poly Mirror Finish                                                                                                                                                | ○    |
| <b>BRAKES*</b>                             |  |                                                                                                                                                                   |      |
| DHA7200                                    |  | Pull Brakes                                                                                                                                                       | £78  |
| DHA7210                                    |  | Standard Push to Lock Brakes, Lever Bent                                                                                                                          | ◆    |
| DHA7220                                    |  | Standard Push to Lock Brakes, Lever Straight                                                                                                                      | ○    |
| DHA7230                                    |  | Performance Brakes, Push to Lock. Brake barrel flips away from tyre (NC 7300, 7810, 7820, 7830)                                                                   | £78  |
| DHA7240                                    |  | Active/Scissor Brakes (NC 7300)                                                                                                                                   | £78  |
| DHA7250                                    |  | Light/Compact Active, Aluminium Brakes (NC 7300)                                                                                                                  | £136 |
| <b>BRAKE LEVER EXTENSION</b>               |  |                                                                                                                                                                   |      |
| DHA7300                                    |  | Brake Lever Extension                                                                                                                                             | £49  |
| <b>SAFETY OPTIONS</b>                      |  |                                                                                                                                                                   |      |
| DHA7500                                    |  | Anti-Tipper Right                                                                                                                                                 | £122 |
| DHA7510                                    |  | Anti-Tipper Left                                                                                                                                                  | £122 |
| DHA7520                                    |  | Anti-Tipper Both sides                                                                                                                                            | £244 |
| DHA7800                                    |  | Standard Pelvic belt                                                                                                                                              | £74  |
| <b>BODYPOINT BELTS</b>                     |  |                                                                                                                                                                   |      |
| DHA7810                                    |  | Bodypoint Posture Belt Evoflex Small                                                                                                                              | £81  |
| DHA7820                                    |  | Bodypoint Posture Belt Evoflex Medium                                                                                                                             | £81  |
| DHA7830                                    |  | Bodypoint Posture Belt Evoflex Large                                                                                                                              | £81  |

| <b>OTHER OPTIONS</b> |  |                                                                                                                            |      |
|----------------------|--|----------------------------------------------------------------------------------------------------------------------------|------|
| DHA8010              |  | Front Frame Protectors                                                                                                     | £37  |
| DHA8100              |  | Tipper-Aid, Mounted Right Side Max. 2 adaptions per side (Adaptations = Anti tipper, Tipper aid and Cane holder)           | £46  |
| DHA8200              |  | Cane Holder, Mounted Right Side Max. 2 adaptions per side (Adaptations = Anti tipper, Tipper aid and cane holder)          | £75  |
| DHA8500              |  | Air Pump                                                                                                                   | £13  |
| DHA8550              |  | Tool Kit                                                                                                                   | ◆    |
| DHA8680              |  | 24" Mountain Bike Wheels with bag. Fitted with Specific Black Tyres, Grey Aluminium Handrims; Quick Release Axles Included | £320 |
| DHA8700              |  | Travel bag for Champion, Wheels can remain on the wheelchair                                                               | £254 |

| <b>BODYPOINT ACCESSORIES</b> |  |                        |     |
|------------------------------|--|------------------------|-----|
| DHA8800                      |  | Bodypoint Calf Strap   | £57 |
| DHA8850                      |  | Bodypoint Mobility Bag | £55 |

|                    |  |
|--------------------|--|
| <b>TOTAL PRICE</b> |  |
|--------------------|--|

**Non Conforming Options:**

|         |                                                                                                                                                |
|---------|------------------------------------------------------------------------------------------------------------------------------------------------|
| DHA1320 | 360 mm                                                                                                                                         |
| DHA1325 | 380 mm                                                                                                                                         |
| DHA1350 | 480 mm                                                                                                                                         |
| DHA1946 | 430 mm                                                                                                                                         |
| DHA1948 | 440 mm                                                                                                                                         |
| DHA1950 | 450 mm                                                                                                                                         |
| DHA1952 | 460 mm                                                                                                                                         |
| DHA1954 | 470 mm                                                                                                                                         |
| DHA1956 | 480 mm                                                                                                                                         |
| DHA1958 | 490 mm                                                                                                                                         |
| DHA1960 | 500 mm                                                                                                                                         |
| DHA2200 | Black Standard Light Upholstery, without Tension Adjustment                                                                                    |
| DHA2240 | Without Backrest Upholstery Select if fitting a Matrix/Fixed backrest, no upholstery required                                                  |
| UHA2241 | Neon Yellow                                                                                                                                    |
| UHA2242 | Neon Pink                                                                                                                                      |
| UHA2243 | Blue                                                                                                                                           |
| UHA2244 | White                                                                                                                                          |
| UHA2245 | Red                                                                                                                                            |
| UHA0002 | Your personal colour stitching. Add RAL colour of your choice in special comments box, for a personalised colour stitching on back upholstery. |
| UHA0003 | Your personal design embroidery. Supply your design with your order, subject to copyright. Design size up to 250mm x 250mm                     |
| UHA0004 | Your personal flock design. Supply your design with your order, subject to copyright. Design size up to 280mm x 400mm                          |
| DHA2400 | Without Push Handles                                                                                                                           |
| DHA2430 | Push Handles Height Adjustable, Rear Set                                                                                                       |
| DHA2710 | Flo-Tech Lite Cushion 78mm high- Comfort/Low Risk                                                                                              |
| DHA2720 | Matrix Libra Cushion 83mm high, Black Smooth Cover - Very High Risk                                                                            |
| DHA2748 | Contour Visco NG Cushion 100mm high - High Risk                                                                                                |
| DCA2750 | Vicair Active O2 lightweight cushion 60mm high - High risk                                                                                     |
| DHA2810 | Fluid supplement sac (Pair)                                                                                                                    |
| DHA2820 | Tapered wedge (Pair)                                                                                                                           |
| DHA2830 | Abductor                                                                                                                                       |
| DHA2840 | Lateral side wedge                                                                                                                             |
| DHA2850 | Anterior front wedge                                                                                                                           |
| DHA2860 | Seat rigidizer                                                                                                                                 |
| DHA3480 | Tubular Armrest. Height Adjustable, Removable, Non-locking, Swing-away and Comfort Padded, no Clothesguard included                            |
| DHA3490 | Armrest, Height Adjustable, Short Armpad (T-Arm), Non-locking, no Clothesguard included                                                        |
| UHA4066 | Ceramic glaze finish on frame, for added durability                                                                                            |
| DHA4190 | Bold Decal (Full coloured transfer on frame)                                                                                                   |
| DHA4195 | Tag (Small Kuschall brand transfer on front of frame)                                                                                          |
| DHA4200 | Black                                                                                                                                          |
| DHA4210 | Dark Grey                                                                                                                                      |
| DHA4220 | Light Grey                                                                                                                                     |
| DHA4230 | Red                                                                                                                                            |
| DHA4240 | Blue                                                                                                                                           |

|         |                                                                                                                                  |
|---------|----------------------------------------------------------------------------------------------------------------------------------|
| DHA4250 | Gold                                                                                                                             |
| DHA1111 | No bold colour, select only with Option DHA4195 Tag or DHA4197 Without Tag                                                       |
| DHA4300 | Black                                                                                                                            |
| DHA4310 | Red                                                                                                                              |
| DHA2222 | No tag colour, select only with Option DHA4190 Bold Decal or DHA4197 Without Tag                                                 |
| DHA5000 | 75mm x 25mm (3" x 1") Supersport                                                                                                 |
| DHA5050 | 100mm x 30mm (4" x 1.2") Everyday Black hub and tyre                                                                             |
| DHA5100 | 100mm x 30mm (4" x 1.2") Comfort Black hub and tyre                                                                              |
| DHA5150 | 100mm x 25mm (4" x 1") Starec Black hub and tyre                                                                                 |
| DHA5200 | 100mm x 36mm (4" x 1.4") Aluminium Softroll Black hub and tyre (Frogs Legs) Use this option to access 4" anodised colour options |
| DHA5250 | 100mm x 36mm (4" x 1.4") Aluminium Softroll Silver hub, Black tyre (Frog Legs)                                                   |
| UHA5502 | Pink                                                                                                                             |
| UHA5504 | Orange                                                                                                                           |
| UHA5506 | Blue                                                                                                                             |
| UHA5507 | Purple                                                                                                                           |
| UHA5508 | Red                                                                                                                              |
| UHA5509 | Gold                                                                                                                             |
| UHA5510 | Ceramic castor bearings for Aluminium soft roll castors                                                                          |
| DHA6180 | 9th from front                                                                                                                   |
| DHA6410 | 24"                                                                                                                              |
| DHA6420 | 25"                                                                                                                              |
| DHA6430 | 26"                                                                                                                              |
| DHA6510 | 24", 26" Kuschall Starec Wheel, 24" & 26" Black Cross Spokes, Rim & Hub.                                                         |
| DHA6520 | 24", 25" High Performance Wheel Black hub and rim                                                                                |
| DHA6550 | 24", 25" Spinergy LXX, Black hub, rim and 18 spokes in Black (Unique coloured spokes available)                                  |
| UHA6551 | 24", 25" Spinergy Lite Extreme LX, Weight 687 grams per wheel (Unique coloured rims, hub and spokes available)                   |
| UHA6552 | 24", 25" Spinergy Lite Extreme X-laced XLX, Weight 821 grams per wheel (Unique coloured rims, hub and spokes available)          |
| UHA6553 | 24", 25" Spinergy X-Laced Sport Lite Extreme XSLX, Weight 757 grams per wheel (Unique coloured rims, hub and spokes available)   |
| UHA6554 | 24", 25" Spinergy Carbon Fibre CLX, Weight 530 grams per wheel (Unique coloured hub and spokes available)                        |
| SWO0907 | 24", 25" Spinergy SPOX, Chrome hub and 18 spokes in Black. (Unique coloured spokes available)                                    |
| UHA6557 | Blue                                                                                                                             |
| UHA6558 | Gold                                                                                                                             |
| UHA6559 | Green                                                                                                                            |
| UHA6560 | Orange                                                                                                                           |
| UHA6561 | Red                                                                                                                              |
| UHA6562 | Purple                                                                                                                           |
| UHA6563 | Blue                                                                                                                             |
| UHA6564 | Gold                                                                                                                             |
| UHA6565 | Silver                                                                                                                           |
| UHA6566 | Green                                                                                                                            |
| UHA6567 | Orange                                                                                                                           |
| UHA6568 | Red                                                                                                                              |
| UHA6569 | Purple                                                                                                                           |

|         |                                                                                                                                                                                                                                 |
|---------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| UHA6570 | Blue                                                                                                                                                                                                                            |
| UHA6571 | Green                                                                                                                                                                                                                           |
| UHA6572 | Orange                                                                                                                                                                                                                          |
| UHA6573 | Pink                                                                                                                                                                                                                            |
| UHA6574 | Red                                                                                                                                                                                                                             |
| UHA6575 | Purple                                                                                                                                                                                                                          |
| UHA6576 | White                                                                                                                                                                                                                           |
| UHA6577 | Yellow                                                                                                                                                                                                                          |
| UHA6579 | Ceramic Bearing for Spinergy rear wheels                                                                                                                                                                                        |
| DHA6600 | 24" x 1" Maxxis Lite guard Pneumatic tyre. Herringbone tread, Silvershield puncture protection                                                                                                                                  |
| DHA6630 | 24" x 1" MBL Speedlite Super Lightweight (Dual layer very high puncture resistance)                                                                                                                                             |
| DHA6660 | 24" x 1 3/8" Solid Profile Black Tyre                                                                                                                                                                                           |
| DHA6700 | 24", 26" Supergripp, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance                                                                                                                                      |
| DHA6710 | 26" Steel, short tab with 1 hole, mounted at 3 cm distance                                                                                                                                                                      |
| DHA6720 | 24" Titanium, long tab with 2 holes for 3cm or 4cm mounting, mounted at 3cm distance                                                                                                                                            |
| DHA6740 | 24", 25" Carbolife Gekko, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3 cm distance. Unique shape for low-effort propulsion, with fully integrated rubber strip for extra grip                                     |
| DHA6750 | 24", 25" Carbolife Gekko S, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion, with fully integrated rubber strip for extra grip, smaller profile for smaller hands |
| DHA6760 | 24" Carbolife Curve L, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion                                                                                            |
| DHA6770 | 24" Carbolife Quadro, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion, with full anti-slip coating over the entire handrim                                        |
| DHA6900 | 24", 26" Transparent Spokeguards                                                                                                                                                                                                |
| DHA7000 | Fixation Strap (Hook and loop strap under the seat holds the frame together when folded for easy transportation, designed for clients with limited hand function)                                                               |
| DHA7010 | Fixation Clip (Plastic clip on either side of the seat holds the frame together when folded for easy transportation)                                                                                                            |
| DHA7230 | Performance Brakes, Push to Lock. Brake barrel flips away from tyre                                                                                                                                                             |
| DHA7300 | Brake Lever Extension                                                                                                                                                                                                           |
| DHA7810 | Bodypoint Posture Belt Evoflex Small                                                                                                                                                                                            |
| DHA7820 | Bodypoint Posture Belt Evoflex Medium                                                                                                                                                                                           |
| DHA7830 | Bodypoint Posture Belt Evoflex Large                                                                                                                                                                                            |



For a more complex customer requirement  
please refer to our Invacare Unique special order form.  
email [specialsuk@invacare.com](mailto:specialsuk@invacare.com)

**Special Instructions**

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_

Description of special request (please include relevant option code or part number):  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_