



**Yes, you can:**

Tel: 01656 776222  
 Email: ordersuk@invacare.com  
 Online Spares available at: www.invacare.co.uk

# Power Retail Prescription Form

JANUARY 2025  
 Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## STORM 4 Modulite

### Standard Features

- Max user weight 150Kg (23.5 stone)
- Crash tested ISO 7176 -Part 19
- Modulite highly adjustable seating system
- Telescopic adjustable Seat Widths from 330mm - 530mm (13" - 21") with an additional 50mm, (2") by arms
- Seat Depths 410mm - 510mm Stepless Adjustment (16" - 20")
- Tension Adjustable Backrest
- Standard T-Armrests, detachable with Short armpad
- Manual settable seat Tilt 0° - 20°
- Manual settable Recliner (0° - 30° in 3.75° steps)
- VARI F with Fixed Footplate
- LinX Control REM110 with standard bracket
- Choice of 6 Colours
- Puncture Proof Castors
- Puncture Proof Drive Wheels
- 73.5 Ah Batteries
- 6 km/h (4mph) Motors
- 8A (200 - 250V) Charger
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)
- Anti-Tippers
- Rear wheel suspension

**Standard Base price: £6,491**



Q'STRAIT  
 QLK-150  
 Compatible



Please make your selection with a cross in the box    ○ Free of charge    ◆ standard feature    POA special request

Customer weight =                      KG

Code	Description	
<b>CHAIR TYPE*</b>		
ATB0000	Storm 4 Modulite	£6,491
ATZ0000	Special Storm 4 Modulite (Please provide details in special section)	POA
<b>SEATING*</b>		
ATB2100	Modulite Seating Unit (Flex3 seating unit)	◆
<b>SEATING OPTION*</b>		
ATB2111	Telescopic Seat Frame with metal seat plate (with stepless adjustment on armrests, legrests and backrest)	◆

<b>SEAT WIDTH*</b>			
ATBXXXX		Seat width 330mm (13") seat width achieved on a 380mm seat frame (2121) with inrigged armrests (2183)	○
ATB2121		Set at 380mm (adjustable from 380mm to 430mm (15" - 17") with stepless adjustment on armrests, legrests and backrest) (NC 2703, 2705)	○
ATB2122		Set at 430mm (adjustable from 430mm to 480mm (17" - 19) with stepless adjustment on armrests, legrests and backrest) (NC 2703, 2705)	○
ATB2123		Set at 480mm (adjustable from 480mm to 530mm (19" - 21") with stepless adjustment on armrests, legrests and backrest) (NC 2704, 2706)	○
ATB2124		Set at 530mm (adjustable from 530mm to 580mm (21" - 23") with stepless adjustment on armrests, legrests and backrest) (NC 2704, 2706)	○
<b>SEAT DEPTH*</b>			
ATB2132		410mm (16") (adjustable 410mm - 510mm stepless adjustment) (NC 2701, 2702, 2703, 2704, 2705, 2706)	○
ATB2134		460mm (18") (adjustable 410mm - 510mm stepless adjustment)	○
ATB2136		510mm (20") (adjustable 410mm - 510mm stepless adjustment)	○
<b>SEAT OPTION</b>			
ATBSWO		"H" Frame Interface Kit	£447
<b>SEAT CUSHION*</b>			
ATB2142		Matrx Flo-tech Lite	£112
ATB2143		Matrx Flo-tech Contour Visco (NC 2155)	£112
ATB2147		Matrx Libra 90mm thick	£391
ASB2140		Without seat cushion	○
<b>BACKREST SUPPORT*</b>			
ATB2162		Backposts 540mm (Tall) must be selected for Flex3 Laguna backrest	◆
ATB2163		Backposts (Short) (NC 2173)	○
<b>PUSH HANDLES FOR BACKPOSTS*</b>			
ATB2107		Backposts WITH push handles	◆
ATB2108		Backposts WITHOUT push handles, Push handles will be cut down on the vertical section. (NC 2219, 2227, 2228, 2245)	£36
<b>BACKREST UPHOLSTERY / CUSHION*</b>			
ATB2171		Basic Tension Adjustable Backrest (540mm or 480mm upholstery height) Straps easily accessible and adjustable from the rear, Black fabric only (NC 2155, 1496, 2230, 2247, 2324, 2325)	◆
ATB2179		Flex3 Tension Adjustable (540mm or 480mm) Black multi stretch PU coated fabric only (NC 2155, 2230, 2247, 2324, 2325)	£138
ATB2172		Flex3 adjustable Backrest Plates (without Backrest Cushion) (NC 2128)	£224
ATB2173		Flex3 adjustable Backrest Plates with Laguna Backrest cushion (NC 2128)	£363
ATB2174		Matrx Posture Back, height 400mm (Mesh fabric)	£496
ATB2175		Matrx Posture Back, height 500mm (Mesh fabric)	£496
ATB2176		Matrx Elite E2 Back , height 500 mm (Mesh fabric) (NC 2121)	£496
ATBS177		Matrx Elite TR Back, height 500mm (Mesh fabric) (NC 2121)	£496
ATB2178		Matrx Elite Deep E2 Back, height 500 mm (Mesh fabric)	£550
ATB2170		Without Backrest Upholstery (NC 2151, 2152, 2155, 2230, 2247)	○

UPHOLSTERY TYPE*				
ATB2151		Black fabric	(NC 2170)	◆
ATB2152		Black multi stretch PU coated fabric	(NC 2170)	○
ATB2155		Black Spacer fabric	(NC 2171, 2179, 2143, 2147)	○
ATB3002		No Upholstry Colour required - No backrest upholstery	(NC 2179, 2173, 2174, 2175, 2176, 2178)	○
ARMREST (PAIR)*				
ATB2191		T-Armrest Standard, Range of adjustment from seat plate 245mm to 315mm		◆
ATB2181		T-Armrest Extended height, Range of adjustment from seat plate 295mm - 365mm.		£36
ATB2183		In-rigged Armrest (reduces seat width by 25mm each side)		○
ATB2193		Flip-up Armrest (height and depth adjustable) (height adjustable from 230mm - 360mm)		£58
ATB2192		Following Armrest (height adjustable from 230mm - 300mm)		£165
ATB2194		Following Armrest (height adjustable form 300mm - 360mm)		£165
ARMPAD (PAIR)*				
ATB2195		Short Armpad (325mm)		◆
ATB2196		Long Armpad (375mm)		○
ATB2197		Long Ergonomic Armpad (397mm)	(NC 2230, 2247)	○
ATB2199		Soft Armpad (REA Azalea Wide Pad)	(NC 2179)	£86
HIP SUPPORTS LH (Not available with 410mm seat depth aNd MaNual CeNtre legrests)				
ATB2184		SMALL 140mm x 100mm, for flip up and following armrests only	(NC 2702, 2703, 2704, 2705, 2706)	£226
ATB2186		MEDIUM 230mm x 100mm, for flip up and following armrests only	(NC 2702, 2703, 2704, 2705, 2706)	£226
ATB2188		LARGE 330mm x 100mm, for flip up and following armrests only	(NC 2702, 2703, 2704, 2705, 2706)	£226
HIP SUPPORTS RH (Not available with 410mm seat depth aNd MaNual CeNtre legrests)				
ATB2185		SMALL 140mm x 100mm, for flip up and following armrests only	(NC 2702, 2703, 2704, 2705, 2706)	£226
ATB2187		MEDIUM 230mm x 100mm, for flip up and following armrests only	(NC 2702, 2703, 2704, 2705, 2706)	£226
ATB2189		LARGE 330mm x 100mm, for flip up and following armrests only	(NC 2702, 2703, 2704, 2705, 2706)	£226
MOUNTING HARDWARE FOR HIP SUPPORTS				
ATB2165		Mounting hardware for hip supports (8.1" / 206mm high)		○
HEADREST / NECKREST - Adequate headrest option must be selected*				
ATB1350		neckrest: height and depth adjustable	(NC 2247)	£72
ATB1432		Headrest, height, depth and sideways adjustable	(NC 2247)	£203
ATB1427		Headrest with cheek support: height and depth adjustable	(NC 2247)	£203
ATB2324		Elan Standard Pad 355mm (10") wide with standard cover complete with Elan Headrest Hardware	(NC 2230, 2247)	£294
ATB2325		Elan Occipital Pad 230mm x 180mm (9" x 7") with standard cover complete with Elan Headrest Hardware	(NC 2230, 2247)	£294
ATB1496		non Invacare Headrest but profile mount fitting Flex3 hard back	(NC 2171)	○
ATB1495		Push bar dia. 22mm in conjunction with telescopic seat frame / backposts, no Headrest - no mounting hardware.		○
LATERAL SUPPORTS				
ATB1445		Multi-adjustable pad (Azalea type, multi stretch PU coated fabric)	(NC 2171)	£150

<b>SEAT HEIGHT AND SEAT RISER*</b>			
ATB1023		Fixed Seat Height 450mm (17.75")	○
ATB1025		Fixed Seat Height 480mm (19")	○
ATB1030		Electrical Seat Riser (400mm - 650mm) (16" - 25.5")	£1,088
ATB1040		Electrical Seat Riser (440mm - 690mm) incl. adapter (17.25" - 27")	£1,156
<b>SEAT TILT*</b>			
ATB1002		Manual Seat Tilt by spindle	◆
ATB1014		Electrical Tilt 0° to 25°	£540
<b>BACKREST RECLINER*</b>			
ATB0975		Manual Settable Recliner (0° - 30°, in 3.75° steps) Foldable	◆
ATB0977		Powered Recliner (0° - 30°) Foldable	£453
<b>LEGREST FIXATION POSITION*</b>			
ATB0590		Position for seat cushion up to 70mm thick	◆
<b>FOOTREST HANGER TYPE (PAIR)*</b>			
ATB0503		Vari F : Angle adjustable with tools 0° - 70°, length 290mm - 460mm	○
ATB0533		Vari A : Angle adjustable with hand wheel 0° - 70°, length 290mm - 460mm	£209
ATB0531		ADE Powered : Angle adjustable powered by joystick, 0° - 80° length 290mm - 460mm	£600
ATB0501		Without Hanger (NC 0504, 0513)	○
<b>FOOTPLATE TYPE*</b>			
ATB0512		Fixed Footplate (NC 0533, 0531)	○
ATB0504		Angle Adjustable Footplate	£65
ATB0513		Angle and depth adjustable Footplate (Adapt)	£106
ATBXXXX		Without footplates for option 0501 without hanger (NC 0533, 0531)	○
<b>CONTROLS*</b>			
ATB2201		REM110 Drive only (NC 0977, 0531, 1437, 2209, 2210, 2230, 2247, 2257, 2233, 2240, 2270, 2219, 2245, 1908, 1909)	◆
ATB2203		REM211 Drive & Actuators (NC 2209, 2210, 2230, 2247, 2257, 2233, 2240, 2270, 1908, 1909)	○
ATB2204		REM216 Drive, Actuators & Lights (NC 2210, 2230, 2247, 2257, 2233, 2240, 2270)	£135
ATB2205		REM400, touch-screen Interface, built in Bluetooth HID (NC 2210)	£306
ATB2206		REM500, without joystick for alternative access, built in Bluetooth HID (NC 2228, 2270, 1709, 1710)	£696
<b>CONTROL BRACKETRY*</b>			
ATB1207		Standard Bracket RH (NC 2206)	○
ATB1208		Standard Bracket LH (NC 2206)	○
ATB1707		Swing Away RH Bracket (NC 2206)	£119
ATB1708		Swing Away LH Bracket (NC 2206)	£119
ATB1709		Matrx Resolve Swing in/out bracket RH (NC 2213, 2235, 1404)	£219
ATB1710		Matrx Resolve Swing in/out bracket LH (NC 2214, 2236, 1405)	£219
ATB2213		Centre mount swing away to RH nucleus Bracket (NC 1405)	£447
ATB2214		Centre mount swing away to LH nucleus Bracket (NC 1405)	£447
ATB2227		Bracket RH mounted for Attendant (NC 1404, 1405, 2128)	£57
ATB2228		Bracket LH mounted for Attendant (NC 1404, 1405, 2128)	£57
ATZ2239		Individual mounting option per (special) - Requires additional information / specification	POA
ATB2209		REM500 Swing away RH multi-adjustable Bracket (NC 2203, 2204, 2205, 1404)	○
ATB2210		REM500 Swing away LH multi-adjustable Bracket (NC 2203, 2204, 2205)	○

<b>REM400 TOGGLE OPTION</b>			
ATB2270		Assignable Toggle Switch Module for REM400, (as default set to left On / OFF, right speed dial)	(NC 2203, 2204) £282
<b>CHIN CONTROL</b>			
UNI9250		Mo-vis Micro joystick (Please specify in Special Comments box mounting request if required)	POA
UNI9251		Mo-vis Multi joystick (Please specify in Special Comments box mounting request if required)	POA
<b>DUAL CONTROL</b>			
ATB2220		REM050 Dual Control (Please select user control and bracket) (Please select dual control bracket below)	(NC 1222, 2227, 2228, 2245, 2128, 2230, 2247, 2257, 2233, 2240) £354
ATB2219		Compact Dual Control (DLX-ACU200) (Please select user control and bracket) (Please select dual control bracket below)	(NC 1222, 2227, 2228, 2245, 2128, 2201) £354
<b>DUAL CONTROL BRACKETRY</b>			
ATB1207		Dual Control Standard Bracket RH	(NC 1222, 2227, 2228, 2245) ○
ATB1208		Dual Control Standard Bracket LH	(NC 1222, 2227, 2228, 2245) ○
<b>INVACARE LINX IDC DUAL CONTROL</b>			
ATB2245		IDC Dual control (Please select IDC dual control bracket below)	(NC 2219, 2227, 2228, 2201) £673
<b>INVACARE LINX IDC DUAL CONTROL BRACKETRY</b>			
ATB1221		IDC Dual control Detachable Bracket RH	(NC 2219, 2227, 2228) ○
ATB1222		IDC Dual control Detachable Bracket LH	(NC 2219, 2227, 2228) ○
<b>LINX COMPACT REMOTE</b>			
ATB2231		Proportional Compact Remote DLX-CR400 (Please also select Proportional Bracket)	(NC 2203, 2204) £741
ATB2232		Proportional Compact Remote Low Force DLX-CR400LF (Please select Proportional Bracketry)	(NC 2203, 2204) £990
<b>LINX COMPACT REMOTE BRACKET</b>			
ATB2207		Proportional Compact RH Mount Bracket	(NC 2203, 2204) ○
ATB2208		Proportional Compact LH Mount Bracket	(NC 2203, 2204) ○
ATB2213		Proportional Compact Centre Mount swing away to RH (nucleus)	(NC 2203, 2204) £447
ATB2214		Proportional Compact Centre Mount swing away to LH (nucleus)	(NC 2203, 2204, 1404, 1405) £447
ATZ2215		Individual mounting option per request (Special) Requires additional information / specification	(NC 2203, 2204) POA
<b>LINX HEAD CONTROL</b>			
ATB2230		Head Control with ultrasonic proximity sensors (ATOM , ASL104) with standard wings M4	(NC 2179, 2197, 2201, 2203, 2204, 2220, 2233, 2324, 2325, 1350, 1432, 1427, 2257, 2240) £3,835
ATB2247		Head Control with ultrasonic proximity sensors (PROTON , ASL104P) with PROTON wings	(NC 2179, 2197, 2201, 2203, 2204, 2220, 2233, 2324, 2325, 1350, 1432, 1427, 2257, 2240) £4,491
<b>FOUR-SWITCH PROXIMITY ARRAY REMOTE</b>			
ATZ2257		Four-Switch Proximity Array (ASL106)	(NC 2203, 2204, 2220, 2230, 2247, 2233, 2240, 2128) £5,939
<b>FOUR-SWITCH PROXIMITY ARRAY BRACKETRY</b>			
ATZ2275		Centre Mount, swing-away to RH, nucleus (includes Eclipse Tray ASL619, 360mm x 180mm)	(NC 1405, 1415) £447
ATZ2276		Centre Mount, swing-away to LH, nucleus (includes Eclipse Tray ASL619, 360mm x 180mm)	(NC 1405, 1339) £447
ATZ2215		Individual mounting option per request (Special) Requires additional information / specification	POA
<b>WIRELESS ACCESSORIES FOR ATOM</b>			
ATB2249		Wireless Triple Switch Receiver (ASL557-3) (Sub D plug)	£306
ATB2250		Wireless Mouse Emulator (ASL558)	£733
<b>LINX MICRO EXTREMITY CONTROL</b>			
ATB2233		MEC Proportional Joystick (Micro Extremity Control) (Please select LinX MEC Bracketry)	(NC 2203, 2204, 2220, 2230, 2247, 2257, 2240, 2128) £3,723

<b>LINX MICRO EXTREMITY CONTROL BRACKET</b>				
ATB2235		MEC Tray Mount RH (ASL603, 100mm x 150mm) (Please also select MEC Bracketry)	(NC 1415)	○
ATB2236		MEC Tray Mount LH (ASL603, 100mm x 150mm)	(NC 1339)	○
ATB2213		Centre Mount, swing away to RH, nucleus	(NC 1405, 1415)	£447
ATB2214		Centre Mount, swing away to LH, nucleus	(NC 1405, 1339)	£447
ATB2241		Centre Mount, swing away to RH, nucleus (Inc Tray ASL603, 100m x 150mm)	(NC 1405, 1415)	£447
ATB2242		Centre Mount, swing away to LH, nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 1405, 1339)	£447
ATZ2215		Individual mounting option per request (Special) Requires additional information / specification		POA
<b>LINX HEAVY DUTY REMOTE</b>				
ATZ2240		Proportional Heavy Duty Remote Joystick with mode switch (ASL133) (Please also select LinX Heavy Duty Bracketry)	(NC 2203, 2204, 2220, 2230, 2247, 2257, 2233, 2128)	£2,076
<b>LINX HEAVY DUTY REMOTE BRACKET</b>				
ATZ2235		Tray Mount RH, (ASL603, 100mm x 150mm)	(NC 1415)	○
ATZ2236		Tray Mount LH, (ASL603, 100mm x 150mm)	(NC 1339)	○
ATZ2213		Centre Mount, swing away to RH, nucleus	(NC 1405, 1415)	£447
ATZ2214		Centre Mount, swing away to LH, nucleus	(NC 1405, 1339)	£447
ATZ2241		Centre Mount, swing away to RH, nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 1405, 1415)	£447
ATZ2242		Centre Mount, swing away to LH, nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 1405, 1339)	£447
ATZ2215		Individual mounting option per request (Special) Requires additional information / specification		POA
<b>LINX 10-WAY SWITCH</b>				
ATB2251		10-Way Switch DLX-SW10-A for direct Powered Seating Access (included loom to connect + additional jack for mode switch)		£453
<b>LINX 10-WAY SWITCH BRACKET</b>				
ATB2217		10-Way Switch Bracket RH, Mounted for User	(NC 2213, 2275, 2235, 2241)	○
ATB2218		10-Way Switch Bracket LH, Mounted for User	(NC 2214, 2276, 2236, 2242)	○
ATB2227		10-Way Switch Bracket RH, Mounted for Attendant		£57
ATB2228		10-Way Switch Bracket LH, Mounted for Attendant		£57
<b>LINX ACCESSORY</b>				
ATB2252		Loom providing jack sockets for Drive and Mode (includes connection for 10-way switch)		£90
<b>LINX USB CHARGER PORT</b>				
ATB2253		USB Charger Port DLX-USB02-A	(NC 2702, 2703, 2704, 2705, 2706)	£143
<b>LINX USB CHARGER PORT BRACKET</b>				
ATB2207		USB Charger Port RH Mount		○
ATB2208		USB Charger Port LH Mount		○
<b>FURTHER LINX COMPONENTS / ACCESSORIES</b>				
ATB2300		To be enclosed only		○
ATB2301		To be connected and programmed (special)		POA
<b>BUDDY BUTTONS</b>				
ATB1685		Red Buddy Button		£79
ATB1686		Yellow Buddy Button		£79
ATB1687		Green Buddy Button		£79
ATB1690		Blue Buddy Button		£79
<b>EGG SWITCH</b>				
ATB2254		Egg Switch, black		£79

PIKO BUTTONS			
ATB1711	<input type="checkbox"/>	Red Piko Button	£142
ATB1712	<input type="checkbox"/>	Yellow Piko Button	£142
ATB1713	<input type="checkbox"/>	Green Piko Button	£142
ATB1714	<input type="checkbox"/>	Orange Piko Button	£142
ATB1715	<input type="checkbox"/>	Blue Piko Button	£142
ATB1716	<input type="checkbox"/>	Black Piko Button	£142
TOGGLES / BUTTONS			
ATZ2261	<input type="checkbox"/>	Stereo Toggle - direct control for Tilt only	£240
ATZ2262	<input type="checkbox"/>	Stereo Button - direct control for Tilt only	£231
ATZ2263	<input type="checkbox"/>	Four-Way Toggle - direct control for Tilt & Recline	£301
ATZ2264	<input type="checkbox"/>	Four-Way Button - direct control for Tilt & Recline	£321
ATZ2207	<input type="checkbox"/>	Toggle / Button RH Mount	○
ATZ2208	<input type="checkbox"/>	Toggle / Button LH Mount	○
ATZ2225	<input type="checkbox"/>	Without Mounting	○
BODYPOINT JOYSTICK KNOBS			
ATB2721	<input type="checkbox"/>	3" U-Shaped, Joystick Handle, Bodypoint PC101A	£64
ATB2722	<input type="checkbox"/>	4" U-Shaped, Joystick Handle, Bodypoint PC102A	£64
ATB2723	<input type="checkbox"/>	Domed Rubber, Joystick Handle, Bodypoint PC107A	£64
ATB2724	<input type="checkbox"/>	Mushroom Grooved, Joystick Handle, Bodypoint PC110A	£64
REMOTE STOP SWITCH			
ATZ2258	<input type="checkbox"/>	Remote Stop Switch (ASL504)	£941
SHROUD COLOUR*			
ATB0202	<input type="checkbox"/>	Blue	○
ATB0228	<input type="checkbox"/>	Black	○
ATB0204	<input type="checkbox"/>	Transparent Cherry Red	○
ATZ0299	<input type="checkbox"/>	Special shroud colour (RAL code required)	POA
FRONT WHEEL*			
ATB0851	<input type="checkbox"/>	3.00 - 6 Grey Pneumatic	○
ATB0853	<input type="checkbox"/>	3.00 - 6 Grey Puncture Proof (sealed)	○
ATB0855	<input type="checkbox"/>	3.00 - 6 Grey Puncture Proof (solid inlay)	£97
REAR WHEEL*			
ATB0771	<input type="checkbox"/>	3.00 - 8 Grey Pneumatic	○
ATB0773	<input type="checkbox"/>	3.00 - 8 Grey Puncture Proof (sealed)	○
ATB0750	<input type="checkbox"/>	3.00 - 8 Grey Puncture Proof (solid inlay)	£97
MOTOR TYPE*			
ATB1277	<input type="checkbox"/>	6 km/h (4mph)	○
ATB1908	<input type="checkbox"/>	10 km/h (6mph) (includes LED lights system with indicators)	£464
ATB1909	<input type="checkbox"/>	12 km/h (8mph) (includes LED light system with indicators)	£464
BATTERY*			
ATB1364	<input type="checkbox"/>	70 Ah Gel	◆

<b>SUSPENSION*</b>			
ATB1240		Suspension low (up to 80kg user weight)	○
ATB1241		Suspension middle (from 70kg up to 100kg user weight)	○
ATB1242		Suspension hard (from 90kg up to 150kg user weight)	○
<b>G TRAC</b>			
ATB1215		G-Trac (ensures the powerchair takes the intended direction as indicated by the user control, without the need for continual correction by the client on uneven or sloping terrain)	£738
<b>CHARGER*</b>			
ATB1343		OptiCharge Charger 8A (200 - 250V)	◆
ATB1313		OptiCharge Charger 10A (200 - 250V) Charging time will be 20% faster	£135
<b>LIGHTS*</b>			
ATB1437		LED Light system with indicators (30 LUX) (NC 1909)	£429
ATB1254		Without Light system	◆
<b>ANTI TIPPER*</b>			
ATB1205		Anti-tipper (standard version)	○
<b>KERB RIDER</b>			
ATB1244		Central Kerbrider (with parking position)	£190
<b>TRANSPORTATION KIT*</b>			
ATB1450		Preparation for DAHL-Docking station, (Adequate headrest is compulsory, see options below)	£366
ATB1494		4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory, see options below)	◆
<b>FIXATION BELT*</b>			
ATB1329		Lap Belt with Auto Lock (NC 2702, 2703, 2704, 2705, 2706)	○
ATB2113		Lap belt with auto lock (flexible mount in side rails) (NC 2702, 2703, 2704, 2705, 2706)	◆
ATB1332		H-Belt	£68
ATB2701		Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size L (padded length 620mm) (NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113)	£85
ATB2702		Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size M (padded length 460mm) (NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113)	£85
ATB2703		Bodypoint 2-Point Hip Belt, Centre pull - push button, Cinch mount size L (padded length 620mm) (NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113)	£78
ATB2704		Bodypoint 2-Point Hip Belt, Centre pull - push button, Cinch mount size M (padded length 460mm) (NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113)	£78
ATB2705		Bodypoint 2-Point Hip Belt, Dual pull Side release button, Cinch mount size L (padded length 620mm) (NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113)	£78
ATB2706		Bodypoint 2-Point Hip Belt, Dual pull Side release button, Cinch mount size M (padded length 460mm) (NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113)	£78
<b>REAR VIEW MIRROR</b>			
ATB1415		Rear view Mirror RH - Must be selected with 10 km/h Motor (NC 2213, 2275, 2235, 2241)	£38
ATB1339		Rear view Mirror LH - Must be selected with 10 km/h Motor (NC 2214, 2276, 2236, 2242)	£38
<b>BAG HOOKS</b>			
ATB2128		Bag Hooks (Max load up to 6kg = impact on max user weight) (NC 2172, 2173, 2227, 2228, 2220, 2219, 2245, 2257, 2233, 2240)	£69
<b>CANE HOLDER</b>			
ATB1312		Cane holder RH (In combination with Attendant bracket or Dual control the cane holder has to be mounted on the opposite side of the remote)	£38
ATB1363		Cane holder LH (In combination with Attendant bracket or Dual control the cane holder has to be mounted on the opposite side of the remote)	£38



TRAY				
ATB1404		Tray Swing Away RH	(NC 1707, 2213, 2214, 2209, 2235, 2236, 2241, 2242, 2207, 2275, 2276)	£314
ATB1405		Tray Swing Away LH	(NC 1708, 2213, 2214, 2210, 2235, 2236, 2241, 2242, 2208, 2275, 2276)	£314
FENDER				
ATB1416		Fender for front Wheels		£108
ATB1417		Fender for rear Wheels		£108

<b>TOTAL PRICE</b>				
--------------------	--	--	--	--

**Non Conforming Options:**

ATB2121	Set at 380mm (adjustable from 380mm to 430mm (15" - 17") with stepless adjustment on armrests, legrests and backrest)
ATB2122	Set at 430mm (adjustable from 430mm to 480mm (17" - 19) with stepless adjustment on armrests, legrests and backrest)
ATB2132	410mm (16") (adjustable 410mm - 510mm stepless adjustment)
ATB2143	Matrx Flo-tech Contour Visco
ATB2147	Matrx Libra 90mm thick
ATB2171	Basic Tension Adjustable Backrest (540mm or 480mm upholstery height) Straps easily accessible and adjustable from the rear, Black fabric only
ATB2179	Flex3 Tension Adjustable (540mm or 480mm) Black multi stretch PU coated fabric only
ATB2172	Flex3 adjustable Backrest Plates (without Backrest Cushion)
ATB2173	Flex3 adjustable Backrest Plates with Laguna Backrest cushion
ATB2174	Matrx Posture Back, height 400mm (Mesh fabric)
ATB2175	Matrx Posture Back, height 500mm (Mesh fabric)
ATB2176	Matrx Elite E2 Back , height 500 mm (Mesh fabric)
ATB2178	Matrx Elite Deep E2 Back, height 500 mm (Mesh fabric)
ATB2170	Without Backrest Upholstery
ATB2151	Black fabric
ATB2152	Black multi stretch PU coated fabric
ATB2155	Black Spacer fabric
ATB2197	Long Ergonomic Armpad (397mm)
ATB2184	SMALL 140mm x 100mm, for flip up and following armrests only
ATB2186	MEDIUM 230mm x 100mm, for flip up and following armrests only
ATB2188	LARGE 330mm x 100mm, for flip up and following armrests only
ATB2185	SMALL 140mm x 100mm, for flip up and following armrests only
ATB2187	MEDIUM 230mm x 100mm, for flip up and following armrests only
ATB2189	LARGE 330mm x 100mm, for flip up and following armrests only
ATB1350	neckrest: height and depth adjustable
ATB1432	Headrest, height, depth and sideways adjustable
ATB1427	Headrest with cheek support: height and depth adjustable
ATB2324	Elan Standard Pad 355mm (10") wide with standard cover complete with Elan Headrest Hardware
ATB2325	Elan Occipital Pad 230mm x 180mm (9" x 7") with standard cover complete with Elan Headrest Hardware
ATB1496	non Invacare Headrest but profile mount fitting Flex3 hard back
ATB0977	Powered Recliner (0° - 30°) Foldable
ATB0533	Vari A : Angle adjustable with hand wheel 0° - 70°, length 290mm - 460mm
ATB0531	ADE Powered : Angle adjustable powered by joystick, 0° - 80° length 290mm - 460mm
ATB0504	Angle Adjustable Footplate
ATB0513	Angle and depth adjustable Footplate (Adapt)
ATB2201	REM110 Drive only
ATB2203	REM211 Drive & Actuators
ATB2204	REM216 Drive, Actuators & Lights
ATB2205	REM400, touch-screen Interface, built in Bluetooth HID
ATB2206	REM500, without joystick for alternative access, built in Bluetooth HID
ATB1707	Swing Away RH Bracket
ATB1708	Swing Away LH Bracket
ATB1709	Matrx Resolve Swing in/out bracket RH

ATB1710	Matrx Resolve Swing in/out bracket LH
ATB2213	Centre mount swing away to RH nucleus Bracket
ATB2214	Centre mount swing away to LH nucleus Bracket
ATB2227	Bracket RH mounted for Attendant
ATB2228	Bracket LH mounted for Attendant
ATB2209	REM500 Swing away RH multi-adjustable Bracket
ATB2210	REM500 Swing away LH multi-adjustable Bracket
ATB2270	Assignable Toggle Switch Module for REM400, (as default set to left On / OFF, right speed dial)
ATB2220	REM050 Dual Control (Please select user control and bracket) (Please select dual control bracket below)
ATB2219	Compact Dual Control (DLX-ACU200) (Please select user control and bracket) (Please select dual control bracket below)
ATB2245	IDC Dual control (Please select IDC dual control bracket below)
ATB1222	IDC Dual control Detachable Bracket LH
ATB2207	Proportional Compact RH Mount Bracket
ATB2208	Proportional Compact LH Mount Bracket
ATB2213	Proportional Compact Centre Mount swing away to RH (nucleus)
ATB2214	Proportional Compact Centre Mount swing away to LH (nucleus)
ATB2230	Head Control with ultrasonic proximity sensors (ATOM , ASL104) with standard wings M4
ATB2247	Head Control with ultrasonic proximity sensors (PROTON , ASL104P) with PROTON wings
ATZ2257	Four-Switch Proximity Array (ASL106)
ATZ2275	Centre Mount, swing-away to RH, nucleus (includes Eclipse Tray ASL619, 360mm x 180mm)
ATZ2276	Centre Mount, swing-away to LH, nucleus (includes Eclipse Tray ASL619, 360mm x 180mm)
ATB2233	MEC Proportional Joystick (Micro Extremity Control) (Please select LinX MEC Bracketry)
ATB2235	MEC Tray Mount RH (ASL603, 100mm x 150mm) (Please also select MEC Bracketry)
ATB2236	MEC Tray Mount LH (ASL603, 100mm x 150mm)
ATB2213	Centre Mount, swing away to RH, nucleus
ATB2214	Centre Mount, swing away to LH, nucleus
ATB2241	Centre Mount, swing away to RH, nucleus (Inc Tray ASL603, 100m x 150mm)
ATB2242	Centre Mount, swing away to LH, nucleus (Inc Tray ASL603, 100mm x 150mm)
ATZ2240	Proportional Heavy Duty Remote Joystick with mode switch (ASL133) (Please also select LinX Heavy Duty Bracketry)
ATZ2235	Tray Mount RH, (ASL603, 100mm x 150mm)
ATZ2236	Tray Mount LH, (ASL603, 100mm x 150mm)
ATZ2213	Centre Mount, swing away to RH, nucleus
ATZ2214	Centre Mount, swing away to LH, nucleus
ATZ2241	Centre Mount, swing away to RH, nucleus (Inc Tray ASL603, 100mm x 150mm)
ATZ2242	Centre Mount, swing away to LH, nucleus (Inc Tray ASL603, 100mm x 150mm)
ATB2227	10-Way Switch Bracket RH, Mounted for Attendant
ATB2228	10-Way Switch Bracket LH, Mounted for Attendant
ATB2253	USB Charger Port DLX-USB02-A
ATB2207	USB Charger Port RH Mount
ATB2208	USB Charger Port LH Mount
ATZ2207	Toggle / Button RH Mount
ATZ2208	Toggle / Button LH Mount
ATB1908	10 km/h (6mph) (includes LED lights system with indicators)
ATB1909	12 km/h (8mph) (includes LED light system with indicators)

ATB1437	LED Light system with indicators (30 LUX)
ATB1329	Lap Belt with Auto Lock
ATB2113	Lap belt with auto lock (flexible mount in side rails)
ATB2701	Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size L (padded length 620mm)
ATB2702	Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size M (padded length 460mm)
ATB2703	Bodypoint 2-Point Hip Belt, Centre pull - push button, Cinch mount size L (padded length 620mm)
ATB2704	Bodypoint 2-Point Hip Belt, Centre pull - push button, Cinch mount size M (padded length 460mm)
ATB2705	Bodypoint 2-Point Hip Belt, Dual pull Side release button, Cinch mount size L (padded length 620mm)
ATB2706	Bodypoint 2-Point Hip Belt, Dual pull Side release button, Cinch mount size M (padded length 460mm)
ATB1415	Rear view Mirror RH - Must be selected with 10 km/h Motor
ATB1339	Rear view Mirror LH - Must be selected with 10 km/h Motor
ATB2128	Bag Hooks (Max load up to 6kg = impact on max user weight)
ATB1404	Tray Swing Away RH
ATB1405	Tray Swing Away LH



For a more complex customer requirement  
please refer to our Invacare Unique special order form.  
email [specialsuk@invacare.com](mailto:specialsuk@invacare.com)

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):