



**Yes, you can:**

Tel: 01656 776222 Fax: 01656 776220  
Email: ordersuk@invacare.com  
Online Spares available at: www.invacare.co.uk

SEPTEMBER 2024  
Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

**TDX SP2 HD**

**Standard Features**

- Max user weight up to 180kg (28 stone)
- Crash tested ISO 7176 -Part 19
- Compact, manoeuvrable, CWD TDX SP2 Heavy Duty
- Telescopic adjustable Seat Widths from 330mm - 530mm (13" - 21") with an additional 50mm, (2") by arms
- Seat Depths 410mm - 560mm Stepless Adjustment (16" - 22")
- Tension Adjustable backrest
- Standard T-Armrests, detachable with Short armpad
- Manual settable seat Tilt 0° - 20°
- Manual settable backrest Recline 0° - 30°, in 3.75° steps
- Standard 80° hangers with fixed footplate
- LiNX Control REM110 with standard bracket
- Selection of 10 Stunning Shroud and Rim colours
- Centre (Driving Wheels) 355mm x 76mm (14" x 3") Black Tyres
- Suspension and Sure Step Technology
- 6 km/h (4mph) 4 Pole Motors
- 8A (200 - 250V) Charger
- 73.5 Ah Batteries
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)

**Standard Base price: £5,382**



Please make your selection with a cross in the box     Free of charge     standard feature    POA special request

Customer weight =                      KG

Code	Description	
<b>CHAIR TYPE*</b>		
ARJ0000	TDX SP2 Heavy Duty	£5,382
ARZ0001	Special TDX SP2 Heavy Duty (Please provide details in special section)	POA

User Weight*			
ARJ1918		User weight less 100 kg	○
ARJ1919		User weight 101 to 125 kg	○
ARJ1921		User weight 126 to 136 kg	○
ARJ1922		User weight 137 to 150 kg (NC 1092, 2052, 2053)	○
ARJ1923		User weight 151 to 160 kg (NC 2192, 2193, 2194, 1092, 2052, 2053)	○
ARJ1924		User weight 161 to 180 kg (NC 2192, 2193, 2194, 1092, 1093, 1094, 1909, 2052, 2053)	○
SEATING TYPE*			
ARJ2100		Modulite Seating unit	◆
SEAT FRAME AND SEAT PLATE*			
ARJ2111		Telescopic Seat Frame with metal seat plate (with stepless adjustment on armrests, legrests and backrest)	◆
SEAT WIDTH*			
ASJXXXX		Seat width 330mm (13") seat width achieved on a 380mm seat frame (2121) with inrigged armrests (2183)	○
ARJ2121		Set at 380mm (15") (adjustable from 380mm to 430mm (15" - 17") with stepless adjustment on armrests, legrests and backrest) (NC 2143, 2701, 2703, 2705)	○
ARJ2122		Set at 430mm (17") (adjustable from 430mm to 480mm (17" - 19") with stepless adjustment on armrests, legrests and backrest) (NC 2703, 2705)	○
ARJ2123		Set at 480mm (19") (adjustable from 480mm to 530mm (19" - 21") with stepless adjustment on armrests, legrests and backrest) (NC 2702, 2704, 2706)	○
ARJ2124		Set at 530mm (21") (adjustable from 530mm to 580mm (21" - 23") with stepless adjustment on armrests, legrests and backrest) (NC 2174, 2175, 2176, 2177, 2178, 2143, 2702, 2704, 2706)	○
Body Fit Kit			
UNI2106		Body fit Kit Outrigging Backrest Brackets and Pivot Point (outrig: 1 inch/2,54 cm per side) (NC 2124, 2192, 2194)	○
SEAT DEPTH RAIL LENGTH*			
ARJ2129		Standard Depth Package (16" - 20", 410mm - 510mm simple stepless adjustment from the side) (NC 2138)	○
ARJ2130		Extended Depth Package (20" - 22", 510mm - 560mm simple stepless adjustment from the side) (NC 2132, 2134)	○
SEAT DEPTH SETTING*			
ARJ2132		410mm (16") (410mm - 510mm simple stepless adjustment from the side) (NC 2196, 2197, 2701, 2702, 2703, 2704, 2705, 2706)	○
ARJ2134		460mm (18") (410mm - 510mm simple stepless adjustment from the side)	○
ARJ2136		510mm (20") (410mm - 510mm simple stepless adjustment from the side)	○
ARJ2138		560mm (22") (510mm - 560mm simple stepless adjustment from the side) (Cushion (option 2142) is only available with Seat Width (option 2124) 530mm) (NC 2143)	○
SEAT OPTION			
ARJSWO		"H" Frame Interface Kit	£446
SEAT CUSHION*			
ARJ2142		Matrx Flo-tech Lite	£112
ARJ2143		Matrx Flo-tech Contour Visco (17x17 only) (NC 2136, 2121, 2124, 2138)	£112
ARJ2147		Matrx Libra 90mm thick	£391
ARJ2140		Without seat cushion	○
BACKREST SUPPORT*			
ARJ2162		Backposts 540mm (Tall) must be selected for Flex3 adjustable backrest plates	◆
ARJ2163		Backposts 480mm (Short) (NC 2173)	○

<b>PUSH HANDLES FOR BACKPOSTS*</b>			
ARJ2107		Backposts WITH push handles	◆
ARJ2108		Backposts WITHOUT push handles, Push handles will be cut down on the vertical section.	(NC 1332, 2220, 2219, 2245, 2227, 2228) £36
<b>BACKREST UPHOLSTERY / CUSHION*</b>			
ARJ2171		Tension Adjustable Backrest (540mm or 480mm upholstery height) Straps easily accessible and adjustable from the rear (Black fabric only)	(NC 1445, 1496, 2155, 2230, 2247) ◆
ARJ2179		Flex3 Tension Adjustable (540mm or 480mm) Black Dartex only	(NC 2151, 2192, 2193, 2194, 2199, 2155, 2230, 2247, 2128) £139
ARJ2172		Flex3 adjustable Backrest Plates (without Backrest Cushion)	(NC 2163, 2128) £224
ARJ2173		Flex3 adjustable Backrest Plates with Laguna Backrest cushion	(NC 2163, 2128) £363
ARJ2174		Matrx Posture Back, height 400mm (Mesh fabric)	(NC 1496, 2124) £496
ARJ2175		Matrx Posture Back, height 500mm (Mesh fabric)	(NC 1496, 2124) £496
ARJ2176		Matrx Elite E2 Back , height 500 mm (Mesh fabric)	(NC 1496, 2124) £496
ARJ2177		Matrx Elite TR Back, height 500mm (Mesh fabric)	(NC 1496, 2124) £496
ARJ2178		Matrx Elite Deep E2 Back, height 500 mm (Mesh fabric)	(NC 1496, 2124) £550
ARJ2170		Without Backrest Upholstery	(NC 1496, 2230, 2247) ○
<b>UPHOLSTERY COLOUR*</b>			
ARJ2151		Black fabric	◆
ARJ2152		Black multi stretch PU coated fabric	(NC 2171) ○
ARJ2155		Black Spacer fabric	(NC 2171, 2179) ○
ARZ1002		Without backrest upholstery (for option 2170)	(NC 2179, 2172, 2173, 2174, 2175, 2176, 2177, 2178) ○
<b>ARMREST*</b>			
ARJ2191		T-Armrest Standard, Range of adjustment from seat plate 245mm to 315mm	◆
ARJ2181		T-Armrest Extended height, Range of adjustment from seat plate 295mm - 365mm.	£36
ARJ2183		In-rigged Armrests (reduces seat width by 25mm each side)	(NC 2185, 2186, 2187, 2188, 2189) £25
ARJ2193		Flip-up Armrest (height and depth adjustable) (height adjustable from 230mm - 360mm) 150 kg max user weight	(NC 2176, 2199) £58
ARJ2192		Following Armrest (height adjustable from 230mm - 300mm) 150 kg max user weight	(NC 1445, 2179) £166
ARJ2194		Following Armrest (height adjustable from 300mm - 360mm) 150 kg max user weight	(NC 1445, 2179) £166
<b>ARMPAD*</b>			
ARJ2195		Short Armpad (325mm)	◆
ARJ2196		Long Armpad (375mm)	○
ARJ2197		Long Ergonomic Armpad (397mm)	(NC 2247) ○
ARJ2199		Soft Armpad (REA Azalea Wide Pad)	£86
<b>HIP SUPPORTS LH (Not available with 410mm seat depth and Manual Centre legrests)</b>			
ARJ2184		SMALL 140mm x 100mm, for flip up and following armrests only	(NC 2702, 2703, 2704, 2705, 2706) £226
ARJ2186		MEDIUM 230mm x 100mm, for flip up and following armrests only	(NC 2702, 2703, 2704, 2705, 2706) £226
ARJ2188		LARGE 330mm x 100mm, for flip up and following armrests only	(NC 2702, 2703, 2704, 2705, 2706) £226
<b>HIP SUPPORTS RH (Not available with 410mm seat depth and Manual Centre legrests)</b>			
ARJ2185		SMALL 140mm x 100mm, for flip up and following armrests only	(NC 2702, 2703, 2704, 2705, 2706) £226
ARJ2187		MEDIUM 230mm x 100mm, for flip up and following armrests only	(NC 2702, 2703, 2704, 2705, 2706) £226
ARJ2189		LARGE 330mm x 100mm, for flip up and following armrests only	(NC 2702, 2703, 2704, 2705, 2706) £226

<b>MOUNTING HARDWARE FOR HIP SUPPORTS</b>			
ARJ2165		Mounting hardware for hip supports (8.1" / 206mm high)	○
<b>HEADREST / NECKREST - Adequate headrest option must be selected*</b>			
ARJ1350		Neckrest height and depth adjustable (NC 2247)	£72
ARJ1432		Headrest, height, depth and sideways adjustable (NC 2247)	◆
ARJ1427		Headrest with cheek support, height and depth adjustable (NC 2247)	£203
ARJ2324		Elan Standard Pad 355mm (10") wide with standard cover complete with Elan Headrest Hardware (NC 2247)	£294
ARJ2325		Elan Occipital Pad 230mm x 180mm (9" x 7") with standard cover complete with Elan Headrest Hardware (NC 2247)	£294
ARJ2327		Elan 4 Point Headrest Pad 280mm x 255mm (11" x 10") with standard cover complete with Elan Headrest Hardware (NC 2247)	£396
ARJ1496		All bracketry as option 1432 without headrest cushion	○
ARJ1495		Push bar dia. 22mm in conjunction with telescopic seat frame/backposts, No Headrest - no mounting hardware.	○
<b>LATERAL SUPPORT</b>			
ARJ1445		Multi-adjustable pad (Azalea type, Dartex) (NC 2171)	£149
<b>SEAT LIFTER / SEAT TILT*</b>			
ARJ1091		Manually settable Seat Tilt (0° - 20°)	◆
ARJ1092		Powered "fixed pivot" Tilt 20° (Max User weight 136kg)	£362
ARJ1095		Powered "fixed pivot" Tilt 20°, heavy duty (Max User weight 180kg)	£362
ARJ1093		Electric Tilt (0° - 30°) with CoG shift Seat Height 460mm (160KG max user weight)	£961
ARJ2052		Electric Tilt (0°-50°) with COG shift Seat height 460mm	£1,322
ARJ1094		Seat lifter with Electrical Tilt (300mm, 30°) with CoG shift Seat Height 460mm - 490mm (160KG max user weight)	£1,696
ARJ2053		Seat lifter with Electrical Tilt (300mm, 45°) with CoG shift Seat Height 460mm	£2,359
<b>SEAT HEIGHT SETTING*</b>			
ARJ0354		Seat height set at 450mm (17.75") (from floor to front of seat plate at 0° tilt) (NC 1094)	○
ARZ1003		Seat height set at 460mm (18") (from floor to front of seat plate at 0° tilt) (NC 1092, 1095)	○
ARJ0364		Seat height set at 475mm (18.75") (from floor to front of seat plate at 0° tilt) (NC 1093, 1094)	○
<b>SEAT TILT OPTION</b>			
ARJ2054		5° Fixed Anterior Tilt Assist mount (will reduce max tilt by 5°) (NC 1092, 1095)	£361
<b>BACKREST RECLINER*</b>			
ARJ0975		Manual Settable Recliner (0° - 30°, in 3.75° steps) Foldable	◆
ARJ0977		Powered Recliner (0° - 30°) Foldable	£453
<b>CENTRE MOUNT LEGREST*</b>			
ARJ0556		Centre Mounted: Manual : Hanger, centre mounted, angle settable (0° - 90°) : knee to heel length 390 - 510mm, (including calf pad and angle adjustable footplate) (NC 0504, 0513)	£334
ARJ0589		Centre Mounted : LNX Powered: Centre mounted legrest, length compensating (including calf pad and angle adjustable footplate) (NC 0504, 0513)	£1,470
UNI0592		LNX Power centre mount legrest with telescoping foot platform (lowers to the ground) inc calf pads: Length 380-419mm (Selectable length depends on STF height) Custom-made, requires customised setting (NC 0512, 0504)	£2,417
ARX9260		Without Centre mount Legrest	○

<b>FOOTREST HANGER TYPE*</b>			
ARJ0526		Side Mounted : Standard 80° hangers	(NC 0559, 0556, 0589) ◆
ARJ0503		Side Mounted : Vari F : Angle adjustable with tools 0° - 70°, length 290mm - 460mm	(NC 0556, 0589) £107
ARJ0533		Side Mounted : Vari A : Angle adjustable with hand wheel 0° - 70°, length 290mm - 460mm	(NC 0559, 0556, 0589) £210
ARJ0531		Side Mounted : ADE Powered : Angle adjustable powered by joystick, 0° - 80° length 290mm - 460mm	(NC 0513, 0559, 0556, 0589) £600
ARJ0501		Without Hanger	(NC 0504, 0513) ○
<b>AMPUTEE SUPPORT</b>			
ARJ0554		Amputee Support RH (this replaces the RH hanger)	£174
ARJ0555		Amputee Support LH (this replaces the LH hanger)	£174
<b>FOOTPLATE TYPE*</b>			
ARJ0512		Fixed Footplate	(NC 0533, 0556, 0589, 0501) ○
ARJ0504		Angle Adjustable Footplate	(NC 0589, 0501) £66
ARJ0513		Angle and depth adjustable Footplate (Adapt)	(NC 0589, 0501) £106
ARJXXXX		Without footplates for option 0501 without hanger	(NC 0503, 0533, 0531) ○
<b>CENTRE MOUNT FOOTPLATES / FOOT PLATFORM*</b>			
ARJ2365		Angle Adjustable Footplates Small 127mm x 190mm (5" x 7.5") for Centre Mounted Hanger only	(NC 0503, 0533, 0531) £93
ARJ0559		Angle Adjustable Footplate Medium 140mm x 240mm (5.5" x 9.5") for Centre Mounted Hanger only	(NC 0503, 0533, 0531) ○
ARJ0516		Angle Adjustable Footplates Large 150mm x 290mm (6" x 11.5") for Centre Mounted Hanger only	(NC 0503, 0533, 0531) £121
ARJ0517		Angle Adjustable foot platform Medium 305mm x 275mm (12" x 10.75") for Centre Mounted Hanger only	(NC 0526, 0503, 0533, 0531) £243
UNI0559		Angle & width adjustable footplate (240mm depth)	(NC 0526, 0503, 0533, 0531, 0556, 0589) ○
UNI0516		Angle & width adjustable footplates (290mm depth)	(NC 0526, 0503, 0533, 0531, 0503, 0556, 0589) £121
UNI0517		Angle adjustable foot platform (240mm depth, 305mm width)	(NC 0526, 0503, 0533, 0531, 0504, 0556, 0589) £243
ARX9259		Without Centre mount Footplate	○
<b>CONTROLS*</b>			
ARJ2201		REM110 Drive only	(NC 1095, 1093, 1094, 0977, 0589, 0531, 2209, 2210, 2213, 2214, 2230, 2247, 2233, 2240, 2270, 2219, 2245, 2258, 1909) ◆
ARJ2203		REM211 Drive & Actuators	(NC 2209, 2210, 2213, 2214, 2231, 2232, 2230, 2247, 2233, 2240, 2270, 2258, 1909) ○
ARJ2204		REM216 Drive, Actuators & Lights	(NC 2210, 2213, 2214, 2231, 2232, 2230, 2247, 2233, 2240, 2270, 2258) £137
ARJ2205		REM400, touch-screen Interface, built in Bluetooth HID	(NC 2210, 2213, 2215) £305
ARJ2206		REM500, without joystick for alternative access, built in Bluetooth HID (Secondary User input required)	(NC 1208, 1707, 1708, 2227, 2228, 2270) £695

<b>CONTROL BRACKETRY*</b>			
ARJ1207		Standard Bracket RH	○
ARJ1208		Standard Bracket LH	○
ARJ1707		Swing Away RH Bracket	£119
ARJ1708		Swing Away LH Bracket	£119
ARJ1709		Matrx Resolve Swing in/out bracket RH	£219
ARJ1710		Matrx Resolve Swing in/out bracket LH	£219
ARJ2227		Bracket RH mounted for Attendant	£57
ARJ2228		Bracket LH mounted for Attendant	£57
ARJ2209		REM500 Swing away RH multi-adjustable Bracket (NC 2203, 2204, 2205)	○
ARJ2210		REM500 Swing away LH multi-adjustable Bracket (NC 2203, 2204, 2205)	○
ARJ2213		REM500 Centre mount swing away to RH Nucleus Bracket (NC 2203, 2204, 2205)	£446
ARJ2214		REM500 Centre mount swing away to LH Nucleus Bracket (NC 2203, 2204, 2205)	£446
ARZ2239		Individual mounting option per (special) - Requires additional information / specification	POA
<b>REM400 TOGGLE OPTION</b>			
ARJ2270		Assignable Toggle Switch Module for REM400, (as default set to left ON / OFF, right speed dial) (NC 2203, 2204)	£284
<b>CHIN CONTROL &amp; BRACKETS (Requires REM400 or REM500)</b>			
UNI9250		Mo-vis Micro joystick	POA
UNI9251		Mo-vis Multi joystick	POA
ARJ1289		Chin Control Swing Away Bracket RH Manual (NC 2203, 2204, 2171, 2179)	○
ARJ1288		Chin Control Swing Away Bracket LH Manual (NC 2203, 2204, 2171, 2179)	○
ARJ1287		Chin Control Swing Away Bracket RH Powered (NC 2203, 2204, 2227, 1404, 1415, 1312, 2171, 2179)	£2,251
ARJ1286		Chin Control Swing Away Bracket LH Powered (NC 2203, 2204, 2228, 1405, 1339, 1363, 2171, 2179)	£2,251
<b>DUAL CONTROL</b>			
ARJ2220		REM050 Dual Control (Please select user control and bracket) (Please select dual control bracket below) (NC 1221, 1222, 2227, 2228, 2205, 2206, 2245, 2128, 2230, 2247, 2233, 2240)	£355
ARJ2219		Compact Dual Control (DLX-ACU200) (Please select user control and bracket) (Please select dual control bracket below) (NC 1222, 2227, 2228, 2245, 2128, 2201)	£355
<b>DUAL CONTROL BRACKETRY</b>			
ARJ1207		Dual Control Standard Bracket RH	○
ARJ1208		Dual Control Standard Bracket LH	○
<b>INVACARE IDC DUAL CONTROL</b>			
ARJ2245		Invacare IDC Dual control (Please select IDC dual control bracket below) (NC 2219, 2227, 2228, 2128, 2201)	£673
<b>INVACARE IDC DUAL CONTROL BRACKETRY</b>			
ARJ1221		IDC Dual control Detachable Bracket RH (NC 2219, 2227, 2228)	○
ARJ1222		IDC Dual control Detachable Bracket LH (NC 2219, 2227, 2228)	○
<b>LINX COMPACT REMOTE</b>			
ARJ2231		Proportional Compact Remote DLX-CR400 (Requires Remote REM400 or REM500) (Please select Proportional Bracket) (NC 2203, 2204)	£740
ARJ2232		Proportional Compact Remote Low Force DLX-CR400LF (Requires Primary Remote REM400 or REM500) (Please select Proportional Bracketry) (NC 2203, 2204)	£990



LINX COMPACT REMOTE BRACKET				
ARJ2207		Proportional Compact RH Mount Bracket	(NC 2203, 2204)	○
ARJ2208		Proportional Compact LH Mount Bracket	(NC 2203, 2204)	○
ARJ2213		Proportional Compact Centre Mount swing away to RH (Nucleus)	(NC 2203, 2204)	£446
ARJ2214		Proportional Compact Centre Mount swing away to LH (Nucleus)	(NC 2203, 2204)	£446
ARZ2215		Individual mounting option per request (Special) Requires additional information / specification	(NC 2203, 2204)	POA
LINX HEAD CONTROL				
ARJ2230		Head Control with ultrasonic proximity sensors (ATOM , ASL104) with standard wings M3 (Requires Primary Remote REM400 or REM500)	(NC 2179, 2197, 2201, 2203, 2204, 2220, 2233, 2324, 2325, 2327)	£3,835
ARJ2247		Head Control with ultrasonic proximity sensors (PROTON , ASL104P) with PROTON wings (Requires Primary Remote REM400 or REM500)	(NC 2179, 2197, 2201, 2203, 2204, 2220, 2233, 2324, 2325, 2327)	£4,491
LINX MICRO EXTREMITY CONTROL				
ARZ2233		MEC Proportional Joystick (Micro Extremity Control) (Please select MEC Bracketry) (Requires Primary Remote REM400 or REM500)	(NC 2203, 2204)	£3,723
LINX MICRO EXTREMITY CONTROL BRACKET				
ARZ2235		MEC Tray Mount RH (ASL603, 100mm x 150mm) (Please select MEC Bracketry)	(NC 2203, 2204)	○
ARZ2236		MEC Tray Mount LH (ASL603, 100mm x 150mm)	(NC 2203, 2204)	○
ARZ2213		Centre Mount, swing away to RH, Nucleus	(NC 2203, 2204)	£446
ARZ2214		Centre Mount, swing away to LH, Nucleus	(NC 2203, 2204)	£446
ARZ2241		Centre Mount, swing away to RH, Nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 2203, 2204)	£446
ARZ2242		Centre Mount, swing away to LH, Nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 2203, 2204)	£446
ARZ2215		Individual mounting option per request (Special) Requires additional information / specification	(NC 2203, 2204)	POA
LINX HEAVY DUTY REMOTE				
ARZ2240		Proportional Heavy Duty Remote Joystick with mode switch (ASL133) (Please select Heavy Duty Bracketry) (Requires Primary Remote REM400 or REM500)	(NC 2203, 2204)	£2,076
LINX HEAVY DUTY REMOTE BRACKET				
ARZ2235		Tray Mount RH, (ASL603, 100mm x 150mm)	(NC 2203, 2204)	○
ARZ2236		Tray Mount LH, (ASL603, 100mm x 150mm)	(NC 2203, 2204)	○
ARZ2213		Centre Mount, swing away to RH, Nucleus	(NC 2203, 2204)	£446
ARZ2214		Centre Mount, swing away to LH, Nucleus	(NC 2203, 2204)	£446
ARZ2241		Centre Mount, swing away to RH, Nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 2203, 2204)	£446
ARZ2242		Centre Mount, swing away to LH, Nucleus (Inc Tray ASL603, 100mm x 150mm)	(NC 2203, 2204)	£446
ARZ2215		Individual mounting option per request (Special) Requires additional information / specification	(NC 2203, 2204)	POA
LINX 10-WAY SWITCH				
ARJ2251		10-Way Switch DLX-SW10-A for direct Powered Seating Access (included loom to connect + additional jack for mode switch) (Please select 10 way switch bracket below)		£453
LINX 10-WAY SWITCH BRACKET				
ARJ2217		10-Way Switch Bracket RH, Mounted for User	(NC 2213)	○
ARJ2218		10-Way Switch Bracket LH, Mounted for User	(NC 2214)	○
ARJ2227		10-Way Switch Bracket RH, Mounted for Attendant		£57
ARJ2228		10-Way Switch Bracket LH, Mounted for Attendant		£57

LINX ACCESSORY			
ARJ2252		Loom providing jack sockets for Drive and Mode (includes connection for 10-way switch)	£90
LINX USB CHARGER PORT			
ARJ2253		USB Charger Port DLX-USB02-A (NC 2701, 2703, 2704, 2705, 2706)	£143
LINX USB CHARGER PORT BRACKET			
ARJ2207		USB Charger Port RH Mount	○
ARJ2208		USB Charger Port LH Mount	○
FURTHER LINX COMPONENTS / ACCESSORIES			
ARJ2300		To be enclosed only (Accessory loom 2252 required if no REM400 or REM500 is selected)	○
ARZ2301		To be connected and programmed (special)	POA
BUDDY BUTTONS			
ARJ1685		Red Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	£79
ARJ1686		Yellow Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	£79
ARJ1687		Green Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	£79
ARJ1690		Blue Buddy Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	£79
EGG SWITCH			
ARJ2254		Egg Switch, black (to be enclosed only (Accessory loom 2252 required if no REM400 or REM500 is selected)	£79
PIKO BUTTONS			
ARJ1711		Red Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	£142
ARJ1712		Yellow Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	£142
ARJ1713		Green Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	£142
ARJ1714		Orange Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	£142
ARJ1715		Blue Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	£142
ARJ1716		Black Piko Button (Accessory loom 2252 required if no REM400 or REM500 is selected)	£142
TOGGLES / BUTTONS			
ARZ2261		Stereo Toggle - direct control for Tilt only	£242
ARZ2262		Stereo Button - direct control for Tilt only	£230
ARZ2263		Four-Way Toggle - direct control for Tilt & Recline	£300
ARZ2264		Four-Way Button - direct control for Tilt & Recline	£321
ARZ2207		Toggle / Button RH Mount	○
ARZ2208		Toggle / Button LH Mount	○
ARZ2225		Without Mounting	○
BODYPOINT JOYSTICK KNOBS			
ARJ2721		3" U-Shaped, Joystick Handle, Bodypoint PC101A	£65
ARJ2722		4" U-Shaped, Joystick Handle, Bodypoint PC102A	£65
ARJ2723		Domed Rubber, Joystick Handle, Bodypoint PC107A	£65
ARJ2724		Mushroom Grooved, Joystick Handle, Bodypoint PC110A	£65
REMOTE STOP SWITCH			
ARZ2258		Remote Stop Switch (ASL504)	£942



<b>LINX INPUT MODULE</b>			
UNI2273		IN200 Module (Jack port for external power button ) (NC 2220)	£242
UNI9252		IN500 Module (for DB9 Drive input)	£785
<b>SHROUD COLOUR*</b>			
ARJ2640		Obsidian Black Frame Shrouds	○
ARJ2641		Ruby Red Frame Shrouds	○
ARJ2642		Purple Passion Frame Shrouds	○
ARJ2643		Peacock Blue Frame Shrouds	○
ARJ2644		Arctic Blue Frame Shrouds	○
<b>RIM COLOUR*</b>			
ARJ2650		Obsidian Black Frame Shrouds	○
ARJ2651		Ruby Red Frame Shrouds	○
ARJ2652		Purple Passion Frame Shrouds	○
ARJ2653		Peacock Blue Frame Shrouds	○
ARJ2654		Arctic Blue Frame Shrouds	○
<b>FORKS*</b>			
ARJ2631		Single sided fork (New) includes Black castor	£316
ARJ2632		Double sided fork	○
<b>CENTRE WHEELS (DRIVING WHEELS)*</b>			
ARJ0704		DW 14" (355 mm) x 3" (75 mm), Pneumatic	○
ARJ0705		DW 14" (355 mm) x 3" (75 mm), Solid inlay, dual density foam insert	○
ARJ0706		DW 14" (355 mm) x 3" (75 mm), Puncture protected, gel-sealed	○
<b>MOTOR TYPE*</b>			
ARJ1294		3 km/h (2mph)	○
ARJ1282		6 km/h (4mph) (4 pole)	◆
ARJ1909		12 km/h (8mph) (includes light system and indicators) G-Trac recommended	£464
<b>BATTERY*</b>			
ARJ1260		73.5 Ah (Gel)	◆
<b>G TRAC</b>			
ARJ1215		G-Trac (ensures the powerchair takes the intended direction as indicated by the user control, without the need for continual correction by the client on uneven or sloping terrain)	£738
<b>CHARGER*</b>			
ARJ1343		OptiCharge Charger 8A (200 - 250V)	◆
ARJ1313		OptiCharge Charger 10A (200 - 250V) Charging time will be 20% faster	£137
<b>LIGHTS*</b>			
ARJ1437		Light system with Indicators & LED Headlights (NC 1909)	£431
ARJ1254		Without light system (NC 1909)	◆
ARZ1006		Lights Included with Motor Package 10 & 12kmh (NC 1282)	○
<b>TRANSPORTATION KIT*</b>			
ARJ1450		Preparation for DAHL Docking station: Docking station for car implementation, please order docking station separately (Adequate headrest is compulsory, see options below)	£524
ARJ1494		4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory, see options below)	◆

<b>FIXATION BELT*</b>			
ARJ2113		Lapbelt (adjustable both sides) adjustable mounting on seat rail extrusion (NC 2702, 2703, 2704, 2705, 2706)	○
ARJ1332		H-Belt	£66
ARJ2701		Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size L (padded length 620mm) (NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 2113)	£84
ARJ2702		Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size M (padded length 460mm) (NC 2123, 2124, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 2113)	£84
ARJ2703		Bodypoint 2-Point Hip Belt, Centre pull - push button, Cinch mount size L (padded length 620mm) (NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 2113)	£78
ARJ2704		Bodypoint 2-Point Hip Belt, Centre pull - push button, Cinch mount size M (padded length 460mm) (NC 2123, 2124, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 2113)	£78
ARJ2705		Bodypoint 2-Point Hip Belt, Dual pull Side release button, Cinch mount size L (padded length 620mm) (NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 2113)	£78
ARJ2706		Bodypoint 2-Point Hip Belt, Dual pull Side release button, Cinch mount size M (padded length 460mm) (NC 2123, 2124, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 2113)	£78
<b>REAR VIEW MIRROR</b>			
ARJ1415		Rear view Mirror RH - Must be selected with 10 km/h Motor (NC 2213, 1287)	£38
ARJ1339		Rear view Mirror LH - Must be selected with 10 km/h Motor (NC 2214, 1286)	£38
<b>BAG HOOKS</b>			
ARJ2128		Bag Hooks (Max load up to 6kg = impact on max user weight)	£69
<b>CANE HOLDER</b>			
ARJ1312		Cane holder RH (In combination with Attendant bracket or Dual control the cane holder has to be mounted on the opposite side of the remote)	£38
ARJ1363		Cane holder LH (In combination with Attendant bracket or Dual control the cane holder has to be mounted on the opposite side of the remote)	£38
<b>TRAY</b>			
ARJ1404		Tray swing away RH Select opposite side to user controls (NC 1707, 2213, 2214, 2209, 1287)	£303
ARJ1405		Tray swing away LH Select opposite side to user controls (NC 1708, 2213, 2214, 2210, 1286)	£303

<b>TOTAL PRICE</b>	
--------------------	--

**Non Conforming Options:**

ARJ2121	Set at 380mm (15") (adjustable from 380mm to 430mm (15" - 17") with stepless adjustment on armrests, legrests and backrest) (NC 2138 + 2142)
ARJ2122	Set at 430mm (17") (adjustable from 430mm to 480mm (17" - 19) with stepless adjustment on armrests, legrests and backrest) (NC 2138 + 2142)
ARJ2123	Set at 480mm (19") (adjustable from 480mm to 530mm (19" - 21") with stepless adjustment on armrests, legrests and backrest) (NC 2138 + 2142)
ARJ2124	Set at 530mm (21") (adjustable from 530mm to 580mm (21" - 23") with stepless adjustment on armrests, legrests and backrest)
ARJ2132	410mm (16") (410mm - 510mm simple stepless adjustment from the side)
ARJ2134	460mm (18") (410mm - 510mm simple stepless adjustment from the side)
ARJ2136	510mm (20") (410mm - 510mm simple stepless adjustment from the side)
ARJ2138	560mm (22") (510mm - 560mm simple stepless adjustment from the side) (Cushion (option 2142) is only available with Seat Width (option 2124) 530mm)
ARJ2143	Matrx Flo-tech Contour Visco (17x17 only)
ARJ2163	Backposts 480mm (Short)
ARJ2171	Tension Adjustable Backrest (540mm or 480mm upholstery height) Straps easily accessible and adjustable from the rear (Black fabric only)
ARJ2179	Flex3 Tension Adjustable (540mm or 480mm) Black Dartex only
ARJ2172	Flex3 adjustable Backrest Plates (without Backrest Cushion)
ARJ2173	Flex3 adjustable Backrest Plates with Laguna Backrest cushion
ARJ2174	Matrx Posture Back, height 400mm (Mesh fabric)
ARJ2175	Matrx Posture Back, height 500mm (Mesh fabric)
ARJ2176	Matrx Elite E2 Back , height 500 mm (Mesh fabric)
ARJ2177	Matrx Elite TR Back, height 500mm (Mesh fabric)
ARJ2178	Matrx Elite Deep E2 Back, height 500 mm (Mesh fabric)
ARJ2151	Black fabric
ARJ2155	Black Spacer fabric
ARJ2193	Flip-up Armrest (height and depth adjustable) (height adjustable from 230mm - 360mm) 150 kg max user weight
ARJ2192	Following Armrest (height adjustable from 230mm - 300mm) 150 kg max user weight
ARJ2194	Following Armrest (height adjustable from 300mm - 360mm) 150 kg max user weight
ARJ2196	Long Armpad (375mm)
ARJ2197	Long Ergonomic Armpad (397mm)
ARJ2199	Soft Armpad (REA Azalea Wide Pad)
ARJ2184	SMALL 140mm x 100mm, for flip up and following armrests only
ARJ2186	MEDIUM 230mm x 100mm, for flip up and following armrests only
ARJ2188	LARGE 330mm x 100mm, for flip up and following armrests only
ARJ2185	SMALL 140mm x 100mm, for flip up and following armrests only
ARJ2187	MEDIUM 230mm x 100mm, for flip up and following armrests only
ARJ2189	LARGE 330mm x 100mm, for flip up and following armrests only
ARJ2324	Elan Standard Pad 355mm (10") wide with standard cover complete with Elan Headrest Hardware
ARJ2325	Elan Occipital Pad 230mm x 180mm (9" x 7") with standard cover complete with Elan Headrest Hardware
ARJ2327	Elan 4 Point Headrest Pad 280mm x 255mm (11" x 10") with standard cover complete with Elan Headrest Hardware
ARJ1496	All bracketry as option 1432 without headrest cushion
ARJ1445	Multi-adjustable pad (Azalea type, Dartex)
ARJ1092	Powered "fixed pivot" Tilt 20° (Max User weight 136kg)
ARJ1095	Powered "fixed pivot" Tilt 20°, heavy duty (Max User weight 180kg)
ARJ1093	Electric Tilt (0° - 30°) with CoG shift Seat Height 460mm (160KG max user weight)

ARJ2052	Electric Tilt (0°-50°) with COG shift Seat height 460mm
ARJ1094	Seat lifter with Electrical Tilt (300mm, 30°) with CoG shift Seat Height 460mm - 490mm (160KG max user weight)
ARJ2053	Seat lifter with Electrical Tilt (300mm, 45°) with CoG shift Seat Height 460mm
ARJ0977	Powered Recliner (0° - 30°) Foldable
ARJ0556	Centre Mounted: Manual : Hanger, centre mounted, angle settable (0° - 90°) : knee to heel length 390 - 510mm, (including calf pad and angle adjustable footplate)
ARJ0589	Centre Mounted : LNX Powered: Centre mounted legrest, length compensating (including calf pad and angle adjustable footplate)
ARJ0526	Side Mounted : Standard 80° hangers
ARJ0503	Side Mounted : Vari F : Angle adjustable with tools 0° - 70°, length 290mm - 460mm
ARJ0533	Side Mounted : Vari A : Angle adjustable with hand wheel 0° - 70°, length 290mm - 460mm
ARJ0531	Side Mounted : ADE Powered : Angle adjustable powered by joystick, 0° - 80° length 290mm - 460mm
ARJ0501	Without Hanger
ARJ0512	Fixed Footplate
ARJ0504	Angle Adjustable Footplate
ARJ0513	Angle and depth adjustable Footplate (Adapt)
ARJ0559	Angle Adjustable Footplate Medium 140mm x 240mm (5.5" x 9.5") for Centre Mounted Hanger only
UNI0559	Angle & width adjustable footplate (240mm depth)
ARJ2201	REM110 Drive only
ARJ2203	REM211 Drive & Actuators
ARJ2204	REM216 Drive, Actuators & Lights
ARJ2205	REM400, touch-screen Interface, built in Bluetooth HID
ARJ2206	REM500, without joystick for alternative access, built in Bluetooth HID (Secondary User input required)
ARJ1208	Standard Bracket LH
ARJ1707	Swing Away RH Bracket
ARJ1708	Swing Away LH Bracket
ARJ2227	Bracket RH mounted for Attendant
ARJ2228	Bracket LH mounted for Attendant
ARJ2209	REM500 Swing away RH multi-adjustable Bracket
ARJ2210	REM500 Swing away LH multi-adjustable Bracket
ARJ2213	REM500 Centre mount swing away to RH Nucleus Bracket
ARJ2214	REM500 Centre mount swing away to LH Nucleus Bracket
ARJ2270	Assignable Toggle Switch Module for REM400, (as default set to left ON / OFF, right speed dial)
ARJ1287	Chin Control Swing Away Bracket RH Powered
ARJ1286	Chin Control Swing Away Bracket LH Powered
ARJ2220	REM050 Dual Control (Please select user control and bracket) (Please select dual control bracket below)
ARJ2219	Compact Dual Control (DLX-ACU200) (Please select user control and bracket) (Please select dual control bracket below)
ARJ1208	Dual Control Standard Bracket LH
ARJ2245	Invacare IDC Dual control (Please select IDC dual control bracket below)
ARJ1221	IDC Dual control Detachable Bracket RH
ARJ1222	IDC Dual control Detachable Bracket LH
ARJ2231	Proportional Compact Remote DLX-CR400 (Requires Remote REM400 or REM500) (Please select Proportional Bracket)
ARJ2232	Proportional Compact Remote Low Force DLX-CR400LF (Requires Primary Remote REM400 or REM500) (Please select Proportional Bracketry)
ARJ2213	Proportional Compact Centre Mount swing away to RH (Nucleus)
ARJ2214	Proportional Compact Centre Mount swing away to LH (Nucleus)

ARZ2215	Individual mounting option per request (Special) Requires additional information / specification
ARJ2230	Head Control with ultrasonic proximity sensors (ATOM , ASL104) with standard wings M3 (Requires Primary Remote REM400 or REM500)
ARJ2247	Head Control with ultrasonic proximity sensors (PROTON , ASL104P) with PROTON wings (Requires Primary Remote REM400 or REM500)
ARZ2233	MEC Proportional Joystick (Micro Extremity Control) (Please select MEC Bracketry) (Requires Primary Remote REM400 or REM500)
ARZ2213	Centre Mount, swing away to RH, Nucleus
ARZ2214	Centre Mount, swing away to LH, Nucleus
ARZ2215	Individual mounting option per request (Special) Requires additional information / specification
ARZ2240	Proportional Heavy Duty Remote Joystick with mode switch (ASL133) (Please select Heavy Duty Bracketry) (Requires Primary Remote REM400 or REM500)
ARZ2213	Centre Mount, swing away to RH, Nucleus
ARZ2214	Centre Mount, swing away to LH, Nucleus
ARZ2215	Individual mounting option per request (Special) Requires additional information / specification
ARJ2227	10-Way Switch Bracket RH, Mounted for Attendant
ARJ2228	10-Way Switch Bracket LH, Mounted for Attendant
ARJ2253	USB Charger Port DLX-USB02-A
ARZ2258	Remote Stop Switch (ASL504)
ARJ1282	6 km/h (4mph) (4 pole)
ARJ1909	12 km/h (8mph) (includes light system and indicators) G-Trac recommended
ARJ2113	Lapbelt (adjustable both sides) adjustable mounting on seat rail extrusion
ARJ1332	H-Belt
ARJ2701	Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size L (padded length 620mm)
ARJ2702	Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size M (padded length 460mm)
ARJ2703	Bodypoint 2-Point Hip Belt, Centre pull - push button, Cinch mount size L (padded length 620mm)
ARJ2704	Bodypoint 2-Point Hip Belt, Centre pull - push button, Cinch mount size M (padded length 460mm)
ARJ2705	Bodypoint 2-Point Hip Belt, Dual pull Side release button, Cinch mount size L (padded length 620mm)
ARJ2706	Bodypoint 2-Point Hip Belt, Dual pull Side release button, Cinch mount size M (padded length 460mm)
ARJ1415	Rear view Mirror RH - Must be selected with 10 km/h Motor
ARJ1339	Rear view Mirror LH - Must be selected with 10 km/h Motor
ARJ2128	Bag Hooks (Max load up to 6kg = impact on max user weight)
ARJ1312	Cane holder RH (In combination with Attendant bracket or Dual control the cane holder has to be mounted on the opposite side of the remote)
ARJ1363	Cane holder LH (In combination with Attendant bracket or Dual control the cane holder has to be mounted on the opposite side of the remote)
ARJ1404	Tray swing away RH Select opposite side to user controls
ARJ1405	Tray swing away LH Select opposite side to user controls



For a more complex customer requirement please refer to our Invacare Unique special order form. email [specialsuk@invacare.com](mailto:specialsuk@invacare.com) or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Description of special request (please include relevant option code or part number):