



**Yes, you can:**

Tel: 01656 776222 Fax: 01656 776220  
Email: ordersuk@invacare.com  
Online Spares available at: www.invacare.co.uk

SEPTEMBER 2024  
Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

**TDX SP2 Low Rider**

**Standard Features**

- Max user weight 160kg (25.2 Stone)
- Crash tested ISO7176-Part 19
- Compact, manoeuvrable, CWD TDX SP2
- Telescopic adjustable Seat Widths from 430mm - 530mm (17" - 21")
- Seat Depths 410mm - 560mm Stepless Adjustment (16" - 22")
- Tension Adjustable backrest
- Standard T-Armrests, detachable with Short Armpad
- Seat lifter with powered Tilt (300mm, 30° ) CoG Shift
- Choice of seat height either 420mm or 403mm
- Manual settable Recliner 0° - 30° , in 3.75° steps
- Standard 80° Hanger with fixed footplate
- LiNX Control REM211 with standard bracket
- Selection of 10 Stunning Shroud and Rim colours
- Centre (Driving Wheels) 355mm x 76mm (14" x 3") Black Tyres
- Suspension and Sure Step Technology
- 6 km/h 4 Pole Motors
- 8A (200 - 25V) Charge
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)

**Standard Base price: £7,328**



Q'STRAIT  
QLK-150  
Compatible

Please make your selection with a cross in the box    ○ Free of charge    ◆ standard feature    POA special request

Customer weight =                      KG

Code	Description	
<b>BASE PRODUCT*</b>		
ARQ0071	TDX SP2 Low Rider 73ah Batteries 420mm seat to ground height	○
ARQ0072	TDX SP2 Low Rider 60ah Batteries 403mm seat to ground height	○
ARZ0075	Special TDX SP2 Low Rider 73ah Batteries 420mm seat to ground height	POA
ARZ0076	Special TDX SP2 Low Rider 60ah Batteries 403mm seat to ground height	POA

<b>User Weight*</b>			
ARQ1918	<input type="checkbox"/>	User weight less 100 kg	<input type="radio"/>
ARQ1919	<input type="checkbox"/>	User weight 101 to 125 kg	<input type="radio"/>
ARQ1921	<input type="checkbox"/>	User weight 126 to 136 kg (NC 1350, 1432, 1427)	<input type="radio"/>
ARQ1922	<input type="checkbox"/>	User weight 137 to 150 kg (NC 1350, 1432, 1427, 2053)	<input type="radio"/>
ARQ1923	<input type="checkbox"/>	User weight 151 to 160 kg (NC 2192, 2193, 2194, 1908, 1350, 1432, 1427, 2053)	<input type="radio"/>
<b>SEATING TYPE*</b>			
ARQ2100	<input type="checkbox"/>	Modulite Seating Unit	£7,328
<b>SEAT FRAME*</b>			
ARQ2111	<input type="checkbox"/>	Telescopic Seat Frame with metal seat plate	◆
<b>SEAT WIDTH*</b>			
ARQ2122	<input type="checkbox"/>	Set at 430mm (17") (adjustable from 430mm to 480mm (17" - 19)) (NC 2703, 2705)	<input type="radio"/>
ARQ2123	<input type="checkbox"/>	Set at 480mm(19") (adjustable from 480mm to 530mm (19" - 21")) (NC 2702, 2704, 2706)	<input type="radio"/>
ARQ2124	<input type="checkbox"/>	Set at 530mm (21") (adjustable from 530mm to 580mm (21" - 23")) (NC 2174, 2175, 2176, 2177, 2178, 2143, 2702, 2704, 2706)	<input type="radio"/>
<b>SEAT DEPTH RAIL LENGTH*</b>			
ARQ2129	<input type="checkbox"/>	Standard Depth Package (16" - 20", 410mm - 510mm ) (NC 2138)	<input type="radio"/>
ARQ2130	<input type="checkbox"/>	Extended Depth Package (20" - 22", 510mm - 560mm) (NC 2132, 2134)	<input type="radio"/>
<b>SEAT DEPTH*</b>			
ARQ2132	<input type="checkbox"/>	410mm (16") (410mm - 510mm) (NC 2196, 2197, 2701, 2702, 2703, 2704, 2705, 2706)	<input type="radio"/>
ARQ2134	<input type="checkbox"/>	460mm (18") (410mm - 510mm)	<input type="radio"/>
ARQ2136	<input type="checkbox"/>	510mm (20") (410mm - 510mm)	<input type="radio"/>
ARQ2138	<input type="checkbox"/>	560mm (22") (510mm - 560mm) (NC 2143)	<input type="radio"/>
<b>SEAT OPTION</b>			
ARQSWO	<input type="checkbox"/>	"H" Frame Interface Kit (NC 2143)	£446
<b>SEAT CUSHION*</b>			
ARQ2142	<input type="checkbox"/>	Matrx Flo-tech Lite	£112
ARQ2143	<input type="checkbox"/>	Matrx Flo-tech Contour Visco (NC 2136, 2124, 2138)	£112
ARQ2147	<input type="checkbox"/>	Matrx Libra 90mm thick	£391
ARQ2140	<input type="checkbox"/>	Without seat cushion	<input type="radio"/>
<b>BACKREST SUPPORT*</b>			
ARQ2162	<input type="checkbox"/>	Backposts 540mm (Tall) must be selected for Flex3 adjustable backrest plates	<input type="radio"/>
ARQ2163	<input type="checkbox"/>	Backposts 480mm (Short) (NC 2173)	<input type="radio"/>
<b>PUSH HANDLES FOR BACKPOSTS*</b>			
ARQ2107	<input type="checkbox"/>	Backposts WITH push handles	<input type="radio"/>
ARQ2108	<input type="checkbox"/>	Backposts WITHOUT push handles (NC 1332, 2220, 2219, 2245, 2227, 2228)	£36

<b>BACKREST UPHOLSTERY / CUSHION*</b>			
ARQ2171		Basic Tension Adjustable Backrest	(NC 1445, 1496, 2155, 2230, 2247) ○
ARQ2179		Flex3 Tension Adjustable (540mm or 480mm) Black multi stretch PU coated fabric only	(NC 2151, 2192, 2193, 2194, 2199, 2155, 2230, 2247, 2128) £139
ARQ2172		Flex3 Adjustable Backrest Plates (without Backrest Cushion)	(NC 2163, 2128) £224
ARQ2173		Flex3 Adjustable Backrest Plates with Laguna Backrest cushion	(NC 2163, 2128) £363
ARQ2174		Matrx Posture Back, height 400mm (Mesh fabric)	(NC 1496, 2124) £496
ARQ2175		Matrx Posture Back, height 500mm (Mesh fabric)	(NC 1496, 2124) £496
ARQ2176		Matrx Elite E2 Back , height 500 mm (Mesh fabric)	(NC 1496, 2124) £496
ARQ2177		Matrx Elite TR Back, height 500mm (Mesh fabric)	(NC 1496, 2124) £496
ARQ2178		Matrx Elite Deep E2 Back, height 500 mm (Mesh fabric)	(NC 1496, 2124) £550
ARQ2170		Without Backrest Upholstery	(NC 1496, 2230, 2247) ○
<b>UPHOLSTERY TYPE*</b>			
ARQ2151		Black fabric	○
ARQ2152		Black multi stretch PU coated fabric	(NC 2171) ○
ARQ2155		Black Spacer fabric	(NC 2171, 2179) ○
<b>ARMREST*</b>			
ARQ2191		T-Armrest Standard, Range of adjustment	○
ARQ2181		T-Armrest Extended height	£36
ARQ2183		In-rigged Armrests	(NC 2185, 2186, 2187, 2188, 2189) £25
ARQ2193		Flip-up Armrest (from 230mm - 360mm)	(NC 2176, 2199) £58
ARQ2192		Following Armrest (from 230mm - 300mm)	(NC 1445, 2179) £166
ARQ2194		Following Armrest (from 300mm - 360mm)	(NC 1445, 2179) £166
<b>ARMPAD*</b>			
ARQ2195		Short Armpad (325mm)	○
ARQ2196		Long Armpad (375mm)	○
ARQ2197		Long Ergonomic Armpad (397mm)	(NC 2247) ○
ARQ2199		Soft Arpad (REA Azalea Wide Pad)	£86
<b>HIP SUPPORTS LH</b>			
ARQ2184		SMALL 140mm x 100mm	(NC 2702, 2703, 2704, 2705, 2706) £226
ARQ2186		MEDIUM 230mm x 100mm	(NC 2702, 2703, 2704, 2705, 2706) £226
ARQ2188		LARGE 330mm x 100mm	(NC 2702, 2703, 2704, 2705, 2706) £226
<b>HIP SUPPORTS RH</b>			
ARQ2185		SMALL 140mm x 100mm	(NC 2702, 2703, 2704, 2705, 2706) £226
ARQ2187		MEDIUM 230mm x 100mm	(NC 2702, 2703, 2704, 2705, 2706) £226
ARQ2189		LARGE 330mm x 100mm	(NC 2702, 2703, 2704, 2705, 2706) £226
<b>MOUNTING HARDWARE FOR HIP SUPPORTS</b>			
ARQ2165		Mounting hardware for hip supports (8.1" / 206mm high)	○

<b>HEADREST / NECKREST*</b>			
ARQ1350		Neckrest height and depth adjustable, Max User weight 136kg (NC 2247)	£72
ARQ1432		Headrest, height, depth and sideways adjustable, Max User weight 136kg (NC 2247)	○
ARQ1427		Headrest with cheek support, height and depth adjustable, Max User weight 136kg (NC 2247)	£203
ARQ2324		Elan Standard Pad 355mm (10") wide (NC 2247)	£294
ARQ2325		Elan Occipital Pad 230mm x 180mm (9" x 7") (NC 2247)	£294
ARQ2327		Elan 4 Point Headrest Pad 280mm x 255mm (11" x 10") (NC 2247)	£396
ARQ1496		All Bracketry as above option without Headrest cushion	○
ARQ1495		Push bar dia. 22mm	○
<b>LATERAL SUPPORT</b>			
ARQ1445		Multi-adjustable pad (Azalea type, multi stretch PU coated fabric) (NC 2171)	£150
<b>SEAT LIFTER / SEAT TILT*</b>			
ARQ1094		Seat lifter with Electrical Tilt	○
ARQ2053		Seat lifter with Electrical Tilt (300mm, 45°) with CoG shift Seat Height 460mm	£663
<b>SEAT TILT OPTION</b>			
ARQ2054		5° Fixed Anterior Tilt Assist mount	£361
<b>BACKREST RECLINER*</b>			
ARQ0975		Manual Settable Recliner	○
ARQ0977		Powered Recliner (0° - 30°) Foldable	£453
<b>FOOTREST HANGER TYPE*</b>			
ARQ0556		Centre Mounted: Manual : Hanger, centre mounted (NC 0504, 0513)	£334
ARQ0589		Centre Mounted : LNX Powered: Centre mounted legrest (NC 0504, 0513)	£1,470
UNI0592		LNX Power centre mount legrest with telescoping foot platform (lowers to the ground) inc calf pads: Length 380-419mm (Selectable length depends on STF height) Custom-made, requires customised setting (NC 0512, 0504)	£2,417
ARQ0526		Side Mounted : Standard 80° hangers (NC 0559)	○
ARQ0503		Side Mounted : Vari F : Angle adjustable	£107
ARQ0533		Side Mounted : Vari A : Angle adjustable (NC 0559)	£210
ARQ0531		Side Mounted : ADE Powered (NC 0513, 0559)	£600
ARQ0501		Without Hanger (NC 0504, 0513, 0559)	○
<b>FOOTPLATE TYPE*</b>			
ARQ0512		Fixed Footplate (NC 0533, 0556, 0589, 0501)	○
ARQ0504		Angle Adjustable Footplate (NC 0589, 0501)	£64
ARQ0513		Angle and depth adjustable Footplate (Adapt) (NC 0589, 0501)	£106
ARQ0559		Angle Adjustable Footplate for Centre mounted legrest (NC 0503, 0533, 0531, 0501)	○
ARQXXXX		Without footplates f (NC 0589, 0526, 0503, 0533, 0531)	POA
<b>CONTROLS*</b>			
ARQ2203		REM211 Drive & Actuators (NC 2209, 2210, 2213, 2214, 2231, 2232, 2230, 2247, 2233, 2240, 2270, 2258, 1908, 1909)	○
ARQ2204		REM216 Drive, Actuators & Lights (NC 2210, 2213, 2214, 2231, 2232, 2230, 2247, 2233, 2240, 2270, 2258)	£137
ARQ2205		REM400, touch-screen Interface, built in Bluetooth HID (NC 2210, 2213, 2215)	£305
ARQ2206		REM500 (NC 1208, 1707, 1708, 2227, 2228, 2270)	£695

<b>CONTROL BRACKETRY*</b>			
ARQ1207		Standard Bracket RH	○
ARQ1208		Standard Bracket LH	○
ARQ1707		Swing Away RH Bracket	£119
ARQ1708		Swing Away LH Bracket	£119
ARQ1709		Matrx Resolve Swing in/out bracket RH	£219
ARQ1710		Matrx Resolve Swing in/out bracket LH	£219
ARQ2213		Centre mount swing away to RH Nucleus	£446
ARQ2214		Centre mount swing away to LH Nucleus	£446
ARQ2227		Bracket RH mounted for Attendant (NC 2203, 2204, 2205)	£57
ARQ2228		Bracket LH mounted for Attendant (NC 2203, 2204, 2205)	£57
ARZ2239		Individual mounting option per (special)	POA
ARQ2209		REM500 Swing away RH multi-adjustable (NC 2203, 2204, 2205)	○
ARQ2210		REM500 Swing away LH multi-adjustable (NC 2203, 2204, 2205)	○
<b>REM400 TOGGLE OPTION</b>			
ARQ2270		Assignable Toggle Switch Module for REM400, (NC 2203, 2204)	£284
<b>CHIN CONTROL &amp; BRACKETS (Requires REM400 or REM500)</b>			
UNI9250		Mo-vis Micro joystick (NC 0072)	POA
UNI9251		Mo-vis Multi joystick (NC 0072)	POA
<b>CHIN CONTROL BRACKETRY</b>			
ARQ1289		Chin Control Swing Away Bracket RH Manual	○
ARQ1288		Chin Control Swing Away Bracket LH Manual	○
ARQ1287		Chin Control Swing Away Bracket RH Powered	£2,251
ARQ1286		Chin Control Swing Away Bracket LH Powered	£2,251
<b>DUAL CONTROL</b>			
ARQ2220		REM050 Dual Control (NC 1221, 1222, 2227, 2228, 2205, 2206, 2245, 2128, 2230, 2247, 2233, 2240)	£254
ARQ2219		Compact Dual Control (DLX-ACU200) (NC 1222, 2227, 2228, 2245, 2128)	£355
<b>DUAL CONTROL BRACKETRY</b>			
ARQ1207		Dual Control Standard Bracket RH	○
ARQ1208		Dual Control Standard Bracket LH	○
<b>INVACARE LINX IDC DUAL CONTROL</b>			
ARQ2245		IDC Dual control (NC 2219, 2227, 2228, 2128)	£673
<b>INVACARE LINX IDC DUAL CONTROL BRACKETRY</b>			
ARQ1221		IDC Dual Control Detachable Bracket RH (NC 2219, 2227, 2228)	○
ARQ1222		IDC Dual Control Detachable Bracket LH (NC 2219, 2227, 2228)	○
<b>LINX COMPACT REMOTE</b>			
ARQ2231		Proportional Compact Remote DLX-CR400 (NC 2203, 2204)	£737
ARQ2232		Proportional Compact Remote Low Force DLX-CR400LF (NC 2203, 2204)	£990
<b>LINX COMPACT REMOTE BRACKET</b>			
ARQ2207		Proportional Compact RH Mount Bracket (NC 2203, 2204)	○
ARQ2208		Proportional Compact LH Mount Bracket (NC 2203, 2204)	○
ARQ2213		Proportional Compact Centre Mount swing away to RH (Nucleus) (NC 2203, 2204)	£446
ARQ2214		Proportional Compact Centre Mount swing away to LH (Nucleus) (NC 2203, 2204)	£446
ARZ2215		Individual mounting option per request (Special) (NC 2203, 2204)	POA
<b>LINX HEAD CONTROL</b>			
ARQ2230		Head Control with ultrasonic proximity sensors (NC 2179, 2197, 2203, 2204, 2220, 2233, 2324, 2325, 2327)	£3,835
ARQ2247		Head Control with ultrasonic proximity sensors (NC 2179, 2197, 2203, 2204, 2220, 2233, 2324, 2325, 2327)	£4,491

<b>LINX MICRO EXTREMITY CONTROL</b>			
ARZ2233		MEC Proportional Joystick	(NC 2203, 2204) £3,723
<b>LINX MICRO EXTREMITY CONTROL BRACKET</b>			
ARZ2235		MEC Tray Mount RH	(NC 2203, 2204) ○
ARZ2236		MEC Tray Mount LH	(NC 2203, 2204) ○
ARZ2213		Centre Mount, swing away to RH, Nucleus	(NC 2203, 2204) £446
ARZ2214		Centre Mount, swing away to LH, Nucleus	(NC 2203, 2204) £446
ARZ2241		Centre Mount, swing away to RH, Nucleus	(NC 2203, 2204) £446
ARZ2242		Centre Mount, swing away to LH, Nucleus	(NC 2203, 2204) £446
ARZ2215		Individual mounting option per request (Special)	(NC 2203, 2204) POA
<b>LINX HEAVY DUTY REMOTE</b>			
ARZ2240		Proportional Heavy Duty Remote Joystick with mode switch	(NC 2203, 2204) £2,076
<b>LINX HEAVY DUTY REMOTE BRACKET</b>			
ARZ2235		Tray Mount RH, (ASL603, 100mm x 150mm)	(NC 2203, 2204) ○
ARZ2236		Tray Mount LH, (ASL603, 100mm x 150mm)	(NC 2203, 2204) ○
ARZ2213		Centre Mount, swing away to RH, Nucleus	(NC 2203, 2204) £446
ARZ2214		Centre Mount, swing away to LH, Nucleus	(NC 2203, 2204) £446
ARZ2241		Centre Mount, swing away to RH, Nucleus	(NC 2203, 2204) £446
ARZ2242		Centre Mount, swing away to LH, Nucleus	(NC 2203, 2204) £446
ARZ2215		Individual mounting option per request (Special)	(NC 2203, 2204) POA
<b>LINX 10-WAY SWITCH</b>			
ARQ2251		10-Way Switch DLX-SW10-A for direct Powered Seating Access	£453
<b>LINX 10-WAY SWITCH BRACKET</b>			
ARQ2217		10-Way Switch Bracket RH, Mounted for User	(NC 2213) ○
ARQ2218		10-Way Switch Bracket LH, Mounted for User	(NC 2214) ○
ARQ2227		10-Way Switch Bracket RH, Mounted for Attendant	£57
ARQ2228		10-Way Switch Bracket LH, Mounted for Attendant	£57
<b>LINX ACCESSORY</b>			
ARQ2252		Loom providing jack sockets for Drive and Mode	£90
<b>LINX USB CHARGER PORT</b>			
ARQ2253		USB Charger Port DLX-USB02-A	(NC 2701, 2703, 2704, 2705, 2706) £143
<b>LINX USB CHARGER PORT BRACKET</b>			
ARQ2207		USB Charger Port RH Mount	○
ARQ2208		USB Charger Port LH Mount	○
<b>LINX INPUT MODULE</b>			
UNI2273		IN200 Module (Jack port for external power button )	(NC 0072, 2220) £242
UNI9252		IN500 Module (for DB9 Drive input)	(NC 0072) £785
<b>FURTHER LINX COMPONENTS / ACCESSORIES</b>			
ARQ2300		To be enclosed only	○
ARZ2301		To be connected and programmed (special)	POA
<b>BUDDY BUTTONS</b>			
ARQ1685		Red Buddy Button	£79
ARQ1686		Yellow Buddy Button	£79
ARQ1687		Green Buddy Button	£79
ARQ1690		Blue Buddy Button	£79
<b>EGG SWITCH</b>			
ARQ2254		Egg Switch, black	£79

PIKO BUTTONS			
ARQ1711		Red Piko Button	£142
ARQ1712		Yellow Piko Button	£142
ARQ1713		Green Piko Button	£142
ARQ1714		Orange Piko Button	£142
ARQ1715		Blue Piko Button	£142
ARQ1716		Black Piko Button	£142
TOGGLES / BUTTONS			
ARZ2261		Stereo Toggle - direct control for Tilt only	£242
ARZ2262		Stereo Button - direct control for Tilt only	£230
ARZ2263		Four-Way Toggle - direct control for Tilt & Recline	£300
ARZ2264		Four-Way Button - direct control for Tilt & Recline	£321
ARZ2207		Toggle / Button RH Mount	○
ARZ2208		Toggle / Button LH Mount	○
ARZ2225		Without Mounting	○
BODYPOINT JOYSTICK KNOBS			
ARQ2721		3" U-Shaped, Joystick Handle, Bodypoint PC101A	£65
ARQ2722		4" U-Shaped, Joystick Handle, Bodypoint PC102A	£65
ARQ2723		Domed Rubber, Joystick Handle, Bodypoint PC107A	£65
ARQ2724		Mushroom Grooved, Joystick Handle, Bodypoint PC110A	£65
REMOTE STOP SWITCH			
ARZ2258		Remote Stop Switch (ASL504)	£942
SHROUD COLOUR*			
ARQ2640		Obsidian Black Frame Shrouds	○
ARQ2641		Ruby Red Frame Shrouds	○
ARQ2642		Purple Passion Frame Shrouds	○
ARQ2643		Peacock Blue Frame Shrouds	○
ARQ2644		Arctic Blue Frame Shrouds	○
RIM COLOUR*			
ARQ2650		Obsidian Black Rim Inserts	○
ARQ2651		Ruby Red Rim Inserts	○
ARQ2652		Purple Passion Rim Inserts	○
ARQ2653		Peacock Blue Rim Inserts	○
ARQ2654		Arctic Blue Rim Inserts	○
FORKS*			
ARQ2631		Single sided fork (New) includes Black castor	£316
ARQ2632		Double sided fork	○
CENTRE WHEELS (DRIVING WHEELS)*			
ARQ0704		DW 14" (355 mm) x 3" (75 mm), Pneumatic	○
ARQ0705		DW 14" (355 mm) x 3" (75 mm), Solid inlay, dual density foam insert	○
ARQ0706		DW 14" (355 mm) x 3" (75 mm), Puncture protected, gel-sealed	○

<b>MOTOR TYPE*</b>			
ARQ1294	<input type="checkbox"/>	3 km/h (2mph)	○
ARQ1282	<input type="checkbox"/>	6 km/h (4mph) (4 pole)	○
ARQ1908	<input type="checkbox"/>	10 km/h (6mph) (includes light system and indicators) G-Trac recommended	£464
ARQ1909	<input type="checkbox"/>	12 km/h (8mph) (includes light system and indicators) G-Trac recommended	£464
<b>G TRACK</b>			
ARQ1215	<input type="checkbox"/>	G-Trac	£738
<b>CHARGER*</b>			
ARQ1343	<input type="checkbox"/>	OptiCharge Charger 8A (200 - 250V)	○
ARQ1313	<input type="checkbox"/>	OptiCharge Charger 10A (200 - 250V)	£137
<b>LIGHTS*</b>			
ARQ1437	<input type="checkbox"/>	Light system with Indicators & LED Headlights	(NC 1908, 1909) £431
ARQ1254	<input type="checkbox"/>	Without Light system	(NC 1909) ○
ARQ1006	<input type="checkbox"/>	Lights Included with Motor Package 10 & 12kmh	(NC 1282) ○
<b>TRANSPORTATION KIT</b>			
ARQ1450	<input type="checkbox"/>	Preparation for DAHL Docking station:	£524
ARQ1494	<input type="checkbox"/>	4-Point tie-down	○
<b>FIXATION BELT*</b>			
ARQ1329	<input type="checkbox"/>	Lapbelt (adjustable both sides)	(NC 2702, 2703, 2704, 2705, 2706) ○
ARQ2113	<input type="checkbox"/>	Lapbelt (adjustable both sides)	○
ARQ1332	<input type="checkbox"/>	H-Belt	£66
ARQ2701	<input type="checkbox"/>	Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size L (padded length 620mm)	(NC 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113) £84
ARQ2702	<input type="checkbox"/>	Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size M (padded length 460mm)	(NC 2123, 2124, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113) £84
ARQ2703	<input type="checkbox"/>	Bodypoint 2-Point Hip Belt, Centre pull - push button	(NC 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113) £78
ARQ2704	<input type="checkbox"/>	Bodypoint 2-Point Hip Belt, Centre pull - push button	(NC 2123, 2124, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113) £78
ARQ2705	<input type="checkbox"/>	Bodypoint 2-Point Hip Belt, Dual pull Side release	(NC 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113) £78
ARQ2706	<input type="checkbox"/>	Bodypoint 2-Point Hip Belt, Dual pull Side release	(NC 2123, 2124, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113) £78
<b>REAR-VIEW MIRROR</b>			
ARQ1415	<input type="checkbox"/>	Rear View Mirror RH	(NC 2213, 1287) £38
ARQ1339	<input type="checkbox"/>	Rear View Mirror LH	(NC 2214, 1286) £38
<b>BAG HOOKS</b>			
ARQ2128	<input type="checkbox"/>	Bag Hooks (Max load up to 6kg = impact on max user weight)	£69
<b>CANE HOLDER</b>			
ARQ1312	<input type="checkbox"/>	Cane holder RH	£38
ARQ1363	<input type="checkbox"/>	Cane holder LH	£38
<b>TRAY</b>			
ARQ1404	<input type="checkbox"/>	Tray Swing Away RH	(NC 1707, 2213, 2214, 2209, 1287) £303
ARQ1405	<input type="checkbox"/>	Tray Swing Away LH	(NC 1708, 2213, 2214, 2210, 1286) £303
<b>TOTAL PRICE</b>			



**Non Conforming Options:**

ARQ0072	TDX SP2 Low Rider 60ah Batteries 403mm seat to ground height
ARQ2122	Set at 430mm (17") (adjustable from 430mm to 480mm (17" - 19))
ARQ2123	Set at 480mm(19") (adjustable from 480mm to 530mm (19" - 21"))
ARQ2124	Set at 530mm (21") (adjustable from 530mm to 580mm (21" - 23"))
ARQ2132	410mm (16") (410mm - 510mm)
ARQ2134	460mm (18") (410mm - 510mm)
ARQ2136	510mm (20") (410mm - 510mm)
ARQ2138	560mm (22") (510mm - 560mm)
ARQ2143	Matrx Flo-tech Contour Visco
ARQ2163	Backposts 480mm (Short)
ARQ2171	Basic Tension Adjustable Backrest
ARQ2179	Flex3 Tension Adjustable (540mm or 480mm) Black multi stretch PU coated fabric only
ARQ2173	Flex3 Adjustable Backrest Plates with Laguna Backrest cushion
ARQ2174	Matrx Posture Back, height 400mm (Mesh fabric)
ARQ2175	Matrx Posture Back, height 500mm (Mesh fabric)
ARQ2176	Matrx Elite E2 Back , height 500 mm (Mesh fabric)
ARQ2177	Matrx Elite TR Back, height 500mm (Mesh fabric)
ARQ2178	Matrx Elite Deep E2 Back, height 500 mm (Mesh fabric)
ARQ2151	Black fabric
ARQ2155	Black Spacer fabric
ARQ2193	Flip-up Armrest (from 230mm - 360mm)
ARQ2192	Following Armrest (from 230mm - 300mm)
ARQ2194	Following Armrest (from 300mm - 360mm)
ARQ2196	Long Armpad (375mm)
ARQ2197	Long Ergonomic Armpad (397mm)
ARQ2199	Soft Arpad (REA Azalea Wide Pad)
ARQ2184	SMALL 140mm x 100mm
ARQ2186	MEDIUM 230mm x 100mm
ARQ2188	LARGE 330mm x 100mm
ARQ2185	SMALL 140mm x 100mm
ARQ2187	MEDIUM 230mm x 100mm
ARQ2189	LARGE 330mm x 100mm
ARQ1350	Neckrest height and depth adjustable, Max User weight 136kg
ARQ1432	Headrest, height, depth and sideways adjustable, Max User weight 136kg
ARQ1427	Headrest with cheek support, height and depth adjustable, Max User weight 136kg
ARQ2324	Elan Standard Pad 355mm (10") wide
ARQ2325	Elan Occipital Pad 230mm x 180mm (9" x 7")
ARQ2327	Elan 4 Point Headrest Pad 280mm x 255mm (11" x 10")
ARQ1496	All Bracketry as above option without Headrest cushion
ARQ1445	Multi-adjustable pad (Azalea type, multi stretch PU coated fabric)
ARQ2053	Seat lifter with Electrical Tilt (300mm, 45°) with CoG shift Seat Height 460mm
ARQ0556	Centre Mounted: Manual : Hanger, centre mounted
ARQ0589	Centre Mounted : LNX Powered: Centre mounted legrest

ARQ0526	Side Mounted : Standard 80° hangers
ARQ0503	Side Mounted : Vari F : Angle adjustable
ARQ0533	Side Mounted : Vari A : Angle adjustable
ARQ0531	Side Mounted : ADE Powered
ARQ0501	Without Hanger
ARQ0512	Fixed Footplate
ARQ0504	Angle Adjustable Footplate
ARQ0513	Angle and depth adjustable Footplate (Adapt)
ARQ0559	Angle Adjustable Footplate for Centre mounted legrest
ARQ2203	REM211 Drive & Actuators
ARQ2204	REM216 Drive, Actuators & Lights
ARQ2205	REM400, touch-screen Interface, built in Bluetooth HID
ARQ2206	REM500
ARQ1208	Standard Bracket LH
ARQ1707	Swing Away RH Bracket
ARQ1708	Swing Away LH Bracket
ARQ2213	Centre mount swing away to RH Nucleus
ARQ2214	Centre mount swing away to LH Nucleus
ARQ2227	Bracket RH mounted for Attendant
ARQ2228	Bracket LH mounted for Attendant
ARQ2209	REM500 Swing away RH multi-adjustable
ARQ2210	REM500 Swing away LH multi-adjustable
ARQ2270	Assignable Toggle Switch Module for REM400,
ARQ1287	Chin Control Swing Away Bracket RH Powered
ARQ1286	Chin Control Swing Away Bracket LH Powered
ARQ2220	REM050 Dual Control
ARQ2219	Compact Dual Control (DLX-ACU200)
ARQ1208	Dual Control Standard Bracket LH
ARQ2245	IDC Dual control
ARQ1221	IDC Dual Control Detachable Bracket RH
ARQ1222	IDC Dual Control Detachable Bracket LH
ARQ2231	Proportional Compact Remote DLX-CR400
ARQ2232	Proportional Compact Remote Low Force DLX-CR400LF
ARQ2213	Proportional Compact Centre Mount swing away to RH (Nucleus)
ARQ2214	Proportional Compact Centre Mount swing away to LH (Nucleus)
ARZ2215	Individual mounting option per request (Special)
ARQ2230	Head Control with ultrasonic proximity sensors
ARQ2247	Head Control with ultrasonic proximity sensors
ARZ2233	MEC Proportional Joystick
ARZ2213	Centre Mount, swing away to RH, Nucleus
ARZ2214	Centre Mount, swing away to LH, Nucleus
ARZ2215	Individual mounting option per request (Special)
ARZ2240	Proportional Heavy Duty Remote Joystick with mode switch
ARZ2213	Centre Mount, swing away to RH, Nucleus

ARZ2214	Centre Mount, swing away to LH, Nucleus
ARZ2215	Individual mounting option per request (Special)
ARQ2227	10-Way Switch Bracket RH, Mounted for Attendant
ARQ2228	10-Way Switch Bracket LH, Mounted for Attendant
ARQ2253	USB Charger Port DLX-USB02-A
ARZ2258	Remote Stop Switch (ASL504)
ARQ1282	6 km/h (4mph) (4 pole)
ARQ1908	10 km/h (6mph) (includes light system and indicators) G-Trac recommended
ARQ1909	12 km/h (8mph) (includes light system and indicators) G-Trac recommended
ARQ1329	Lapbelt (adjustable both sides)
ARQ2113	Lapbelt (adjustable both sides)
ARQ1332	H-Belt
ARQ2701	Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size L (padded length 620mm)
ARQ2702	Bodypoint Evoflex Pelvic, Stabilizer, Rehab Latch size M (padded length 460mm)
ARQ2703	Bodypoint 2-Point Hip Belt, Centre pull - push button
ARQ2704	Bodypoint 2-Point Hip Belt, Centre pull - push button
ARQ2705	Bodypoint 2-Point Hip Belt, Dual pull Side release
ARQ2706	Bodypoint 2-Point Hip Belt, Dual pull Side release
ARQ2128	Bag Hooks (Max load up to 6kg = impact on max user weight)



For a more complex customer requirement please refer to our Invacare Unique special order form. email [specialsuk@invacare.com](mailto:specialsuk@invacare.com) or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: \_\_\_\_\_ Email: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Description of special request (please include relevant option code or part number):