



Yes, you can:

Tel: 01656 776222 Fax: 01656 776220
 Email: ordersuk@invacare.com
 Online Spares available at: www.invacare.co.uk

Power Retail Prescription Form

SEPTEMBER 2024
 Customer Ref:

| | |
|--------------------------|--------------------|
| Customer Ref* | Account No* |
| Contact No* | Order No* |
| Delivery Address* | Order Qty* |
| | Order date* |

TDX SP2 NB

Standard Features

- Max User Weight 150kg (23.5 Stone)
- Crash test ISO7176 Part-19
- Compact, manoeuvrable, CWD TDX SP2 Narrow Base
- Telescopic adjustable seat widths from 330mm - 530mm (13" - 21") with an additional 50mm, (2") by arms
- Seat depths 410mm - 560mm stepless adjustment (16" - 22")
- Tension adjustable backrest
- Standard T-Armrests, detachable with short armpad
- Manual settable seat Tilt 0° - 20°
- Manual settable backrest Recline 0° - 30°, in 3.75° steps
- Standard 80° hangers with fixed footplate
- LiNX control REM110 with standard bracket
- Selection of 10 stunning shroud and Rim colours
- Centre (driving wheels) 355mm x 76mm (14" x 3") Black tyres with Gel inserts (Solid)
- 6 km/h (4mph) 4 Pole motors and 50Ah batteries
- 8A (200 - 250V) charger
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)

Standard Base price: £5,370



Please make your selection with a cross in the box ○ Free of charge ◆ standard feature POA special request

Customer weight = KG

| Code | Description | |
|----------------------|---|--------|
| BASE PRODUCT* | | |
| ARM2100 | <input type="checkbox"/> TDX SP2 NARROW BASE | £5,370 |
| USER WEIGHT* | | |
| ARM1920 | <input type="checkbox"/> User weight up to 125kg | ○ |
| ARM1921 | <input type="checkbox"/> User weight 126kg to 136kg | ○ |
| ARM1922 | <input type="checkbox"/> User weight 137kg to 150kg (NC 1092, 2052, 2053) | ○ |
| BASE MODEL* | | |
| ARM0000 | <input type="checkbox"/> TDX SP2 Narrow Base | ◆ |
| ARZ0001 | <input type="checkbox"/> Special TDX SP2 Narrow Base | POA |

| SEATING TYPE* | | | |
|--------------------------------|--------------------------|---|--|
| ARM2100 | <input type="checkbox"/> | Modulite Seating Unit | ◆ |
| SEAT FRAME* | | | |
| ARM2111 | <input type="checkbox"/> | Telescopic Seat Frame | ◆ |
| SEAT WIDTH* | | | |
| ARZ1001 | <input type="checkbox"/> | Seat width 330mm (13") seat width achieved on a 380mm seat frame (2121) with inrigged armrests (2183) | (NC 2181, 2193, 2192, 2194) ○ |
| ARM2121 | <input type="checkbox"/> | 380mm (15") | (NC 2143, 2701, 2703, 2705) ○ |
| ARM2122 | <input type="checkbox"/> | 430mm (17") | (NC 2703, 2705) ○ |
| ARM2123 | <input type="checkbox"/> | 480mm (19") | (NC 2702, 2704, 2706) ○ |
| ARM2124 | <input type="checkbox"/> | 530mm (21") | (NC 2174, 2175, 2176, 2177, 2178, 2143, 2702, 2704, 2706) ○ |
| Body Fit Kit | | | |
| UNI2106 | <input type="checkbox"/> | Body fit Kit Outrigging Backrest Brackets and Pivot Point (outrig: 1 inch/2,54 cm per side) | (NC 2124, 2192, 2194) ○ |
| SEAT DEPTH RAIL LENGTH* | | | |
| ARM2129 | <input type="checkbox"/> | Standard 410mm - 510mm | (NC 2138) ○ |
| ARM2190 | <input type="checkbox"/> | Extended 510mm - 560mm | (NC 2132, 2134) ○ |
| SEAT DEPTH* | | | |
| ARM2132 | <input type="checkbox"/> | 410mm (16") | (NC 2196, 2197, 2701, 2702, 2703, 2704, 2705, 2706) ○ |
| ARM2134 | <input type="checkbox"/> | 460mm (18") | ○ |
| ARM2136 | <input type="checkbox"/> | 510mm (20") | ○ |
| ARM2138 | <input type="checkbox"/> | 560mm (22") | (NC 2143) ○ |
| SEAT OPTION | | | |
| ARZSWO | <input type="checkbox"/> | H Frame Interface Kit | £446 |
| SEAT CUSHION* | | | |
| ARM2142 | <input type="checkbox"/> | Matrx Flo-Tech Lite | £112 |
| ARM2143 | <input type="checkbox"/> | Matrx Flo-tech Contour Visco (430mm x 430mm only - up to 125kg) | (NC 2136, 2121, 2124, 2138) £112 |
| ARM2147 | <input type="checkbox"/> | Matrx Libra | £391 |
| ARM2140 | <input type="checkbox"/> | Without seat cushion | ○ |
| BACKREST SUPPORT* | | | |
| ARM2163 | <input type="checkbox"/> | Short 480mm | (NC 2172, 2173) ○ |
| ARM2162 | <input type="checkbox"/> | Tall 540mm | ◆ |
| PUSH HANDLES FOR BACKPOSTS* | | | |
| ARM2107 | <input type="checkbox"/> | WITH push handles | ◆ |
| ARM2108 | <input type="checkbox"/> | WITHOUT push handles | (NC 1332, 2219, 2220, 2227, 2228, 2245) £36 |
| BACKREST UPHOLSTERY / CUSHION* | | | |
| ARM2171 | <input type="checkbox"/> | Tension Adjustable Backrest | (NC 1445, 1496, 2155, 2230, 2247, 1288, 1289, 2052, 2053) ◆ |
| ARM2179 | <input type="checkbox"/> | Flex3 Tension Adjustable | (NC 2151, 2192, 2193, 2194, 2199, 2155, 2230, 2247, 2128) £139 |
| ARM2172 | <input type="checkbox"/> | Flex3 Backrest plates without cushion | (NC 2163, 2128) £224 |
| ARM2173 | <input type="checkbox"/> | Flex3 Backrest plates with Laguna cushion | (NC 2163, 2128) £363 |
| ARM2174 | <input type="checkbox"/> | Matrx PB 400mm | (NC 1496, 2124) £496 |
| ARM2175 | <input type="checkbox"/> | Matrx PB 500mm | (NC 1496, 2124) £496 |
| ARM2176 | <input type="checkbox"/> | Matrx Elite E2 Back , height 500 mm (Mesh fabric) | (NC 1496, 2124) £496 |
| ARM2177 | <input type="checkbox"/> | Matrx Elite TR | (NC 1496, 2124) £496 |
| ARM2178 | <input type="checkbox"/> | Matrx Elite Deep E2 Back, height 500 mm (Mesh fabric) | (NC 1496, 2124) £550 |
| ARM2170 | <input type="checkbox"/> | Without Backrest | (NC 1496, 2230, 2247) ○ |

| UPHOLSTERY TYPE* | | | |
|--|--|--|--------|
| ARM2151 | | Black fabric | ◆ |
| ARM2152 | | Black multi stretch PU coated fabric (NC 2171) | ○ |
| ARM2155 | | Black Spacer fabric (NC 2171, 2179) | ○ |
| ARZ1002 | | Without Backrest Upholstery (NC 2179, 2172, 2173, 2174, 2175, 2176, 2177, 2178) | ○ |
| ARMREST* | | | |
| ARM2191 | | T-Armrest standard | ◆ |
| ARM2181 | | T-Armrest extended height | £36 |
| ARM2183 | | In-rigged armrests (NC 2185, 2186, 2187, 2188, 2189) | £25 |
| ARM2193 | | Flip-up armrests (NC 1445, 2179) | £58 |
| ARM2192 | | Following 230mm - 300mm (NC 1445, 2179) | £166 |
| ARM2194 | | Following 300mm - 360mm (NC 1445, 2179) | £166 |
| ARMPAD* | | | |
| ARM2195 | | Short 325mm Armpad | ◆ |
| ARM2196 | | Long 375mm Armpad | ○ |
| ARM2197 | | Long ergonomic 397mm Armpad (NC 2230, 2247) | ○ |
| ARM2199 | | Soft Armpad (NC 2179) | £86 |
| HIP SUPPORTS LH | | | |
| ARM2184 | | Small 140mm x 100mm (NC 2183, 2701, 2702, 2703, 2704, 2705, 2706) | £226 |
| ARM2186 | | Medium 230mm x 100mm (NC 2183, 2701, 2702, 2703, 2704, 2705, 2706) | £226 |
| ARM2188 | | Large 330mm x 100mm (NC 2183, 2701, 2702, 2703, 2704, 2705, 2706) | £226 |
| HIP SUPPORTS RH | | | |
| ARM2185 | | Small 140mm x 100mm (NC 2183, 2701, 2702, 2703, 2704, 2705, 2706) | £226 |
| ARM2187 | | Medium 230mm x 100mm (NC 2183, 2701, 2702, 2703, 2704, 2705, 2706) | £226 |
| ARM2189 | | Large 330mm x 100mm (NC 2183, 2701, 2702, 2703, 2704, 2705, 2706) | £226 |
| MOUNTING HARDWARE FOR HIP SUPPORTS | | | |
| ARM2165 | | Extended height (+2") | ○ |
| HEADREST / NECKREST - Adequate headrest option must be selected* | | | |
| ARM1350 | | Neckrest height and depth adjustable (NC 2247) | £72 |
| ARM1432 | | Headrest, height, depth and sideways adjustable (NC 2247) | £203 |
| ARM1427 | | Headrest with cheek support, height and depth adjustable (NC 2247) | £203 |
| ARM2324 | | Elan standard pad 255mm (10") (NC 2247) | £294 |
| ARM2325 | | Elan occipital pad 230mm x 180mm (9" X 7") (NC 2247) | £294 |
| ARM2327 | | Elan 4 point headrest pad 280mm x 255mm (11" X 10") (NC 2247) | £396 |
| ARM1496 | | All bracketry as per option 1432 (NC 2171, 2174, 2175, 2176, 2177, 2178) | ○ |
| ARM1495 | | Push bar 22mm | ○ |
| LATERAL SUPPORT | | | |
| ARM1445 | | Multi-adjustable pad (NC 2171, 2174, 2175, 2176, 2177, 2178, 2192, 2193, 2194) | £150 |
| SEAT LIFTER / SEAT TILT* | | | |
| ARM1091 | | Manually settable Seat Tilt (0° - 20°) Seat Height 450mm | ◆ |
| ARM1092 | | Electric Tilt (0° - 20°) (NC 2201) | £362 |
| ARM1095 | | Electric Tilt (0° - 20°) Max User weight 150kg (NC 2201) | £362 |
| ARM1093 | | Electric Tilt (0° - 30°) with CoG shift Seat Height 460mm (NC 0354, 0364, 2201) | £961 |
| ARM2052 | | Electric Tilt (0°-50°) with COG shift Seat height 460mm | £1,322 |
| ARM1094 | | Seat lifter with Electrical Tilt (300mm, 30°) with CoG shift Seat Height 460mm - 490mm (NC 0354, 0364, 2201) | £1,696 |
| ARM2053 | | Seat lifter with Electrical Tilt (300mm, 45°) with CoG shift Seat Height 460mm | £2,359 |

| SEAT HEIGHT SETTING* | | | |
|---|--|---|--------|
| ARM0354 | | Seat height 450mm (17.75") | ○ |
| ARZ1003 | | Seat height set at 460mm (18") (from floor to front of seat plate at 0° tilt) For 30° tilt and Lifter (NC 1092, 1095) | ○ |
| ARM0364 | | Seat height 475mm (18.75") | ○ |
| SEAT TILT OPTION | | | |
| ARM2054 | | 5° fixed anterior Tilt assist mount (NC 1092, 1095) | £290 |
| BACKREST RECLINER* | | | |
| ARM0975 | | Manual settable Recliner | ◆ |
| ARM0977 | | Powered Recliner (NC 2201) | £453 |
| CENTRE MOUNT LEGREST* | | | |
| ARM0556 | | Centre mounted : Manual (NC 0517, 0512, 0504, 0513) | £334 |
| ARM0589 | | Centre mounted : LNX Powered (NC 2201, 0504, 0512, 0504, 0513) | £1,470 |
| UNI0592 | | LNX Power centre mount legrest with telescoping foot platform (lowers to the ground) inc calf pads: Length 380-419mm (Selectable length depends on STF height) Custom-made, requires customised setting (NC 2201, 0504, 0512, 0504) | £2,417 |
| ARX9260 | | Without Centre mount Legrest | ○ |
| FOOTREST HANGER TYPE* | | | |
| ARM0526 | | Side mounted : Standard 80° hanger | ◆ |
| ARM0503 | | Side mounted : VARI F | £107 |
| ARM0533 | | Side mounted : VARI A | £210 |
| ARM0531 | | Side mounted : ADE Powered (NC 2201) | £600 |
| ARM0501 | | Without Hanger (NC 0504, 0513) | ○ |
| AMPUTEE SUPPORT | | | |
| ARM0555 | | LH Amputee Support | £174 |
| ARM0554 | | RH Amputee Support | £174 |
| CENTRE MOUNT FOOTPLATES / FOOT PLATFORM* | | | |
| ARM2365 | | Angle Adjustable Footplates Small 127mm x 190mm (5" x 7.5") for Centre Mounted Hanger only (NC 0503, 0533, 0531) | £93 |
| ARM0559 | | Angle Adjustable Footplate Medium 140mm x 240mm (5.5" x 9.5") for Centre Mounted Hanger only (NC 0503, 0533, 0531) | ○ |
| ARM0516 | | Angle Adjustable Footplates Large 150mm x 290mm (6" x 11.5") for Centre Mounted Hanger only (NC 0503, 0533, 0531) | £121 |
| ARM0517 | | Angle Adjustable foot platform Medium 305mm x 275mm (12" x 10.75") for Centre Mounted Hanger only (NC 0526, 0503, 0533, 0531) | £243 |
| UNI0559 | | Angle & width adjustable footplate (240mm depth) (NC 0526, 0503, 0533, 0531, 0556, 0589) | ○ |
| UNI0516 | | Angle & width adjustable footplates (290mm depth) (NC 0526, 0503, 0533, 0531, 0503, 0556, 0589) | £121 |
| UNI0517 | | Angle adjustable foot platform (240mm depth, 305mm width) (NC 0526, 0503, 0533, 0531, 0504, 0556, 0589) | £243 |
| ARX9259 | | Without Centre mount Footplate | ○ |
| FOOTPLATE TYPE Side Mounted* | | | |
| ARM0512 | | Fixed Footplate (NC 0533, 0531) | ◆ |
| ARM0504 | | Angle adjustable Footplate (NC 0526) | £64 |
| ARM0513 | | Angle and Depth adjustable | ○ |
| ARM1004 | | Without Footplates (NC 0503, 0533, 0531) | ○ |
| CONTROLS* | | | |
| ARM2201 | | REM110 (NC 1437, 1092, 1095, 1093, 1094, 0977, 0589, 0531, 2231, 2232, 2230, 2247, 2257, 2233, 2240, 2219, 2245, 2270, 1908, 2209, 2210) | ◆ |
| ARM2203 | | REM211 (NC 1437, 2231, 2232, 2230, 2247, 2257, 2233, 2240, 2270, 1908, 2209, 2210) | ○ |
| ARM2204 | | REM216 (NC 2231, 2232, 2230, 2247, 2257, 2233, 2240, 2270, 2209, 2210) | £137 |
| ARM2205 | | REM400 (NC 2209, 2210) | £305 |
| ARM2206 | | REM500 (NC 1708, 2213, 2214, 2227, 2228, 2220, 2270) | £695 |

| CONTROL BRACKETRY* | | | |
|---|--|---|--------|
| ARM1207 | | Standard Bracket RH - Not Compatible with REM500 | ○ |
| ARM1208 | | Standard Bracket LH - Not Compatible with REM500 | ○ |
| ARM1707 | | Swing away RH Bracket | £119 |
| ARM1708 | | Swing away LH Bracket | £119 |
| ARM1709 | | Matrx Resolve Swing in/out bracket RH (NC 2213, 2275, 2235, 2213, 2241, 1404) | £219 |
| ARM1710 | | Matrx Resolve Swing in/out bracket LH (NC 2214, 2276, 2236, 2214, 2242, 1405) | £219 |
| ARM2213 | | Centre mount, swing away to RH Nucleus (NC 1405, 1415) | £446 |
| ARM2214 | | Centre mount, swing away to LH Nucleus (NC 1405, 1339) | £446 |
| ARM2227 | | RH Bracket mounted for Attendant (NC 2220, 2108, 1286, 2128) | £57 |
| ARM2228 | | LH Bracket mounted for Attendant (NC 2220, 2108, 1286, 2128) | £57 |
| ARM2209 | | REM500 Swing away RH multi-adjustable Bracket (NC 1415) | ○ |
| ARM2210 | | REM500 Swing away LH multi-adjustable Bracket (NC 1339) | ○ |
| ARZ2239 | | Individual mounting option per request (Special) | POA |
| REM400 TOGGLE OPTION | | | |
| ARM2270 | | Assignable toggle switch module (NC 2203, 2204, 2206) | £284 |
| DUAL CONTROL | | | |
| ARM2220 | | Dual control REM050 (NC 2108, 2205, 2206, 1221, 1222, 2230, 2247, 2257, 2233, 2240) | £355 |
| ARM2219 | | Compact (ACU200) (NC 1222, 2108, 2201) | £355 |
| ARM2245 | | IDC Dual control (NC 2201) | £673 |
| DUAL CONTROL BRACKETRY | | | |
| ARM1207 | | Standard Bracket RH | ○ |
| ARM1208 | | Standard Bracket LH | ○ |
| ARM1221 | | IDC Dual control detachable Bracket RH (NC 2220, 2227) | ○ |
| ARM1222 | | IDC Dual control detachable Bracket LH (NC 2220, 2228) | ○ |
| LINX COMPACT REMOTE & BRACKETS (Requires REM400 or REM500) | | | |
| ARM2231 | | Proportional Compact Remote (NC 2203, 2204) | £740 |
| ARM2232 | | Proportional Compact Remote low force (NC 2203, 2204) | £990 |
| ARM2207 | | RH mount Bracket (NC 2201, 2203, 2204) | ○ |
| ARM2208 | | LH mount Bracket (NC 2201, 2203, 2204) | ○ |
| ARM2213 | | Centre mount, swing away to RH Nucleus Bracket (NC 1405, 1415, 2201, 2203, 2204) | £446 |
| ARM2214 | | Centre mount, swing away to LH Nucleus Bracket (NC 1405, 1339, 2210, 2203, 2204) | £446 |
| ARZ2215 | | Individual mounting option per request (Special) | POA |
| CHIN CONTROL & BRACKETS (Requires REM400 or REM500) | | | |
| UNI9250 | | Mo-vis Micro joystick (NC 9251) | POA |
| UNI9251 | | Mo-vis Multi joystick | POA |
| ARM1289 | | Swing away RH Manual Bracket (NC 1432, 1427, 1496, 1495, 2171, 2201, 2203, 2204, 1288, 1287, 1286) | ○ |
| ARM1288 | | Swing away LH Manual Bracket (NC 1432, 1427, 1496, 1495, 2171, 2201, 2203, 2204, 1289, 1287, 1286) | ○ |
| ARM1287 | | Swing away RH Powered Bracket (NC 1404, 1415, 1312, 2201, 2203, 2204, 1289, 1288, 1286) | £2,251 |
| ARM1286 | | Swing away LH Powered Bracket (NC 1339, 1363, 2201, 2203, 2204, 1289) | £2,251 |
| LINX HEAD CONTROL ATOM TECHNOLOGY (Requires REM400 or REM500) | | | |
| ARM2230 | | ATOM with ultrasonic proximity sensor (NC 2171, 2179, 2197, 2201, 2203, 2204, 2220, 2233, 2324, 2325, 2327, 1350, 1432, 1427, 2257, 2240) | £3,835 |
| ARM2247 | | PROTON with ultrasonic proximity sensor (NC 2171, 2179, 2197, 2201, 2203, 2204, 2220, 2233, 2324, 2325, 2327, 1350, 1432, 1427, 2257, 2240) | £4,491 |
| FOUR-SWITCH PROXIMITY ARRAY (with ATOM TECHNOLOGY) | | | |
| ARM2257 | | Four-switch Array (ASL106) (NC 2203, 2204, 2220, 2230, 2247, 2233, 2240, 2128) | £5,938 |
| ARM2275 | | Centre mount, swing-away to RH nucleus (NC 2203, 2204, 1404, 1405, 1415) | £446 |
| ARM2276 | | Centre mount, swing-away to LH nucleus (NC 2203, 2204, 1404, 1405, 1339) | £446 |
| ARZ2215 | | Individual mounting option per request (Special) | POA |

| WIRELESS ACCESSORIES FOR ATOM | | | |
|--|--|---|--|
| ARM2249 | | Wireless triple switch receiver (ASL557-3) | £305 |
| ARM2250 | | Wireless mouse emulator (ASL558) | £733 |
| LINX PROPORTIONAL REMOTE / INPUTS & BRACKETS (Requires REM400 or REM500) | | | |
| ARM2233 | | MEC proportional joystick | (NC 2203, 2204, 2220, 2230, 2247, 2257, 2240, 2128) £3,723 |
| ARM2240 | | Proportional heavy duty with mode switch (ASL133) | (NC 2203, 2204, 2220, 2230, 2247, 2257, 2233, 2128) £2,076 |
| ARM2235 | | MEC tray mount RH (ASL603) | (NC 1415, 2201, 2203, 2204) ○ |
| ARM2236 | | MEC tray mount LH (ASL603) | (NC 1339, 2201, 2203, 2204) ○ |
| ARM2213 | | Centre mount, swing away to RH nucleus | (NC 1405, 1415, 2201, 2203, 2204) £446 |
| ARM2214 | | Centre mount, swing away to LH nucleus | (NC 1405, 1339, 2201, 2203, 2204) £446 |
| ARM2241 | | Centre mount, swing away to RH nucleus | (NC 1405, 1415, 2201, 2203, 2204) £446 |
| ARM2242 | | Centre mount, swing away to LH nucleus | (NC 1405, 1339, 2201, 2203, 2204) £446 |
| ARZ2215 | | Individual mounting option per request (Special) | POA |
| LINX 10-WAY SWITCH | | | |
| ARM2251 | | DLX-SW10-A | £453 |
| ARM2217 | | 10-Way Switch Bracket RH, Mounted for User | (NC 2213, 2275, 2235, 2241, 1415) ○ |
| ARM2218 | | 10-Way Switch Bracket LH, Mounted for User | (NC 2214, 2276, 2236, 2242, 1339) ○ |
| ARM2227 | | 10-Way Switch Bracket RH, Mounted for Attendant | (NC 2128, 1287, 1221) £57 |
| ARM2228 | | 10-Way Switch Bracket LH, Mounted for Attendant | (NC 2128, 1286, 1222) £57 |
| LINX ACCESSORY | | | |
| ARM2252 | | Loom providing jack sockets for drive and mode | £90 |
| LINX USB CHARGER PORT | | | |
| ARM2253 | | USB charger port DLX-USB02-A | (NC 2701, 2702, 2703, 2704, 2705, 2706) £143 |
| ARM2207 | | RH mount | ○ |
| ARM2208 | | LH mount | ○ |
| REMOTE STOP SWITCH | | | |
| ARM2258 | | Remote Stop Switch (ASL504) | (NC 2203, 2204) £942 |
| LINX INPUT MODULE | | | |
| UNI2273 | | IN200 Module (Jack port for external power button) | (NC 2220) £242 |
| UNI9252 | | IN500 Module (for DB9 Drive input) | £785 |
| FURTHER LINX COMPONENTS / ACCESSORIES | | | |
| ARM2300 | | To be enclosed only | ○ |
| ARZ2301 | | To be connected and programmed (special) | POA |
| BUDDY BUTTONS | | | |
| ARM1685 | | Red | £79 |
| ARM1686 | | Yellow | £79 |
| ARM1687 | | Green | £79 |
| ARM1690 | | Blue | £79 |
| EGG SWITCH | | | |
| ARM2254 | | Egg Switch, Black | £79 |
| PIKO BUTTONS | | | |
| ARM1711 | | Red | £142 |
| ARM1712 | | Yellow | £142 |
| ARM1713 | | Green | £142 |
| ARM1714 | | Orange | £142 |
| ARM1715 | | Blue | £142 |
| ARM1716 | | Black | £142 |

| TOGGLES / BUTTONS | | | |
|---------------------------------|--|---|----------------------------|
| ARZ2261 | | Stereo Toggle | £242 |
| ARZ2262 | | Stereo Button | £230 |
| ARZ2263 | | Four-way Toggle | £300 |
| ARZ2264 | | Four-way Button | £321 |
| ARZ2207 | | RH mount | ○ |
| ARZ2208 | | LH mount | ○ |
| ARZ2225 | | Without mounting | ○ |
| BODYPOINT JOYSTICK KNOBS | | | |
| ARM2721 | | 3" U-shaped, joystick handle | £65 |
| ARM2722 | | 4" U-shaped, joystick handle | £65 |
| ARM2723 | | Domed rubber, joystick handle | £65 |
| ARM2724 | | Mushroom grooved, joystick handle | £65 |
| SHROUD COLOUR* | | | |
| ARM2640 | | Obsidian Black Frame Shrouds | ○ |
| ARM2641 | | Ruby Red Frame Shrouds | ○ |
| ARM2642 | | Purple Passion Frame Shrouds | ○ |
| ARM2643 | | Peacock Blue Frame Shrouds | ○ |
| ARM2644 | | Arctic Blue Frame Shrouds | ○ |
| RIM COLOUR* | | | |
| ARM2650 | | Obsidian Black Rim Inserts | ○ |
| ARM2651 | | Ruby Red Rim Inserts | ○ |
| ARM2652 | | Purple Passion Rim Inserts | ○ |
| ARM2653 | | Peacock Blue Rim Inserts | ○ |
| ARM2654 | | Arctic Blue Rim Inserts | ○ |
| FORKS* | | | |
| ARM2631 | | Single sided | £316 |
| ARM2632 | | Double sided | ◆ |
| CENTRE WHEELS (DRIVING WHEELS)* | | | |
| ARM0704 | | DW 14" (355 mm) x 3" (75 mm), Pneumatic | ◆ |
| ARM0705 | | DW 14" (355 mm) x 3" (75 mm), Solid inlay, dual density foam insert | ○ |
| MOTOR TYPE* | | | |
| ARM1294 | | 3 km/h (2mph) | ○ |
| ARM1282 | | 6 km/h (4mph) | ◆ |
| ARM1908 | | 10 km/h (6mph) | (NC 2203, 1437, 1254) £464 |
| BATTERY* | | | |
| ARM1257 | | 50 Ah Gel | ◆ |
| PERFORMANCE ACCESSORIES | | | |
| ARM1215 | | G-Trac | £738 |
| CHARGER* | | | |
| ARM1343 | | OptiCharge 8A (200 - 250V) | ◆ |
| ARM1313 | | OptiCharge 10A (200 - 250V) | £137 |
| LIGHTS* | | | |
| ARM1437 | | LED headlights & indicators | (NC 2203, 1908) £407 |
| ARZ1006 | | Lights Included with Motor Package 10 & 12kmh | (NC 1282) ○ |
| ARM1254 | | Without light system | ◆ |

| TRANSPORTATION KIT* | | | |
|---------------------|--|---|--|
| ARM1450 | | Preparation for DAHL docking station | £524 |
| ARM1494 | | 4-Point tie-down | ◆ |
| FIXATION BELT* | | | |
| ARM1329 | | Lap Belt fitted to base of backrest | (NC 2702, 2703, 2704, 2705, 2706) ◆ |
| ARM2113 | | Lap Belt adjustable mounting on seat rail extrusion | (NC 2701, 2702, 2703, 2704, 2705, 2706) ○ |
| ARM1332 | | H-Belt | (NC 2108) £66 |
| ARM2701 | | Evoflex pelvic stabilizer, size L | (NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113) £84 |
| ARM2702 | | Evoflex pelvic stabilizer, size M | (NC 2123, 2124, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113) £84 |
| ARM2703 | | 2-Point hip belt centre pull, size L | (NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113) £78 |
| ARM2704 | | 2-Point hip belt centre pull, size M | (NC 2123, 2124, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113) £78 |
| ARM2705 | | 2-Point hip belt dual pull, size L | (NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113) £78 |
| ARM2706 | | 2-Point Hip Belt Dual pull, size M | (NC 2123, 2124, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113) £78 |
| REAR VIEW MIRROR | | | |
| ARM1415 | | RH Mirror | (NC 2213, 1287, 2217, 2275, 2235, 2241) £38 |
| ARM1339 | | LH Mirror | (NC 2214, 1286, 2218, 2276, 2236, 2242) £38 |
| OTHER OPTIONS | | | |
| ARM2128 | | Bag hooks | (NC 2173, 2227, 2228, 2257, 2233, 2240) £69 |
| ARM1312 | | RH Cane holder | £38 |
| ARM1363 | | LH Cane holder | £38 |
| ARM1404 | | Swing away RH Tray | (NC 1707, 2209, 2213, 2214, 2235, 2236, 2241, 2242, 1287, 2207, 2275, 2276) £174 |
| ARM1405 | | Swing away LH Tray | (NC 1708, 2210, 2213, 2214, 2235, 2236, 2241, 2242, 1286, 2208, 2275, 2276) £174 |

| | |
|--------------------|--|
| TOTAL PRICE | |
|--------------------|--|

Non Conforming Options:

| | |
|---------|--|
| ARM2121 | 380mm (15") |
| ARM2122 | 430mm (17") |
| ARM2123 | 480mm (19") |
| ARM2124 | 530mm (21") |
| ARM2132 | 410mm (16") |
| ARM2134 | 460mm (18") |
| ARM2136 | 510mm (20") |
| ARM2138 | 560mm (22") |
| ARM2143 | Matrix Flo-tech Contour Visco (430mm x 430mm only - up to 125kg) |
| ARM2163 | Short 480mm |
| ARM2108 | WITHOUT push handles |
| ARM2171 | Tension Adjustable Backrest |
| ARM2179 | Flex3 Tension Adjustable |
| ARM2172 | Flex3 Backrest plates without cushion |
| ARM2173 | Flex3 Backrest plates with Laguna cushion |
| ARM2174 | Matrix PB 400mm |
| ARM2175 | Matrix PB 500mm |
| ARM2176 | Matrix Elite E2 Back , height 500 mm (Mesh fabric) |
| ARM2177 | Matrix Elite TR |
| ARM2178 | Matrix Elite Deep E2 Back, height 500 mm (Mesh fabric) |
| ARM2151 | Black fabric |
| ARM2155 | Black Spacer fabric |
| ARM2181 | T-Armrest extended height |
| ARM2183 | In-rigged armrests |
| ARM2193 | Flip-up armrests |
| ARM2192 | Following 230mm - 300mm |
| ARM2194 | Following 300mm - 360mm |
| ARM2196 | Long 375mm Armpad |
| ARM2197 | Long ergonomic 397mm Armpad |
| ARM2199 | Soft Armpad |
| ARM2184 | Small 140mm x 100mm |
| ARM2186 | Medium 230mm x 100mm |
| ARM2188 | Large 330mm x 100mm |
| ARM2185 | Small 140mm x 100mm |
| ARM2187 | Medium 230mm x 100mm |
| ARM2189 | Large 330mm x 100mm |
| ARM1350 | Neckrest height and depth adjustable |
| ARM1432 | Headrest, height, depth and sideways adjustable |
| ARM1427 | Headrest with cheek support, height and depth adjustable |
| ARM2324 | Elan standard pad 255mm (10") |
| ARM2325 | Elan occipital pad 230mm x 180mm (9" X 7") |
| ARM2327 | Elan 4 point headrest pad 280mm x 255mm (11" X 10") |
| ARM1496 | All bracketry as per option 1432 |

| | |
|---------|---|
| ARM1495 | Push bar 22mm |
| ARM1445 | Multi-adjustable pad |
| ARM1092 | Electric Tilt (0°- 20°) |
| ARM1095 | Electric Tilt (0°- 20°) Max User weight 150kg |
| ARM1093 | Electric Tilt (0° - 30°) with CoG shift Seat Height 460mm |
| ARM2052 | Electric Tilt (0°-50°) with COG shift Seat height 460mm |
| ARM1094 | Seat lifter with Electrical Tilt (300mm, 30°) with CoG shift Seat Height 460mm - 490mm |
| ARM2053 | Seat lifter with Electrical Tilt (300mm, 45°) with CoG shift Seat Height 460mm |
| ARM0354 | Seat height 450mm (17.75") |
| ARM0364 | Seat height 475mm (18.75") |
| ARM0977 | Powered Recliner |
| ARM0556 | Centre mounted : Manual |
| ARM0589 | Centre mounted : LNX Powered |
| ARM0526 | Side mounted : Standard 80° hanger |
| ARM0503 | Side mounted : VARI F |
| ARM0533 | Side mounted : VARI A |
| ARM0531 | Side mounted : ADE Powered |
| ARM0517 | Angle Adjustable foot platform Medium 305mm x 275mm (12" x 10.75") for Centre Mounted Hanger only |
| UNI0517 | Angle adjustable foot platform (240mm depth, 305mm width) |
| ARM0512 | Fixed Footplate |
| ARM0504 | Angle adjustable Footplate |
| ARM0513 | Angle and Depth adjustable |
| ARM2201 | REM110 |
| ARM2203 | REM211 |
| ARM2204 | REM216 |
| ARM2205 | REM400 |
| ARM2206 | REM500 |
| ARM1707 | Swing away RH Bracket |
| ARM1708 | Swing away LH Bracket |
| ARM2213 | Centre mount, swing away to RH Nucleus |
| ARM2214 | Centre mount, swing away to LH Nucleus |
| ARM2227 | RH Bracket mounted for Attendant |
| ARM2228 | LH Bracket mounted for Attendant |
| ARM2209 | REM500 Swing away RH multi-adjustable Bracket |
| ARM2210 | REM500 Swing away LH multi-adjustable Bracket |
| ARM2270 | Assignable toggle switch module |
| ARM2220 | Dual control REM050 |
| ARM2219 | Compact (ACU200) |
| ARM2245 | IDC Dual control |
| ARM1221 | IDC Dual control detachable Bracket RH |
| ARM1222 | IDC Dual control detachable Bracket LH |
| ARM2231 | Proportional Compact Remote |
| ARM2232 | Proportional Compact Remote low force |
| ARM2207 | RH mount Bracket |

| | |
|---------|---|
| ARM2208 | LH mount Bracket |
| ARM2213 | Centre mount, swing away to RH Nucleus Bracket |
| ARM2214 | Centre mount, swing away to LH Nucleus Bracket |
| UNI9251 | Mo-vis Multi joystick |
| ARM1289 | Swing away RH Manual Bracket |
| ARM1288 | Swing away LH Manual Bracket |
| ARM1287 | Swing away RH Powered Bracket |
| ARM1286 | Swing away LH Powered Bracket |
| ARM2230 | ATOM with ultrasonic proximity sensor |
| ARM2247 | PROTON with ultrasonic proximity sensor |
| ARM2257 | Four-switch Array (ASL106) |
| ARM2275 | Centre mount, swing-away to RH nucleus |
| ARM2276 | Centre mount, swing-away to LH nucleus |
| ARM2233 | MEC proportional joystick |
| ARM2240 | Proportional heavy duty with mode switch (ASL133) |
| ARM2235 | MEC tray mount RH (ASL603) |
| ARM2236 | MEC tray mount LH (ASL603) |
| ARM2213 | Centre mount, swing away to RH nucleus |
| ARM2214 | Centre mount, swing away to LH nucleus |
| ARM2241 | Centre mount, swing away to RH nucleus |
| ARM2242 | Centre mount, swing away to LH nucleus |
| ARM2217 | 10-Way Switch Bracket RH, Mounted for User |
| ARM2218 | 10-Way Switch Bracket LH, Mounted for User |
| ARM2227 | 10-Way Switch Bracket RH, Mounted for Attendant |
| ARM2228 | 10-Way Switch Bracket LH, Mounted for Attendant |
| ARM2253 | USB charger port DLX-USB02-A |
| ARM2207 | RH mount |
| ARM2208 | LH mount |
| ARZ2207 | RH mount |
| ARZ2208 | LH mount |
| ARM1282 | 6 km/h (4mph) |
| ARM1908 | 10 km/h (6mph) |
| ARM1437 | LED headlights & indicators |
| ARM1254 | Without light system |
| ARM1329 | Lap Belt fitted to base of backrest |
| ARM2113 | Lap Belt adjustable mounting on seat rail extrusion |
| ARM1332 | H-Belt |
| ARM2701 | Evoflex pelvic stabilizer, size L |
| ARM2702 | Evoflex pelvic stabilizer, size M |
| ARM2703 | 2-Point hip belt centre pull, size L |
| ARM2704 | 2-Point hip belt centre pull, size M |
| ARM2705 | 2-Point hip belt dual pull, size L |
| ARM2706 | 2-Point Hip Belt Dual pull, size M |
| ARM1415 | RH Mirror |

| | |
|---------|--------------------|
| ARM1339 | LH Mirror |
| ARM2128 | Bag hooks |
| ARM1312 | RH Cane holder |
| ARM1363 | LH Cane holder |
| ARM1404 | Swing away RH Tray |
| ARM1405 | Swing away LH Tray |



For a more complex customer requirement please refer to our Invacare Unique special order form. email specialsuk@invacare.com or Fax: 01656 776220

Special Instructions

Please supply the following information:

Name: _____ Email: _____

Phone: _____ Fax: _____

Description of special request (please include relevant option code or part number):