



**Yes, you can:**

Tel: 01656 776222  
 Email: ordersuk@invacare.com  
 Online Spares available at: www.invacare.co.uk

# Power Retail Prescription Form

JANUARY 2025  
 Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## TDX SP2 NB

### Standard Features

- Max User Weight 150kg (23.5 Stone)
- Crash test ISO7176 Part-19
- Compact, manoeuvrable, CWD TDX SP2 Narrow Base
- Telescopic adjustable seat widths from 330mm - 530mm (13" - 21") with an additional 50mm, (2") by arms
- Seat depths 410mm - 560mm stepless adjustment (16" - 22")
- Tension adjustable backrest
- Standard T-Armrests, detachable with short armpad
- Manual settable seat Tilt 0° - 20°
- Manual settable backrest Recline 0° - 30°, in 3.75° steps
- Standard 80° hangers with fixed footplate
- LiNX control REM110 with standard bracket
- Selection of 10 stunning shroud and Rim colours
- Centre (driving wheels) 355mm x 76mm (14" x 3") Black tyres with Gel inserts (Solid)
- 6 km/h (4mph) 4 Pole motors and 50Ah batteries
- 8A (200 - 250V) charger
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)

**Standard Base price: £5,370**



Please make your selection with a cross in the box    ○ Free of charge    ◆ standard feature    POA special request

Customer weight =                      KG

Code	Description	
<b>BASE PRODUCT*</b>		
ARM2100	<input type="checkbox"/> TDX SP2 NARROW BASE	£5,370
<b>USER WEIGHT*</b>		
ARM1920	<input type="checkbox"/> User weight up to 125kg	○
ARM1921	<input type="checkbox"/> User weight 126kg to 136kg	○
ARM1922	<input type="checkbox"/> User weight 137kg to 150kg (NC 1092, 2052, 2053)	○
<b>BASE MODEL*</b>		
ARM0000	<input type="checkbox"/> TDX SP2 Narrow Base	◆
ARZ0001	<input type="checkbox"/> Special TDX SP2 Narrow Base	POA

SEATING TYPE*			
ARM2100	<input type="checkbox"/>	Modulite Seating Unit	◆
SEAT FRAME*			
ARM2111	<input type="checkbox"/>	Telescopic Seat Frame	◆
SEAT WIDTH*			
ARZ1001	<input type="checkbox"/>	Seat width 330mm (13") seat width achieved on a 380mm seat frame (2121) with inrigged armrests (2183)	(NC 2181, 2193, 2192, 2194) ○
ARM2121	<input type="checkbox"/>	380mm (15")	(NC 2143, 2701, 2703, 2705) ○
ARM2122	<input type="checkbox"/>	430mm (17")	(NC 2703, 2705) ○
ARM2123	<input type="checkbox"/>	480mm (19")	(NC 2702, 2704, 2706) ○
ARM2124	<input type="checkbox"/>	530mm (21")	(NC 2174, 2175, 2176, 2177, 2178, 2143, 2702, 2704, 2706) ○
Body Fit Kit			
UNI2106	<input type="checkbox"/>	Body fit Kit Outrigging Backrest Brackets and Pivot Point (outrig: 1 inch/2,54 cm per side)	(NC 2124, 2192, 2194) ○
SEAT DEPTH RAIL LENGTH*			
ARM2129	<input type="checkbox"/>	Standard 410mm - 510mm	(NC 2138) ○
ARM2190	<input type="checkbox"/>	Extended 510mm - 560mm	(NC 2132, 2134) ○
SEAT DEPTH*			
ARM2132	<input type="checkbox"/>	410mm (16")	(NC 2196, 2197, 2701, 2702, 2703, 2704, 2705, 2706) ○
ARM2134	<input type="checkbox"/>	460mm (18")	○
ARM2136	<input type="checkbox"/>	510mm (20")	○
ARM2138	<input type="checkbox"/>	560mm (22")	(NC 2143) ○
SEAT OPTION			
ARZSWO	<input type="checkbox"/>	H Frame Interface Kit	£446
SEAT CUSHION*			
ARM2142	<input type="checkbox"/>	Matrx Flo-Tech Lite	£112
ARM2143	<input type="checkbox"/>	Matrx Flo-tech Contour Visco (430mm x 430mm only - up to 125kg)	(NC 2136, 2121, 2124, 2138) £112
ARM2147	<input type="checkbox"/>	Matrx Libra	£391
ARM2140	<input type="checkbox"/>	Without seat cushion	○
BACKREST SUPPORT*			
ARM2163	<input type="checkbox"/>	Short 480mm	(NC 2172, 2173) ○
ARM2162	<input type="checkbox"/>	Tall 540mm	◆
PUSH HANDLES FOR BACKPOSTS*			
ARM2107	<input type="checkbox"/>	WITH push handles	◆
ARM2108	<input type="checkbox"/>	WITHOUT push handles	(NC 1332, 2219, 2220, 2227, 2228, 2245) £36
BACKREST UPHOLSTERY / CUSHION*			
ARM2171	<input type="checkbox"/>	Tension Adjustable Backrest	(NC 1445, 1496, 2155, 2230, 2247, 2052, 2053) ◆
ARM2179	<input type="checkbox"/>	Flex3 Tension Adjustable	(NC 2151, 2192, 2193, 2194, 2199, 2155, 2230, 2247, 2128) £139
ARM2172	<input type="checkbox"/>	Flex3 Backrest plates without cushion	(NC 2163, 2128) £224
ARM2173	<input type="checkbox"/>	Flex3 Backrest plates with Laguna cushion	(NC 2163, 2128) £363
ARM2174	<input type="checkbox"/>	Matrx PB 400mm	(NC 1496, 2124) £496
ARM2175	<input type="checkbox"/>	Matrx PB 500mm	(NC 1496, 2124) £496
ARM2176	<input type="checkbox"/>	Matrx Elite E2 Back , height 500 mm (Mesh fabric)	(NC 1496, 2124) £496
ARM2177	<input type="checkbox"/>	Matrx Elite TR	(NC 1496, 2124) £496
ARM2178	<input type="checkbox"/>	Matrx Elite Deep E2 Back, height 500 mm (Mesh fabric)	(NC 1496, 2124) £550
ARM2170	<input type="checkbox"/>	Without Backrest	(NC 1496, 2230, 2247) ○

<b>UPHOLSTERY TYPE*</b>			
ARM2151		Black fabric	◆
ARM2152		Black multi stretch PU coated fabric (NC 2171)	○
ARM2155		Black Spacer fabric (NC 2171, 2179)	○
ARZ1002		Without Backrest Upholstery (NC 2179, 2172, 2173, 2174, 2175, 2176, 2177, 2178)	○
<b>ARMREST*</b>			
ARM2191		T-Armrest standard	◆
ARM2181		T-Armrest extended height	£36
ARM2183		In-rigged armrests (NC 2185, 2186, 2187, 2188, 2189)	£25
ARM2193		Flip-up armrests (NC 1445, 2179)	£58
ARM2192		Following 230mm - 300mm (NC 1445, 2179)	£166
ARM2194		Following 300mm - 360mm (NC 1445, 2179)	£166
<b>ARMPAD*</b>			
ARM2195		Short 325mm Armpad	◆
ARM2196		Long 375mm Armpad	○
ARM2197		Long ergonomic 397mm Armpad (NC 2230, 2247)	○
ARM2199		Soft Armpad (NC 2179)	£86
<b>HIP SUPPORTS LH</b>			
ARM2184		Small 140mm x 100mm (NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	£226
ARM2186		Medium 230mm x 100mm (NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	£226
ARM2188		Large 330mm x 100mm (NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	£226
<b>HIP SUPPORTS RH</b>			
ARM2185		Small 140mm x 100mm (NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	£226
ARM2187		Medium 230mm x 100mm (NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	£226
ARM2189		Large 330mm x 100mm (NC 2183, 2701, 2702, 2703, 2704, 2705, 2706)	£226
<b>MOUNTING HARDWARE FOR HIP SUPPORTS</b>			
ARM2164		Standard height	○
ARM2165		Extended height (+2")	○
<b>HEADREST / NECKREST - Adequate headrest option must be selected*</b>			
ARM1350		Neckrest height and depth adjustable (NC 2247)	£72
ARM1432		Headrest, height, depth and sideways adjustable (NC 2247)	£203
ARM1427		Headrest with cheek support, height and depth adjustable (NC 2247)	£203
ARM2324		Elan standard pad 255mm (10") (NC 2247)	£294
ARM2325		Elan occipital pad 230mm x 180mm (9" X 7") (NC 2247)	£294
ARM2327		Elan 4 point headrest pad 280mm x 255mm (11" X 10") (NC 2247)	£396
ARM1496		All bracketry as per option 1432 (NC 2171, 2174, 2175, 2176, 2177, 2178)	○
ARM1495		Push bar 22mm	○
<b>LATERAL SUPPORT</b>			
ARM1445		Multi-adjustable pad (NC 2171, 2174, 2175, 2176, 2177, 2178, 2192, 2193, 2194)	£150

<b>SEAT LIFTER / SEAT TILT*</b>			
ARM1091		Manually settable Seat Tilt (0° - 20°) Seat Height 450mm	◆
ARM1092		Electric Tilt (0° - 20°)	(NC 2201) £362
ARM1095		Electric Tilt (0° - 20°) Max User weight 150kg	(NC 2201) £362
ARM1093		Electric Tilt (0° - 30°) with CoG shift Seat Height 460mm	(NC 0354, 0364, 2201) £961
ARM2052		Electric Tilt (0°-50°) with COG shift Seat height 460mm (136kg Max User Weight)	(NC 0364, 2201) £1,322
ARM1094		Seat lifter with Electrical Tilt (300mm, 30°) with CoG shift Seat Height 460mm - 490mm	(NC 0354, 0364, 2201) £1,696
ARM2053		Seat lifter with Electrical Tilt (300mm, 45°) with CoG shift Seat Height 460mm (136kg Max User Weight)	(NC 0364, 2201) £2,359
<b>SEAT HEIGHT SETTING*</b>			
ARM0354		Seat height 450mm (17.75")	○
ARZ1003		Seat height set at 460mm (18") (from floor to front of seat plate at 0° tilt) For 30° tilt and Lifter	(NC 1092, 1095) ○
ARM0364		Seat height 475mm (18.75")	○
<b>SEAT TILT OPTION</b>			
ARM2054		5° fixed anterior Tilt assist mount	(NC 1092, 1095) £290
<b>BACKREST RECLINER*</b>			
ARM0975		Manual settable Recliner	◆
ARM0977		Powered Recliner	(NC 2201) £453
<b>CENTRE MOUNT LEGREST*</b>			
ARM0556		Centre mounted : Manual	(NC 0517, 0512, 0504, 0513, 0526, 0503, 0533, 0531, 0555, 0554) £334
ARM0589		Centre mounted : LNX Powered	(NC 2201, 0504, 0512, 0504, 0513, 0526, 0503, 0533, 0531, 0555, 0554) £1,470
UNI0592		LNX Power centre mount legrest with telescoping foot platform (lowers to the ground) inc calf pads: Length 380-419mm (Selectable length depends on STF height) Custom-made, requires customised setting	(NC 2201, 0504, 0512, 0504, 0526, 0503, 0533, 0531, 0555, 0554) POA
ARX9260		Without Centre mount Legrest	○
<b>FOOTREST HANGER TYPE*</b>			
ARM0526		Side mounted : Standard 80° hanger	◆
ARM0503		Side mounted : VARI F	£107
ARM0533		Side mounted : VARI A	£210
ARM0531		Side mounted : ADE Powered	(NC 2201) £600
ARM0501		Without Hanger	(NC 0504, 0513) ○
<b>AMPUTEE SUPPORT</b>			
ARM0555		LH Amputee Support	(NC 0556, 0589, 0592) £174
ARM0554		RH Amputee Support	(NC 0556, 0589, 0592) £174
<b>CENTRE MOUNT FOOTPLATES / FOOT PLATFORM*</b>			
ARM2365		Angle Adjustable Footplates Small 127mm x 190mm (5" x 7.5") for Centre Mounted Hanger only	(NC 0503, 0533, 0531) £93
ARM0559		Angle Adjustable Footplate Medium 140mm x 240mm (5.5" x 9.5") for Centre Mounted Hanger only	(NC 0503, 0533, 0531) ○
ARM0516		Angle Adjustable Footplates Large 150mm x 290mm (6" x 11.5") for Centre Mounted Hanger only	(NC 0503, 0533, 0531) £121
ARM0517		Angle Adjustable foot platform Medium 305mm x 275mm (12" x 10.75") for Centre Mounted Hanger only	(NC 0526, 0503, 0533, 0531) £243
UNI0559		Angle & width adjustable footplate (240mm depth)	(NC 0526, 0503, 0533, 0531, 0556, 0589) ○
UNI0516		Angle & width adjustable footplates (290mm depth)	(NC 0526, 0503, 0533, 0531, 0503, 0556, 0589) £121
UNI0517		Angle adjustable foot platform (240mm depth, 305mm width)	(NC 0526, 0503, 0533, 0531, 0504, 0556, 0589) £243
ARX9259		Without Centre mount Footplate	○

FOOTPLATE TYPE Side Mounted*				
ARM0512		Fixed Footplate	(NC 0533, 0531)	◆
ARM0504		Angle adjustable Footplate	(NC 0526)	£64
ARM0513		Angle and Depth adjustable		○
ARM1004		Without Footplates	(NC 0503, 0533, 0531)	○
CONTROLS*				
ARM2201		REM110	(NC 1437, 1092, 1095, 1093, 1094, 0977, 0589, 0531, 2231, 2232, 2230, 2247, 2257, 2233, 2240, 2219, 2245, 2270, 1908, 2209, 2210)	◆
ARM2203		REM211	(NC 1437, 2231, 2232, 2230, 2247, 2257, 2233, 2240, 2270, 1908, 2209, 2210)	○
ARM2204		REM216	(NC 2231, 2232, 2230, 2247, 2257, 2233, 2240, 2270, 2209, 2210)	£137
ARM2205		REM400	(NC 2209, 2210)	£305
ARM2206		REM500	(NC 1708, 1709, 1710, 2213, 2214, 2227, 2228, 2220, 2270)	£695
CONTROL BRACKETRY*				
ARM1207		Standard Bracket RH - Not Compatible with REM500		○
ARM1208		Standard Bracket LH - Not Compatible with REM500		○
ARM1707		Swing away RH Bracket		£119
ARM1708		Swing away LH Bracket		£119
ARM1709		Matrx Resolve Swing in/out bracket RH	(NC 2213, 2275, 2235, 2213, 2241, 1404)	£219
ARM1710		Matrx Resolve Swing in/out bracket LH	(NC 2214, 2276, 2236, 2214, 2242, 1405)	£219
ARM2213		Centre mount, swing away to RH Nucleus	(NC 1405, 1415)	£446
ARM2214		Centre mount, swing away to LH Nucleus	(NC 1405, 1339)	£446
ARM2227		RH Bracket mounted for Attendant	(NC 2220, 2108, 2128)	£57
ARM2228		LH Bracket mounted for Attendant	(NC 2220, 2108, 2128)	£57
ARM2209		REM500 Swing away RH multi-adjustable Bracket	(NC 1415)	○
ARM2210		REM500 Swing away LH multi-adjustable Bracket	(NC 1339)	○
ARZ2239		Individual mounting option per request (Special)		POA
REM400 TOGGLE OPTION				
ARM2270		Assignable toggle switch module	(NC 2203, 2204, 2206)	£284
DUAL CONTROL				
ARM2220		Dual control REM050	(NC 2108, 2205, 2206, 1221, 1222, 2230, 2247, 2257, 2233, 2240)	£355
ARM2219		Compact (ACU200)	(NC 1222, 2108, 2201)	£355
ARM2245		IDC Dual control	(NC 2201)	£673
DUAL CONTROL BRACKETRY				
ARM1207		Standard Bracket RH		○
ARM1208		Standard Bracket LH		○
ARM1221		IDC Dual control detachable Bracket RH	(NC 2220, 2227)	○
ARM1222		IDC Dual control detachable Bracket LH	(NC 2220, 2228)	○
LINX COMPACT REMOTE & BRACKETS (Requires REM400 or REM500)				
ARM2231		Proportional Compact Remote	(NC 2203, 2204)	£740
ARM2232		Proportional Compact Remote low force	(NC 2203, 2204)	£990
ARM2207		RH mount Bracket	(NC 2201, 2203, 2204)	○
ARM2208		LH mount Bracket	(NC 2201, 2203, 2204)	○
ARM2213		Centre mount, swing away to RH Nucleus Bracket	(NC 1405, 1415, 2201, 2203, 2204)	£446
ARM2214		Centre mount, swing away to LH Nucleus Bracket	(NC 1405, 1339, 2210, 2203, 2204)	£446
ARZ2215		Individual mounting option per request (Special)		POA
CHIN CONTROL				
UNI9250		Mo-vis Micro joystick (Please specify in Special Comments box mounting request if required)	(NC 9251)	POA
UNI9251		Mo-vis Multi joystick (Please specify in Special Comments box mounting request if required)		POA

<b>LINX HEAD CONTROL ATOM TECHNOLOGY (Requires REM400 or REM500)</b>			
ARM2230		ATOM with ultrasonic proximity sensor	(NC 2171, 2179, 2197, 2201, 2203, 2204, 2220, 2233, 2324, 2325, 2327, 1350, 1432, 1427, 2257, 2240) £3,835
ARM2247		PROTON with ultrasonic proximity sensor	(NC 2171, 2179, 2197, 2201, 2203, 2204, 2220, 2233, 2324, 2325, 2327, 1350, 1432, 1427, 2257, 2240) £4,491
<b>FOUR-SWITCH PROXIMITY ARRAY (with ATOM TECHNOLOGY)</b>			
ARM2257		Four-switch Array (ASL106)	(NC 2203, 2204, 2220, 2230, 2247, 2233, 2240, 2128) £5,938
ARM2275		Centre mount, swing-away to RH nucleus	(NC 2203, 2204, 1404, 1405, 1415) £446
ARM2276		Centre mount, swing-away to LH nucleus	(NC 2203, 2204, 1404, 1405, 1339) £446
ARZ2215		Individual mounting option per request (Special)	POA
<b>WIRELESS ACCESSORIES FOR ATOM</b>			
ARM2249		Wireless triple switch receiver (ASL557-3)	£305
ARM2250		Wireless mouse emulator (ASL558)	£733
<b>LINX PROPORTIONAL REMOTE / INPUTS &amp; BRACKETS (Requires REM400 or REM500)</b>			
ARM2233		MEC proportional joystick	(NC 2203, 2204, 2220, 2230, 2247, 2257, 2240, 2128) £3,723
ARM2240		Proportional heavy duty with mode switch (ASL133)	(NC 2203, 2204, 2220, 2230, 2247, 2257, 2233, 2128) £2,076
ARM2235		MEC tray mount RH (ASL603)	(NC 1415, 2201, 2203, 2204) ○
ARM2236		MEC tray mount LH (ASL603)	(NC 1339, 2201, 2203, 2204) ○
ARM2213		Centre mount, swing away to RH nucleus	(NC 1405, 1415, 2201, 2203, 2204) £446
ARM2214		Centre mount, swing away to LH nucleus	(NC 1405, 1339, 2201, 2203, 2204) £446
ARM2241		Centre mount, swing away to RH nucleus	(NC 1405, 1415, 2201, 2203, 2204) £446
ARM2242		Centre mount, swing away to LH nucleus	(NC 1405, 1339, 2201, 2203, 2204) £446
ARZ2215		Individual mounting option per request (Special)	POA
<b>LINX 10-WAY SWITCH</b>			
ARM2251		DLX-SW10-A	£453
ARM2217		10-Way Switch Bracket RH, Mounted for User	(NC 2213, 2275, 2235, 2241, 1415) ○
ARM2218		10-Way Switch Bracket LH, Mounted for User	(NC 2214, 2276, 2236, 2242, 1339) ○
ARM2227		10-Way Switch Bracket RH, Mounted for Attendant	(NC 2128, 1221) £57
ARM2228		10-Way Switch Bracket LH, Mounted for Attendant	(NC 2128, 1222) £57
<b>LINX ACCESSORY</b>			
ARM2252		Loom providing jack sockets for drive and mode	£90
<b>LINX USB CHARGER PORT</b>			
ARM2253		USB charger port DLX-USB02-A	(NC 2701, 2702, 2703, 2704, 2705, 2706) £143
ARM2207		RH mount	○
ARM2208		LH mount	○
<b>REMOTE STOP SWITCH</b>			
ARM2258		Remote Stop Switch (ASL504)	(NC 2203, 2204) £942
<b>LINX INPUT MODULE</b>			
UNI2273		IN200 Module (Jack port for external power button )	(NC 2220) £242
UNI9252		IN500 Module (for DB9 Drive input)	£785
<b>FURTHER LINX COMPONENTS / ACCESSORIES</b>			
ARM2300		To be enclosed only	○
ARZ2301		To be connected and programmed (special)	POA
<b>BUDDY BUTTONS</b>			
ARM1685		Red	£79
ARM1686		Yellow	£79
ARM1687		Green	£79
ARM1690		Blue	£79
<b>EGG SWITCH</b>			
ARM2254		Egg Switch, Black	£79

PIKO BUTTONS			
ARM1711		Red	£142
ARM1712		Yellow	£142
ARM1713		Green	£142
ARM1714		Orange	£142
ARM1715		Blue	£142
ARM1716		Black	£142
TOGGLES / BUTTONS			
ARZ2261		Stereo Toggle	£242
ARZ2262		Stereo Button	£230
ARZ2263		Four-way Toggle	£300
ARZ2264		Four-way Button	£321
ARZ2207		RH mount	○
ARZ2208		LH mount	○
ARZ2225		Without mounting	○
BODYPOINT JOYSTICK KNOBS			
ARM2721		3" U-shaped, joystick handle	£65
ARM2722		4" U-shaped, joystick handle	£65
ARM2723		Domed rubber, joystick handle	£65
ARM2724		Mushroom grooved, joystick handle	£65
SHROUD COLOUR*			
ARM2640		Obsidian Black Frame Shrouds	○
ARM2641		Ruby Red Frame Shrouds	○
ARM2642		Purple Passion Frame Shrouds	○
ARM2643		Peacock Blue Frame Shrouds	○
ARM2644		Arctic Blue Frame Shrouds	○
RIM COLOUR*			
ARM2650		Obsidian Black Rim Inserts	○
ARM2651		Ruby Red Rim Inserts	○
ARM2652		Purple Passion Rim Inserts	○
ARM2653		Peacock Blue Rim Inserts	○
ARM2654		Arctic Blue Rim Inserts	○
FORKS*			
ARM2631		Single sided	£316
ARM2632		Double sided	◆
CENTRE WHEELS (DRIVING WHEELS)*			
ARM0704		DW 14" (355 mm) x 3" (75 mm), Pneumatic	◆
ARM0705		DW 14" (355 mm) x 3" (75 mm), Solid inlay, dual density foam insert	○
MOTOR TYPE*			
ARM1294		3 km/h (2mph)	○
ARM1282		6 km/h (4mph)	◆
ARM1908		10 km/h (6mph)	(NC 2203, 1437, 1254) £464
BATTERY*			
ARM1257		50 Ah Gel	◆

PERFORMANCE ACCESSORIES			
ARM1215		G-Trac	£738
CHARGER*			
ARM1343		OptiCharge 8A (200 - 250V)	◆
ARM1313		OptiCharge 10A (200 - 250V)	£137
LIGHTS*			
ARM1437		LED Light System with Indicators (30 LUX)	(NC 2203, 1908) £407
ARZ1006		Lights Included with Motor Package 10 & 12kmh	(NC 1282) ○
ARM1254		Without light system	◆
TRANSPORTATION KIT*			
ARM1450		Preparation for DAHL docking station	£524
ARM1494		4-Point tie-down	◆
FIXATION BELT*			
ARM1329		Lap Belt fitted to base of backrest	(NC 2702, 2703, 2704, 2705, 2706) ◆
ARM2113		Lap Belt adjustable mounting on seat rail extrusion	(NC 2701, 2702, 2703, 2704, 2705, 2706) ○
ARM1332		H-Belt	(NC 2108) £66
ARM2701		Evoflex pelvic stabilizer, size L	(NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113) £84
ARM2702		Evoflex pelvic stabilizer, size M	(NC 2123, 2124, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113) £84
ARM2703		2-Point hip belt centre pull, size L	(NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113) £78
ARM2704		2-Point hip belt centre pull, size M	(NC 2123, 2124, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113) £78
ARM2705		2-Point hip belt dual pull, size L	(NC 2121, 2122, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113) £78
ARM2706		2-Point Hip Belt Dual pull, size M	(NC 2123, 2124, 2132, 2184, 2185, 2186, 2187, 2188, 2189, 2253, 1329, 2113) £78
REAR VIEW MIRROR			
ARM1415		RH Mirror	(NC 2213, 2217, 2275, 2235, 2241) £38
ARM1339		LH Mirror	(NC 2214, 2218, 2276, 2236, 2242) £38
OTHER OPTIONS			
ARM2128		Bag hooks	(NC 2173, 2227, 2228, 2257, 2233, 2240) £69
ARM1312		RH Cane holder	£38
ARM1363		LH Cane holder	£38
ARM1404		Swing away RH Tray	(NC 1707, 2209, 2213, 2214, 2235, 2236, 2241, 2242, 2207, 2275, 2276) £174
ARM1405		Swing away LH Tray	(NC 1708, 2210, 2213, 2214, 2235, 2236, 2241, 2242, 2208, 2275, 2276) £174

<b>TOTAL PRICE</b>	
--------------------	--



**Non Conforming Options:**

ARM2121	380mm (15")
ARM2122	430mm (17")
ARM2123	480mm (19")
ARM2124	530mm (21")
ARM2132	410mm (16")
ARM2134	460mm (18")
ARM2136	510mm (20")
ARM2138	560mm (22")
ARM2143	Matrix Flo-tech Contour Visco (430mm x 430mm only - up to 125kg)
ARM2163	Short 480mm
ARM2108	WITHOUT push handles
ARM2171	Tension Adjustable Backrest
ARM2179	Flex3 Tension Adjustable
ARM2172	Flex3 Backrest plates without cushion
ARM2173	Flex3 Backrest plates with Laguna cushion
ARM2174	Matrix PB 400mm
ARM2175	Matrix PB 500mm
ARM2176	Matrix Elite E2 Back , height 500 mm (Mesh fabric)
ARM2177	Matrix Elite TR
ARM2178	Matrix Elite Deep E2 Back, height 500 mm (Mesh fabric)
ARM2151	Black fabric
ARM2155	Black Spacer fabric
ARM2181	T-Armrest extended height
ARM2183	In-rigged armrests
ARM2193	Flip-up armrests
ARM2192	Following 230mm - 300mm
ARM2194	Following 300mm - 360mm
ARM2196	Long 375mm Armpad
ARM2197	Long ergonomic 397mm Armpad
ARM2199	Soft Armpad
ARM2184	Small 140mm x 100mm
ARM2186	Medium 230mm x 100mm
ARM2188	Large 330mm x 100mm
ARM2185	Small 140mm x 100mm
ARM2187	Medium 230mm x 100mm
ARM2189	Large 330mm x 100mm
ARM1350	Neckrest height and depth adjustable
ARM1432	Headrest, height, depth and sideways adjustable
ARM1427	Headrest with cheek support, height and depth adjustable
ARM2324	Elan standard pad 255mm (10")
ARM2325	Elan occipital pad 230mm x 180mm (9" X 7")
ARM2327	Elan 4 point headrest pad 280mm x 255mm (11" X 10")
ARM1496	All bracketry as per option 1432

ARM1445	Multi-adjustable pad
ARM1092	Electric Tilt (0°- 20°)
ARM1095	Electric Tilt (0°- 20°) Max User weight 150kg
ARM1093	Electric Tilt (0° - 30°) with CoG shift Seat Height 460mm
ARM2052	Electric Tilt (0°-50°) with COG shift Seat height 460mm (136kg Max User Weight)
ARM1094	Seat lifter with Electrical Tilt (300mm, 30°) with CoG shift Seat Height 460mm - 490mm
ARM2053	Seat lifter with Electrical Tilt (300mm, 45°) with CoG shift Seat Height 460mm (136kg Max User Weight)
ARM0354	Seat height 450mm (17.75")
ARM0364	Seat height 475mm (18.75")
ARM0977	Powered Recliner
ARM0556	Centre mounted : Manual
ARM0589	Centre mounted : LNX Powered
UNI0592	LNX Power centre mount legrest with telescoping foot platform (lowers to the ground) inc calf pads: Length 380-419mm (Selectable length depends on STF height) Custom-made, requires customised setting
ARM0526	Side mounted : Standard 80° hanger
ARM0503	Side mounted : VARI F
ARM0533	Side mounted : VARI A
ARM0531	Side mounted : ADE Powered
ARM0555	LH Amputee Support
ARM0554	RH Amputee Support
ARM0517	Angle Adjustable foot platform Medium 305mm x 275mm (12" x 10.75") for Centre Mounted Hanger only
UNI0517	Angle adjustable foot platform (240mm depth, 305mm width)
ARM0512	Fixed Footplate
ARM0504	Angle adjustable Footplate
ARM0513	Angle and Depth adjustable
ARM2201	REM110
ARM2203	REM211
ARM2204	REM216
ARM2205	REM400
ARM2206	REM500
ARM1707	Swing away RH Bracket
ARM1708	Swing away LH Bracket
ARM1709	Matrx Resolve Swing in/out bracket RH
ARM1710	Matrx Resolve Swing in/out bracket LH
ARM2213	Centre mount, swing away to RH Nucleus
ARM2214	Centre mount, swing away to LH Nucleus
ARM2227	RH Bracket mounted for Attendant
ARM2228	LH Bracket mounted for Attendant
ARM2209	REM500 Swing away RH multi-adjustable Bracket
ARM2210	REM500 Swing away LH multi-adjustable Bracket
ARM2270	Assignable toggle switch module
ARM2220	Dual control REM050
ARM2219	Compact (ACU200)
ARM2245	IDC Dual control
ARM1221	IDC Dual control detachable Bracket RH

ARM1222	IDC Dual control detachable Bracket LH
ARM2231	Proportional Compact Remote
ARM2232	Proportional Compact Remote low force
ARM2207	RH mount Bracket
ARM2208	LH mount Bracket
ARM2213	Centre mount, swing away to RH Nucleus Bracket
ARM2214	Centre mount, swing away to LH Nucleus Bracket
UNI9251	Mo-vis Multi joystick (Please specify in Special Comments box mounting request if required)
ARM2230	ATOM with ultrasonic proximity sensor
ARM2247	PROTON with ultrasonic proximity sensor
ARM2257	Four-switch Array (ASL106)
ARM2275	Centre mount, swing-away to RH nucleus
ARM2276	Centre mount, swing-away to LH nucleus
ARM2233	MEC proportional joystick
ARM2240	Proportional heavy duty with mode switch (ASL133)
ARM2235	MEC tray mount RH (ASL603)
ARM2236	MEC tray mount LH (ASL603)
ARM2213	Centre mount, swing away to RH nucleus
ARM2214	Centre mount, swing away to LH nucleus
ARM2241	Centre mount, swing away to RH nucleus
ARM2242	Centre mount, swing away to LH nucleus
ARM2217	10-Way Switch Bracket RH, Mounted for User
ARM2218	10-Way Switch Bracket LH, Mounted for User
ARM2227	10-Way Switch Bracket RH, Mounted for Attendant
ARM2228	10-Way Switch Bracket LH, Mounted for Attendant
ARM2253	USB charger port DLX-USB02-A
ARM2207	RH mount
ARM2208	LH mount
ARZ2207	RH mount
ARZ2208	LH mount
ARM1282	6 km/h (4mph)
ARM1908	10 km/h (6mph)
ARM1437	LED Light System with Indicators (30 LUX)
ARM1254	Without light system
ARM1329	Lap Belt fitted to base of backrest
ARM2113	Lap Belt adjustable mounting on seat rail extrusion
ARM1332	H-Belt
ARM2701	Evoflex pelvic stabilizer, size L
ARM2702	Evoflex pelvic stabilizer, size M
ARM2703	2-Point hip belt centre pull, size L
ARM2704	2-Point hip belt centre pull, size M
ARM2705	2-Point hip belt dual pull, size L
ARM2706	2-Point Hip Belt Dual pull, size M
ARM1415	RH Mirror

---

ARM1339	LH Mirror
ARM2128	Bag hooks
ARM1404	Swing away RH Tray
ARM1405	Swing away LH Tray

---



For a more complex customer requirement  
please refer to our Invacare Unique special order form.  
email [specialsuk@invacare.com](mailto:specialsuk@invacare.com)

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):