



Yes, you can:

Tel: 01656 776222
 Email: ordersuk@invacare.com
 Online Spares available at: www.invacare.co.uk

Manual Retail Prescription Form

JANUARY 2025
 Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Rea Azalea Base

Standard Features

- Max User Weight 135KG/21 Stone
- Crash Tested ISO 7176-Part 19
- Transit Propulsion, Attendant hub brakes and Solid black tyres
- 4 Seat Widths 340mm to 490mm (+50mm, +100mm with Transit)
- 2 Seat Heights 400mm or 450mm
- 2 seat depths 380-450mm (with 340mm width) or 430-500mm
- Seat plate or Trapeze Interface (without Interface option)
- Manually settable backrest recline
- Assistant operated gas piston -1° to +23° seat tilt
- Tilt lock - Pre-settable angle
- Fixed castor attachment
- Pearl Silver frame
- User manoeuvred tyre brakes
- Anti Tippers, Front tiedown brackets

Standard Base price: £1,302



Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

Code	Description	
BASE PRODUCT*		
BQA0010	Rea Azalea Base	£1,302
BQA8000	SPECIAL REQUEST - PLEASE INCLUDE QUOTATION REFERENCE IN SPECIAL DESCRIPTION AT THE END OF THE FORM WHERE APPLICABLE	POA
PROPULSION*		
BQA0020	Transit 16" wheels with hub brakes (NC 0705, 0715, 0724, 0810, 0850, 0863, 0867)	◆
BQA0030	Self Propel 24"/22" wheels with user brakes (NC 0100, 0430, 0480, 0702, 0726)	£31
BQA0040	Self Propel 24"/22" wheels with assistant hub brakes (NC 0100, 0430, 0480, 0724)	£140

SEAT WIDTH*			
BQA0050	<input type="checkbox"/>	340 mm/13" (Extra small chassis)	(NC 0120, 0164, 0168, 0169, 0180, 0186, 0210, 0260, 0600, 0665, 0675, 0863, 0867, 0900, 0902) ○
BQA0060	<input type="checkbox"/>	390 mm/15" (Small chassis)	(NC 0110) ○
BQA0070	<input type="checkbox"/>	440 mm/17" (Medium chassis)	(NC 0110) ○
BQA0080	<input type="checkbox"/>	490 mm/19" (Large chassis)	(NC 0110) ○
SEAT WIDTH EXTENSION*			
BQA0090	<input type="checkbox"/>	+ 50 mm seat width pre-adjustment on armrests and legrests	○
BQA0100	<input type="checkbox"/>	+ 100 mm seat width pre-adjustment on armrests and legrests	(NC 0863, 0867) ○
BQA0105	<input type="checkbox"/>	Without seat width extension, set at 0mm adjustment but can be adjusted wider if required	○
SEAT DEPTH*			
BQA0110	<input type="checkbox"/>	380-450 mm (15" - 18") only possible with seat width 340mm	○
BQA0120	<input type="checkbox"/>	430-500 mm (17" - 20")	○
SEAT HEIGHT*			
BQA0130	<input type="checkbox"/>	400 mm / 16"	(NC 0180, 0340, 0665, 0675, 0715) ○
BQA0140	<input type="checkbox"/>	450 mm / 18"	(NC 0705) ◆
SEAT ANGLE			
BQA0145	<input type="checkbox"/>	Negative Seat Angle -1 To +23° (For attendant operated gas piston in negative tilt)	○
SEAT INTERFACE*			
BQA0150	<input type="checkbox"/>	Seat Plate standard	◆
BQA0155	<input type="checkbox"/>	Seat Plate with wedges (Seat Angle +/- 5°)	£212
BQA0157	<input type="checkbox"/>	H-Bracket	£77
BQA0170	<input type="checkbox"/>	Without interface (allows mounting of trapeze interface)	○
BACKREST REGULATION*			
BQA0160	<input type="checkbox"/>	Assistant operated gas piston Backrest Recline	£85
BQA0161	<input type="checkbox"/>	Backrest regulation with rigging screw	◆
TILT LOCK*			
BQA0932	<input type="checkbox"/>	Tilt Lock (locking of pre-set angle)	◆
BACKREST*			
BQA0164	<input type="checkbox"/>	Matrx PB Elite Deep (includes hardware mounted to frame, backrest supplied in box)	£546
BQA0168	<input type="checkbox"/>	Matrx Elite (includes hardware mounted to frame, backrest supplied in box)	(NC 0567) £494
BQA0169	<input type="checkbox"/>	Matrx Posture Back (includes hardware mounted to frame, backrest supplied in box)	£353
BQA0167	<input type="checkbox"/>	Without Backrest	(NC 0172, 0173, 0567) ○
MATRX BACKREST HEIGHT			
BQA0172	<input type="checkbox"/>	Matrx Backrest Height 410 mm / 16"	○
BQA0173	<input type="checkbox"/>	Matrx Backrest Height 510 mm / 20"	○
PUSH HANDLE*			
BQA0174	<input type="checkbox"/>	Push Handles connected	◆
CENTRAL LEGREST			
BQA0180	<input type="checkbox"/>	Central Legrest, includes calfpad and footrest. Also select options BQA0230, BQA0280, BQA1111 and BQA2222	(NC 0186, 0188, 0190, 0200, 0210, 0220, 0240, 0250, 0260, 0270, 0290, 0295, 0310, 0340, 0370, 0400, 0580, 0590, 0960) £303

LEGREST HANGER ATTACHMENT*			
BQA0185		Standard position for legrest mounting	◆
SWO0189		Standard position with 0-30° Legrest Abduction function (NC 0220, 0230, 0270, 0280, 0290, 0295)	£100
BQA0186		High position 50mm higher for leg length compensation	○
BQA0188		90° Hanger attachment to achieve a 90° legrest position on the Elevating Legrest (NC 0200, 0210, 0220, 0250, 0260, 0270)	○
LEGREST HANGER LH*			
BQA0190		Elevating Legrest, angle adjustable 0 to 75° swing in/swing out with knee pad, calf support with locking clip LH (NC 0310, 0595, 3333)	£97
BQA0200		Legrest Hanger fixed 80 deg, swing in/swing out LH (NC 0576, 0577, 0580)	£57
BQA0210		Legrest Hanger fixed 90 deg, swing in/swing out LH (NC 0340, 0576, 0577, 0580, 0600, 0665, 0675)	£57
BQA0220		Legrest Amputee LH (NC 0290, 0295, 0310, 0340, 0580, 0576, 0577, 0595)	£137
BQA0230		Without Legrest, LH (NC 0290, 0295, 0310, 0340, 0576, 0577, 0580, 0595)	○
LEGREST HANGER RH*			
BQA0240		Elevating Legrest, angle adjustable 0 to 75° swing In/swing out with knee pad and calf support with locking clip RH (NC 0370, 0595, 3333)	£97
BQA0250		Legrest Hanger fixed 80 deg, swing in/swing out RH (NC 0586, 0587, 0590)	£57
BQA0260		Legrest Hanger fixed 90 deg, swing in/swing out RH (NC 0400, 0586, 0587, 0590, 0600, 0665, 0675)	£57
BQA0270		Legrest Amputee RH (NC 0290, 0295, 0370, 0400, 0586, 0587, 0590, 0595)	£137
BQA0280		Without Legrest, RH (NC 0290, 0295, 0370, 0400, 0586, 0587, 0590, 0595)	○
FOOTREST LH*			
BQA0290		Footrest one piece, depth and angle adjustable (NC 0295, 0310, 0340, 0370, 0400, 0580, 0590)	£200
BQA0295		Footbox, angle adjustable, padded with short tube (NC 0310, 0340, 0370, 0400, 0580, 0590, 0595)	£543
BQA0310		Footrest lower part, depth and angle adjustable with short tube and heel strap, LH	£35
BQA0340		Footrest lower part, depth and angle adjustable with long tube and heel strap, LH	£35
BQA1111		Without footrest. Only select with options BQA0180, Central legrest, BQA0220 Amputee legrest and BQA0230 Without legrest	○
FOOTREST RH*			
BQA0370		Footrest lower part, depth and angle adjustable with short tube and heel strap, RH	£35
BQA0400		Footrest lower part, depth and angle adjustable with long tube and heel strap, RH	£35
BQA2222		Without footrest. Only select with options BQA0180, Central legrest, BQA0270 Amputee legrest, BQA0280 Without legrest, BQA0290 One piece footrest and BQA0295 Footbox	○
ARMREST LH*			
BQA0420		Armrest Ducktail Pad LH, height adjustable with autolock	£83
BQA0430		Armrest Wide Pad LH, height adjustable with autolock Compatible with 16" wheels only	£83
BQA0435		Armrests Reinforced Wide Pad, LH, height adjustable with autolock (NC 0900, 0902, 0930)	£89
BQA0440		Armrest Padded LH, height adjustable with autolock (NC 0520, 0900, 0902, 0920, 0930)	£64
BQA0450		Armrest Hemiplegic LH, height adjustable with autolock (NC 0900, 0902, 0910, 0915, 0920, 0930)	£200
BQA0460		Without Armrest LH (NC 0520, 0900, 0902, 0910, 0915, 0920, 0930)	○

ARMREST RH*			
BQA0470		Armrest Ducktail Pad RH, height adjustable with autolock	£83
BQA0480		Armrest Wide Pad RH, height adjustable with autolock compatible with 16" wheels only	£83
BQA0485		Armrest Reinforced Wide Pad, RH, height adjustable with autolock	(NC 0900, 0902, 0930) £89
BQA0490		Armrest Padded, RH, height adjustable with autolock.	(NC 0520, 0900, 0902, 0920, 0930) £64
BQA0500		Armrest Hemiplegic, RH, height adjustable with autolock	(NC 0900, 0902, 0910, 0915, 0920, 0930) £200
BQA0510		Without armrest RH	(NC 0520, 0900, 0902, 0910, 0915, 0920, 0930) ○
ARMREST OPTIONS			
BQA0520		Side Support Cover, removable with tools (wipeable cover for the sideguard)	(NC 9999) £51
BQA9999		Side Support Cover with foam cushion. 2cm thick cushion supplied complete with PU fabric cover. Reduces seat width by 2cm for a comfortable fit. Consists of options 0520 and 0530 (Removable with tools)	£51
BACKREST OPTIONS			
BQA0540		Headrest, Sideway Adjustable	(NC 0560, 0567) £108
BQA0550		Headrest with cheek support, sideway adjustable	(NC 0545, 0560, 0567) £137
BQA0545		Cover for headrest	£31
BQA0560		Neckrest	(NC 0567) £108
BQA0567		Headrest Matrix, 250 mm Soft Pad	(NC 0570) £185
BQA0570		Handles Clinic Kit for Head-/Neckrest (Levers which replace screws for easy adjustment)	£42
ELEVATING LEGREST OPTIONS*			
BQA0576		Calf pad cover, for LH elevating legrest	○
BQA0577		Knee pad for LH elevating legrest	◆
BQA0586		Calf pad cover, for RH elevating legrest	○
BQA0587		Knee pad for RH elevating legrest	◆
BQA3333		Without knee pad, select only with Fixed legrests, Amputee legrests and without legrests	○
LEGREST OPTIONS			
BQA0580		Footplate Extender LH	£56
BQA0590		Footplate Extender RH	£56
BQA0595		Calf Strap (for One piece footrest)	○
CHASSIS COLOUR*			
BQA0591		Pearl Grey	◆
BQA0592		Azurite	○
BQA0593		Sand	○
BQA0594		Electric Green	○
BQA0596		Happy Red	○
CHASSIS SAFETY			
BQA0600		Chassis extender, castors extend 50mm forward for added stability	○
CASTOR ATTACHMENT*			
BQA0610		Fixed castor attachment	◆

FRONT CASTORS*			
BQA0638		140x45 mm (5 1/2" x 1 3/4") Medium Hard Black tyre	£22
BQA0645		150x30 mm (6" x 1 1/4") Black Medium Hard tyre	£22
BQA0665		200x32 mm (8" x 1 1/4") Black Solid Tyre	○
BQA0675		200x45 mm (8" X 1 3/4") Black Solid tyre	◆
REAR WHEEL SIZE*			
BQA0702		16" / 405 mm Transit rear wheel with fixed attachment (NC 0724, 0745, 0790, 0810, 0850, 0863, 0867)	◆
BQA0705		22" / 560 mm Self Propel rear wheel with QR attachment	○
BQA0715		24" / 600 mm Self Propel rear wheel with QR attachment	○
REAR WHEEL TYPE*			
BQA0724		Standard self propel with colourless hub, clear anodized rim and Stainless Steel spokes	○
BQA0726		Hub Brake Self Propel or Transit 16" Transit with Black Plastic Hub and Rim 22" and 24" / 560 mm and 600 mm Self Propel with colourless hub, clear anodized rim and Stainless Steel spokes	◆
REAR WHEEL TYRE*			
BQA0745		Pneumatic Profile Black 1" 3/8 tyres	○
BQA0775		Solid Black tyres	◆
HANDRIM			
BQA0790		Aluminium handrim	○
BQA0810		Plastic coated handrim	£48
SPOKEGUARD			
BQA0850		Frost spokeguard	£166
BRAKE*			
BQA0860		User manoeuvred tyre brakes	◆
BQA0863		One Arm Brake LH	£126
BQA0867		One Arm Brake RH	£126
SAFETY OPTIONS			
BQA0870		Pelvic belt	○
TRANSPORTATION BRACKETS*			
BQA0895		Transportation Brackets, mounted to the front of the wheelchair. (Important note selection of this option will require that a suitable risk assessment is carried out on the final seat up of the chair. Please refer to the Instructions for Use for further information regarding using the Rea Azalea Base as a seat in transit)	◆
OTHER OPTIONS			
BQA0900		Table Tray without elbow pads includes mounting kit (NC 0930)	£138
BQA0902		Table Tray with elbow Pads includes mounting kit (NC 0930)	£138
BQA0910		Swing away table without elbow pads (NC 0920, 0930)	£195
BQA0915		Swing away table with elbow pads (NC 0920, 0930)	£195
BQA0905		Table tray padded cushion (To be ordered in combination with the fixed & swingaway table tray)	£51
BQA0920		Locking device for table tray (NC 0930)	£15
BQA0930		Mounting kit for table (if ordered separately)	£9
BQA0960		Clinic kit	£55
BQA0970		Harness attachment , must be selected in combination with Head/Neckrest	£28

TOTAL PRICE	
--------------------	--

Non Conforming Options:

BQA0100	+ 100 mm seat width pre-adjustment on armrests and legrests
BQA0110	380-450 mm (15" - 18") only possible with seat width 340mm
BQA0120	430-500 mm (17" - 20")
BQA0164	Matrx PB Elite Deep (includes hardware mounted to frame, backrest supplied in box)
BQA0168	Matrx Elite (includes hardware mounted to frame, backrest supplied in box)
BQA0169	Matrx Posture Back (includes hardware mounted to frame, backrest supplied in box)
BQA0172	Matrx Backrest Height 410 mm / 16"
BQA0173	Matrx Backrest Height 510 mm / 20"
BQA0180	Central Legrest, includes calfpad and footrest. Also select options BQA0230, BQA0280, BQA1111 and BQA2222
BQA0186	High position 50mm higher for leg length compensation
BQA0188	90° Hanger attachment to achieve a 90° legrest position on the Elevating Legrest
BQA0190	Elevating Legrest, angle adjustable 0 to 75° swing in/swing out with knee pad, calf support with locking clip LH
BQA0200	Legrest Hanger fixed 80 deg, swing in/swing out LH
BQA0210	Legrest Hanger fixed 90 deg, swing in/swing out LH
BQA0220	Legrest Amputee LH
BQA0230	Without Legrest, LH
BQA0240	Elevating Legrest, angle adjustable 0 to 75° swing in/swing out with knee pad and calf support with locking clip RH
BQA0250	Legrest Hanger fixed 80 deg, swing in/swing out RH
BQA0260	Legrest Hanger fixed 90 deg, swing in/swing out RH
BQA0270	Legrest Amputee RH
BQA0280	Without Legrest, RH
BQA0290	Footrest one piece, depth and angle adjustable
BQA0295	Footbox, angle adjustable, padded with short tube
BQA0310	Footrest lower part, depth and angle adjustable with short tube and heel strap, LH
BQA0340	Footrest lower part, depth and angle adjustable with long tube and heel strap, LH
BQA0370	Footrest lower part, depth and angle adjustable with short tube and heel strap, RH
BQA0400	Footrest lower part, depth and angle adjustable with long tube and heel strap, RH
BQA0430	Armrest Wide Pad LH, height adjustable with autolock Compatible with 16" wheels only
BQA0480	Armrest Wide Pad RH, height adjustable with autolock compatible with 16" wheels only
BQA0520	Side Support Cover, removable with tools (wipeable cover for the sideguard)
BQA9999	Side Support Cover with foam cushion. 2cm thick cushion supplied complete with PU fabric cover. Reduces seat width by 2cm for a comfortable fit. Consists of options 0520 and 0530 (Removable with tools)
BQA0545	Cover for headrest
BQA0560	Neckrest
BQA0567	Headrest Matrx, 250 mm Soft Pad
BQA0570	Handles Clinic Kit for Head-/Neckrest (Levers which replace screws for easy adjustment)
BQA0576	Calf pad cover, for LH elevating legrest
BQA0577	Knee pad for LH elevating legrest
BQA0586	Calf pad cover, for RH elevating legrest
BQA0587	Knee pad for RH elevating legrest
BQA3333	Without knee pad, select only with Fixed legrests, Amputee legrests and without legrests
BQA0580	Footplate Extender LH
BQA0590	Footplate Extender RH
BQA0595	Calf Strap (for One piece footrest)

BQA0600	Chassis extender, castors extend 50mm forward for added stability
BQA0665	200x32 mm (8" x 1 1/4") Black Solid Tyre
BQA0675	200x45 mm (8" X 1 3/4") Black Solid tyre
BQA0702	16" / 405 mm Transit rear wheel with fixed attachment
BQA0705	22" / 560 mm Self Propel rear wheel with QR attachment
BQA0715	24" / 600 mm Self Propel rear wheel with QR attachment
BQA0724	Standard self propel with colourless hub, clear anodized rim and Stainless Steel spokes
BQA0726	Hub Brake Self Propel or Transit 16" Transit with Black Plastic Hub and Rim 22" and 24" / 560 mm and 600 mm Self Propel with colourless hub, clear anodized rim and Stainless Steel spokes
BQA0745	Pneumatic Profile Black 1" 3/8 tyres
BQA0790	Aluminium handrim
BQA0810	Plastic coated handrim
BQA0850	Frost spokeguard
BQA0863	One Arm Brake LH
BQA0867	One Arm Brake RH
BQA0900	Table Tray without elbow pads includes mounting kit
BQA0902	Table Tray with elbow Pads includes mounting kit
BQA0910	Swing away table without elbow pads
BQA0915	Swing away table with elbow pads
BQA0920	Locking device for table tray
BQA0930	Mounting kit for table (if ordered separately)
BQA0960	Clinic kit



For a more complex customer requirement
please refer to our Invacare Unique special order form.
email specialsuk@invacare.com

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):