



**Yes, you can:**

Tel: 01656 776222  
 Email: ordersuk@invacare.com  
 Online Spares available at: www.invacare.co.uk

OCTOBER 2024  
 Customer Ref:

<b>Customer Ref*</b>	<b>Account No*</b>
<b>Contact No*</b>	<b>Order No*</b>
<b>Delivery Address*</b>	<b>Order Qty*</b>
	<b>Order date*</b>

## Rea Azalea Transit

### Standard Features

- Max User Weight 135KG/21 Stone
- Crash Tested ISO 7176-Part 19
- 16" Transit Propulsion with Attendant Hub Brake and Solid Black Tyres
- Choice of Seat Widths 390mm, 440mm or 490mm (+50mm or +100mm Width Adj)
- Choice of 2 Seat Heights 400mm or 450mm
- Depth Adjustable Seat 430mm-500mm
- Without Seat Cushion, Supplied with Seat Board
- Flex 3 Backrest Plate, Width and Height Adjustable
- Passad 2 Backrest Cushion with Multi-stretch Black PU coated fabric Cover and Push Handles Connected
- Assistant Operated Gas Piston Seat Tilt +1° - +24°
- Backrest Recline 30° with Tilt Lock
- Angle Adjustable Legrests with Calf Support, Height and Depth Adjustable Footplates
- Height Adjustable Wide Pad Armrests with Autolock
- Pearl Silver Frame Colour, Solid Black Front Castors
- User Manoeuvred Tyre Brakes
- Pelvic Belt, Adjustable both sides
- Anti Tippers
- Front Transportation Brackets

**Standard Base price: £2,393**



Please make your selection with a cross in the box    ○ Free of charge    ◆ standard feature    POA special request

Customer weight =                      KG

Code	Description	
<b>BASE PRODUCT*</b>		
BNE9002	Rea Azalea Transit	£2,393
BNE8000	SPECIAL REQUEST	POA
<b>PROPULSION*</b>		
BNE0020	Transit 16" wheels with hub brakes	◆
<b>SEAT WIDTH*</b>		
BNE0040	390 mm/15" (Small chassis)	(NC 0244, 0245, 0246, 0247) ○
BNE0050	440 mm/17" (Medium chassis)	(NC 0241, 0242, 0246, 0247) ○
BNE0060	490 mm/19" (Large chassis)	(NC 0206, 0241, 0242, 0243, 0244) ○

SEAT WIDTH EXTENSION*			
BNE0065		+ 50 mm pre-adjustment on armrests and legrests	○
BNE0075		+ 100 mm pre-adjustment on armrests and legrests (NC 0110)	○
BNE0085		Without seat width extension, set at 0mm adjustment but can be adjusted wider if required (NC 0206)	○
SEAT DEPTH*			
BNE0100		430-500 mm / 17" - 20"	◆
SEAT HEIGHT*			
BNE0110		400 mm / 16" Seat Height (NC 1015, 1025)	○
BNE0120		450 mm / 18" Seat Height	○
SEAT ANGLE			
BNE0127		Negative Seat Angle -1 To +23° (For attendant operated gas piston in negative tilt) (NC 0360)	○
SEAT TYPE*			
BNE0130		Depth Adjustable	◆
SEAT CUSHION*			
BNE0165		Cushion Flosshape seat upholstery (NC 1111)	£57
BNE0167		Cushion Vicair Adjuster O2 10cm (Black) (NC 0201, 0202)	£547
SWO0190		Cushion Flosshape seat upholstery Shorter depth 450mm deep. Only use with seat depth in shortest position (NC 0202, 1111)	£57
BNE0180		Without Seat Cushion, Seat Plate only (NC 0201, 0202)	○
SEAT UPHOLSTERY*			
BNE0201		Multi-stretch Black PU coated fabric	◆
BNE0202		Black Spacer TR36	○
BNE1111		No upholstery, select with options 0167 Vicair Adjuster O2 or 0180 Without seat cushion	○
BACKREST WIDTH EXTENSION*			
BNE0205		+ 50 mm Backrest width pre-adjustment, set at +50mm but can be adjusted if required (NC 0221, 0223, 0225)	○
BNE0206		+ 100 mm Backrest width pre-adjustment, set at +100mm but can be adjusted if required (NC 0220, 0221, 0223, 0225)	○
BNE0207		Without Backrest width pre-adjustment, set at 0mm but can be adjusted if required	○
BACKREST*			
BNE0210		Flex 3 Backrest, Fixed Plate Width Adjustable. Headrest attachment built into backrest plate (NC 0233, 0234, 0235, 0240, 0241, 0242, 0243, 0244, 0245, 0246, 0247, 0343)	◆
BNE0220		Flex 3 Backrest Tension (Width Adjustable) (NC 0233, 0234, 0235, 0241, 0242, 0243, 0244, 0245, 0246, 0247, 0343)	£102
BNE0221		Matrx Elite E2 (includes hardware mounted to frame, backrest supplied in box) (NC 0240, 0250, 0260, 0270, 0305, 0321, 0322, 0800, 0810, 0820, 0830, 0880, 0890, 0895, 0896, 0900, 0920, 7777)	£494
BNE0223		Matrx Elite E2 Deep (includes hardware mounted to frame, backrest supplied in box) (NC 0240, 0250, 0260, 0270, 0305, 0321, 0322, 0800, 0810, 0820, 0830, 0880, 0890, 0895, 0896, 0900, 0920, 7777)	£546
BNE0225		Matrx Posture Back (includes hardware mounted to frame, backrest supplied in box) (NC 0235, 0240, 0241, 0242, 0243, 0244, 0245, 0246, 0247, 0250, 0260, 0270, 0305, 0321, 0322, 0335, 0800, 0810, 0820, 0830, 0880, 0890, 0895, 0896, 0900, 0920)	£353
BNE0230		Backrest Tubes Without Backrest (NC 0233, 0234, 0235, 0240, 0241, 0242, 0243, 0244, 0245, 0246, 0247, 0250, 0260, 0270, 0305, 0321, 0322, 0345, 0740, 0800, 0810, 0820, 0830, 0880, 0890, 0895, 0896, 0900, 0920)	○
MATRX BACKREST HEIGHT			
BNE0233		Matrx Backrest Height 410 mm / 16"	○
BNE0235		Matrx Backrest Height 460 mm / 18"	○
BNE0234		Matrx Backrest Height 510 mm / 20"	○

<b>MATRIX BACKREST WIDTH</b>			
BNE0241		Matrix Backrest Width 360mm / 14" (Chest width 265mm / 10 1/2")	○
BNE0242		Matrix Backrest Width 380mm / 15" (Chest width 290mm / 11 1/2")	○
BNE0243		Matrix Backrest Width 410mm / 16" (Chest width 315mm / 12 1/2")	○
BNE0244		Matrix Backrest Width 430mm / 17" (Chest width 340mm / 13")	○
BNE0245		Matrix Backrest Width 460mm / 18" (Chest width 370mm / 14 1/2")	○
BNE0246		Matrix Backrest Width 480mm / 19" (Chest width 390mm / 15")	○
BNE0247		Matrix Backrest Width 510mm / 20" (Chest width 420mm / 16 1/2")	○
BNE7777		No backrest width, select with Matrix Posture back only	○
<b>BACKREST CUSHION*</b>			
BNE0240		Cover for tension adjustable upholstery (NC 2222)	○
BNE0250		Cushion, Laguna (Lateral Support) (NC 2222)	○
BNE0260		Cushion, Mistral 2 (Waist Support) (NC 0322, 2222)	£144
BNE0270		Cushion, Passad 2 (Shoulder Support) (NC 2222)	◆
BNE0300		Without Backrest Cover/ Cushion (NC 0305, 0321, 0322)	○
<b>WIDER BACKREST COVER</b>			
BNE0305		Extra large backrest cover or cushion. One size wider backrest cushion or tension cover fitted if selected	○
<b>BACKREST UPHOLSTERY*</b>			
BNE0321		Multi-stretch Black PU coated fabric	◆
BNE0322		Black Spacer TR36	○
BNE2222		No upholstery, select with Matrix, and Without backrest options	○
<b>PUSH HANDLES*</b>			
BNE0335		Push Handles connected	◆
BNE0343		Push Handles connected,with Rea head/neckrest attachment	£72
BNE0345		Push Bar, Angle Adjustable	£102
<b>SEAT TILT AND BACKREST REGULATION*</b>			
BNE0350		Assistant operated gas piston seat tilt and backrest recline (NC 5555)	◆
BNE0360		Electrical operated seat tilt and backrest recline (Tilt range is -1° To +23°) (NC 1263)	£616
<b>TILT LOCK*</b>			
BNE1263		Tilt Lock (locking of pre-set angle)	◆
BNE5555		No Tilt Lock, Select with 0360 Electric operated seat and backrest	○
<b>LEGREST HANGER ATTACHMENT*</b>			
BNE0375		Standard position for legrest mounting	◆
SWO0189		Standard position with 0-30° Legrest Abduction function (NC 0475, 0476, 0477)	£100
BNE0376		High position 50mm higher for leg length compensation	○
BNE0378		90° Hanger attachment to achieve a 90° legrest position on the Elevating Legrest only (NC 0390, 0400, 0410, 0440, 0450, 0460)	○

<b>LEGREST HANGER LH*</b>			
BNE0380		Elevating Legrest, angle adjustable 0 to 75° swing in/swing out with knee pad, calf support with locking clip LH	(NC 0475, 0490, 0942, 0945, 6666) ◆
BNE0390		Legrest Hanger fixed 80 deg, swing in/swing out LH	(NC 0928, 0929, 0930) ○
BNE0400		Legrest Hanger fixed 90 deg, swing in/swing out LH	(NC 0520, 0928, 0929, 0930) ○
BNE0410		Legrest Amputee LH	(NC 0475, 0476, 0477, 0490, 0520, 0928, 0929, 0930, 0941, 0942, 0945, 0952, 0954, 0956) £137
BNE0420		Without Legrest, LH	(NC 0475, 0476, 0477, 0490, 0520, 0928, 0929, 0930, 0941, 0942, 0945, 0952, 0954, 0956) ○
<b>LEGREST HANGER RH*</b>			
BNE0430		Elevating Legrest, angle adjustable 0 to 75° swing In/swing out with knee pad and calf support with locking clip RH	(NC 0475, 0550, 0942, 0945, 6666) ◆
BNE0440		Legrest Hanger fixed 80 deg, swing in/swing out RH	(NC 0931, 0932, 0940) ○
BNE0450		Legrest Hanger fixed 90 deg, swing in/swing out RH	(NC 0580, 0931, 0932, 0940) ○
BNE0460		Legrest Amputee RH	(NC 0475, 0476, 0477, 0550, 0580, 0931, 0932, 0940, 0941, 0942, 0945, 0952, 0954, 0956) £137
BNE0470		Without Legrest, RH	(NC 0475, 0476, 0477, 0550, 0580, 0931, 0932, 0940, 0941, 0942, 0945, 0952, 0954, 0956) ○
<b>FOOTREST LH*</b>			
BNE0475		Footrest one piece, depth and angle adjustable short tube	(NC 0477, 0490, 0520, 0550, 0580, 0930, 0940, 0952, 0954, 0956) £118
BNE0476		Footrest one piece, depth and angle adjustable long tube	(NC 0477, 0490, 0520, 0550, 0580, 0930, 0940, 0952, 0954, 0956) £118
BNE0477		Footbox, angle adjustable, padded with short tube	(NC 0490, 0520, 0550, 0580, 0930, 0940, 0941, 0942, 0945, 0952, 0954, 0956) £543
BNE0490		Footrest lower part, depth and angle adjustable with short tube and heel strap, LH. Not compatible with Elevating legrest.	○
BNE0520		Footrest lower part, depth and angle adjustable with long tube and heel strap, LH. All knee to heel ranges can be achieved when used with Elevating legrests.	○
BNE3333		Without footrest, only select with option 0370 Central legrest, 0410 Amputee legrest and 0420 Without legrest	○
<b>FOOTREST RH*</b>			
BNE0550		Footrest lower part, depth and angle adjustable with short tube and heel strap, RH. Not compatible with Elevating legrest.	○
BNE0580		Footrest lower part, depth and angle adjustable with long tube and heel strap, RH. All knee to heel ranges can be achieved when used with Elevating legrests.	○
BNE4444		Without footrest, only select with options 0370 Central legrest, 0460 Amputee legrest, 0470 Without legrest, 0475 and 0476 One piece footrests or 0477 Footbox	○
<b>ARMREST LH*</b>			
BNE0600		Armrest Ducktail Pad LH, height adjustable with autolock	(NC 1232) ○
BNE0605		Armrest Wide Pad LH, height adjustable with autolock. Compatible with 16" wheels only	(NC 1232) ◆
BNE0607		Armrests Reinforced Wide Pad, LH, height adjustable with autolock	(NC 1230, 1231, 1232, 1250, 1260) £89
BNE0610		Armrest Padded LH, height adjustable with autolock	(NC 0680, 0695, 1230, 1231, 1232, 1250, 1260, 9999) £64
BNE0615		Armrest Hemiplegic LH, height adjustable with autolock	(NC 0695, 1230, 1231, 1232, 1240, 1245, 1250, 1260) £179
BNE0620		Without Armrest LH Select when choosing BNE1232 Half Tray LH, Narrow Armrest is included	(NC 0695, 1230, 1231, 1240, 1245, 1250, 1260) ○

ARMREST RH*			
BNE0630		Armrest Ducktail Pad RH, height adjustable with autolock (NC 1233)	○
BNE0635		Armrest Wide Pad RH, height adjustable with autolock. Compatible with 16" wheels only (NC 1233)	◆
BNE0637		Armrest Reinforced Wide Pad, RH, height adjustable with autolock (NC 1230, 1231, 1233, 1250, 1260)	£89
BNE0640		Armrest Padded RH, height adjustable with autolock. (NC 0680, 0695, 1230, 1231, 1233, 1250, 1260, 9999)	£64
BNE0645		Armrest Hemiplegic, RH, height adjustable with autolock (NC 0695, 1230, 1231, 1233, 1240, 1245, 1250, 1260)	£179
BNE0650		Without Armrest RH Select when choosing BNE1233 Half Tray RH, Narrow Armrest is included (NC 0695, 1230, 1231, 1240, 1245, 1250, 1260)	○
ARMREST OPTIONS			
BNE0680		Side Support Cover, PU fabric cover for the sideguards (Removable with tools) (NC 0695, 9999)	£51
BNE9999		Side Support Cover with foam cushion. 2cm thick cushion supplied complete with PU fabric cover. Reduces seat width by 2cm for a comfortable fit. Consists of options 0680 and 0690 (Removable with tools) (NC 0695)	£51
BNE0695		Pad for sideguard, removable with hook and loop fastening	£51
SEAT OPTIONS			
BNE0710		Incontinence Cover	£40
BNE0720		Pommel Abduction Cushion	£118
BNE0725		Tilt Scale for seat tilt	£44
HEADREST OPTIONS*			
BNE0730		Headrest, Sideway Adjustable Fixed hardware included (NC 0742, 0744, 0746, 0750, 0770, 0782)	◆
BNE0750		Headrest with cheek support, sideway adjustable Fixed hardware included (NC 0740, 0742, 0744, 0746, 0770, 0782)	£107
BNE0731		Matrx Elan Standard Pad (254 mm / 10") (NC 0740, 0770, 0782, 0790)	£124
BNE0732		Matrx Elan Standard Pad (355 mm / 14") (NC 0740, 0770, 0782, 0790)	£124
BNE0733		Matrx Elan Occipital Pad (229 mm / 9") (NC 0740, 0770, 0782, 0790)	£124
BNE0734		Matrx Elan Occipital Pad (305 mm / 12") (NC 0740, 0770, 0782, 0790)	£124
BNE0735		Matrx Elan 4-Point Pad (229W x 177H mm / 9"W x 7"H) (NC 0740, 0770, 0782, 0790)	£232
BNE0736		Matrx Elan 4-Point Pad (280W x 254H mm / 11"W x 10"H) (NC 0740, 0770, 0782, 0790)	£232
BNE9991		Without headrest (NC 0742, 0744, 0746)	○
HEADREST HARDWARE*			
BNE0742		Matrx Elan Hardware (NC 0740, 0770, 0790)	£184
BNE0744		Matrx Elan Extended Hardware (Extra link) (NC 0740, 0770, 0790)	£284
BNE0746		Matrx Elan LOXX Hardware (NC 0740, 0770, 0790)	£166
BNE8888		Without Elan headrest hardware (Only select with option 0730 & 0750 Headrests) (NC 0732, 0733, 0734, 0735, 0736)	○
OTHER HEADREST OPTIONS			
BNE0740		Cover for headrest	£31
BNE0770		Neckrest	○
BNE9992		Without neckrest	○
BNE0782		Alternative headrest/neckrest mounting bracket	£23
BNE0790		Handles Clinic Kit for Head/Neckrest (Lever which replace screws for easy adjustment)	○

TRUNK SUPPORT LH				
BNE0800		Trunk Support Multi Adjustable with standard arm LH	(NC 0900)	○
BNE0820		Trunk Support Multi Adjustable with extended height arm LH	(NC 0900)	○
BNE0880		Trunk Support Swingaway LH	(NC 0900)	£105
BNE0895		Trunk Support Multi Adjustable Swingaway LH	(NC 0900)	£113
TRUNK SUPPORT RH				
BNE0810		Trunk Support Multi Adjustable with standard arm RH	(NC 0900)	○
BNE0830		Trunk Support Multi Adjustable with extended height arm RH	(NC 0900)	○
BNE0890		Trunk Support Swingaway RH	(NC 0900)	£105
BNE0896		Trunk Support Multi Adjustable Swingaway RH	(NC 0900)	£113
BACKREST OPTIONS				
BNE0900		Trunk Support Attachment to mount laterals later. If trunk supports are chosen, the attachment bracket is already included		○
BNE0915		Tilt Scale for backrest recline		£44
BNE0920		Wedges (One Pair)		£11
ELEVATING LEGREST OPTIONS*				
BNE0928		Calf pad cover, for LH elevating legrest		○
BNE0929		Knee pad for LH elevating legrest		◆
BNE0931		Calf pad cover, for RH elevating legrest		○
BNE0932		Knee pad for RH elevating legrest		◆
BNE6666		Without Elevating Legrest Options (Select with Fixed legrests, Amputee legrests and Without legrests only)		○
LEGREST OPTIONS				
BNE0930		Footplate Extender LH	(NC 0941)	£34
BNE0940		Footplate Extender RH	(NC 0941)	£34
BNE0941		Cover Footplate, Black Plush Fabric		£54
BNE0942		Calf Strap Padded		£73
BNE0945		Calf Strap (for One piece footrest)		○
BODYPOINT ANKLE HUGGERS				
BNE0952		Bodypoint Ankle Huggers Large (Supplied with fixing kit, unmounted)		£114
BNE0954		Bodypoint Ankle Huggers Medium (Supplied with fixing kit, unmounted)		£114
BNE0956		Bodypoint Ankle Huggers Small (Supplied with fixing kit, unmounted)		£114
CHASSIS COLOUR*				
BNE0944		Silver Pearl		◆
BNE0946		Azurite (Dark Blue)		○
BNE0947		Sand (Beige)		○
BNE0948		Electric Green		○
BNE0949		Happy Red		○
CHASSIS SAFETY				
BNE0950		Chassis extender, castors extend 50mm forward for added stability		○
CASTOR ATTACHMENT*				
BNE0960		Fixed castor attachment		◆



<b>FRONT CASTORS*</b>			
BNE0992		150x30 mm (6" x 1 1/4") Medium Hard Black tyre	○
BNE1015		200x32 mm (8" x 1 1/4") Black Solid Tyre	○
BNE1025		200x45 mm (8" X 1 3/4") Black Solid tyre	◆
<b>REAR WHEEL SIZE*</b>			
BNE1047		16" / 405 mm Transit rear wheel with fixed attachment	◆
<b>REAR WHEEL TYPE*</b>			
BNE1076		Hub Brake 16" Transit with Black Plastic Hub and Rim	◆
<b>REAR WHEEL TYRE*</b>			
BNE1092		Solid Black tyres	◆
<b>BRAKE*</b>			
BNE1170		User manoeuvred tyre brakes, includes extended handles	◆
<b>PELVIC BELT*</b>			
BNE1200		Pelvic belt	◆
BNE1201		Bodypoint 2-Point Centre-pull Push Button Medium (Mounted to chair)	£83
BNE1202		Bodypoint 2-Point Centre-pull Push Button Large (Mounted to chair)	£83
BNE1203		Bodypoint 4-Point Centre-pull Push Button Medium (Mounted to chair)	£98
BNE1204		Bodypoint 4-Point Centre-pull Push Button Large (Mounted to chair)	£98
<b>BODYPOINT HARNESS</b>			
BNE1206		Bodypoint Pivotfit Harness Medium (Supplied with fixing kit, unmounted)	£109
BNE1207		Bodypoint Pivotfit Harness Large (Supplied with fixing kit, unmounted)	£109
<b>TRANSPORTATION BRACKETS*</b>			
BNE1225		Transportation Brackets, mounted to the front of the wheelchair	◆
<b>OTHER OPTIONS</b>			
BNE1230		Fixed Table Tray without elbow pads includes mounting kit (NC 1237, 1260)	£138
BNE1231		Fixed Table Tray with elbow Pads includes mounting kit (NC 1237, 1260)	£138
BNE1232		Half table tray, clear plastic, flip up, depth adjustable, left side (armrest with straight narrow pad included) (NC 1235, 1260)	£103
BNE1233		Half table tray, clear, plastic, flip up, depth adjustable, right side (armrest with straight narrow pad included) (NC 1235, 1260)	£103
BNE1240		Swing away table without elbow pads (NC 1237, 1250, 1260)	£152
BNE1245		Swing away table with elbow pads (NC 1237, 1250, 1260)	£152
BNE1235		Table tray padded cushion (To be ordered in combination with the fixed & swingaway table tray)	£51
BNE1237		Half table tray padded cushion	£46
BNE1250		Locking device for table tray	£15
BNE1260		Mounting kit for table (if ordered separately)	£9
BNE1290		Handles Clinic Kit for backrest, armrest, legrest and seat frame (Levers which replace screws for easy adjustment)	£55
BNE1300		Harness attachment , must be selected in combination with Head/Neckrest	£85
BNE1305		Front castor locks, includes front tiedown brackets. Allows castors to remain parallel for easy access into confined spaces such as lifts or vehicles	£185
BNE1310		Attachment for positioning belt	£32

<b>TOTAL PRICE</b>	
--------------------	--



**Non Conforming Options:**

BNE0110	400 mm / 16" Seat Height
BNE0201	Multi-stretch Black PU coated fabric
BNE0202	Black Spacer TR36
BNE1111	No upholstery, select with options 0167 Vicair Adjuster O2 or 0180 Without seat cushion
BNE0206	+ 100 mm Backrest width pre-adjustment, set at +100mm but can be adjusted if required
BNE0220	Flex 3 Backrest Tension (Width Adjustable)
BNE0221	Matrx Elite E2 (includes hardware mounted to frame, backrest supplied in box)
BNE0223	Matrx Elite E2 Deep (includes hardware mounted to frame, backrest supplied in box)
BNE0225	Matrx Posture Back (includes hardware mounted to frame, backrest supplied in box)
BNE0233	Matrx Backrest Height 410 mm / 16"
BNE0235	Matrx Backrest Height 460 mm / 18"
BNE0234	Matrx Backrest Height 510 mm / 20"
BNE0241	Matrx Backrest Width 360mm / 14" (Chest width 265mm / 10 1/2")
BNE0242	Matrx Backrest Width 380mm / 15" (Chest width 290mm / 11 1/2")
BNE0243	Matrx Backrest Width 410mm / 16" (Chest width 315mm / 12 1/2")
BNE0244	Matrx Backrest Width 430mm / 17" (Chest width 340mm / 13")
BNE0245	Matrx Backrest Width 460mm / 18" (Chest width 370mm / 14 1/2")
BNE0246	Matrx Backrest Width 480mm / 19" (Chest width 390mm / 15")
BNE0247	Matrx Backrest Width 510mm / 20" (Chest width 420mm / 16 1/2")
BNE7777	No backrest width, select with Matrx Posture back only
BNE0240	Cover for tension adjustable upholstery
BNE0250	Cushion, Laguna (Lateral Support)
BNE0260	Cushion, Mistral 2 (Waist Support)
BNE0270	Cushion, Passad 2 (Shoulder Support)
BNE0305	Extra large backrest cover or cushion. One size wider backrest cushion or tension cover fitted if selected
BNE0321	Multi-stretch Black PU coated fabric
BNE0322	Black Spacer TR36
BNE2222	No upholstery, select with Matrx, and Without backrest options
BNE0335	Push Handles connected
BNE0343	Push Handles connected,with Rea head/neckrest attachment
BNE0345	Push Bar, Angle Adjustable
BNE0360	Electrical operated seat tilt and backrest recline (Tilt range is -1° To +23°)
BNE1263	Tilt Lock (locking of pre-set angle)
BNE5555	No Tilt Lock, Select with 0360 Electric operated seat and backrest
BNE0390	Legrest Hanger fixed 80 deg, swing in/swing out LH
BNE0400	Legrest Hanger fixed 90 deg, swing in/swing out LH
BNE0410	Legrest Amputee LH
BNE0440	Legrest Hanger fixed 80 deg, swing in/swing out RH
BNE0450	Legrest Hanger fixed 90 deg, swing in/swing out RH
BNE0460	Legrest Amputee RH
BNE0475	Footrest one piece, depth and angle adjustable short tube
BNE0476	Footrest one piece, depth and angle adjustable long tube
BNE0477	Footbox, angle adjustable, padded with short tube

BNE0490	Footrest lower part, depth and angle adjustable with short tube and heel strap, LH. Not compatible with Elevating legrest.
BNE0520	Footrest lower part, depth and angle adjustable with long tube and heel strap, LH. All knee to heel ranges can be achieved when used with Elevating legrests.
BNE0550	Footrest lower part, depth and angle adjustable with short tube and heel strap, RH. Not compatible with Elevating legrest.
BNE0580	Footrest lower part, depth and angle adjustable with long tube and heel strap, RH. All knee to heel ranges can be achieved when used with Elevating legrests.
BNE0680	Side Support Cover, PU fabric cover for the sideguards (Removable with tools)
BNE9999	Side Support Cover with foam cushion. 2cm thick cushion supplied complete with PU fabric cover. Reduces seat width by 2cm for a comfortable fit. Consists of options 0680 and 0690 (Removable with tools)
BNE0695	Pad for sideguard, removable with hook and loop fastening
BNE0750	Headrest with cheek support, sideways adjustable Fixed hardware included
BNE0732	Matrx Elan Standard Pad (355 mm / 14")
BNE0733	Matrx Elan Occipital Pad (229 mm / 9")
BNE0734	Matrx Elan Occipital Pad (305 mm / 12")
BNE0735	Matrx Elan 4-Point Pad (229W x 177H mm / 9"W x 7"H)
BNE0736	Matrx Elan 4-Point Pad (280W x 254H mm / 11"W x 10"H)
BNE0742	Matrx Elan Hardware
BNE0744	Matrx Elan Extended Hardware (Extra link)
BNE0746	Matrx Elan LOXX Hardware
BNE0740	Cover for headrest
BNE0770	Neckrest
BNE0782	Alternative headrest/neckrest mounting bracket
BNE0790	Handles Clinic Kit for Head/Neckrest (Levers which replace screws for easy adjustment)
BNE0800	Trunk Support Multi Adjustable with standard arm LH
BNE0820	Trunk Support Multi Adjustable with extended height arm LH
BNE0880	Trunk Support Swingaway LH
BNE0895	Trunk Support Multi Adjustable Swingaway LH
BNE0810	Trunk Support Multi Adjustable with standard arm RH
BNE0830	Trunk Support Multi Adjustable with extended height arm RH
BNE0890	Trunk Support Swingaway RH
BNE0896	Trunk Support Multi Adjustable Swingaway RH
BNE0900	Trunk Support Attachment to mount laterals later. If trunk supports are chosen, the attachment bracket is already included
BNE0920	Wedges (One Pair)
BNE0928	Calf pad cover, for LH elevating legrest
BNE0929	Knee pad for LH elevating legrest
BNE0931	Calf pad cover, for RH elevating legrest
BNE0932	Knee pad for RH elevating legrest
BNE6666	Without Elevating Legrest Options (Select with Fixed legrests, Amputee legrests and Without legrests only)
BNE0930	Footplate Extender LH
BNE0940	Footplate Extender RH
BNE0941	Cover Footplate, Black Plush Fabric
BNE0942	Calf Strap Padded
BNE0945	Calf Strap (for One piece footrest)
BNE0952	Bodypoint Ankle Huggers Large (Supplied with fixing kit, unmounted)
BNE0954	Bodypoint Ankle Huggers Medium (Supplied with fixing kit, unmounted)
BNE0956	Bodypoint Ankle Huggers Small (Supplied with fixing kit, unmounted)

BNE1015	200x32 mm (8" x 1 1/4") Black Solid Tyre
BNE1025	200x45 mm (8" X 1 3/4") Black Solid tyre
BNE1230	Fixed Table Tray without elbow pads includes mounting kit
BNE1231	Fixed Table Tray with elbow Pads includes mounting kit
BNE1232	Half table tray, clear plastic, flip up, depth adjustable, left side (armrest with straight narrow pad included)
BNE1233	Half table tray, clear, plastic, flip up, depth adjustable, right side (armrest with straight narrow pad included)
BNE1240	Swing away table without elbow pads
BNE1245	Swing away table with elbow pads
BNE1235	Table tray padded cushion (To be ordered in combination with the fixed & swingaway table tray)
BNE1237	Half table tray padded cushion
BNE1250	Locking device for table tray
BNE1260	Mounting kit for table (if ordered separately)



For a more complex customer requirement  
please refer to our Invacare Unique special order form.  
email [specialsuk@invacare.com](mailto:specialsuk@invacare.com)

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):