



Yes, you can:

Tel: 01656 776222
 Email: ordersuk@invacare.com
 Online Spares available at: www.invacare.co.uk

Manual Retail Prescription Form

JANUARY 2025
 Customer Ref:

Customer Ref*	Account No*
Contact No*	Order No*
Delivery Address*	Order Qty*
	Order date*

Rea Clematis Pro

Standard Features

- Max user weight 135kg/21 Stone
- Crash tested ISO 7176 Part 19
- Self propulsion, 24" rear wheels
- Choice of 3 seat widths 390mm, 440mm and 490mm (+25mm width adj)
- Seat height 450mm
- Depth adjustable seat 420mm - 480mm
- Backrest plate, Laguna 2 cushion with height adjustable push handles
- Flo-shape seat cushion, Multi stretch PU coated upholstery TR26
- Assistant operated gas piston 25 Deg seat tilt and 30 Deg backrest recline
- Angle adjustable legrests with calf support
- Height adjustable standard pad armrests (including autolock)
- Pearl beige frame
- Quick release self propel wheels, solid tyres and aluminium handrims
- Anti tip devices
- Headrest, Pelvic belt, User manoeuvred tyre brakes

Standard Base price: £1,813



Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

Code	Description	
BASE PRODUCT*		
SNB0001	CLEMATIS PRO	£1,813
SNBSPEC	SPECIAL REQUEST - PLEASE INCLUDE QUOTATION REFERENCE IN SPECIAL DESCRIPTION AT THE END OF THE FORM WHERE APPLICABLE	<input type="radio"/>

PROPULSION*		
SNB0100	<input type="checkbox"/>	SELF PROPULSION (24" REAR WHEELS) ◆
FRAME*		
SNB0200	<input type="checkbox"/>	FIXED FRAME WITH 25MM WIDTH ADJUSTMENT EACH SIDE ◆
SEAT WIDTH*		
SNB0300	<input type="checkbox"/>	390 MM ○
SNB0310	<input type="checkbox"/>	440 MM ○
SNB0320	<input type="checkbox"/>	490 MM ○
SEAT DEPTH*		
SNB0360	<input type="checkbox"/>	420 - 480 MM ◆
SEAT HEIGHT*		
SNB0390	<input type="checkbox"/>	450 MM ◆
SEAT CUSHION*		
SNB0430	<input type="checkbox"/>	CUSHION FLO-SHAPE BLACK MULTI STRETCH PU COATED TR26 ◆
BACKREST*		
SNB0500	<input type="checkbox"/>	RIGID BACKREST PLATE ◆
BACKREST CUSHION*		
SNB0530	<input type="checkbox"/>	CUSHION LAGUNA 2 BLACK MULTI STRETCH PU COATED FABRIC ◆
PUSH HANDLES*		
SNB0560	<input type="checkbox"/>	HEIGHT ADJUSTABLE PUSH HANDLES ◆
SNB0570	<input type="checkbox"/>	PUSH BAR ANGLE ADJUSTABLE £102
SEAT TILT AND BACKREST RECLINE*		
SNB0600	<input type="checkbox"/>	ASSISTANT OPERATED GAS PISTON 30 DEG BACKREST RECLINE AND 25 DEG SEAT TILT ◆
TILT LOCK*		
SNB0630	<input type="checkbox"/>	TILT LOCK (LOCKING OF PRE-SET ANGLE) ◆
LEGREST HANGER ATTACHMENT*		
SNB0700	<input type="checkbox"/>	STANDARD LEGREST ATTACHMENT ◆
LEGREST HANGER LH*		
SNB0720	<input type="checkbox"/>	LEGREST ANGLE ADJUSTABLE, INCL CALF SUPPORT, LH ◆
LEGREST HANGER RH*		
SNB0740	<input type="checkbox"/>	LEGREST ANGLE ADJUSTABLE, INCL CALF SUPPORT, RH ◆
FOOTREST LH*		
SNB0770	<input type="checkbox"/>	FOOTREST LOWER PART, DEPTH, WIDTH AND ANGLE ADJUSTABLE WITH HEEL STRAP, LH ◆
FOOTREST RH*		
SNB0790	<input type="checkbox"/>	FOOTREST LOWER PART, DEPTH, WIDTH AND ANGLE ADJUSTABLE WITH HEEL STRAP, RH ◆
ARMREST LH*		
SNB0800	<input type="checkbox"/>	DUCKTAIL PAD, LH DEPTH ADJUSTABLE ○
SNB0810	<input type="checkbox"/>	CONCAVE ARMREST PAD, LH DEPTH ADJUSTABLE ◆
ARMREST RH*		
SNB0850	<input type="checkbox"/>	DUCKTAIL PAD, RH DEPTH ADJUSTABLE ○
SNB0860	<input type="checkbox"/>	CONCAVE ARMREST PAD, RH DEPTH ADJUSTABLE ◆

SEAT OPTION			
SNB0900	<input type="checkbox"/>	SIDE SUPPORT COVER	£51
SNB0910	<input type="checkbox"/>	POMMEL ABDUCTION CUSHION	£118
ARMREST SPACER			
SNB0920	<input type="checkbox"/>	ARMREST SPACER, INCREASES WIDTH TO A TOTAL +25 MM PER SIDE (ONLY SELECT WITH OPTION SNA1510, ALWAYS INCLUDED WITH OPTION SNA1525)	◆
BACKREST OPTIONS			
SNB1100	<input type="checkbox"/>	HEADREST ADJUSTABLE WITH HANDWHEEL	◆
SNB1110	<input type="checkbox"/>	HEADREST WITH CHEEK SUPPORT, SIDEWISE ADJUSTABLE	(NC 1130) £107
SNB1120	<input type="checkbox"/>	NECKREST FIXED	(NC 1130) £31
SNB1130	<input type="checkbox"/>	COVER FOR HEADREST	£44
SNB1140	<input type="checkbox"/>	COVER FOR NECKREST	£22
TRUNK SUPPORT LH			
SNB1200	<input type="checkbox"/>	TRUNK SUPPORT SMALL, LH	£44
SNB1210	<input type="checkbox"/>	TRUNK SUPPORT SWINGAWAY (LH)	£48
TRUNK SUPPORT RH			
SNB1220	<input type="checkbox"/>	TRUNK SUPPORT SMALL, RH	£44
SNB1230	<input type="checkbox"/>	TRUNK SUPPORT SWINGAWAY (RH)	£48
FOOTPLATE PAD			
SNB1300	<input type="checkbox"/>	FOOTPLATE PADDING	£54
CHASSIS COLOUR*			
SNB1400	<input type="checkbox"/>	PEARL BEIGE	◆
CASTOR ATTACHMENT*			
SNB1410	<input type="checkbox"/>	CASTOR ATTACHMENT QUICK RELEASE	◆
FRONT CASTORS*			
SNB1440	<input type="checkbox"/>	200X45 MM BLACK SOLID TYRE	◆
REAR WHEEL TYPE*			
SNB1510	<input type="checkbox"/>	24" (610 MM) REAR WHEELS QR ATTACHMENT	◆
SNB1525	<input type="checkbox"/>	24" (610 MM) DRUM BRAKE WHEELS QR ATTACHMENT (INCLUDES SNA0920 ARMREST SPACER)	£133
REAR WHEEL TYRE*			
SNB1540	<input type="checkbox"/>	BLACK TYRES SOLID	◆
HANDRIM*			
SNB1550	<input type="checkbox"/>	HANDRIM ALUMINIUM	◆
SPOKEGUARD			
SNB1560	<input type="checkbox"/>	FROST	£166
BRAKE*			
SNB1600	<input type="checkbox"/>	OCCUPANT OPERATED BRAKES	◆
SAFETY OPTIONS*			
SNB1650	<input type="checkbox"/>	POSTURE BELT	◆
SNB1660	<input type="checkbox"/>	ANTI-TIP DEVICE INCLUDED IN THE CHASSIS	◆
TRANSPORTATION KIT*			
SNB1680	<input type="checkbox"/>	TRANSPORTATION KIT (WITHOUT POSTURE BELT)	◆

OTHER OPTIONS

SNB1700		TABLE TRAY INCL MOUNTING KIT, WITH ELBOW PADS (NC 1800)	£138
SNB1710		TABLE TRAY INCL MOUNTING KIT (NC 1800)	£138
SNB1740		TABLE TRAY CUSHION (TO BE ORDERED IN COMBINATION WITH THE FIXED & SWINGAWAY TABLE TRAY)	£51
SNB1790		LOCKING DEVICE FOR TABLE TRAY	£15
SNB1800		MOUNTING KIT TABLE (IF TABLE ORDERED SEPARATELY)	£9

TOTAL PRICE

Non Conforming Options:

SNB1130 COVER FOR HEADREST

SNB1800 MOUNTING KIT TABLE (IF TABLE ORDERED SEPARATELY)



For a more complex customer requirement
please refer to our Invacare Unique special order form.
email specialsuk@invacare.com

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):