



Yes, you can:

Tel: 01656 776222
 Email: ordersuk@invacare.com
 Online Spares available at: www.invacare.co.uk

Manual Retail Prescription Form

FEBRUARY 2025
 Customer Ref:

| | |
|--------------------------|--------------------|
| Customer Ref* | Account No* |
| Contact No* | Order No* |
| Delivery Address* | Order Qty* |
| | Order date* |

Küschall K-Series One Solution

Standard Features

- Max User Weight 130kg (20.5 stone)
- Crash Tested ISO 7176-Part 19
- Active Rigid Wheelchair with Height Adjustable Rigid Footrest
- Multiple adjustments, Front and Rear Seat Heights and Centre of Gravity
- Full Hydroformed Aluminium Front Frame 75° Front Angle
- 24" Self Propel Q/R Wheels with Grey Aluminium Handrims and Pneumatic Tyres
- Choice of 10 Seat Widths 360mm to 500mm (14" - 20")
- Choice of 8 Seat Depths 350mm to 525mm (13.8" - 20.5")
- Range of Front Seat Heights from 450mm to 520mm
- Range of Rear Seat Heights from 390mm to 490mm
- Lockdown Folding, Angle Adjustable Fixed Height Backrest
- Without Push Handles, Carbon Clothesguards
- Black Aluminium Axle Tube with 3° Camber
- Black Comfort Front Castor
- Standard Brake Lever Bent
- Tool Kit

Standard Base price: £2,540



Please make your selection with a cross in the box Free of charge standard feature POA special request

Customer weight = KG

| Code | Description | |
|---|--|-----------------------|
| BASE PRODUCT* | | |
| DKA0067 | <input type="checkbox"/> KÜSCHALL K-SERIES | £2,540 |
| UNIQUE SPECIALS (For special modifications only) Additional lead time of 4 weeks | | |
| DKA0000 | <input type="checkbox"/> UNIQUE Special Request - For use when requesting any special options or builds | POA |
| UNIQUE CUSTOMISATION AND PERSONALISATION OPTIONS Additional 4 weeks lead time for Unique options | | |
| UKA0001 | <input type="checkbox"/> UNIQUE Personalisation - For use when selecting UKA customisation and personalisation options | <input type="radio"/> |

| FRAME* | | | |
|---|--|---|--|
| DKA1040 | | V Style Standard Frame Angle 75° Aluminium (Footrest width = seat width minus 120mm) | (NC 1475, 1490, 1495, 7240) ◆ |
| DKA1050 | | Tapered Style Standard Frame Angle 75° Aluminium (Footrest width = seat width minus 100mm) | (NC 1475, 1490, 1495, 1740, 1745, 1750, 3230, 7240) ○ |
| DKA1060 | | V Style +5cm Frame Angle 75° Aluminium. Frame to accommodate a longer seat depth and seat height for taller clients (Footrest width = seat width minus 120mm) | (NC 1410, 1420, 1435, 1740, 1745, 1750, 1755, 1760) ○ |
| DKA1070 | | V Style Standard Frame Angle 90° Aluminium (Footrest width = seat width minus 120mm) | (NC 1475, 1490, 1495, 7240) ○ |
| DKA1080 | | Tapered Style Standard Frame Angle 90° Aluminium. Taper before the frame bend creating a parallel style front frame with a wider footrest space (Footrest width = seat width minus 100mm) | (NC 1475, 1490, 1495, 1740, 1745, 1750, 3230, 5300, 5350, 5400, 5450, 5500, 7240) ○ |
| DKA1090 | | V Style +5cm Frame Angle 90° Aluminium. Frame to accommodate a longer seat depth and seat height for taller clients (Footrest width = seat width minus 120mm) | (NC 1410, 1420, 1435, 1740, 1745, 1750, 1755, 1760) ○ |
| SUSPENSION | | | |
| DKA1200 | | Seat Suspension, Weight Range 40-65 KG | (NC 1350, 1355, 1410, 1520, 1815, 1820, 1825, 1830, 1835, 6030, 6040, 6420, 6430) £434 |
| DKA1210 | | Seat Suspension, Weight Range 65-85 KG | (NC 1350, 1355, 1410, 1520, 1815, 1820, 1825, 1830, 1835, 6030, 6040, 6420, 6430) £434 |
| DKA1220 | | Seat Suspension, Weight Range 85-100 KG | (NC 1350, 1355, 1410, 1520, 1815, 1820, 1825, 1830, 1835, 6030, 6040, 6420, 6430) £434 |
| SEAT WIDTH* | | | |
| DKA1320 | | 360 mm | (NC 2364, 2366, 2368, 2370, 2372, 2710, 2748, 2750) ○ |
| DKA1325 | | 380 mm | (NC 2366, 2368, 2370, 2372, 2750) ○ |
| DKA1330 | | 400 mm | (NC 2360, 2366, 2368, 2370, 2372) ◆ |
| DKA1335 | | 420 mm | (NC 2360, 2362, 2368, 2370, 2372) ○ |
| DKA1340 | | 440 mm | (NC 2360, 2362, 2364, 2370, 2372, 2630) ○ |
| DKA1345 | | 460 mm | (NC 2360, 2362, 2364, 2372, 2630) ○ |
| DKA1350 | | 480 mm | (NC 2300, 2305, 2310, 2315, 2345, 2360, 2362, 2364, 2366, 2630, 3230) ○ |
| DKA1355 | | 500 mm | (NC 2300, 2305, 2310, 2315, 2345, 2360, 2362, 2364, 2366, 2368, 2630, 3210, 3230) ○ |
| SEAT DEPTH* | | | |
| DKA1410 | | 350 mm | (NC 1500, 1510, 2630, 2710, 2720, 2748, 2860, 7810, 7820, 7830) ○ |
| DKA1420 | | 375 mm | (NC 2710, 2720, 2748, 7810, 7820, 7830) ○ |
| DKA1435 | | 400 mm | ◆ |
| DKA1450 | | 425 mm | ○ |
| DKA1465 | | 450 mm | (NC 2630) ○ |
| DKA1475 | | 475 mm | (NC 2630, 2748) ○ |
| DKA1490 | | 500 mm | (NC 2630, 2710, 2748) ○ |
| DKA1495 | | 525 mm | (NC 2630, 2610, 2710, 2748) ○ |
| SEAT POSITION. Seat module can be moved forward and backward along the frame to achieve a more active or stable rear wheel position* | | | |
| DKA1500 | | Position 1 = 140 mm MOST ACTIVE POSITION | ○ |
| DKA1510 | | Position 2 = 125 mm | ○ |
| DKA1520 | | Position 3 = 110 mm | ◆ |
| DKA1530 | | Position 4 = 95 mm | ○ |
| DKA1540 | | Position 5 = 80 mm LEAST ACTIVE POSITION | ○ |

| BACK HEIGHT. Transportation kit is only possible with backrest heights of 405mm and above* | | | | |
|---|--|--------|---|---|
| DKA1605 | | 270 mm | (NC 2325, 2330, 2331, 2350, 2355, 2357, 7600) | ○ |
| DKA1610 | | 285 mm | (NC 2325, 2330, 2331, 2350, 2355, 2357, 7600) | ○ |
| DKA1615 | | 300 mm | (NC 2325, 2330, 2331, 2350, 2355, 2357, 7600) | ○ |
| DKA1620 | | 315 mm | (NC 7600) | ○ |
| DKA1625 | | 330 mm | (NC 7600) | ○ |
| DKA1630 | | 345 mm | (NC 7600) | ○ |
| DKA1635 | | 360 mm | (NC 7600) | ○ |
| DKA1640 | | 375 mm | (NC 7600) | ○ |
| DKA1645 | | 390 mm | (NC 7600) | ○ |
| DKA1650 | | 405 mm | | ○ |
| DKA1655 | | 420 mm | | ○ |
| FRONT SEAT HEIGHT. For 75° frames, the maximum lower leg must be front seat height minus 30mm and minus 40mm for 90° frames* | | | | |
| DKA1740 | | 450 mm | (NC 5300, 5350, 5400, 5450, 5500) | ○ |
| DKA1745 | | 460 mm | (NC 5300, 5350, 5400, 5450, 5500) | ○ |
| DKA1750 | | 470 mm | | ○ |
| DKA1755 | | 480 mm | | ○ |
| DKA1760 | | 490 mm | | ○ |
| DKA1765 | | 500 mm | | ◆ |
| DKA1770 | | 510 mm | | ○ |
| DKA1775 | | 520 mm | | ○ |
| REAR SEAT HEIGHT* | | | | |
| DKA1815 | | 390 mm | (NC 3490, 3495, 6030, 6040, 6430) | ○ |
| DKA1820 | | 400 mm | (NC 3490, 3495, 6030, 6040, 6430) | ○ |
| DKA1825 | | 410 mm | (NC 3490, 3495) | ○ |
| DKA1830 | | 420 mm | | ○ |
| DKA1835 | | 430 mm | | ○ |
| DKA1840 | | 440 mm | | ◆ |
| DKA1845 | | 450 mm | | ○ |
| DKA1850 | | 460 mm | | ○ |
| DKA1855 | | 470 mm | | ○ |
| DKA1860 | | 480 mm | | ○ |
| DKA1865 | | 490 mm | (NC 6410) | ○ |

| LOWER LEG LENGTH Seat canvas to footplate measurement, cushion height must be deducted.* | | | |
|---|--------------------------|--|---|
| DKA1918 | <input type="checkbox"/> | 290 mm | ○ |
| DKA1920 | <input type="checkbox"/> | 300 mm | ○ |
| DKA1922 | <input type="checkbox"/> | 310 mm | ○ |
| DKA1924 | <input type="checkbox"/> | 320 mm | ○ |
| DKA1926 | <input type="checkbox"/> | 330 mm | ○ |
| DKA1928 | <input type="checkbox"/> | 340 mm | ○ |
| DKA1930 | <input type="checkbox"/> | 350 mm | ○ |
| DKA1932 | <input type="checkbox"/> | 360 mm | ○ |
| DKA1934 | <input type="checkbox"/> | 370 mm | ○ |
| DKA1936 | <input type="checkbox"/> | 380 mm | ○ |
| DKA1938 | <input type="checkbox"/> | 390 mm | ○ |
| DKA1940 | <input type="checkbox"/> | 400 mm | ○ |
| DKA1942 | <input type="checkbox"/> | 410 mm | ○ |
| DKA1944 | <input type="checkbox"/> | 420 mm | ○ |
| DKA1946 | <input type="checkbox"/> | 430 mm | ◆ |
| DKA1948 | <input type="checkbox"/> | 440 mm | ○ |
| DKA1950 | <input type="checkbox"/> | 450 mm | ○ |
| DKA1952 | <input type="checkbox"/> | 460 mm | ○ |
| DKA1954 | <input type="checkbox"/> | 470 mm | ○ |
| DKA1956 | <input type="checkbox"/> | 480 mm | ○ |
| BACKREST MATERIAL* | | | |
| DKA2000 | <input type="checkbox"/> | Aluminium | ◆ |
| BACKREST TYPE* | | | |
| DKA2150 | <input type="checkbox"/> | Backrest with Lower Backbrace for use with Matrix fixed back. Lockdown folding for easy transportation | ◆ |
| BACKREST UPHOLSTERY* | | | |
| DKA2240 | <input type="checkbox"/> | Without Backrest Upholstery | ◆ |

MATRIX BACKREST: Please ensure the top of the Matrix backrest is not positioned higher than 480mm, which is the highest possible setting on the K-Series.*

| | | | | |
|---------|--|---|-----------------------|------|
| DKA2300 | | Matrix MX2, height 150 mm with quick release mini set hardware | (NC 2370, 2372, 7600) | £751 |
| DKA2305 | | Matrix MX2, height 230 mm with quick release mini set hardware | (NC 2370, 2372, 7600) | £751 |
| DKA2310 | | Matrix MX2, height 305 mm with quick release mini set hardware | (NC 2370, 2372) | £751 |
| DKA2315 | | Matrix MX2, height 405 mm with quick release mini set hardware | (NC 2370, 2372) | £751 |
| DKA2316 | | Matrix Elite Deep E2, height 160 mm with Mini-Set Mounting Hardware | (NC 7600) | £546 |
| DKA2317 | | Matrix Elite Deep E2, Height 210 mm with Mini-Set Mounting Hardware | (NC 7600) | £546 |
| DKA2318 | | Matrix Elite Deep E2, Height 260 mm with Mini-Set Mounting Hardware | (NC 7600) | £546 |
| DKA2320 | | Matrix Elite Deep E2, Height 310 mm with Mini-Set Mounting Hardware | | £546 |
| DKA2322 | | Matrix Elite Deep E2, Height 360 mm with Mini-Set Mounting Hardware | | £546 |
| DKA2325 | | Matrix Elite Deep E2, Height 410 mm with Mini-Set Mounting Hardware | | £546 |
| DKA2330 | | Matrix Elite Deep E2, Height 460 mm with Easy-Set Mounting Hardware | | £546 |
| DKA2331 | | Matrix Elite Deep E2, Height 510 mm with Easy-Set Mounting Hardware | | £546 |
| DKA2332 | | Matrix Elite E2, Height 160 mm with Mini-Set Mounting Hardware | (NC 7600) | £494 |
| DKA2334 | | Matrix Elite E2, Height 210 mm with Mini-Set Mounting Hardware | (NC 7600) | £494 |
| DKA2335 | | Matrix Elite E2, Height 260 mm with Mini-Set Mounting Hardware | (NC 7600) | £494 |
| DKA2340 | | Matrix Elite E2, Height 310 mm with Mini-Set Mounting Hardware | | £494 |
| DKA2345 | | Matrix Elite E2, Height 360 mm with Mini-Set Mounting Hardware | | £494 |
| DKA2350 | | Matrix Elite E2, Height 410 mm with Mini-Set Mounting Hardware | | £494 |
| DKA2355 | | Matrix Elite E2, Height 460 mm with Easy-Set Mounting Hardware | | £494 |
| DKA2357 | | Matrix Elite E2, Height 510 mm with Easy-Set Mounting Hardware | | £494 |

MATRIX BACKREST WIDTH CHOICE*

| | | | | |
|---------|--|--|--|---|
| DCA2360 | | 360mm Backrest width, to fit chest width 265mm | | ○ |
| DCA2362 | | 380mm Backrest width, to fit chest width 290mm | | ○ |
| DCA2364 | | 410mm Backrest width, to fit chest width 315mm | | ○ |
| DCA2366 | | 430mm Backrest width, to fit chest width 340mm | | ○ |
| DCA2368 | | 460mm Backrest width, to fit chest width 370mm | | ○ |
| DCA2370 | | 480mm Backrest width, to fit chest width 390mm | | ○ |
| DCA2372 | | 510mm Backrest width, to fit chest width 420mm | | ○ |

PUSH HANDLES*

| | | | | |
|---------|--|--|--|------|
| DKA2400 | | Without Push Handles | | ◆ |
| DKA2410 | | Push Handles Fixed Short/Mini Hand Grip | | £133 |
| DKA2420 | | Push Handles Fixed Long Hand Grip | | £133 |
| DKA2440 | | Push Handles Height Adjustable, Integrated in backrest (Height range with Matrix backs 340mm to 460mm from the seat) | | £219 |
| DKA2450 | | Push Handles Foldable | | £198 |

| BACKREST ANGLE* | | | |
|--|--|--|------|
| DKA2520 | | 78° angle from the seat | ○ |
| DKA2540 | | 82° angle from the seat | ○ |
| DKA2560 | | 86° angle from the seat | ○ |
| DKA2580 | | 90° angle from the seat | ◆ |
| DKA2590 | | 94° angle from the seat | ○ |
| SEAT UPHOLSTERY* | | | |
| DKA2600 | | Nylon Seat Upholstery with one 5cm central band (Loop) | ○ |
| DKA2601 | | Standard Seat Upholstery with 1 central loop and 2 side hook | ◆ |
| DKA2610 | | Tension Adjustable Seat Upholstery with one 5cm central band (Loop) | £45 |
| DKA2630 | | Light Seat Upholstery, mesh material to reduce weight | £69 |
| SEAT CUSHION | | | |
| DKA2700 | | Standard Seat Cushion, 60mm high Double Foam with Visco Top Layer - Comfort/Low Risk (NC 2810, 2820, 2830, 2840, 2850) | £66 |
| DKA2710 | | Flo-Tech Lite Cushion 78mm high - Comfort/Low Risk (NC 2810, 2820, 2830, 2840, 2850) | £121 |
| DKA2720 | | Matrx Libra Cushion, Black Smooth Cover 70mm High - Very High Risk | £410 |
| DKA2748 | | Contour Visco NG Cushion 100mm High - High Risk (NC 2810, 2820, 2830, 2840, 2850) | £169 |
| DKA2750 | | Vicair Active O2 lightweight cushion 60mm High - High risk (NC 2810, 2820, 2830, 2840, 2850, 2860) | £521 |
| SEAT BAG | | | |
| DKA2800 | | Seat Bag, Click in and removable | £74 |
| MATRX LIBRA CUSHION ACCESSORIES | | | |
| DKA2810 | | Fluid supplement sac (Pair) | £44 |
| DKA2820 | | Tapered Wedge Pair | £30 |
| DKA2830 | | Abductor | £23 |
| DKA2840 | | Lateral Side Wedge | £20 |
| DKA2850 | | Anterior Front Wedge | £27 |
| SEAT RIGIDISER | | | |
| DKA2860 | | Seat Rigidiser | £39 |
| FOOTREST* | | | |
| DKA3200 | | Footrest Tube, Rigid - Aluminium Hydroformed profile. Weight approx. 260g | ◆ |
| DKA3210 | | Footrest Tube, Rigid - Titanium Rounded profile tube. Weight approx. 210g (NC 6962) | £185 |
| DKA3220 | | Carbon Footplate, Angle Adjustable (NC 3290, 6962) | £225 |
| DKA3230 | | Footplate Flip-to Side, Angle Adjustable (NC 3290, 6962, 6963, 6964, 6965, 6966, 6967, 6968, 6969, 6970, 6971, 6972) | £143 |
| FOOTREST COVER | | | |
| DKA3290 | | Footrest Cover | £37 |
| CALF STRAP* | | | |
| DKA3300 | | Calf strap mounted on the frame with hook and loop | ◆ |
| DKA3305 | | Calf strap mounted with bracket support | ○ |
| CLOTHESGUARDS/MUDGUARDS* | | | |
| DKA3405 | | Composite Clothesguards, Flip up, Fold in | ○ |
| DKA3420 | | Carbon clothes guards, Flip up, Fold in. | ◆ |
| DKA3440 | | Carbon Mudguards, Removable (NC 3490, 3495) | £307 |
| DKA3450 | | Composite mudguards, removable (NC 3490, 3495) | £65 |

| ARMREST | | | |
|----------------------------|--|---|------|
| DKA3480 | | Tubular Armrest. Height Adjustable, Removable, Non-locking, Swing-away and Comfort Padded, no Clothesguard included | £241 |
| DKA3490 | | Armrest, Height Adjustable, Short Armpad (T-Arm), Non-locking, no Clothesguard included | £164 |
| DKA3495 | | Armrest, Height Adjustable, Long Armpad (T-Arm), Non-locking, no Clothesguard included | £164 |
| FRAME COLOUR* | | | |
| DKA4000 | | Space Black (RAL 9005) (NC 4066) | ○ |
| DKA4010 | | Anthracite Grey (RAL 7016) (NC 4066) | ○ |
| DKA4020 | | Pure White (RAL 9016) (NC 4066, 4220) | ○ |
| DKA4030 | | Carmine Red (RAL 3002) (NC 4066) | ○ |
| DKA4040 | | Azure Blue (RAL 5017) (NC 4066) | ○ |
| DKA4060 | | Flashy Yellow (RAL 1516) (Additional 15 day lead time) (NC 4066) | £31 |
| DKA4100 | | Black Glossy (RAL 9005) (Additional 15 day lead time) (NC 4066) | £31 |
| DKA4130 | | Silver Grey (RAL 9006) (Additional 15 day lead time) (NC 4066) | £31 |
| DKA4150 | | Ruby Red (RAL 1503) (Additional 15 day lead time) (NC 4066) | £31 |
| DKA4170 | | Dynamite Orange (Additional 15 day lead time) (NC 4066) | £31 |
| DKA4180 | | Forest Green (RAL 6029) (Additional 15 day lead time) (NC 4066) | £31 |
| UKA4101 | | Mirror Polish with standard varnish finish (Premium ceramic glaze finish available, for this option select UKA4066) | £399 |
| UKA0008 | | Your personal frame colour. Add RAL colour of your choice in special comments box, for a personalised colour frame | £299 |
| CERAMIC GLAZE FRAME FINISH | | | |
| UKA4066 | | Ceramic glaze finish on frame, for added durability | £336 |
| FRAME DECAL TYPE* | | | |
| DKA4190 | | Bold Decal (Full coloured transfer on frame) (NC 4300, 4310, 1111) | £31 |
| DKA4195 | | Tag (Small Kuschall brand transfer on front of frame) (NC 4200, 4210, 4220, 4230, 4240, 4250, 2222) | ◆ |
| DKA4197 | | Without Tag (Kuschall branding only) (NC 4200, 4210, 4220, 4230, 4240, 4250, 4300, 4310) | ○ |
| BOLD DECAL COLOUR* | | | |
| DKA4200 | | Black | ○ |
| DKA4210 | | Dark Grey | ○ |
| DKA4220 | | Light Grey | ○ |
| DKA4230 | | Red | ○ |
| DKA4240 | | Blue | ○ |
| DKA4250 | | Gold | ○ |
| DKA1111 | | No bold colour, select only with Option DKA4195 Tag or DKA4197 Without Tag (NC 4190) | ○ |
| TAG COLOUR* | | | |
| DKA4300 | | Black | ○ |
| DKA4310 | | Red | ○ |
| DKA2222 | | No tag colour, select only with Option DKA4190 Bold Decal or DKA4197 Without Tag (NC 4195) | ○ |

| FRONT CASTORS* | | | | |
|----------------------------------|--|--|---|------|
| DKA5000 | | 75mm x 25mm (3" x 1") Supersport | (NC 5502, 5504, 5506, 5507, 5508, 5509, 5510) | ○ |
| DKA5050 | | 100mm x 30mm (4" x 1.2") Everyday Black hub and tyre | (NC 5502, 5504, 5506, 5507, 5508, 5509, 5510) | ○ |
| DKA5100 | | 100mm x 30mm (4" x 1.2") Comfort Black hub and tyre | (NC 5502, 5504, 5506, 5507, 5508, 5509, 5510) | ○ |
| DKA5150 | | 100mm x 25mm (4" x 1") Starec Black hub and tyre | (NC 5502, 5504, 5506, 5507, 5508, 5509, 5510) | £83 |
| DKA5200 | | 100mm x 36mm (4" x 1.4") Aluminium Softroll Black hub and tyre (Frogs Legs) Use this option to access 4" anodised colour options | | £93 |
| DKA5250 | | 100mm x 36mm (4" x 1.4") Aluminium Softroll Silver hub, Black tyre (Frog Legs) | (NC 5502, 5504, 5506, 5507, 5508, 5509) | £93 |
| DKA5300 | | 125mm x 30mm (5" x 1.2") Everyday Black hub and tyre | (NC 5502, 5504, 5506, 5507, 5508, 5509, 5510) | ○ |
| DKA5350 | | 125mm x 30mm (5" x 1.2") Comfort Black hub and tyre | (NC 5502, 5504, 5506, 5507, 5508, 5509, 5510) | ○ |
| DKA5400 | | 125mm x 25mm (5" x 1") Starec Black hub and tyre | (NC 5502, 5504, 5506, 5507, 5508, 5509, 5510) | £83 |
| DKA5450 | | 125mm x 36mm (5" x 1.4") Aluminium Softroll Black hub and tyre (Frog Legs) Use this option to access 5" anodised colour options | | £93 |
| DKA5500 | | 125mm x 36mm (5" x 1.4") (Frog Legs) Aluminium Softroll Silver hub, Black tyre (Frogs Legs) | (NC 5502, 5504, 5506, 5507, 5508, 5509) | £93 |
| ANODISED COLOUR SOFTROLL CASTORS | | | | |
| UKA5502 | | Pink | | £50 |
| UKA5504 | | Orange | | £50 |
| UKA5506 | | Blue | | £50 |
| UKA5507 | | Purple | | £50 |
| UKA5508 | | Red | | £50 |
| UKA5509 | | Gold | | £50 |
| CERAMIC CASTOR BEARINGS | | | | |
| UKA5510 | | Ceramic castor bearings for Aluminium soft roll castors | | £58 |
| REAR WHEEL POSITION* | | | | |
| DKA6020 | | Standard, Fixed Rear Wheel Position | (NC 6100, 6110, 6120, 6130, 6140, 6230, 6240) | ◆ |
| DKA6030 | | Light Rear Wheel Extension (3 Positions) | (NC 6070, 6130, 6140, 6230, 6240) | £113 |
| DKA6040 | | Rear Wheel Extension (5 Positions) | (NC 6070) | £131 |
| CAMBER AXLE TYPE* | | | | |
| DKA6050 | | Vario-Ax, Adjustable Camber can achieve 1°, 3°, 7° and 10°. Only active/scissor brake will lock 10° camber position | (NC 6200, 6210, 6220, 6270, 7230) | £287 |
| DKA6060 | | Standard Aluminium Axle | (NC 4444) | ◆ |
| DKA6070 | | Carbon Axle | (NC 4444) | £212 |
| POSITION OF AXLE | | | | |
| DKA6100 | | 1st Position from the front. If 2nd axle is selected, 1st axle must be placed in 1st position and 2nd in the 5th position | | ○ |
| DKA6110 | | 2nd Position from the front | (NC 6230, 6240) | ○ |
| DKA6120 | | 3rd Position from the front | (NC 6230, 6240) | ○ |
| DKA6130 | | 4th Position from the front | (NC 6230, 6240) | ○ |
| DKA6140 | | 5th Position from the front | (NC 6230, 6240) | ○ |
| CAMBER ANGLE* | | | | |
| DKA6200 | | 0° Camber | (NC 6270) | ○ |
| DKA6210 | | 3° Camber | | ◆ |
| DKA6220 | | 6° Camber | (NC 6240, 6250) | ○ |
| DKA4444 | | Without camber angle, select only with DKA6050 Vario-Axle | | ○ |

| 2ND CAMBER AXLE TYPE. If a 2nd axle is selected, 1st axle must be placed in 1st position and the 2nd axle in the 5th position | | | |
|--|--|--|--------|
| DKA6230 | | Standard Aluminium Axle | £232 |
| DKA6240 | | Vario-Ax, Adjustable Camber can achieve 1°, 3°, 7° and 10° (NC 6250, 6260, 6270, 7230) | £287 |
| CAMBER ANGLE 2ND AXLE. Max of 3° difference in camber from 1st axle. If only 1 of the axles is a Vario-Axle, 2nd axle is only possible in 0° or 3° camber | | | |
| DKA6250 | | 0° Camber | ○ |
| DKA6260 | | 3° Camber | ○ |
| DKA6270 | | 6° Camber | ○ |
| DKA5555 | | Without camber angle, select only with DKA6240 Vario-Axle | ○ |
| QUAD QUICK RELEASE AXLE | | | |
| DKA6300 | | Quad Quick Release/Tetra Clip (NC 6550, 6551, 6552, 6553, 6554, 0907) | £85 |
| REAR WHEEL SIZE* | | | |
| DKA6410 | | 24" | ◆ |
| DKA6420 | | 25" (NC 6510, 6600, 6630, 6660, 6700, 6710, 6720, 6760, 6770, 6900, 6910) | ○ |
| DKA6430 | | 26" (NC 6520, 6550, 6551, 6552, 6553, 6554, 6600, 6630, 6660, 6720, 6740, 6750, 6760, 6770) | ○ |
| REAR WHEEL TYPE* | | | |
| DKA6510 | | 24", 26" Kuschall Starec Wheel, 24" Black Cross Spokes, Rim & Hub. 26" Silver Cross-Spokes, Colourless hub and Black rim (NC 6557, 6558, 6559, 6560, 6561, 6562, 6563, 6564, 6565, 6566, 6567, 6568, 6569, 6570, 6571, 6572, 6573, 6574, 6575, 6576, 6577, 6579) | ◆ |
| DKA6520 | | 24", 25" High Performance Wheel Black hub and rim (NC 6557, 6558, 6559, 6560, 6561, 6562, 6563, 6564, 6565, 6566, 6567, 6568, 6569, 6570, 6571, 6572, 6573, 6574, 6575, 6576, 6577, 6579, 6660) | £182 |
| DKA6550 | | 24", 25" Spinergy LXK, Black hub, rim and 18 spokes in Black (NC 6557, 6558, 6559, 6560, 6561, 6562, 6563, 6564, 6565, 6566, 6567, 6568, 6569, 6660, 6910) | £611 |
| UKA6551 | | 24", 25" Spinergy Lite Extreme LX, Weight 687 grams per wheel (Unique coloured rims, hub and spokes available, supplied in black as standard) (NC 6660, 6910) | £785 |
| UKA6552 | | 24", 25" Spinergy Lite Extreme X-laced XLX, Weight 821 grams per wheel (Unique coloured rims, hub and spokes available, supplied in black as standard) (NC 6660, 6910) | £844 |
| UKA6553 | | 24", 25" Spinergy X-Laced Sport Lite Extreme XSLX, Weight 757 grams per wheel (Unique coloured rims, hub and spokes available, supplied in black as standard) (NC 6660, 6910) | £893 |
| UKA6554 | | 24", 25" Spinergy Carbon Fibre CLX, Weight 530 grams per wheel (Unique coloured hub and spokes available, supplied in black as standard) (NC 6557, 6558, 6559, 6560, 6561, 6562, 6660, 6910) | £1,084 |
| SWO0907 | | 24", 25" Spinergy SPOX, Chrome hub and 18 spokes in Black. (Unique coloured spokes available) (NC 6557, 6558, 6559, 6560, 6561, 6562, 6563, 6564, 6565, 6566, 6567, 6568, 6569, 6660, 6910) | £486 |
| ANODISED COLOUR SPINERGY RIM (DESIGN LAB) | | | |
| UKA6557 | | Blue | £179 |
| UKA6558 | | Gold | £179 |
| UKA6559 | | Green | £179 |
| UKA6560 | | Orange | £179 |
| UKA6561 | | Red | £179 |
| UKA6562 | | Purple | £179 |

| ANODISED COLOUR SPINERGY HUB (DESIGN LAB) | | | |
|---|--|---|------|
| UKA6563 | | Blue | £179 |
| UKA6564 | | Gold | £179 |
| UKA6565 | | Silver | £179 |
| UKA6566 | | Green | £179 |
| UKA6567 | | Orange | £179 |
| UKA6568 | | Red | £179 |
| UKA6569 | | Purple | £179 |
| SPINERGY COLOURED SPOKES | | | |
| UKA6570 | | Blue | £84 |
| UKA6571 | | Green | £84 |
| UKA6572 | | Orange | £84 |
| UKA6573 | | Pink | £84 |
| UKA6574 | | Red | £84 |
| UKA6575 | | Purple | £84 |
| UKA6576 | | White | £84 |
| UKA6577 | | Yellow | £84 |
| CERAMIC BEARING REAR WHEEL | | | |
| UKA6579 | | Ceramic Bearing for Spinergy rear wheels | £170 |
| REAR WHEEL TYRE* | | | |
| DKA6600 | | 24" x 1" Maxxis Lite guard Pneumatic tyre. Herringbone tread, Silvershield puncture protection | ◆ |
| DKA6625 | | 24", 25", 26" x 1" Maxxis Alpha Protect. Lightweight, Pneumatic and highly puncture resistant. | £100 |
| DKA6630 | | 24" x 1" MBL Speedlite Super Lightweight (Dual layer very high puncture resistance) | £97 |
| DKA6660 | | 24" x 1 3/8" Solid Profile Black Tyre (NC 7230) | £66 |
| UNIQUE CUSTOMISATION COLOUR PACKAGE* | | | |
| UKA6960 | | UNIQUE KIT 1 - Includes custom colour forks and bearing cover. Select colour from 10 choices below | £111 |
| UKA6961 | | UNIQUE KIT 2 - Includes custom colour forks, bearing cover, backrest bracket, seat end caps and brake clamps. Select colour from the 10 choices below | £285 |
| UKA6962 | | UNIQUE KIT 3 - Includes custom colour forks, bearing cover, backrest bracket, seat end caps, brake clamps, axle clamps, calf strap bracket, rear wheel extension. Select colour from the 10 choices below | £329 |
| UKA1234 | | Without colour package | ○ |
| UNIQUE CUSTOMISATION COLOURS | | | |
| UKA6963 | | Red | ○ |
| UKA6964 | | Purple | ○ |
| UKA6965 | | Copper | ○ |
| UKA6966 | | Green | ○ |
| UKA6967 | | Orange | ○ |
| UKA6968 | | Gold | ○ |
| UKA6969 | | Champagne | ○ |
| UKA6970 | | Pink | ○ |
| UKA6971 | | Blue | ○ |
| UKA6972 | | Mirror Finish | ○ |

| HANDRIM* | | | |
|------------------------------|--|---|------|
| DKA6700 | | 24", 26" Supergripp, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance | £223 |
| DKA6710 | | 24", 26" Steel, short tab with 1 hole, mounted at 3cm distance | £62 |
| DKA6720 | | 24" Titanium, long tab with 2 holes for 3 cm or 4 cm mountings, mounted at 3cm distance | £192 |
| DKA6740 | | 24" , 25" Carbolife Gekko, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3 cm distance. Unique shape for low-effort propulsion, with fully integrated rubber strip for extra grip. 32mm diameter | £320 |
| DKA6750 | | 24", 25" Carbolife Gekko S, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion, with fully integrated rubber strip for extra grip, smaller profile for smaller hands, 24mm diameter. | £320 |
| DKA6760 | | 24" Carbolife Curve L, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion | £284 |
| DKA6770 | | 24" Carbolife Quadro, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion, with full anti-slip coating over the entire handrim | £356 |
| DKA6780 | | 24", 25", 26" Light Aluminium Grey, long tab with 2 holes for 3 cm or 4cm mountings, mounted at 3cm distance | ◆ |
| DKA6790 | | 24", 25", 26" Light Aluminium Black, short tab with 1 hole for 3cm mounting | £42 |
| SPOKEGUARD | | | |
| DKA6900 | | 24", 26" Transparent Spokeguard | £126 |
| DKA6910 | | 24", 26" Black Spokeguard | £126 |
| BRAKES* | | | |
| DKA7210 | | Standard Push to Lock Brakes, Lever Bent | ◆ |
| DKA7230 | | Performance Brakes, Push to Lock. Brake barrel flips away from tyre (NC 7300) | £78 |
| DKA7240 | | Active/Scissor Brakes (NC 7300) | £78 |
| DKA7250 | | Light/Compact Active, Aluminium Brakes (NC 7300) | £136 |
| BRAKE LEVER EXTENSION | | | |
| DKA7300 | | Brake Lever Extension | £49 |
| SAFETY OPTIONS | | | |
| DKA7500 | | Anti-Tipper Right (NC 7520) | £122 |
| DKA7510 | | Anti-Tipper Left (NC 7520) | £122 |
| DKA7520 | | Anti-Tipper Both sides (NC 8300) | £244 |
| DKA7600 | | Transportation Kit Includes Tiedown Labels | ○ |
| DKA7800 | | Standard lap belt. We recommend the fitment of a seat belt if the Transportation kit is selected | £5 |
| BODYPOINT BELTS | | | |
| DKA7810 | | Bodypoint Posture Belt Evoflex Small | £81 |
| DKA7820 | | Bodypoint Posture Belt Evoflex Medium | £81 |
| DKA7830 | | Bodypoint Posture Belt Evoflex Large | £81 |
| OTHER OPTIONS | | | |
| DKA8010 | | Front Frame Protectors; 2 Piece. Removable with hook and loop | £37 |
| DKA8500 | | Air Pump | £13 |
| DKA8680 | | 24" Mountain Bike Wheels with bag. Fitted with Specific Black Tyres, Grey Aluminium Handrims; Quick Release Axles Included | £320 |

| CANE HOLDER | | |
|------------------------------|--|-------------------------------------|
| DKA8200 | | Cane Holder, Mounted Right Side £75 |
| DKA8210 | | Cane Holder, Mounted Left Side £75 |
| TRANSIT WHEELS | | |
| DKA8300 | | Transit Wheels 3" £230 |
| BODYPOINT ACCESSORIES | | |
| DKA8800 | | Bodypoint Calf Strap £57 |
| DKA8850 | | Bodypoint Mobility Bag £55 |

| | |
|--------------------|--|
| TOTAL PRICE | |
|--------------------|--|

Non Conforming Options:

| | |
|---------|---|
| DKA1350 | 480 mm |
| DKA1355 | 500 mm |
| DKA1410 | 350 mm |
| DKA1420 | 375 mm |
| DKA1435 | 400 mm |
| DKA1475 | 475 mm |
| DKA1490 | 500 mm |
| DKA1495 | 525 mm |
| DKA1500 | Position 1 = 140 mm MOST ACTIVE POSITION |
| DKA1510 | Position 2 = 125 mm |
| DKA1520 | Position 3 = 110 mm |
| DKA1740 | 450 mm |
| DKA1745 | 460 mm |
| DKA1750 | 470 mm |
| DKA1755 | 480 mm |
| DKA1760 | 490 mm |
| DKA1815 | 390 mm |
| DKA1820 | 400 mm |
| DKA1825 | 410 mm |
| DKA1830 | 420 mm |
| DKA1835 | 430 mm |
| DKA2300 | Matrx MX2, height 150 mm with quick release mini set hardware |
| DKA2305 | Matrx MX2, height 230 mm with quick release mini set hardware |
| DKA2310 | Matrx MX2, height 305 mm with quick release mini set hardware |
| DKA2315 | Matrx MX2, height 405 mm with quick release mini set hardware |
| DKA2325 | Matrx Elite Deep E2, Height 410 mm with Mini-Set Mounting Hardware |
| DKA2330 | Matrx Elite Deep E2, Height 460 mm with Easy-Set Mounting Hardware |
| DKA2331 | Matrx Elite Deep E2, Height 510 mm with Easy-Set Mounting Hardware |
| DKA2345 | Matrx Elite E2, Height 360 mm with Mini-Set Mounting Hardware |
| DKA2350 | Matrx Elite E2, Height 410 mm with Mini-Set Mounting Hardware |
| DKA2355 | Matrx Elite E2, Height 460 mm with Easy-Set Mounting Hardware |
| DKA2357 | Matrx Elite E2, Height 510 mm with Easy-Set Mounting Hardware |
| DCA2360 | 360mm Backrest width, to fit chest width 265mm |
| DCA2362 | 380mm Backrest width, to fit chest width 290mm |
| DCA2364 | 410mm Backrest width, to fit chest width 315mm |
| DCA2366 | 430mm Backrest width, to fit chest width 340mm |
| DCA2368 | 460mm Backrest width, to fit chest width 370mm |
| DCA2370 | 480mm Backrest width, to fit chest width 390mm |
| DCA2372 | 510mm Backrest width, to fit chest width 420mm |
| DKA2610 | Tension Adjustable Seat Upholstery with one 5cm central band (Loop) |
| DKA2630 | Light Seat Upholstery, mesh material to reduce weight |
| DKA2710 | Flo-Tech Lite Cushion 78mm high - Comfort/Low Risk |
| DKA2720 | Matrx Libra Cushion, Black Smooth Cover 70mm High - Very High Risk |

| | |
|---------|---|
| DKA2748 | Contour Visco NG Cushion 100mm High - High Risk |
| DKA2750 | Vicair Active O2 lightweight cushion 60mm High - High risk |
| DKA2810 | Fluid supplement sac (Pair) |
| DKA2820 | Tapered Wedge Pair |
| DKA2830 | Abductor |
| DKA2840 | Lateral Side Wedge |
| DKA2850 | Anterior Front Wedge |
| DKA2860 | Seat Rigidiser |
| DKA3210 | Footrest Tube, Rigid - Titanium Rounded profile tube. Weight approx. 210g |
| DKA3230 | Footplate Flip-to Side, Angle Adjustable |
| DKA3290 | Footrest Cover |
| DKA3490 | Armrest, Height Adjustable, Short Armpad (T-Arm), Non-locking, no Clothesguard included |
| DKA3495 | Armrest, Height Adjustable, Long Armpad (T-Arm), Non-locking, no Clothesguard included |
| UKA4066 | Ceramic glaze finish on frame, for added durability |
| DKA4190 | Bold Decal (Full coloured transfer on frame) |
| DKA4195 | Tag (Small Kuschall brand transfer on front of frame) |
| DKA4200 | Black |
| DKA4210 | Dark Grey |
| DKA4220 | Light Grey |
| DKA4230 | Red |
| DKA4240 | Blue |
| DKA4250 | Gold |
| DKA1111 | No bold colour, select only with Option DKA4195 Tag or DKA4197 Without Tag |
| DKA4300 | Black |
| DKA4310 | Red |
| DKA2222 | No tag colour, select only with Option DKA4190 Bold Decal or DKA4197 Without Tag |
| DKA5300 | 125mm x 30mm (5" x 1.2") Everyday Black hub and tyre |
| DKA5350 | 125mm x 30mm (5" x 1.2") Comfort Black hub and tyre |
| DKA5400 | 125mm x 25mm (5" x 1") Starec Black hub and tyre |
| DKA5450 | 125mm x 36mm (5" x 1.4") Aluminium Softroll Black hub and tyre (Frog Legs) Use this option to access 5" anodised colour options |
| DKA5500 | 125mm x 36mm (5" x 1.4") (Frog Legs) Aluminium Softroll Silver hub, Black tyre (Frogs Legs) |
| UKA5502 | Pink |
| UKA5504 | Orange |
| UKA5506 | Blue |
| UKA5507 | Purple |
| UKA5508 | Red |
| UKA5509 | Gold |
| UKA5510 | Ceramic castor bearings for Aluminium soft roll castors |
| DKA6030 | Light Rear Wheel Extension (3 Positions) |
| DKA6040 | Rear Wheel Extension (5 Positions) |
| DKA6070 | Carbon Axle |
| DKA6100 | 1st Position from the front. If 2nd axle is selected, 1st axle must be placed in 1st position and 2nd in the 5th position |
| DKA6110 | 2nd Position from the front |
| DKA6120 | 3rd Position from the front |

| | |
|---------|---|
| DKA6130 | 4th Position from the front |
| DKA6140 | 5th Position from the front |
| DKA6200 | 0° Camber |
| DKA6210 | 3° Camber |
| DKA6220 | 6° Camber |
| DKA4444 | Without camber angle, select only with DKA6050 Vario-Axle |
| DKA6230 | Standard Aluminium Axle |
| DKA6240 | Vario-Ax, Adjustable Camber can achieve 1°, 3°, 7° and 10° |
| DKA6250 | 0° Camber |
| DKA6260 | 3° Camber |
| DKA6270 | 6° Camber |
| DKA6410 | 24" |
| DKA6420 | 25" |
| DKA6430 | 26" |
| DKA6510 | 24", 26" Kuschall Starec Wheel, 24" Black Cross Spokes, Rim & Hub. 26" Silver Cross-Spokes, Colourless hub and Black rim |
| DKA6520 | 24", 25" High Performance Wheel Black hub and rim |
| DKA6550 | 24", 25" Spinergy L XK, Black hub, rim and 18 spokes in Black |
| UKA6551 | 24", 25" Spinergy Lite Extreme LX, Weight 687 grams per wheel (Unique coloured rims, hub and spokes available, supplied in black as standard) |
| UKA6552 | 24", 25" Spinergy Lite Extreme X-laced XLX, Weight 821 grams per wheel (Unique coloured rims, hub and spokes available, supplied in black as standard) |
| UKA6553 | 24", 25" Spinergy X-Laced Sport Lite Extreme XSLX, Weight 757 grams per wheel (Unique coloured rims, hub and spokes available, supplied in black as standard) |
| UKA6554 | 24", 25" Spinergy Carbon Fibre CLX, Weight 530 grams per wheel (Unique coloured hub and spokes available, supplied in black as standard) |
| SWO0907 | 24", 25" Spinergy SPOX, Chrome hub and 18 spokes in Black. (Unique coloured spokes available) |
| UKA6557 | Blue |
| UKA6558 | Gold |
| UKA6559 | Green |
| UKA6560 | Orange |
| UKA6561 | Red |
| UKA6562 | Purple |
| UKA6563 | Blue |
| UKA6564 | Gold |
| UKA6565 | Silver |
| UKA6566 | Green |
| UKA6567 | Orange |
| UKA6568 | Red |
| UKA6569 | Purple |
| UKA6570 | Blue |
| UKA6571 | Green |
| UKA6572 | Orange |
| UKA6573 | Pink |
| UKA6574 | Red |
| UKA6575 | Purple |
| UKA6576 | White |

| | |
|---------|---|
| UKA6577 | Yellow |
| UKA6579 | Ceramic Bearing for Spinergy rear wheels |
| DKA6600 | 24" x 1" Maxxis Lite guard Pneumatic tyre. Herringbone tread, Silvershield puncture protection |
| DKA6630 | 24" x 1" MBL Speedlite Super Lightweight (Dual layer very high puncture resistance) |
| DKA6660 | 24" x 1 3/8" Solid Profile Black Tyre |
| UKA6962 | UNIQUE KIT 3 - Includes custom colour forks, bearing cover, backrest bracket, seat end caps, brake clamps, axle clamps, calf strap bracket, rear wheel extension. Select colour from the 10 choices below |
| UKA6963 | Red |
| UKA6964 | Purple |
| UKA6965 | Copper |
| UKA6966 | Green |
| UKA6967 | Orange |
| UKA6968 | Gold |
| UKA6969 | Champagne |
| UKA6970 | Pink |
| UKA6971 | Blue |
| UKA6972 | Mirror Finish |
| DKA6700 | 24", 26" Supergripp, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance |
| DKA6710 | 24", 26" Steel, short tab with 1 hole, mounted at 3cm distance |
| DKA6720 | 24" Titanium, long tab with 2 holes for 3 cm or 4 cm mountings, mounted at 3cm distance |
| DKA6740 | 24" , 25" Carbolife Gekko, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3 cm distance. Unique shape for low-effort propulsion, with fully integrated rubber strip for extra grip. 32mm diameter |
| DKA6750 | 24", 25" Carbolife Gekko S, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion, with fully integrated rubber strip for extra grip, smaller profile for smaller hands, 24mm diameter. |
| DKA6760 | 24" Carbolife Curve L, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion |
| DKA6770 | 24" Carbolife Quadro, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion, with full anti-slip coating over the entire handrim |
| DKA6900 | 24", 26" Transparent Spokeguard |
| DKA6910 | 24", 26" Black Spokeguard |
| DKA7230 | Performance Brakes, Push to Lock. Brake barrel flips away from tyre |
| DKA7240 | Active/Scissor Brakes |
| DKA7300 | Brake Lever Extension |
| DKA7520 | Anti-Tipper Both sides |
| DKA7600 | Transportation Kit Includes Tiedown Labels |
| DKA7810 | Bodypoint Posture Belt Evoflex Small |
| DKA7820 | Bodypoint Posture Belt Evoflex Medium |
| DKA7830 | Bodypoint Posture Belt Evoflex Large |
| DKA8300 | Transit Wheels 3" |



For a more complex customer requirement
please refer to our Invacare Unique special order form.
email specialsuk@invacare.com

Special Instructions

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):