



**Yes, you can:**

Tel: 01656 776222  
 Email: ordersuk@invacare.com  
 Online Spares available at: www.invacare.co.uk

JANUARY 2025  
 Customer Ref:

|                          |                    |
|--------------------------|--------------------|
| <b>Customer Ref*</b>     | <b>Account No*</b> |
| <b>Contact No*</b>       | <b>Order No*</b>   |
| <b>Delivery Address*</b> | <b>Order Qty*</b>  |
|                          | <b>Order date*</b> |

## Küschall The KSL One Solution

### Standard Features

- Max User Weight 100kg (16 stone)
- Active, Fully Welded, Lightweight, Rigid Wheelchair with Height Adjustable Rigid Footrest
- Fully Hydroformed Aluminium Front Frame 75° or 90° Front Angle
- 24" Self Propel Q/R Wheels with Grey Aluminium Handrims & Pneumatic Tyres
- Choice of 6 Seat Widths 340mm to 440mm (13" - 17")
- Choice of 5 Seat Depths 385mm to 485mm (15" - 19")
- Range of Front Seat Heights from 450mm to 530mm (17" - 21")
- Range of Rear Seat Heights from 400mm to 480mm (15" - 19")
- Lockdown Folding, Angle & Height Adjustable Backrest
- Black Nylon Upholstery
- Without Push Handles, Carbon Clothesguards
- Carbon Fibre Axle Tube with 3° Camber
- Black Comfort Front Castor
- Standard Brake Lever Bent
- Tool Kit

**Standard Base price: £3,150**



Please make your selection with a cross in the box     Free of charge     standard feature    POA special request

Customer weight =                      KG

| Code  | Description  |                       |
|---|--|-----------------------|
| <b>BASE PRODUCT*</b>  |  |                       |
| DSA0085   | <input type="checkbox"/> KÜSCHALL THE KSL 2.0  | £3,150                |
| <b>UNIQUE SPECIALS (For special modifications only) Additional lead time of 4 weeks*</b>                |  |                       |
| DSA0000   | <input type="checkbox"/> UNIQUE Special Request - For use when requesting any special options or builds                | POA                   |
| <b>UNIQUE CUSTOMISATION AND PERSONALISATION OPTIONS Additional 4 weeks lead time for Unique options</b> |  |                       |
| USA0001   | <input type="checkbox"/> UNIQUE Personalisation - For use when selecting USA customisation and personalisation options | <input type="radio"/> |

| <b>FRAME*</b>  |  |   |      |
|--|--|---|------|
| DSA1000  |  | Aluminium Hydroformed V Style Frame Angle 75° to the ground                                   | ◆    |
| DSA1010  |  | Aluminium Hydroformed Tapered Frame Angle 75° to the ground                                   | ○    |
| DSA1020  |  | Aluminium Hydroformed V Style Frame Angle 90° to the ground                                   | ○    |
| DSA1030  |  | Aluminium Hydroformed Tapered Frame Angle 90° to the ground (NC 5300, 5350, 5400, 5450, 5500) | ○    |
| <b>WELDED UNDERSEAT RIGID BAR*</b>   |  |   |      |
| SWO3333  |  | Under seat rigidiser bar (Welded)   | £100 |
| USA5555  |  | Without under seat rigidiser bar (Welded)   | ○    |
| <b>SEAT WIDTH*</b>   |  |   |      |
| DSA1315  |  | 340 mm (NC 2300, 2305, 2310, 2315, 2710, 2720, 2750, 2748)                                    | ○    |
| DSA1320  |  | 360 mm (NC 2362, 2364, 2366, 2368, 2748)  | ○    |
| DSA1325  |  | 380 mm (NC 2364, 2366, 2368)  | ○    |
| DSA1330  |  | 400 mm (NC 2360, 2366, 2368)  | ◆    |
| DSA1335  |  | 420 mm (NC 2360, 2362, 2368)  | ○    |
| DSA1340  |  | 440 mm (NC 2360, 2362, 2364)  | ○    |
| <b>SEAT DEPTH*</b>   |  |   |      |
| DSA1430  |  | 385 mm (NC 2720, 7810, 7820, 7830)  | ○    |
| DSA1440  |  | 410 mm  | ◆    |
| DSA1455  |  | 435 mm  | ○    |
| DSA1470  |  | 460 mm  | ○    |
| DSA1485  |  | 485 mm (NC 2710, 2748)  | ○    |
| <b>BACK HEIGHT*</b>  |  |   |      |
| DSA1605  |  | 270 mm  | ○    |
| DSA1610  |  | 285 mm  | ○    |
| DSA1615  |  | 300 mm  | ○    |
| DSA1620  |  | 315 mm  | ○    |
| DSA1625  |  | 330 mm  | ○    |
| DSA1630  |  | 345 mm  | ○    |
| DSA1635  |  | 360 mm  | ○    |
| DSA1640  |  | 375 mm  | ○    |
| DSA1645  |  | 390 mm  | ○    |
| DSA1650  |  | 405 mm  | ○    |
| DSA1655  |  | 420 mm  | ○    |
| <b>FRONT SEAT HEIGHT. Lower leg must be front height -100mm To -30mm on 75° frames and -80mm To -40mm on 90° frames*</b> |  |   |      |
| DSA1740  |  | 450 mm (NC 1840, 1845, 1850, 1855, 1860, 6420)  | ○    |
| DSA1745  |  | 460 mm (NC 1845, 1850, 1855, 1860)  | ○    |
| DSA1750  |  | 470 mm (NC 1850, 1855, 1860)  | ○    |
| DSA1755  |  | 480 mm (NC 1855, 1860)  | ○    |
| DSA1760  |  | 490 mm (NC 1860)  | ○    |
| DSA1765  |  | 500 mm  | ◆    |
| DSA1770  |  | 510 mm (NC 1820)  | ○    |
| DSA1775  |  | 520 mm (NC 1820, 1825)  | ○    |
| DSA1780  |  | 530 mm (NC 1820, 1825, 1830, 6410)  | ○    |

| <b>REAR SEAT HEIGHT. Rear height = "Front height -100 MM" TO "Front height -20 MM"*</b>         |  |   |             |
|---|--|---|-------------|
| DSA1820   |  | 400 mm  | (NC 6420) ○ |
| DSA1825   |  | 410 mm  | ○           |
| DSA1830   |  | 420 mm  | ○           |
| DSA1835   |  | 430 mm  | ◆           |
| DSA1840   |  | 440 mm  | ○           |
| DSA1845   |  | 450 mm  | ○           |
| DSA1850   |  | 460 mm  | ○           |
| DSA1855   |  | 470 mm  | ○           |
| DSA1860   |  | 480 mm  | ○           |
| <b>LOWER LEG LENGTH Seat canvas to footplate measurement, cushion height must be deducted.*</b> |  |   |             |
| DSA1930   |  | 350 mm  | ○           |
| DSA1932   |  | 360 mm  | ○           |
| DSA1934   |  | 370 mm  | ○           |
| DSA1936   |  | 380 mm  | ○           |
| DSA1938   |  | 390 mm  | ○           |
| DSA1940   |  | 400 mm  | ○           |
| DSA1942   |  | 410 mm  | ○           |
| DSA1944   |  | 420 mm  | ○           |
| DSA1946   |  | 430 mm  | ◆           |
| DSA1948   |  | 440 mm  | ○           |
| DSA1950   |  | 450 mm  | ○           |
| DSA1952   |  | 460 mm  | ○           |
| DSA1954   |  | 470 mm  | ○           |
| DSA1956   |  | 480 mm  | ○           |
| DSA1958   |  | 490 mm  | ○           |
| DSA1960   |  | 500 mm  | ○           |
| <b>BACKREST TYPE*</b>   |  |   |             |
| DSA2110   |  | Fixed Aluminium Hydroformed backrest in frame colour  | ◆           |
| SWO4444   |  | K-Series style lockdown folding backrest, supplied with carbon fibre flip up, fold in clothesguards (NC 3410, 3430, 9998) | £80         |
| SWO9999   |  | K-Series style lockdown folding backrest, supplied with carbon fibre, removable mudguards (NC 3410, 3430, 6666)           | £307        |
| <b>BACKREST UPHOLSTERY*</b>   |  |   |             |
| DSA2240   |  | Without Backrest Upholstery Select if fitting a Matrix/Fixed backrest, no upholstery required                             | ○           |
| <b>MATRIX BACKREST*</b>   |  |   |             |
| DSA2300   |  | Matrix MX2, height 150 mm with quick release mini set hardware  | £715        |
| DSA2305   |  | Matrix MX2, height 230 mm with quick release mini set hardware  | £715        |
| DSA2310   |  | Matrix MX2, height 305 mm with quick release mini set hardware  | £715        |
| DSA2315   |  | Matrix MX2, height 405 mm with quick release mini set hardware  | £715        |

| <b>MATRX BACKREST WIDTH CHOICE*</b>    |  |   |      |
|--|--|---|------|
| DSA2360                                |  | 360mm Backrest width, to fit chest width 265mm  | ○    |
| DSA2362                                |  | 380mm Backrest width, to fit chest width 290mm  | ○    |
| DSA2364                                |  | 410mm Backrest width, to fit chest width 315mm  | ○    |
| DSA2366                                |  | 430mm Backrest width, to fit chest width 340mm  | ○    |
| DSA2368                                |  | 460mm Backrest width, to fit chest width 370mm  | ○    |
| <b>PUSH HANDLES*</b>                   |  |   |      |
| DSA2400                                |  | Without Push Handles  | ◆    |
| DSA2410                                |  | Push Handles Fixed Short/Mini Hand Grip   | £133 |
| DSA2420                                |  | Push Handles Fixed Long Hand Grip   | £133 |
| DSA2440                                |  | Push Handles Height Adjustable, Integrated in backrest  | £219 |
| DSA2450                                |  | Push Handles Foldable   | £198 |
| <b>BACKREST ANGLE*</b>                 |  |   |      |
| DSA2500                                |  | 74° angle from the seat   | ○    |
| DSA2520                                |  | 78° angle from the seat   | ○    |
| DSA2540                                |  | 82° angle from the seat   | ◆    |
| DSA2560                                |  | 86° angle from the seat   | ○    |
| DSA2580                                |  | 90° angle from the seat   | ○    |
| <b>SEAT UPHOLSTERY*</b>                |  |   |      |
| DSA2600                                |  | Standard Nylon Seat Upholstery with one 5cm central band (Loop)                                 | ◆    |
| DSA2630                                |  | Light Seat Upholstery, mesh material to reduce weight   | £69  |
| <b>SEAT CUSHION</b>                    |  |   |      |
| DSA2700                                |  | Standard Seat Cushion, 60 mm Double Foam with Visco Top Layer (NC 2810, 2820, 2830, 2840, 2850) | £66  |
| DSA2710                                |  | Flo-Tech Lite Cushion - Low Risk (NC 2810, 2820, 2830, 2840, 2850)                              | £121 |
| DSA2720                                |  | Matrx Libra Cushion, Black Smooth Cover - Very High Risk  | £411 |
| DSA2748                                |  | Contour Visco NG Cushion 120mm deep - High Risk (NC 2810, 2820, 2830, 2840, 2850)               | £169 |
| DSA2750                                |  | Vicair Active O2 lightweight cushion 6cm deep (NC 2810, 2820, 2830, 2840, 2850, 2860)           | £521 |
| <b>MATRX LIBRA CUSHION ACCESSORIES</b> |  |   |      |
| DSA2810                                |  | Fluid supplement sac (Pair)   | £44  |
| DSA2820                                |  | Tapered Wedge Pair  | £31  |
| DSA2830                                |  | Abductor  | £23  |
| DSA2840                                |  | Lateral Side Wedge  | £20  |
| DSA2850                                |  | Anterior Front Wedge  | £27  |
| <b>SEAT RIGIDISER</b>                  |  |   |      |
| DSA2860                                |  | Seat Rigidiser  | £39  |
| <b>FOOTREST*</b>                       |  |   |      |
| DSA3200                                |  | Footrest Tube, Rigid - Aluminium Hydroformed profile. Weight approx. 260g (NC 6962)             | ◆    |
| DSA3210                                |  | Footrest Tube, Rigid - Titanium Rounded profile tube. Weight approx. 210g (NC 3290, 6962)       | £185 |
| DSA3220                                |  | Carbon Footplate, Angle Adjustable (NC 3290)  | £225 |
| <b>FOOTREST COVER</b>                  |  |   |      |
| DSA3290                                |  | Footrest Cover  | £37  |
| <b>CALF STRAP*</b>                     |  |   |      |
| DSA3300                                |  | Calf strap Mounted on the frame with hook and loop  | ◆    |
| DSA3305                                |  | Calf strap Mounted on the supporter with bracket  | ○    |

| <b>CLOTHESGUARDS/MUDGUARDS*</b>   |  |   |      |
|-----------------------------------|--|---|------|
| DSA3410                           |  | Carbon Clothesguards, fixed between the back rest and seat rails  | ◆    |
| DSA3430                           |  | Carbon Mudguards, fixed between the backrest and seat rails (NC 7230)   | £255 |
| SWO6666                           |  | Carbon fibre flip up clothesguards, for use with option SWO4444 K-Series style folding backrest                     | ○    |
| SWO9998                           |  | Carbon fibre removable mudguards, for use with option SWO9999 K-Series style folding backrest                       | ○    |
| <b>FRAME COLOUR*</b>              |  |   |      |
| DSA4000                           |  | Space Black (RAL 9005) (NC 4066)  | ○    |
| DSA4010                           |  | Anthracite Grey (RAL 7016) (NC 4066)  | ○    |
| DSA4020                           |  | Pure White (RAL 9016) (NC 4066, 4220)   | ○    |
| DSA4030                           |  | Carmine Red (RAL 3002) (NC 4066)  | ○    |
| DSA4040                           |  | Azure Blue (RAL 5017) (NC 4066)   | ○    |
| DSA4060                           |  | Flashy Yellow (RAL 1516) (Additional 15 day lead time) (NC 4066)  | £31  |
| DSA4100                           |  | Black Glossy (RAL 9005) (Additional 15 day lead time) (NC 4066)   | £31  |
| DSA4130                           |  | Silver Grey (RAL 9006) (Additional 15 day lead time) (NC 4066)  | £31  |
| DSA4145                           |  | Pearl Night Blue (Additional 15 day lead time) (NC 4066)  | £31  |
| DSA4150                           |  | Ruby Red (RAL 1503) (Additional 15 day lead time) (NC 4066)   | £31  |
| DSA4170                           |  | Dynamite Orange (Additional 15 day lead time) (NC 4066)   | £31  |
| DSA4180                           |  | Forest Green (RAL 6029) (Additional 15 day lead time) (NC 4066)   | £31  |
| USA4101                           |  | Mirror Polish with standard varnish finish (Premium ceramic glaze finish available, for this option select USA4066) | £399 |
| USA0008                           |  | Your personal frame colour. Add RAL colour of your choice in special comments box, for a personalised colour frame  | £299 |
| <b>CERAMIC GLAZE FRAME FINISH</b> |  |   |      |
| USA4066                           |  | Ceramic glaze finish on frame, for added durability   | £336 |
| <b>FRAME DECAL TYPE*</b>          |  |   |      |
| DSA4190                           |  | Bold Decal (Full coloured transfer on frame) (NC 4300, 4310, 1111)  | £31  |
| DSA4195                           |  | Tag (Small Kuschall brand transfer on front of frame) (NC 4200, 4210, 4220, 4230, 4240, 4250, 2222)                 | ◆    |
| DKA4197                           |  | Without Tag (Kuschall branding only) (NC 4200, 4210, 4220, 4230, 4240, 4250, 4300, 4310)                            | ○    |
| <b>BOLD DECAL COLOUR*</b>         |  |   |      |
| DSA4200                           |  | Black   | ○    |
| DSA4210                           |  | Dark Grey   | ○    |
| DSA4220                           |  | Light Grey  | ○    |
| DSA4230                           |  | Red   | ○    |
| DSA4240                           |  | Blue  | ○    |
| DSA4250                           |  | Gold  | ○    |
| DSA1111                           |  | No bold colour, select only with Option DSA4195 Tag or DSA4197 Without Tag (NC 4190)                                | ○    |
| <b>TAG COLOUR*</b>                |  |   |      |
| DSA4300                           |  | Black   | ○    |
| DSA4310                           |  | Red   | ○    |
| DSA2222                           |  | No tag colour, select only with Option DSA4190 Bold Decal or DSA4197 Without Tag (NC 4195)                          | ○    |

| FRONT CASTORS*   |  |  |   |
|--|--|--|---|
| DSA5000  |  | 75mm x 25mm (3" x 1") Supersport   | (NC 5502, 5504, 5506, 5507, 5508, 5509, 5510) ○                   |
| DSA5050  |  | 100mm x 30mm (4" x 1.2") Everyday Black hub and tyre   | (NC 5502, 5504, 5506, 5507, 5508, 5509, 5510) ○                   |
| DSA5100  |  | 100mm x 30mm (4" x 1.2") Comfort Black hub and tyre  | (NC 5502, 5504, 5506, 5507, 5508, 5509, 5510) ○                   |
| DSA5150  |  | 100mm x 25mm (4" x 1") Starec Black hub and tyre   | (NC 5502, 5504, 5506, 5507, 5508, 5509, 5510) £83                 |
| DSA5200  |  | 100mm x 36mm (4" x 1.4") Aluminium Softroll Black hub and tyre (Frogs Legs) Use this option to access 4" anodised colour options | £93   |
| DSA5250  |  | 100mm x 36mm (4" x 1.4") Aluminium Softroll Silver hub, Black tyre (Frog Legs)   | (NC 5502, 5504, 5506, 5507, 5508, 5509) £93                       |
| DSA5300  |  | 125mm x 30mm (5" x 1.2") Everyday Black hub and tyre   | (NC 5502, 5504, 5506, 5507, 5508, 5509, 5510) ○                   |
| DSA5350  |  | 125mm x 30mm (5" x 1.2") Comfort Black hub and tyre  | (NC 5502, 5504, 5506, 5507, 5508, 5509, 5510) ○                   |
| DSA5400  |  | 125mm x 25mm (5" x 1") Starec Black hub and tyre   | (NC 5502, 5504, 5506, 5507, 5508, 5509, 5510) £83                 |
| DSA5450  |  | 125mm x 36mm (5" x 1.4") Aluminium Softroll Black hub and tyre (Frog Legs) Use this option to access 5" anodised colour options  | £93   |
| DSA5500  |  | 125mm x 36mm (5" x 1.4") (Frog Legs) Aluminium Softroll Silver hub, Black tyre (Frogs Legs)                                      | (NC 5502, 5504, 5506, 5507, 5508, 5509) £93                       |
| ANODISED COLOUR SOFTROLL CASTORS   |  |  |   |
| USA5502  |  | Pink   | £50   |
| USA5504  |  | Orange   | £50   |
| USA5506  |  | Blue   | £50   |
| USA5507  |  | Purple   | £50   |
| USA5508  |  | Red  | £50   |
| USA5509  |  | Gold   | £50   |
| CERAMIC CASTOR BEARINGS  |  |  |   |
| USA5510  |  | Ceramic castor bearings  | £58   |
| ADJUSTABLE REAR WHEEL POSITION (SELECT POSITION FROM THE 5 OPTIONS BELOW)  |  |  |   |
| USA6050  |  | Adjustable centre of gravity (adjustable +/- 1 position)<br>Select from positions 1-5 below (15mm increments adjustment)         | £445  |
| SEAT POSITION WELDED (Measured from centre of backpost to centre of axle)* |  |  |   |
| DSA6100  |  | Position 1 = 140 mm MOST ACTIVE POSITION   | ○   |
| DSA6110  |  | Position 2 = 125 mm  | ○   |
| DSA6120  |  | Position 3 = 110 mm  | ○   |
| DSA6130  |  | Position 4 = 95 mm   | ◆   |
| DSA6140  |  | Position 5 = 80 mm LEAST ACTIVE POSITION   | ○   |
| CAMBER OF CARBON AXLE*   |  |  |   |
| DSA6200  |  | 0° Camber  | ○   |
| DSA6210  |  | 3° Camber  | ◆   |
| DSA6220  |  | 6° Camber  | ○   |
| QUAD QUICK RELEASE AXLE  |  |  |   |
| DSA6300  |  | Quad Quick Release/Tetra Clip  | (NC 6550, 6551, 6552, 6553, 6554, 0907) £86                       |
| REAR WHEEL SIZE*   |  |  |   |
| DSA6410  |  | 24"  | ◆   |
| DSA6420  |  | 25"  | (NC 6510, 6600, 6630, 6700, 6710, 6720, 6760, 6770, 6900, 6910) ○ |

| <b>REAR WHEEL TYPE*</b>                          |  |  |  |
|--|--|--|--|
| DSA6510  |  | 24" KÜSCHALL Starec Wheel, Black Cross Spokes, Rim & Hub.  | (NC 6557, 6558, 6559, 6560, 6561, 6562, 6563, 6564, 6565, 6566, 6567, 6568, 6569, 6570, 6571, 6572, 6573, 6574, 6575, 6576, 6577, 6579) ◆    |
| DSA6520  |  | 24", 25" High Performance Wheel Black hub and rim  | (NC 6557, 6558, 6559, 6560, 6561, 6562, 6563, 6564, 6565, 6566, 6567, 6568, 6569, 6570, 6571, 6572, 6573, 6574, 6575, 6576, 6577, 6579) £182 |
| DSA6550  |  | 24", 25" Spinergy LXX, Black hub, rim and 18 spokes in Black   | (NC 6557, 6558, 6559, 6560, 6561, 6562, 6563, 6564, 6565, 6566, 6567, 6568, 6569, 6910) £611   |
| USA6551  |  | 24", 25" Spinergy Lite Extreme LX, Weight 687 grams per wheel (Unique coloured rims, hub and spokes available)                 | (NC 6910) £785   |
| USA6552  |  | 24", 25" Spinergy Lite Extreme X-laced XLX, Weight 821 grams per wheel (Unique coloured rims, hub and spokes available)        | (NC 6910) £844   |
| USA6553  |  | 24", 25" Spinergy X-Laced Sport Lite Extreme XSLX, Weight 757 grams per wheel (Unique coloured rims, hub and spokes available) | (NC 6910) £893   |
| USA6554  |  | 24", 25" Spinergy Carbon Fibre CLX, Weight 530 grams per wheel (Unique coloured hub and spokes available)                      | (NC 6557, 6558, 6559, 6560, 6561, 6562, 6910) £1,084   |
| SWO0907  |  | 24", 25" Spinergy SPOX, Chrome hub and 18 spokes in Black. Coloured spokes are available via MyKÜSCHALL                        | (NC 6557, 6558, 6559, 6560, 6561, 6562, 6571, 6572, 6573, 6575) £486   |
| <b>ANODISED COLOUR SPINERGY RIM (DESIGN LAB)</b> |  |  |  |
| USA6557  |  | Blue   | £179   |
| USA6558  |  | Gold   | £179   |
| USA6559  |  | Green  | £179   |
| USA6560  |  | Orange   | £179   |
| USA6561  |  | Red  | £179   |
| USA6562  |  | Purple   | £179   |
| <b>ANODISED COLOUR SPINERGY HUB (DESIGN LAB)</b> |  |  |  |
| USA6563  |  | Blue   | £179   |
| USA6564  |  | Gold   | £179   |
| USA6565  |  | Silver   | £179   |
| USA6566  |  | Green  | £179   |
| USA6567  |  | Orange   | £179   |
| USA6568  |  | Red  | £179   |
| USA6569  |  | Purple   | £179   |
| <b>SPINERGY COLOURED SPOKES</b>                  |  |  |  |
| USA6570  |  | Blue   | £84  |
| USA6571  |  | Green  | £84  |
| USA6572  |  | Orange   | £84  |
| USA6573  |  | Pink   | £84  |
| USA6574  |  | Red  | £84  |
| USA6575  |  | Purple   | £84  |
| USA6576  |  | White  | £84  |
| USA6577  |  | Yellow   | £84  |
| <b>CERAMIC BEARING REAR WHEEL</b>                |  |  |  |
| USA6579  |  | Ceramic Bearing for Spinergy rear wheels   | £170   |



| <b>REAR WHEEL TYRE*</b>                    |  |   |      |
|--|--|---|------|
| DSA6600                                    |  | 24" Maxxis Lite guard Pneumatic tyre. Herringbone tread, Silvershield puncture protection   | ◆    |
| DSA6625                                    |  | 24" Maxxis Alpha Protect. Pneumatic and puncture resistant. Only for use with 24" High performance wheels   | £100 |
| DSA6630                                    |  | 24" MBL Speedlite Super Lightweight (Dual layer very high puncture resistance)  | £97  |
| <b>UNIQUE CUSTOMISATION COLOUR PACKAGE</b> |  |   |      |
| USA6960                                    |  | UNIQUE KIT 1 - Includes custom colour forks and bearing cover. Select colour from 10 choices below  | £111 |
| USA6961                                    |  | UNIQUE KIT 2 - Includes custom colour forks, bearing cover, backrest bracket with K-Series backrest, brake clamps. Select colour from the 10 choices below  | £285 |
| USA6962                                    |  | UNIQUE KIT 3 - Includes custom colour forks, bearing cover, backrest bracket with K-Series backrest, brake clamps, rear wheel axle clamps, calf strap bracket. Select colour from the 10 choices below                          | £329 |
| <b>UNIQUE CUSTOMISATION COLOURS</b>        |  |   |      |
| USA6963                                    |  | Red   | ○    |
| USA6964                                    |  | Purple  | ○    |
| USA6965                                    |  | Copper  | ○    |
| USA6966                                    |  | Green   | ○    |
| USA6967                                    |  | Orange  | ○    |
| USA6968                                    |  | Gold  | ○    |
| USA6969                                    |  | Champagne   | ○    |
| USA6970                                    |  | Pink  | ○    |
| USA6971                                    |  | Blue  | ○    |
| USA6972                                    |  | Mirror Finish   | ○    |
| <b>HANDRIM*</b>                            |  |   |      |
| DSA6700                                    |  | 24" Supergripp, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance   | £223 |
| DSA6710                                    |  | 24" Steel, mounted at 3 cm distance   | £62  |
| DSA6720                                    |  | 24" Titanium, long tab with 2 holes for 3 cm or 4 cm mountings, mounted at 3cm distance   | £192 |
| DSA6740                                    |  | 24" , 25" Carbolife Gekko, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3 cm distance. Unique shape for low-effort propulsion, with fully integrated rubber strip for extra grip                                    | £320 |
| DSA6750                                    |  | 24", 25" Carbolife Gekko S, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion, with fully integrated rubber strip for extra grip, smaller profile for smaller hands | £320 |
| DSA6760                                    |  | 24" Carbolife Curve L, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion  | £284 |
| DSA6770                                    |  | 24" Carbolife Quadro, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion, with full anti-slip coating over the entire handrim  | £356 |
| DSA6780                                    |  | 24", 25" Light Aluminium Grey, long tab with 2 holes for 3 cm or 4cm mountings, mounted at 3cm distance   | ◆    |
| DSA6790                                    |  | 25", 26" Light Aluminium Black, short tab with 1 hole for 3cm mounting  | £42  |
| <b>SPOKEGUARD</b>                          |  |   |      |
| DSA6900                                    |  | 24" Transparent Spokeguard  | £126 |
| DSA6910                                    |  | 24" Black Spokeguard  | £126 |



| <b>BRAKES*</b>               |  |  |      |
|------------------------------|--|--|------|
| DSA7210                      |  | Standard Push to Lock Brakes, Lever Bent   | ◆    |
| DSA7230                      |  | Performance Brakes, Push to Lock. Brake barrel flips away from tyre  | £78  |
| DSA7240                      |  | Active/Scissor Brakes  | £78  |
| DSA7250                      |  | Light/Compact Active, Aluminium Brakes   | £136 |
| <b>SAFETY OPTIONS</b>        |  |  |      |
| DSA7500                      |  | Anti-Tipper Right  | £122 |
| DSA7510                      |  | Anti-Tipper Left   | £122 |
| DSA7520                      |  | Anti-Tipper Both sides   | £244 |
| DSA7810                      |  | Bodypoint Posture Belt Evoflex 2 points, Small, Flat-Mount with Push Button. User Weight: 25-50KG                          | £81  |
| DSA7820                      |  | Bodypoint Posture Belt Evoflex 2 points, Medium, Flat-Mount with Rehab Latch. User Weight: 50-75KG                         | £81  |
| DSA7830                      |  | Bodypoint Posture Belt Evoflex, 2 points, Large, Flat-Mount with Rehab Latch. User Weight: 75-100KG                        | £81  |
| <b>OTHER OPTIONS</b>         |  |  |      |
| DSA8010                      |  | Front Frame Protectors; 2 Piece. Removable with hook and loop  | £37  |
| DSA8500                      |  | Air Pump   | £13  |
| DSA8550                      |  | Tool Kit   | ◆    |
| DSA8680                      |  | 24" Mountain Bike Wheels with bag. Fitted with Specific Black Tyres, Grey Aluminium Handrims; Quick Release Axles Included | £320 |
| <b>BODYPOINT ACCESSORIES</b> |  |  |      |
| DSA8800                      |  | Bodypoint Calf Strap   | £57  |
| DSA8850                      |  | Bodypoint Mobility Bag   | £55  |

|                    |  |
|--------------------|--|
| <b>TOTAL PRICE</b> |  |
|--------------------|--|

**Non Conforming Options:**

|         |   |
|---------|---|
| DSA1820 | 400 mm  |
| DSA1825 | 410 mm  |
| DSA1830 | 420 mm  |
| DSA1840 | 440 mm  |
| DSA1845 | 450 mm  |
| DSA1850 | 460 mm  |
| DSA1855 | 470 mm  |
| DSA1860 | 480 mm  |
| DSA2300 | Matrx MX2, height 150 mm with quick release mini set hardware                                   |
| DSA2305 | Matrx MX2, height 230 mm with quick release mini set hardware                                   |
| DSA2310 | Matrx MX2, height 305 mm with quick release mini set hardware                                   |
| DSA2315 | Matrx MX2, height 405 mm with quick release mini set hardware                                   |
| DSA2360 | 360mm Backrest width, to fit chest width 265mm  |
| DSA2362 | 380mm Backrest width, to fit chest width 290mm  |
| DSA2364 | 410mm Backrest width, to fit chest width 315mm  |
| DSA2366 | 430mm Backrest width, to fit chest width 340mm  |
| DSA2368 | 460mm Backrest width, to fit chest width 370mm  |
| DSA2710 | Flo-Tech Lite Cushion - Low Risk  |
| DSA2720 | Matrx Libra Cushion, Black Smooth Cover - Very High Risk  |
| DSA2748 | Contour Visco NG Cushion 120mm deep - High Risk   |
| DSA2750 | Vicair Active O2 lightweight cushion 6cm deep   |
| DSA2810 | Fluid supplement sac (Pair)   |
| DSA2820 | Tapered Wedge Pair  |
| DSA2830 | Abductor  |
| DSA2840 | Lateral Side Wedge  |
| DSA2850 | Anterior Front Wedge  |
| DSA2860 | Seat Rigidiser  |
| DSA3290 | Footrest Cover  |
| DSA3410 | Carbon Clothesguards, fixed between the back rest and seat rails                                |
| DSA3430 | Carbon Mudguards, fixed between the backrest and seat rails                                     |
| SWO6666 | Carbon fibre flip up clothesguards, for use with option SWO4444 K-Series style folding backrest |
| SWO9998 | Carbon fibre removable mudguards, for use with option SWO9999 K-Series style folding backrest   |
| USA4066 | Ceramic glaze finish on frame, for added durability   |
| DSA4190 | Bold Decal (Full coloured transfer on frame)  |
| DSA4195 | Tag (Small Kuschall brand transfer on front of frame)   |
| DSA4200 | Black   |
| DSA4210 | Dark Grey   |
| DSA4220 | Light Grey  |
| DSA4230 | Red   |
| DSA4240 | Blue  |
| DSA4250 | Gold  |
| DSA1111 | No bold colour, select only with Option DSA4195 Tag or DSA4197 Without Tag                      |
| DSA4300 | Black   |

|         |   |
|---------|---|
| DSA4310 | Red   |
| DSA2222 | No tag colour, select only with Option DSA4190 Bold Decal or DSA4197 Without Tag  |
| DSA5300 | 125mm x 30mm (5" x 1.2") Everyday Black hub and tyre  |
| DSA5350 | 125mm x 30mm (5" x 1.2") Comfort Black hub and tyre   |
| DSA5400 | 125mm x 25mm (5" x 1") Starec Black hub and tyre  |
| DSA5450 | 125mm x 36mm (5" x 1.4") Aluminium Softroll Black hub and tyre (Frog Legs) Use this option to access 5" anodised colour options |
| DSA5500 | 125mm x 36mm (5" x 1.4") (Frog Legs) Aluminium Softroll Silver hub, Black tyre (Frogs Legs)                                     |
| USA5502 | Pink  |
| USA5504 | Orange  |
| USA5506 | Blue  |
| USA5507 | Purple  |
| USA5508 | Red   |
| USA5509 | Gold  |
| USA5510 | Ceramic castor bearings   |
| DSA6410 | 24"   |
| DSA6420 | 25"   |
| DSA6510 | 24" KÜSCHALL Starec Wheel, Black Cross Spokes, Rim & Hub.   |
| DSA6550 | 24", 25" Spinergy LXX, Black hub, rim and 18 spokes in Black  |
| USA6551 | 24", 25" Spinergy Lite Extreme LX, Weight 687 grams per wheel (Unique coloured rims, hub and spokes available)                  |
| USA6552 | 24", 25" Spinergy Lite Extreme X-laced XLX, Weight 821 grams per wheel (Unique coloured rims, hub and spokes available)         |
| USA6553 | 24", 25" Spinergy X-Laced Sport Lite Extreme XSLX, Weight 757 grams per wheel (Unique coloured rims, hub and spokes available)  |
| USA6554 | 24", 25" Spinergy Carbon Fibre CLX, Weight 530 grams per wheel (Unique coloured hub and spokes available)                       |
| SWO0907 | 24", 25" Spinergy SPOX, Chrome hub and 18 spokes in Black. Coloured spokes are available via MyKÜSCHALL                         |
| USA6557 | Blue  |
| USA6558 | Gold  |
| USA6559 | Green   |
| USA6560 | Orange  |
| USA6561 | Red   |
| USA6562 | Purple  |
| USA6563 | Blue  |
| USA6564 | Gold  |
| USA6565 | Silver  |
| USA6566 | Green   |
| USA6567 | Orange  |
| USA6568 | Red   |
| USA6569 | Purple  |
| USA6570 | Blue  |
| USA6571 | Green   |
| USA6572 | Orange  |
| USA6573 | Pink  |
| USA6574 | Red   |
| USA6575 | Purple  |
| USA6576 | White   |

|         |  |
|---------|--|
| USA6577 | Yellow   |
| USA6579 | Ceramic Bearing for Spinergy rear wheels   |
| DSA6600 | 24" Maxxis Lite guard Pneumatic tyre. Herringbone tread, Silvershield puncture protection  |
| DSA6630 | 24" MBL Speedlite Super Lightweight (Dual layer very high puncture resistance)   |
| USA6962 | UNIQUE KIT 3 - Includes custom colour forks, bearing cover, backrest bracket with K-Series backrest, brake clamps, rear wheel axle clamps, calf strap bracket. Select colour from the 10 choices below |
| DSA6700 | 24" Supergripp, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance  |
| DSA6710 | 24" Steel, mounted at 3 cm distance  |
| DSA6720 | 24" Titanium, long tab with 2 holes for 3 cm or 4 cm mountings, mounted at 3cm distance  |
| DSA6760 | 24" Carbolife Curve L, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion   |
| DSA6770 | 24" Carbolife Quadro, long tab with 2 holes for 3cm & 4cm mounting, mounted at 3cm distance. Unique shape for low-effort propulsion, with full anti-slip coating over the entire handrim               |
| DSA6900 | 24" Transparent Spokeguard   |
| DSA6910 | 24" Black Spokeguard   |
| DSA7230 | Performance Brakes, Push to Lock. Brake barrel flips away from tyre  |
| DSA7810 | Bodypoint Posture Belt Evoflex 2 points, Small, Flat-Mount with Push Button. User Weight: 25-50KG  |
| DSA7820 | Bodypoint Posture Belt Evoflex 2 points, Medium, Flat-Mount with Rehab Latch. User Weight: 50-75KG   |
| DSA7830 | Bodypoint Posture Belt Evoflex, 2 points, Large, Flat-Mount with Rehab Latch. User Weight: 75-100KG  |



For a more complex customer requirement  
please refer to our Invacare Unique special order form.  
email [specialsuk@invacare.com](mailto:specialsuk@invacare.com)

**Special Instructions**

Please supply the following information:

Name:

Email:

Phone:

Description of special request (please include relevant option code or part number):